Gebr. Brasseler GmbH & Co. KG  
 Trophagener Weg 25  
 32657 Lemgo - Germany  
 Telefon +49 (0) 5261 701-400  
 Telefax +49 (0) 5261 701-580  
 info@kometmedical.de  
 www.kometmedical.de

(doesn't apply to  
 Kirschner/Steinmann US.)

**Guide Pins and Screws | Kirschner/Steinmann US.**  
 © Gebr. Brasseler GmbH & Co. KG  
**309869 - Rev 2018-11**

Do not use in case of damaged packaging  
 No utilizar si el embalaje está dañado

## Rx Only

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.  
 Precaución: La ley federal de EE.UU. limita la venta de este dispositivo a médicos o bajo prescripción facultativa.

Do not re-use! (Single use only)  
 No reutilizar (un solo uso)

Instructions for use and safety recommendations for the single use of pins, wires and screws in the medical sector

The below instructions and recommendations apply to:

- Kirschner wires
- Steinmann pins
- Guide pins for sawing templates
- Guide screws for sawing templates

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

The pins and screws are made of high-grade materials. It is essential that the below described instructions for use and safety recommendations be observed. Improper use can lead to damage to the tissue, premature wear, destruction of the pins, wires and screws and injury to the operator, the patient or third persons.

### Note:

The Kirschner wires, Steinmann pins, guide pins and guide screws are made of stainless steel conforming to ASTM F138 or ASTM F899.

The instruments are supplied in sterile or non sterile packaging, as required. This piece of information is clearly indicated on the label of the instruments.

### 1. Safety Notice

Please read the instructions for use prior to using the pins, wires and screws for orthopedic surgery. The pins, wires and screws may only be used for their intended purpose according to the general rules for orthopaedic surgical treatment, occupational safety and accident prevention. The pins, wires and screws should only be used in accordance with the instructions for use. The pins, wires and screws should only be used by surgeons with training and experience in orthopaedic surgery as well as in diagnosis and pre-operative planning. Store out of children's reach. Intended for surgical use only. The product has been developed solely for use in surgery. Descriptions and data constitute no warranty of attributes and are not binding. In case of contributory negligence by the user, Gebr. Brasseler partly or totally declines liability for all resulting damages, particularly due to the non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

To ensure traceability of the MD by means of a UDF, please make sure to keep the packaging of articles without direct marking.

### 2. Product description

The pins, wires and screws are designed for fixing bones and for temporarily aligning and securing cutting templates and guides on bones.

### 3. Instructions for use

Pins, wires and screws may only be used according to the directions described herein and the instructions for use enclosed with surgical power systems.

In case of use, make sure that the instruments are sterile and have not been damaged. Do not reuse pins, wires or screws.

Insert pins, wires and screws into the hand-piece until the stop position is reached.

Prior to start-up, check that the instrument is securely fitted in its proper position.

Pins, wires and screws may not be used for work on metal or used for leverage.

When using Kirschner wires in cannulated instruments, make sure that the diameter of the wire does not exceed that of the cannulation. Ideally, the diameter of the wire should be slightly smaller than the cannulation.

When using saw templates, make sure that the diameter of the guide pins and screws does not exceed that of the intended holes. Ideally, the diameter of the guide pins and guide screws should be slightly smaller than that of the hole in the template.

The physical properties of the wires, pins and screws are impaired after one use. Therefore the Kirschner wires, Steinmann pins, and guide pins and screws are single use only. Bent, blunt or otherwise damaged pins, wires and screws must be discarded immediately.

### 4. Indications

The Kirschner Wires and Steinmann Pins are indicated for use in fixation of bone fractures, for bone reconstructions, as Guide pins and Guide Screws for insertion of other implants or implantation through the skin so that traction may be applied to the skeletal system.

### 5. Contraindications

All known contraindications for surgery must be generally taken into account. These include, but are not limited to:

- Diseases that inhibit the healing process
- Insufficient bone quality
- Allergies to the implant material. In case of a suspected metal allergy, tests have to be carried out before the start of the treatment.
- Active infection

### 6. Cleaning, disinfection, sterilization and storage

All pins, wires and screws have to be discarded after one use.

Components supplied in a sterile condition:

Do not clean, disinfect or sterilize instruments that will be supplied in a sterile condition. Sterilized instruments and instruments in sterile packages must be stored away from dust, moisture and any sources of contamination. Components supplied in a non-sterile condition:

- The instruments have to be cleaned, disinfected and sterilized prior to use. These devices should be stored in their original packaging at room temperature, in a dry, dust-free environment.

### Instructions for use

Please read carefully and retain for future reference.

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- The instruments are suitable for cleaning in an ultrasonic bath. To avoid damage to the instruments, take care that they do not come in contact with one another.
- They can be cleaned in the gentle cycle (e.g. Vario TD) of a thermo disinfectant. Only use detergents and disinfectants with proven material compatibility. We recommend Neodisher MediClean forte (co. Dr. Weigert).
- For recommendations for use of the cleaning agents please refer to the instructions provided by the manufacturer of the agents used.
- All contact with H<sub>2</sub>O<sub>2</sub> (Hydrogen peroxide) as well as detergents and disinfectants containing chlorine or oxalic acid must be avoided, as these may cause pitting and corrosion.
- Use suitable anti-corrosive disinfecting and cleaning agents at all times.
- Rinse off any disinfecting and cleaning agent immediately after use under running water, as residues of these agents may set off a chemical reaction which could cause corrosion.
- For further information on the use, immersion time and suitability of the disinfecting and cleaning agents please refer to the instructions provided by the manufacturers of these agents.
- During disinfection and sterilization, make sure that the chosen methods are suitable for the instrument to be treated. Returned instruments are only accepted if they are in clean condition and have been disinfected and sterilised.
- The following cycle will achieve a sterility assurance level (SAL) of 10<sup>-6</sup> in accordance with ISO 11737-2

Sterilization of medical devices – Microbiological methods – Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process:

Method:	Steam
Cycle:	Gravity Displacement
Temperature:	132 °C
Exposure Time:	10 minutes (Full cycle)
Drying time:	0-1 minute
Alternatively Method:	Steam
Cycle:	Pre-Vacuum
Temperature:	132 °C
Exposure Time:	4 minutes (Full cycle)

### 7. Warnings

The instruments are supplied in sterile or non-sterile packaging, as required. The sterility information is clearly indicated on the device labels.

K-wires, Steinmann pins, Guide Pins and Screws supplied in non-sterile condition must be cleaned, disinfected and sterilized prior to use on the patient.

Improper use of the pins, wires and screws can lead to damage to the tissue, destruction of the device components and injury to the operator, patient or third parties.

Avoid excessive contact pressure. This can result in damage of the device. At the same time, extreme heat is generated which increases the risk of thermal necrosis. Excessive contact pressure can also result in breakage. Avoid jamming of the device during operation. Do not use as a lever.

Do not use the product if the sterile package is damaged.

Single use devices must be disposed of in compliance with all applicable local, state, and federal laws and regulations concerning medical waste.

### 8. Precautions

Use aseptic technique to open package for delivery of device into sterile field.

To avoid injury and infection, wear gloves or finger protection when inserting or removing the instruments from the hand piece. The tips of the instruments are sharp.

The use of a device in any manner or medical procedure other than those for which it is designed and indicated may result in damage or breakage.

Metallic implants can loosen, fracture, corrode, migrate, or cause pain and stress shielding of bone, even after a fracture has healed, particularly in young, active patients. The fixation devices must be removed once their service as an aid to healing is accomplished. Implant removal should be followed by adequate postoperative management to avoid re-fracture.

The implant system has not been evaluated for safety and compatibility in the MR environment, nor has it been tested for heating or migration in the MR environment.

### 9. Special patient population

Patients must, as a matter of principle, be suitable for surgical bone treatment (see contraindications).

### 10. Adverse reactions

All known reactions on orthopedic treatment have to be considered carefully. Adverse events may include but are not limited to:

- Bone resorption and poor bone formation
- Device failure
- Infection
- Tissue damage as a result of surgical trauma

### 11. Disclaimer of Warranties

ALL WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE EXPRESSLY DISCLAIMED. ALL INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING, BUT NOT LIMITED TO, LOST PROFITS, INJURY TO OTHER PROPERTY, LOSS OF USE OR OTHER COMMERCIAL LOSSES, ARE EXPRESSLY DISCLAIMED. If, notwithstanding the foregoing, consequential and incidental damages cannot be excluded due to operation of law, such damages are expressly limited in amount to the purchase price of any defective products.

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