| | |
|-----------|--|
| fr | Instruments non stériles - Instruments chirurgicaux à usage orthopédique. |
| en | Non-sterile instruments - Surgical instruments for orthopaedic use. |
| de | Nichtsterile Instrumente - Chirurgische Instrumente zum orthopädischen Gebrauch. |
| nl | Niet-steriele instrumenten - Chirurgische instrumenten voor orthopedisch gebruik. |
| es | Instrumentos no estériles - Instrumentos quirúrgicos para uso ortopédico. |
| it | Strumenti non sterili - Strumenti chirurgici per uso ortopedico. |
| pt | Instrumentos não-estéreis - Instrumentos cirúrgicos para utilização ortopédica. |
| sv | Ikkesterile instrument - Kirurgiska instrument för ortopedi. |
| da | Ikke-sterile instrument - Kirurgiske instrument til ortopædisk brug. |
| fi | Epiteerilliset instrumentit - Kirurgiset instrumentit ortopediseen käyttöön. |
| no | Usterile instrumenter - Kirurgiske instrumenter for ortopedisk bruk. |
| pl | Instrumenty niesterylizowane - Instrumenty chirurgiczne do użytku ortopedycznego |
| el | Μη αποστειρωμένα εργαλεία - Χειρουργικά εργαλεία για ορθοπαιδική χρήση. |

| |
|--|
|  |
| novastep ® S.A.S |
| 2, Allée Jacques Frimot |
| 35000 Rennes - France |
| Tel. : +33 (0)2-99 33 86 50 |
| Fax : +33 (0)9 70 29 18 95 |
| contact@novastep-ortho.com |

(CE): 2014-06

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en Novastep instruments are designed to be used in surgical procedures in association with Novastep devices. For more specific instructions, please refer to the instructions for use dedicated to each range of Novastep implants.

PRESENTATION

Non-sterile instrument is delivered in plastic bag or directly in a tray adapted for steam sterilisation.

COMPOSITION

Stainless steel | Titanium alloy | CoCr L605 | Propylux®

TARGET GROUP(S)

Target population

Adult patients

Intended use

For surgery; orthopaedic surgeons

For cleaning and sterilization: Sterilization officer

LIFETIME

Single use device: 60 minutes

INDICATIONS FOR USE

Surgical instruments are indicated:

- to prepare the bone site for the insertion of the implant.
- to hold the bone fragments and instruments and to guide the instruments.
- to define the adapted implant reference.
- to insert the implant in the bone site.

CONTRA-INDICATIONS

- CrCo wire contains a CMR, do not use on pregnant or breastfeeding women.
- Hypersensitivity to one or more components.

DIRECTIONS FOR USE

WARNING: The instrument is supplied non-sterile and must be sterilised before first use and between uses (for reusable device) to avoid cross-contamination. Read and apply the paragraph "Sterilization and reprocessing instructions".

• Before each use, all the instruments of the set should be removed then visually checked to ensure that they are in perfect technical condition. It is possible to refer to **IDM** available on the Website (<https://int.novastep.life/ifu/>) for examples of verification:

- Check **drilling and cutting instruments** for the following indications to avoid any risk of injury to the patient and the user:
 - o blunt (rounded) blades and notches
 - o damage to the shank,
 - o bent or "out of true" instruments,
- Check that the instruments are working properly, including the moving parts.
- Check the legibility of the various markings on the instrument to avoid loss of device traceability (absence of corrosion, etc.).
- Check that the instrument cavities are clean before, during and after the procedure to ensure that no bone debris accumulates and so reduce any risk of obstruction.

- If an instrument is dull, damaged, corroded or its markings are no longer legible, do not use it and return it to Novastep for replacement.
- At the beginning of the procedure check the interface between instruments designed to be connected together.
- The blunt or damaged devices should be handled carefully and should be disposed of according to hospital protocol or returned to Novastep.
- Single use instruments are in direct contact with biological tissue and should be disposed of after use in a biohazardous waste bin according to hospital protocol.
- If used in combination with guide, the cutting device must be inserted into the guide before it is used.
- During the use of rotating instruments, ensure sufficient cooling, otherwise this leads to clogging of the blades which results in increased heat generation and may cause irreversible damage to the bone (thermal necrosis) and the life of the instruments will be reduced.
- During use, take care to prevent the instruments from becoming bent, jammed or levered.

CAUTION: In the case of breakage during use, it is at the practitioner's discretion to take appropriate measures depending on the clinical situation.

CAUTION: Do not use the surgical instruments supplied by Novastep with other instruments or implants from other manufacturers, unless otherwise specified. Please refer to the operating techniques for each device implant on the following website: <https://int.novastep.life/products/>

WARNINGS AND PRECAUTIONS FOR USE

- **Single use device only**: Do not reuse device on several patients to avoid cross contamination.
- Not authorized for work on metal. Do not allow the instrument to come in conflict during use, as this might cause nicks in the blades.
- Do not exert excessive contact pressure and avoid prying on the instrument during operation as this may lead to damage to the functional part of the instruments as well as to blade breakage.

SIDE EFFECTS

No known side effects.

Please report any serious incident that occurred in relation to the device to Novastep (regulatory@novastep-ortho.com) and to the local competent authority for medical devices.

DISPOSAL

- **CAUTION**: The disposal of sharp or cutting instruments must be carried out with the utmost care to prevent all risks of physical harm to users.
- The blunt or damaged devices should be handled carefully and should be disposed of according to hospital protocol or returned to Novastep.
- Single use instruments are in direct contact with biological tissue and should be disposed of after use in a biohazardous waste bin according to hospital protocol.

STERILIZATION AND REPROCESSING INSTRUCTIONS

The information provided relates exclusively to reusable instruments or unused and unopened single use instruments.

Warnings

- Sterilise the instrument before the first use and between each use to avoid cross contamination.

Limitations: processing

- Reusable device: Refer to the **IDM** available on the Website.

- Any single-use instruments with scratches, nicks, residues or debris should be discarded.

Initial treatment at the point of use

- Before first cleaning, remove and discard all original packaging: sealed pouch.
- For assembled instruments, disassemble before cleaning.
- Before steam sterilisation, place the instrument in the Novastep sterilisation tray.

CLEANING PROCEDURE

In accordance with AAMI TIR12, AAMI TIR30, ANSI/AAMI S181 and ISO 17664 the following cleaning procedures have been validated:

Manual cleaning

NOTICE: Due to its significantly lower efficiency and reproducibility, it is preferable to use automatic cleaning rather than manual cleaning

Equipment required : Enzymatic solution (neutral pH)

For a cannulated device: using the cleaning pin (provided with the device), remove gross soil from the device lumen by inserting the wire and moving it back and forth a minimum of 3 times.

- Rinse soiled device under running, cold tap water (target temperature between 22°C and 30°C) for a minimum of 1min. During rinsing, remove gross soil from the outer surfaces using a soft-bristled brush or clean, soft, lint-free cloth.
- For a cannulated device**: Using a syringe, pipette, or water jet, flush the hole with a minimum of 10ml of water.
- 2) Prepare an enzymatic solution (neutral pH) following the detergent manufacturer's instructions (target temperature between 33°C and 36°C).
- 3) Immerse the device completely in the enzymatic solution for at least 5min.
- 4) Manually clean the device for at least 1min in the freshly prepared enzymatic solution. Use a soft-bristled brush to remove soil and debris from the outer surfaces of the device.

For a cannulated device: Using a syringe, pipette or water jet, flush the lumen with a minimum of 10ml of solution.

- Prepare an enzymatic solution (neutral pH) following the detergent manufacturer's instructions (target temperature between 33°C and 36°C).
- Immerse the device in an ultrasonic bath filled with fresh enzyme solution and vibrate it for at least 10min.
- 7) Rinse device under cold running tap water for a minimum of 30sec (target temperature between 22°C and 30°C).

For a cannulated device: Using a syringe, pipette or water jet, flush the lumen with a minimum of 10ml of water.

Automated cleaning

Equipment required : Enzymatic solution (neutral pH)

Perform a manual pre-wash before automated cleaning procedure as follows :

For a cannulated device: using the cleaning pin (provided with the device), remove gross soil from the device lumen by inserting the wire and moving it back and forth a minimum of 3 times.

- Rinse soiled device under running, cold tap water (target temperature between 22°C and 30°C) for a minimum of 1min. During rinsing, remove gross soil from the outer surfaces using a soft-bristled brush or clean, soft, lint-free cloth.
- For a cannulated device**: Using a syringe, pipette, or water jet, flush the lumen with a minimum of 10ml of water.
- 2) Prepare an enzymatic solution (neutral pH) following the detergent manufacturer's instructions (target temperature between 33°C and 36°C).
- 3) Immerse the device completely in the enzymatic solution for at least 5min.
- 4) Manually clean the device for at least 1min in the freshly prepared enzymatic solution. Use a soft-bristled brush to remove soil and debris from the outer surfaces of the device.

For a cannulated device: Using a syringe, pipette or water jet, flush the lumen with a minimum of 10ml of solution.

- Prepare an enzymatic solution (neutral pH) following the detergent manufacturer's instructions (target temperature between 33°C and 36°C).
- Immerse the device in an ultrasonic bath filled with fresh enzyme solution and vibrate it for at least 10min.
- 7) Rinse device under cold running tap water for a minimum of 30sec (target temperature between 22°C and 30°C).

For a cannulated device: Using a syringe, pipette, or water jet to flush the lumen. Remove residual water from the lumen using filtered compressed air or a syringe.

Perform an automated cleaning procedure as follow

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