

Pro Night Splint

PATIENT APPLICATION



BEFORE USING THIS DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USE/INDICATIONS:

Helps reduce pain associated with Plantar Fasciitis.

CONTRAINDICATIONS:

Not recommended for everyday walking, or for patients allergic to neoprene.

WARNINGS AND PRECAUTIONS:

If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, discontinue use and consult your medical professional immediately.

APPLICATION INFORMATION:

1. Apply splint while seated.
2. Loosen all straps.
3. Slide foot into splint. Ensure toes and ball of the foot are comfortably positioned in the forefoot section of the splint.
4. Secure the forefoot straps first, then the calf strap. Make sure all straps are secure but not too tight.

CLEANING INSTRUCTIONS:

Wipe clean with warm water. Line dry. Avoid heat or sun.

MATERIAL CONTENTS:

Polyurethane, Polyester, Nylon, PVC

WARRANTY:

DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

FOR SINGLE PATIENT USE ONLY.

NOT MADE WITH NATURAL RUBBER LATEX.

NOTICE:

WHILE EVERY EFFORT HAS BEEN MADE IN STATE OF THE ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THIS DEVICE IS ONLY ONE ELEMENT IN THE OVERALL TREATMENT PROGRAM ADMINISTERED BY A MEDICAL PROFESSIONAL. THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.



DJO, LLC

T 800.336.6569 D 760.727.1280 F 800.936.6569
2900 Lake Vista Drive | Lewisville, TX 75067 | U.S.A.
DJOglobal.com