Title: INSTRUCTIONS FOR USE – DJO SURGICAL SHOULDER SYSTEM

Revision	ECO Date	ECO	Summary of Changes
J	01/18/11	14318	MULTIPLE CHANGES TO SECTIONS 1, 6 AND 9.
K	10/13/11	16148	MULTIPLE CHANGES TO SECTION 2, 3 AND 9. ADD INDICATIONS FOR THE HUMERAL STEM ADAPTERS. CHANGE THE DJO SURGICAL COMPANY LOGO.
L	11/26/12	18880	REMOVE THE TURKISH REPRESENTATIVE BECAUSE THEY ARE NO LONGER OUR REPRESENTATIVE. UPDATE SECTION 2 (PRODUCT DESCRIPTION) TO ADD THE TURON e+ SHOULDER GLENOID. UPDATE SECTION 9 (STERILIZATION) TO ADD THE TURON e+ GLENOID STERILIZATION REQUIREMENTS. REMOVE THE "USER RESTERILIZATION GUIDELINES FOR IMPLANTS" INFORMATION PER FDA REQUEST.
M	XX/XX/14	DC-00918	UPDATE SECTION 2.0 (PRODUCT DESCRIPTION AND IMPLANT MATERIALS) TO ADD THE BILOX DELTA CERAMIC TURON SHOULDER HUMERAL HEADS.

0400-0061 Rev. M 2014-08

EN

1. Product Handling

Implants are provided sterile and should always be stored unopened in their respective protective containers. Prior to use, inspect package for damage that may compromise sterility. If packaging has been opened or damaged upon receipt, contact the manufacturer's representative. Also inspect the labeling to verify that the expiration date has not passed. If the product is expired, contact Customer Service and do not use the implant. When unpacking the implant, verify the labeling for correct Reference No. and size. When removing the implant from its packaging, the relevant aseptic instructions must be observed. Protect prosthesis from contact with objects that may damage the surface finish. Inspect each implant prior to use for visual damage. This implant is part of a system and should be used only in combination with other original DJO Surgical product belonging to the same shoulder system, unless otherwise specified.

2. Product Description and Implant Materials

Component	Fixation Method	Material	Applicable ASTM Standard	Applicable ISO Standard
FOUNDATION® Humeral Stem ¹	Cemented	Ti6Al4V alloy or Cast Titanium	ASTM F136 / ASTM F1472	ISO 5832-3
		Ti6Al4V Titanium Plasma Spray	ASTM 1580	ISO 5832-3
FOUNDATION® Humeral Head ¹	Cementless or Cemented	CoCrMo	ASTM F799	ISO 5832-4
FOUNDATION® Glenoid ¹	Cemented	Medical grade UltraHigh Molecular Weight Polyethylene	ASTM F648	ISO 5834-1 & 2
TURON™ Shoulder Humeral Stem ²	Cementless	Ti6Al4V alloy or Cast Titanium	ASTM F136 / ASTM F1472	ISO 5832-3
		Ti6Al4V Titanium Plasma Spray	ASTM 1580	ISO 5832-3
TURON™ Shoulder Humeral Head ²	Cemented or Cementless	CoCrMo	ASTM F799	ISO 5832-4
		BioloX® delta Ceramic ³		ISO 6474
TURON™ Shoulder Humeral Neck ²	Cemented or Cementless	CoCrMo	ASTM F799	ISO 5832-4
		Ti6Al4V alloy	ASTM F136/F1472	ISO 5832-3
TURON™ Shoulder Glenoid ²	Cemented	Medical grade UltraHigh Molecular Weight Polyethylene	ASTM F648	ISO 5834-1 & 2
		Highly Cross-Linked Vitamin E UHMWPE	ASTM F2695 ASTM F2565	
Humeral Stem Adapter ⁴	Cementless	Ti6Al4V alloy	ASTM F136 / ASTM F1472	ISO 5832-3
TURON™ Shoulder Metal-Backed Revision Glenoid ²	Cementless	Ti-6Al-4V	ASTM F136 / ASTM F1472	ISO 5832-3
		Hydroxylapatite (HA) Coating	ASTM F1185	ISO 13779-1
		Medical grade UltraHigh Molecular Weight Polyethylene	ASTM F648	ISO 5834-1 & 2

¹ These implants are part of a system and should be used only in combination with other original FOUNDATION® product belonging to the same system.

² These implants are part of a system and should be used only in combination with other original TURON™ SHOULDER product belonging to the same system.

³ These implants are part of the TURON™ SHOULDER system and should be used only in combination with other original TURON™ SHOULDER Ti6Al4V alloy Humeral Heads.

⁴ This implant should be used with the TURON™ SHOULDER Humeral Head.

3. Indications

Joint replacement is indicated for patients suffering from disability due to:

- noninflammatory degenerative joint disease including osteoarthritis and avascular necrosis of the natural humeral head and/or glenoid;
- rheumatoid arthritis;
- correction of functional deformity;
- humeral fracture.

This device may also be indicated in the salvage of previously failed surgical attempts.

Indications for RSP Humeral Stem Adapters:

The Reverse® Shoulder Prosthesis (RSP) is indicated for treatment of patients with a grossly rotator cuff deficient shoulder joint with severe arthropathy or a previously failed joint replacement with a grossly rotator cuff deficient shoulder joint. The patient's joint must be anatomically and structurally suited to receive the selected implant(s), and a functional deltoid muscle is necessary to use the device. The glenoid baseplate is intended for cementless application with the addition of screws for fixation. The humeral stem is intended for cemented use only.

During primary surgery, after the humerus is prepared for the RSP humeral stem (modular and monoblock), if purchase to the glenoid bone is insufficient to bear the load of the glenoid baseplate and alternative glenoid bone reconstruction and/or repair is inadequate, the corresponding RSP humeral stem adapter can be used to convert the RSP humeral stem to hemiarthroplasty prosthesis as a salvage procedure. During revision surgery of an RSP (modular or monoblock), if the glenoid bone stock appears to be insufficient to bear the load of the glenoid baseplate and alternative glenoid bone reconstruction and/or repair is inadequate, the corresponding RSP humeral stem adapter can be used to convert the RSP device to hemiarthroplasty prosthesis as a salvage procedure. For modular RSP stems, the Foundation Shoulder humeral head should be used. For the monoblock stem, the Turon humeral head should be used.

This stem/adapter construct is not approved for use as a surrogate for traditional hemiarthroplasty or anatomic replacement indications.

Indications for use for the TURON™ Shoulder Metal-Backed Glenoid:

The TURON™ Revision Glenoid is indicated for cementless use for the following indications:

- noninflammatory degenerative joint disease including osteoarthritis and avascular necrosis;
- rheumatoid arthritis;
- correction of functional deformity;
- other difficult management problems where arthrodesis or resectional arthroplasty are not acceptable (e.g. revision of a failed primary component);
- cuff tear arthropathy.

4. Intended Use

DJO Surgical shoulder devices are intended for treatment of patients who are candidates for shoulder arthroplasty per the indications for use. While shoulder replacements are not intended to withstand activity levels and loads of normal healthy bone, they are a means of restoring mobility and reducing pain for many patients.

5. Contraindications

Joint replacement is contraindicated where there is:

- infection or sepsis;
- insufficient bone quality which may affect the stability of the implant;
- muscular, neurological or vascular deficiencies, which compromise the affected extremity;
- alcoholism or other addictions;
- materials sensitivity;
- loss of ligamentous structures;
- high levels of physical activity (e.g. competitive sports, heavy physical labor).

6. Precautions and Warnings

An implant should never be reused. Although the implant may appear undamaged, previous stresses could create imperfections that may lead to mechanical failure. It is advised to utilize new prostheses of current design.

Familiarity with, and attention to the surgical technique recommended for this device is imperative for best results. The correct selection as well as the correct seating/placement of the prosthetic implant is extremely important. Only DJO Surgical Shoulder System instruments and trial prostheses should be used.

Care must be taken to protect mating surfaces (i.e. tapers) and polished bearing surfaces from nicks and scratches which could become the focal point for failure. Contouring or bending of the implant may reduce its service life and may cause immediate or eventual failure under load. An implant must not be tampered with, as tampering will adversely affect the performance of the implant.

DJO Surgical (Hip, Knee, and Shoulder) systems have not been evaluated for safety and compatibility in the Magnetic resonance environment. The (Hip, Knee, Shoulder) stems have not been tested for heating or migration in the MR environment.

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

7. Preoperative Planning and Postoperative Care

Preoperative planning provides essential information regarding the appropriate prosthesis and likely combinations of components. Use instrument trial components for fit verification (where applicable) and extra implant components for backup. X-ray templates for all sizes of the FOUNDATION® and TURON™ SHOULDER systems are available upon request.

Accepted surgical practices should be followed for postoperative care. The patient should be made aware of the limitation of total joint reconstruction. Excessive physical activity and trauma affecting the replaced joint have been implicated in premature failure by loosening, fracture, and/or wear of the prosthetic implants. The patient should be cautioned to govern his/her activities accordingly as the risk of implant failure increases with weight and activity levels of the patient.

8. Adverse Effects

- Accelerated wear of the polyethylene articulating surfaces have been reported following total shoulder replacement. Such wear may be initiated by particles of cement, metal, or other debris which can cause abrasion of the articulating surfaces. Accelerated wear shortens the useful life of the prosthesis, and leads to early revision surgery to replace the worn prosthetic components.
- Metallosis and osteolysis may be implicated from wear debris associated with the use of orthopedic implants.
- Peripheral neuropathies have been reported following total joint surgery. Subclinical nerve damage occurs more frequently, possibly the result of surgical trauma.
- Metal sensitivity reactions in patients following joint replacement have been rarely reported. Implantation of foreign material in tissues can result in histological reactions involving macrophages and fibroblasts. The clinical significance of this effect is uncertain, as similar changes may occur as a precursor to, or during the healing process. In some cases, wear debris can initiate the process of histiocytic granuloma formation and consequent osteolysis and loosening of the implant.
- Dislocation and subluxation of implant components can result from improper positioning of the components. Muscle and fibrous tissue laxity can also contribute to these conditions.
- Implants can loosen or migrate due to trauma or loss of fixation.
- Infection can lead to failure of the joint replacement.
- While rare, fatigue fracture of the implant can occur as a result of strenuous activity, improper alignment, or duration of service.
- Fracture of the humerus can occur while press-fitting (seating) the humeral stem into the prepared humeral canal.
- Allergic reactions.

Intraoperative and early postoperative complications can include:

- humeral perforation, or fracture;
- humeral fracture can occur while seating the device;
- damage to blood vessels;
- temporary or permanent nerve damage resulting in pain or numbness of the affected limb;
- undesirable shortening or lengthening of the limb;
- traumatic arthrosis of the shoulder from intraoperative positioning of the extremity;
- cardiovascular disorders including venous thrombosis, pulmonary embolism, or myocardial infarction;
- hematoma;
- delayed wound healing; and,
- infection.

Late postoperative complications can include:

- avulsion as a result of excess muscular weakening;
- non-union due to inadequate reattachment and/or early weight bearing;
- aggravated problems of other joints of the affected limb or muscle deficiencies;
- humeral fracture by trauma or excessive loading, particularly in the presence of poor bone stock;
- periarticular calcification or ossification, with or without impediment to joint mobility;
- inadequate range of motion due to improper selection or positioning of components, by impingement, and calcification.

9. Sterilization

Unless opened or damaged, DJO Surgical implants are supplied sterile in multiple pouches or barrier blister trays. Upon receipt, check all packaging for punctures or other damage. If packaging is opened or damaged, contact the manufacturer or manufacturer's representative for instructions.

Sterilization of implants, other than the Turon e+ glenoid, is by gamma radiation at the minimum dose of 25 kGy to achieve a Sterility Assurance Level (SAL) of 10⁻⁶. Sterilization of the Turon e+ glenoid is performed by hydrogen peroxide gas plasma to achieve a Sterility Assurance Level (SAL) of 10⁻⁶. These glenoids are single-use devices and CANNOT be resterilized by a healthcare facility. Implants are single-use devices. Trials and other instruments are used to determine sizing before the sterile package needs to be opened. Should the original sterile package be inadvertently opened or compromised before implantation, the device cannot be implanted. Contact manufacturer or manufacturer's representative for instructions. Do not resterilize an implant or component that has been in contact with or contaminated by blood or other substances. Do not try to clean an implant since standard procedures cannot be relied upon to remove contamination from porous coating and storage of any opened implant or component should be avoided.

Instruments are provided nonsterile and should be stored in their original packaging until cleaned and sterilized according to the recommended guidelines found in the DJO Surgical Instrumentation Instructions for Use.

WARNING: DO NOT resterilize any shoulder prosthesis distributed by DJO Surgical (Encore Medical, L.P.) if sterile packaging is opened or damaged upon receipt. Return the implant with respective packaging to DJO Surgical for inspection and disposition.

WARNING: Protect all porous coated and polished surfaces. Standard cleaning procedures cannot be relied upon to remove contamination from porous coating.

WARNING: DO NOT resterilize UHMWPE (ultra-high molecular weight polyethylene) or HA (hydroxylapatite) coated implants.

DJO Surgical has validated the above steam sterilization cycles and has data on file. Other sterilization cycles may also be suitable; however individuals or hospitals not using the recommended method are advised to validate any alternate method using appropriate laboratory techniques. Proper validation of the autoclave is essential to ensure proper sterilization temperatures and cycle times.

NOTE: DJO Surgical does not recommend Flash or Chemical Sterilization.

For further information regarding the use of the DJO Surgical Shoulder Systems, contact your DJO Surgical representative or distributor.

DJO Surgical Shoulder Systems are manufactured by ENCORE MEDICAL, L.P.
9800 Metric Blvd., Austin, TX 78758 USA (Made in the USA)

10. Trademarks and patents

FOUNDATION® is a trademark of DJO Surgical.

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