
Enovis™ Knee Systems  
 Enovis™ Knie-Systeme  
 Systèmes de genou Enovis™  
 Sistemas de rodilla Enovis™  
 Sistemi protesici per ginocchio Enovis™  
 Συστήματα γονάτου της Enovis™  
 Enovis™ Diz Sistemleri

Encore Medical, L.P.  
 9800 Metric Blvd.  
 Austin, TX 78758

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0400-0294 Rev. B 2022-10

A printable copy of the IFU for this device can be located at: [www.djosurgicalifus.com](http://www.djosurgicalifus.com). A paper copy can be requested via phone at +1-800-520-8976.

EN

1. Product Handling

Implants are provided sterile and should always be stored unopened in their respective protective containers. Prior to use, inspect package for damage that may compromise sterility. If packaging has been opened or damaged upon receipt, contact the manufacturer's representative. Also inspect the labeling to verify that the expiration date has not passed. If the product is expired, contact Customer Service and do not use the implant. When unpacking the implant, verify the labeling for correct Reference No. and size. When removing the implant from its packaging, the relevant aseptic instructions must be observed. Protect prosthesis from contact with objects that may damage the surface finish. Inspect each implant prior to use for visual damage. This implant is part of a system and should be used only in combination with other original Surgical division of Enovis™ product belonging to the same knee system, unless otherwise specified.

2. Product Description and Implant Materials

Component	Fixation Method	Material	Applicable ASTM Standard	Applicable ISO Standard
<b>3DKnee™ System</b>				
3DKNEE™ Femur	Cemented	CoCrMo	ASTM F75	ISO 5832-4
3DKNEE™ Tibial Baseplate FK/3DKNEE™ Downsize Tibial Baseplate	Cemented	CoCrMo	ASTM F75	ISO 5832-4
3DKNEE™ Insert, e+™ with Attachment Screw	Cementless	Medical grade UltraHigh Molecular Weight Polyethylene (Moderately Cross-Linked)	ASTM F648 / ASTM F2565	ISO 5834-1 / ISO 5834-2
		Vitamin E UHMWPE (a-tocopheral)	ASTM F2695	
		Ti6Al4V alloy	ASTM F136 / ASTM F1472	ISO 5832-3
1. 3DKNEE™ femurs are to be used with 3DKNEE™ tibial inserts and FK/3DKNEE™, or 3DKNEE™ tibial baseplates. If using FK/3DKNEE™ or tibial baseplates with 3DKNEE™ femurs, use baseplates that are one size smaller. 2. Size 3, 5, and 7 3DKNEE™ femurs and inserts are to be used with 3DKNEE™ Tibial Components that are one size larger or FK/3DKNEE™ downsize components that are one size smaller. 3. The 3DKNEE™ femur and the 3DKNEE™ insert size must match. The insert will mate with the same size 3DKNEE™ baseplate or one size smaller or FK/3DKNEE™ baseplate. 4. Cemented items to be used with PMMA bone cement only.				
<b>EMPOWR KNEE System®</b>				
EMPOWR 3D™ Femur <sup>1,4</sup> EMPOWR PS™ Femur <sup>2,3</sup>	Cemented	CoCrMo	ASTM F75	ISO 5832-4
EMPOWR™ Tibial Baseplate	Cemented	CoCrMo	ASTM F75	ISO 5832-4
EMPOWR™ Tibial Insert: - EMPOWR 3D™ Tibial Insert - e+™. <sup>1</sup> - EMPOWR PS™ Tibial Insert - e+™. <sup>2</sup>	Cementless	Medical grade UltraHigh Molecular Weight Polyethylene (Moderately Cross-Linked)	ASTM F648 / ASTM F2565	ISO 5834-1 / ISO 5834-2
		Vitamin E UHMWPE (a-tocopheral)	ASTM F2695	
1. Size interchangeability between the EMPOWR 3D™ femur to EMPOWR 3D™ tibial inserts are compatible with the same size component, or one size smaller femur on a larger insert. 2. Size interchangeability between the EMPOWR PS™ femur and EMPOWR PS™ tibial inserts are compatible with the same size component, one size smaller femur on a larger insert, or one size larger femur on a smaller insert. 3. Cemented items to be used with PMMA bone cement only.				

Patellae				
Patella – All Poly, Domed <sup>1</sup>	Cemented	Medical grade UltraHigh Molecular Weight Polyethylene	ASTM F648	ISO 5834-1 / ISO 5834-2
Patella – All Poly e+™, Domed <sup>2</sup>	Cemented	Medical grade UltraHigh Molecular Weight Polyethylene (Moderately Cross-Linked)	ASTM F648 / ASTM F2565	ISO 5834-1 / ISO 5834-2
		Vitamin E UHMWPE (a-tocopheral)	ASTM F2695	
1. The 3D Knee System is intended for use with the All Poly, Domed Patella only. 2. The EMPOWR Knee System® is intended for use with the e+™ Domed Patella only. 3. Cemented items to be used with PMMA bone cement only.				

Depending on the Enovis™ Knee System, femoral prostheses may be available in left and right configurations. Depending on the Enovis™ Knee System, the stemmed baseplate is available in left and right configurations.

### 3. Indications

In the EU, the 3DKnee is only approved for primary arthroplasty. The 3DKnee is indicated for use as a total knee replacement for the relief of pain and significant disability resulting from degenerative joint diseases such as primary or secondary osteoarthritis or rheumatoid arthritis

The EMPOWR Knee is indicated for patients suffering from disability due to:

- degenerative, post-traumatic or rheumatoid arthritis;
- avascular necrosis of the femoral condyle;
- post-traumatic loss of joint configuration, particularly when there is patellofemoral erosion, dysfunction or prior patellectomy;
- moderate valgus, varus or flexion deformities;
- treatment of fractures that are unmanageable using other techniques.

This device may also be indicated in the salvage of previously failed surgical attempts. All devices are intended for cemented applications.

### 4. Intended Use

Enovis™ knee devices are intended for treatment of patients who are candidates for knee arthroplasty per the Indication for use. While total knee replacements are not intended to withstand activity levels and loads of normal healthy bone, they are a means of restoring mobility and reducing pain for many patients.

### 5. Contraindications

Joint replacement is contraindicated where there is:

- infection (or a history of infection), acute or chronic, local or systemic;
- insufficient bone quality which may affect the stability of the implant;
- muscular, neurological or vascular deficiencies, which compromise the affected extremity;
- obesity;
- alcoholism or other addictions;
- materials sensitivity;
- untreated medial and lateral ligamentous instability
- high levels of physical activity (e.g. competitive sports, heavy physical labor).

### 6. Precautions and Warnings

An implant should never be reused. Although the implant may appear undamaged, previous stresses could create imperfections that may lead to mechanical failure. It is advised to utilize new prostheses of current design.

Familiarity with, and attention to the surgical technique recommended for this device is imperative for best results. The correct selection as well as the correct seating/placement of the prosthetic implant is extremely important. Malposition may predispose the device to excess wear and early failure. Only Enovis™ Knee System implants, instruments and trial prostheses should be used.

Care must be taken to protect mating surfaces (i.e. tapers) and polished bearing surfaces from nicks and scratches which could become the focal point for failure. Contouring or bending of the implant may reduce its service life and may cause immediate or eventual failure under load. An implant must not be tampered with, as tampering will adversely affect the performance of the implant.

The manufacturer has not evaluated the use of the device in children under 18.

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

### 7. Preoperative Planning and Postoperative Care

Preoperative planning provides essential information regarding the appropriate prosthesis and likely combinations of components. Use instrument trial components for fit verification (where applicable) and extra implant components for backup. X-ray templates for all sizes of the Enovis™ Knee systems are available upon request.

Accepted surgical practices should be followed for postoperative care. The patient should be made aware of the limitation of total joint reconstruction. Excessive physical activity and trauma affecting the replaced joint have been implicated in premature failure by loosening, fracture, and/or wear of the prosthetic implants. The patient should be cautioned to govern his activities accordingly as the risk of implant failure increases with weight and activity levels of the patient.

### 8. MRI Compatibility

EMPOWR Knee systems listed in Section 2 of this document have not been evaluated for safety and compatibility in the Magnetic resonance environment. These devices have not been tested for heating, migration, or image artifact in the MR environment. The safety of these Enovis™ components in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Non-clinical testing has demonstrated that the 3DKnee systems contained in Section 2 listed above are MR Conditional. Patients can be scanned safely under the following conditions:

- Static magnetic field of 1.5-Tesla (1.5T) or 3.0-Tesla (3.0T).
- Spatial gradient field of up to:
  - 3,930 G/cm (39.3 T/m) for 1.5T systems.
  - 1,960 G/cm (19.6 T/m) for 3.0T systems.
- Maximum whole body averaged specific absorption rate (SAR) of:
  - 1.0 W/kg for 15 minutes of scanning in Normal Operating Mode at 1.5T.
  - 2.0 W/kg for 15 minutes of scanning in Normal Operating Mode at 3.0T.

### 3.0T RF heating

In non-clinical testing with body coil excitation, representative devices produced a temperature rise of less than 2.0°C at a maximum whole body averaged specific absorption rate (SAR) of 2.0 W/kg, as assessed by calorimetry for 15 minutes of scanning in a 3.0T Siemens Trio (MRC20587) MR scanner with SYNGO MR A30 4VA30A software.

### 1.5T RF heating

In non-clinical testing with body coil excitation, representative devices produced a temperature rise of less than 5.0°C at a maximum whole body averaged specific absorption rate (SAR) of 1.0 W/kg, as assessed by calorimetry for 15 minutes of scanning in a 1.5T Siemens Espree (MRC30732) MR scanner with SYNGO MR B17 software.

Caution: The RF heating behavior does not scale with static field strength. Devices which do not exhibit detectable heating at one field strength may exhibit high values of localized heating at another field strength.

### MR Artifact

In testing using a 3.0T system with gradient-echo sequencing, the shape of the image artifact follows the approximate contour of the device and extends radially up to 7.7 cm from the implant

Note: Patients receiving MRI should be made aware of risks associated with this procedure. This could include the following:

- "The strong, static magnetic field of the MRI scanner will pull on magnetic materials and may cause unwanted movement of the medical device."
- "The radiofrequency energy and magnetic fields that change with time may cause heating of the implanted medical device and the surrounding tissue, which could lead to burns."
- "The presence of the medical device will degrade the quality of the MR image, which may make the MRI scan uninformative or may lead to an inaccurate clinical diagnosis, potentially resulting in inappropriate medical treatment."

### 9. Adverse Effects

Some of the adverse effects that could occur related to total knee arthroplasty are:

- fracture of the tibia or femur;
- transient peroneal palsy secondary to surgical manipulation;
- patellar subluxation or dislocation;
- patella femoral impingement;
- instability, changes in position, or loosening of components from bone interfaces;
- ligamentous laxity;
- dissociation of components;
- infection;
- poor range of motion;
- shortening of limbs;
- lengthening of limbs if severe deformity is present
- stiffness post-operatively
- metal sensitivity reactions.

### 10. Sterilization

Unless opened or damaged, Enovis™ implants are supplied sterile in multiple pouches or barrier blister trays. Upon receipt, check all packaging for punctures or other damage. If during inspection, packaging is found opened or damaged, contact manufacturer or manufacturer's representative for instructions.

Sterilization of implants other than Tibial Inserts and Patellae manufactured from Moderately Cross-Linked Polyethylene with Vitamin E (e+™) is performed by gamma radiation at the minimum dose of 25 kGy to achieve a Sterility Assurance Level (SAL) of 10<sup>-6</sup>.

Sterilization of the Tibial Inserts and Patellae manufactured from Moderately Cross-Linked Polyethylene with Vitamin E (e+™) is performed by hydrogen peroxide gas plasma to achieve a Sterility Assurance Level (SAL) of 10<sup>-6</sup>.

Implants are single-use devices. Trials and other instruments are used to determine sizing before the sterile package needs to be opened. Should the integrity of the original sterile package be lost by being opened, punctured, or torn before implantation in the surgical field, contact manufacturer or manufacturer's representative for instructions. These inserts are single-use devices and CANNOT be resterilized by a healthcare facility. Contact manufacturer or manufacturer's representative for instructions.

Do not resterilize an implant or component that has been opened outside of the surgical field or in contact with or contaminated by blood or other substances. Do not try to clean an implant since standard procedures cannot be relied upon to remove contamination from the implant or component and storage of the opened implant or component should be avoided.

Instruments are provided nonsterile and should be stored in their original packaging until cleaned and sterilized according to the recommended guidelines found in the Enovis™ Instrumentation Instructions for Use.

**WARNING: DO NOT resterilize any knee prosthesis distributed by Enovis™ (Encore Medical, L.P.) if sterile packaging is opened or damaged upon receipt. Return the implant with respective packaging to Enovis™ for inspection and disposition.**

**WARNING: Protect all porous coated and polished surfaces. Standard cleaning procedures cannot be relied upon to remove contamination from porous coating.**

**WARNING: DO NOT resterilize UHMWPE (ultra-high molecular weight polyethylene) and Moderately Cross-Linked Polyethylene Vitamin E (e+™) implants.**

Enovis™ has validated sterilization cycle data on file.

**NOTE:** Enovis™ does not recommend Flash or Chemical Sterilization.

For further information regarding the use of the Enovis™ Knee Systems, contact your Enovis™ representative or distributor.

Enovis™ Knee Systems are manufactured by ENCORE MEDICAL, L.P.  
9800 Metric Blvd., Austin, TX 78758 USA (Made in the USA)

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### 11. Ticari markalar ve patentler

FK/3DKNEE™, 3DKNEE™, EMPower Knee System™, EMPower 3D™, EMPower 3D KNEE®, EMPower PS™, EMPower PS KNEE®, DJO EMPower KNEE®, e+™ ; Enovis™'in tescilli markalarıdır. ABD patentleri: 5,413,604

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Bone Cement Usage – The following legends are displayed on the product labeling to indicate bone cement usage:

Usage	Legend
Implants intended to be used with bone cement	CEMENTED
Implants intended to be used without bone cement	CEMENTLESS
Implants intended to be used optionally	NO LEGEND

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### 11. Trademarks and patents
FK/3DKNEE™, 3DKNEE™, EMPOWR Knee System™, EMPOWR 3D™, EMPOWR 3D KNEE®, EMPOWR PS™, EMPOWR PS KNEE®, DJO EMPOWR KNEE®, e+™ are trademarks of Enovis™.