Reverse® Shoulder Prosthesis System  
 Reverse® Schulterprothesen-System  
 Système de prothèse d'épaule Reverse®  
 Sistema de prótesis de hombro Reverse®  
 Sistema protesico per spalla Reverse®

Encore® Medical, L.P.  
 9800 Metric Blvd.  
 Austin, TX 78758-5445 USA

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0400-0297 REV. C 2024-01

EN

A printable copy of the IFU for this device can be located at: [www.djosurgicalifus.com](http://www.djosurgicalifus.com). A paper copy can be requested via phone at +1-800-520-8976. The legal entity Encore Medical L.P. does business under the company name DJO Surgical®\* and shall be referenced as such throughout this IFU.

### 1. Product Handling

Devices not returned to DJO® should be treated as biohazardous material and disposed of in accordance with local laws and regulations.

Implants are provided sterile and should always be stored unopened in their respective protective containers. Prior to use, inspect package for damage that may compromise sterility. If packaging has been opened or damaged upon receipt, do not use the product and contact the manufacturer's representative. Also inspect the labeling to verify that the expiration date has not passed. If the product is expired, contact Customer Service and do not use the implant. When unpacking the implant, verify the labeling for correct Reference No. and size. When removing the implant from its packaging, the relevant aseptic instructions must be observed. Protect prosthesis from contact with objects that may damage the surface finish. Inspect each implant prior to use for visual damage. This implant is part of a system and should be used only in combination with other original DJO Surgical® product belonging to the same shoulder system, unless otherwise specified.

### 2. Product Description and Implant Materials<sup>1</sup>

Component	Fixation Method	Material	Applicable ASTM Standard	Applicable ISO Standard
<b>Reverse® Shoulder Prosthesis (RSP®) and Altivate Reverse® Humeral Stems</b>				
RSP® Monoblock Humeral Stem (Primary and Revision)	Cemented or Cementless	Ti6Al4V alloy Ti Plasma Coating	ASTM F136 / ASTM F1472 ASTM F1580	ISO 5832-3
Altivate Reverse® Standard Shell Humeral Stem (Primary and Revision)	Cemented or Cementless	Ti6Al4V alloy CP Ti Porous Coating	ASTM F1472 ASTM F67	ISO 5832-3 ISO 5832-2
Altivate Reverse® Small Shell Humeral Stem (Primary and Revision)	Cemented or Cementless	Ti6Al4V alloy CP Ti Porous Coating	ASTM F1472 ASTM F67	ISO 5832-3 ISO 5832-2
Cemented items to be used with PMMA bone cement only.				
<b>Reverse® Shoulder Prosthesis (RSP®) and Altivate Reverse® Humeral Inserts</b>				
RSP® Humeral Socket Insert (Standard and Semi-Constrained)	Cementless	Medical Grade Ultra High Molecular Weight Polyethylene	ASTM F648	ISO 5834-1 / ISO 5834-2
RSP® Humeral Socket Insert e+™* (Standard and Semi-Constrained)	Cementless	Medical Grade Ultra High Molecular Weight Polyethylene (Highly Cross-Linked) Vitamin E UHMWPE (a-tocopheral)	ASTM F648 / ASTM F2565 ASTM F2695	ISO 5834-1 / ISO 5834-2
Altivate Reverse® Small Socket Insert (Neutral and Semi-Constrained)	Cementless	Medical Grade Ultra High Molecular Weight Polyethylene	ASTM F648	ISO 5834-1 / ISO 5834-2
Altivate Reverse® Small Socket Insert e+™ (Neutral and Semi-Constrained)	Cementless	Medical Grade Ultra High Molecular Weight Polyethylene (Highly Cross-Linked) with Vitamin E (a-tocopheral)	ASTM F648 / ASTM F2565 ASTM F2695	ISO 5834-1 / ISO 5834-2
<b>Reverse® Shoulder Prosthesis (RSP®) and Altivate Reverse® Humeral Adapters and Spacers Turon® Humeral Heads</b>				
RSP® Monoblock Hemi-Adapter with Retaining Screw <sup>1</sup>	Cementless	CoCrMo Ti6Al4V alloy	ASTM F799 ASTM F136 / ASTM F1472	ISO 5832-4 ISO 5832-3
RSP® Monoblock Hemi-Adapter Retaining Screw <sup>5</sup>	Cementless	Ti6Al4V alloy	ASTM F136 / ASTM F1472	ISO 5832-3
RSP® Monoblock Spacer with Retaining Screw <sup>2</sup>	Cementless	Ti6Al4V alloy	ASTM F136 / ASTM F1472	ISO 5832-3
RSP® Monoblock Spacer Retaining Screw <sup>6</sup>	Cementless	Ti6Al4V alloy	ASTM F136 / ASTM F1472	ISO 5832-3
Altivate Reverse® Small Spacer with Retaining Screw <sup>3</sup>	Cementless	Ti6Al4V alloy	ASTM F136 / ASTM F1472	ISO 5832-3
Altivate Reverse® Small Hemi-Adapter with Retaining Screw <sup>3</sup>	Cementless	CoCrMo	ASTM F1537	ISO 5832-4

<sup>1</sup> Refer to Section 18 for implant composition.

		Ti6Al4V alloy	ASTM ASTM F1472	ISO 5832-3
AltiVate® Anatomic to Reverse® Conversion Module <sup>7</sup>	Cementless	Ti6Al4V alloy	ASTM F1472	ISO 5832-3
Turon® Shoulder Humeral Head (Neutral and Offset) <sup>4</sup>	Cementless	CoCrMo	ASTM F799	ISO 5832-4
AltiVate® Anatomic Humeral Head (Neutral and Offset) <sup>4</sup>	Cementless	CoCrMo	ASTM F1537	ISO 5832-4
1. The RSP® Monoblock Hemi-Adapter is compatible with the RSP® Monoblock and AltiVate Reverse® Standard Shell Stems above. 2. The RSP® Monoblock Spacer is compatible with the RSP® Monoblock and AltiVate Reverse® Standard Shell Stems above. 3. The AltiVate Reverse® Small Spacer and Small Hemi-Adapter are only compatible with the AltiVate Reverse® Small Shell Stems above. 4. The RSP® Monoblock and AltiVate Reverse® Small Hemi-Adapters are compatible with the Turon® Humeral Heads listed above or AltiVate® Anatomic Humeral Heads <sup>4</sup> 5. The RSP® Monoblock Hemi-Adapter Retaining Screw is also compatible with the AltiVate Reverse® Small Hemi-Adapter. 6. The RSP® Monoblock Spacer Retaining Screw is also compatible with the AltiVate Reverse® Small Spacer. 7. The RSP® Monoblock Socket Inserts and Glenoid Systems are compatible with AltiVate® Anatomic to Reverse® Conversion Module <sup>4</sup>				
<b>Reverse® Shoulder Prosthesis (RSP®) Glenoid System</b>				
RSP® Glenoid Head with Retaining Screw	Cementless	CoCrMo Ti6Al4V alloy	ASTM F799 ASTM F136 / ASTM F1472	ISO 5832-4 ISO 5832-3
RSP® Glenoid Head Retaining Screw	Cementless	Ti6Al4V alloy	ASTM F136 / ASTM F1472	ISO 5832-3
RSP® Glenoid Baseplate with HA Coating	Cementless	Ti6Al4V alloy CP Ti Porous Coating Hydroxyapatite Coating	ASTM F136 / ASTM F1472 ASTM F67 ASTM F1185	ISO 5832-3 ISO 5832-2 ISO 13779-1
RSP® Glenoid Baseplate Screws (Locking and Non-Locking)	Cementless	Ti6Al4V alloy	ASTM F136 / ASTM F1472	ISO 5832-3
The RSP® Glenoid system is compatible with all RSP® and AltiVate Reverse® Shoulder Humeral systems listed above				

\*Reference DJO IFU 0400-0257 for Instructions for Use for the AltiVate® Anatomic components

### 3. Intended Purpose:

DJO Surgical® Reverse Shoulder Systems are intended to replace all or a portion of a shoulder joint in patients meeting the indications for use.

### 4. Intended Use:

DJO Surgical® shoulder devices are intended for treatment of patients who are candidates for shoulder arthroplasty per the Indication for use. While total shoulder replacements are not intended to withstand activity levels and loads of normal healthy bone, they are a means of restoring mobility and reducing pain for many patients.

### 5. Indications for Use:

#### Indications for RSP® Monoblock Stem:

The Reverse® Shoulder Prosthesis Monoblock is indicated for patients with a functional deltoid muscle with a grossly deficient rotator cuff shoulder joint with severe arthropathy or a previously failed joint replacement with a grossly deficient rotator cuff shoulder joint:

- In cases of fracture of glenohumeral joint from trauma or pathologic conditions of the shoulder, including humeral head fracture or displaced 3- or 4-part fractures of proximal humerus.
- In cases of bone defect in proximal humerus.

The patient's joint must be anatomically and structurally suited to receive the selected implant(s).

The glenoid baseplate is intended for cementless application with the addition of screws for fixation. The humeral stem is intended for cemented or Cementless use (Cementless use not cleared in the EU).

The RSP® Monoblock Socket Inserts and Glenoid systems are indicated to be used with the AltiVate® Anatomic to Reverse® Conversion Module and previously implanted AltiVate® Anatomic stem for revision surgeries in patients with a grossly rotator cuff deficient shoulder joint with severe arthropathy or a previously failed joint replacement with a grossly rotator cuff deficient shoulder joint. The patient's joint must be anatomically and structurally suited to receive the selected implant(s), and a functional deltoid muscle is necessary to use the device.

#### Additional Indications for RSP® Monoblock Hemi-Adapter used with RSP® Monoblock:

During primary surgery, after the humerus is prepared for the RSP® Monoblock humeral stem, if the glenoid bone stock appears insufficient to bear the load of the glenoid baseplate, a RSP® Monoblock Hemi-Adapter can be used to convert the RSP® Monoblock humeral stem to hemiarthroplasty prosthesis. During revision surgery of an RSP® Monoblock, if the glenoid bone stock appears to be insufficient to bear the load of the glenoid baseplate, a RSP® Monoblock Hemi-Adapter is used to convert the RSP® Monoblock device to hemiarthroplasty prosthesis.

#### Indications for AltiVate Reverse® Standard Shell and Small Shell Humeral Stems (Primary and Revision)

#### Hemi Shoulder Indications:

When used with a Hemi-Adapter, the AltiVate Reverse® Shoulder Prosthesis Stem is indicated as a hemi shoulder joint replacement for patients suffering from pain and dysfunction due to:

- Noninflammatory degenerative joint disease including osteoarthritis;
- Inflammatory arthritis of the glenohumeral joint including rheumatoid arthritis;
- Post-traumatic arthritis of the glenohumeral joint;
- Avascular necrosis of the humeral head with and without involvement of the glenoid;
- Correction of functional deformity;
- Rotator cuff tear arthropathy;
- Humeral fracture.
- Failed previous shoulder surgery

#### Reverse Total Shoulder Indications:

The AltiVate Reverse® Shoulder Prosthesis is indicated for patients with a functional deltoid muscle and a grossly deficient rotator cuff joint suffering from pain and dysfunction due to:

- Severe arthropathy with a grossly deficient rotator cuff;
- Previously failed joint replacement with a grossly deficient rotator cuff;
- Fracture of glenohumeral joint from trauma or pathologic conditions of the shoulder including humeral head fracture, displaced 3- or 4-part fractures of proximal humerus, or reconstruction after tumor resection;
- Bone defect in proximal humerus;
- Non-inflammatory degenerative disease including osteoarthritis and avascular necrosis of the natural humeral head and/or glenoid;
- Inflammatory arthritis including rheumatoid arthritis;
- Correction of functional deformity.

The glenoid baseplate is intended for cementless application with addition of screws for fixation. This device may also be indicated in the salvage of previously failed surgical attempts for anatomic and hemi procedures

All RSP® Monoblock and AltiVate Reverse® Standard Shell and Small Shell Humeral Stems are intended for cemented or cementless use.

#### Indications for Turon® Humeral Heads:

The Turon® Humeral Head is indicated to be used in a hemiarthroplasty configuration with a hemi-adapter and reverse stem to:

- Convert to a hemiarthroplasty either during a primary or revision reverse arthroplasty procedure if the glenoid bone stock appears insufficient to bear the load of the glenoid baseplate.

The Turon® Humeral Head is indicated to be used in a hemiarthroplasty configuration with either an AltıVate® Anatomic stem or with one of the reverse stems listed in the table above and a hemi-adapter to treat patients with a functional deltoid muscle suffering from pain and dysfunction due to:

- Noninflammatory degenerative joint disease including osteoarthritis;
- Inflammatory arthritis of the glenohumeral joint including rheumatoid arthritis;
- Post-traumatic arthritis of the glenohumeral joint;
- Avascular necrosis of the humeral head with and without involvement of the glenoid;
- Correction of functional deformity, including fracture malunions
- Rotator cuff tear arthropathy;
- Humeral/ humeral head fracture.
- Failed previous shoulder surgery (if sufficient bone stock remains)

**6. Intended Users**

Users of the DJO Surgical® Reverse® Shoulder Systems should be licensed medical professionals.

**7. Performance Characteristics**

The DJO Surgical Reverse Shoulder Systems replace the worn cartilage of the humeral head and glenoid surfaces with artificial surfaces to restore motion in the shoulder joint. The glenohumeral articulation surfaces provide multi-directional rotation and translation.

**8. Contraindications**

Total joint replacement is contraindicated where there is:

- Infection or Sepsis;
- Insufficient bone quality which may affect the stability of the implant;
- Muscular, neurological, or vascular deficiencies, which compromise the affected extremity;
- Alcoholism or other addictions;
- Materials (metals, etc.) sensitivity;
- Loss of ligamentous structures;
- High levels of physical activity (e.g. competitive sports, heavy physical labor);
- Non-functional deltoid muscle.
- Rotator cuff insufficiency for anatomic shoulder arthroplasty with the AltıVate Reverse® humeral stem;
- Intraoperative conversion from a reverse to an anatomic shoulder

Contraindications for the Hemi-Adapters:

Hemi-arthroplasty is contraindicated where there is:

- Non-functional deltoid muscle;
- Active sepsis;
- Excessive glenoid bone loss;
- Pregnancy;
- Muscular, neurological or vascular deficiencies, which compromise the affected extremity;
- Conditions that place excessive demand on the implant (i.e. Charcot's joints, muscle deficiencies, refusal to modify postoperative physical activities, skeletal immaturity);
- Known metal allergy (i.e., jewelry).

**9. Clinical Benefits**

The AltıVate Reverse and RSP Monoblock Shoulder Systems are intended to replace the shoulder joint in patients who are candidates for shoulder arthroplasty per the indications for use in order to increase quality of life for patients suffering from a painful shoulder joint by reducing pain and restoring shoulder motion.

**10. Precautions and Warnings**

An implant should never be reused. Although the implant may appear undamaged, previous stresses could create imperfections that may lead to mechanical failure. It is advised to utilize new prostheses of current design.

Familiarity with, and attention to the surgical technique recommended for this device is imperative for best results. The correct selection as well as the correct seating/placement of the prosthetic implant is extremely important. Malposition may predispose the device to excess wear and early failure. Only DJO Surgical® Shoulder System instruments and trial prostheses should be used.

Care must be taken to protect mating surfaces (i.e. tapers) and polished bearing surfaces from nicks and scratches which could become the focal point for failure. Contouring or bending of the implant may reduce its service life and may cause immediate or eventual failure under load. An implant must not be tampered with, as tampering will adversely affect the performance of the implant.

The manufacturer has not evaluated the use of the device in children under 18.

This shoulder is a semi-constrained device designed to address irreparable soft tissue, irreparable rotator cuffs, musculature and bony deficiencies. Due to the constraints built into the design, there may be limits to the patient's achievable range of motion. In addition, because of the limit to the range of motion, there may be the possibility of impingement and/or additional wear.

**11. Preoperative Planning and Postoperative Care**

Preoperative planning provides essential information regarding the appropriate prosthesis and likely combinations of components. Use instrument trial components for fit verification (where applicable) and extra implant components for backup. X-ray templates for all sizes of the Reverse® Prosthesis Shoulder and AltıVate Reverse® Shoulder Systems are available upon request.

Accepted surgical practices should be followed for postoperative care. The patient should be made aware of the limitation of total joint reconstruction. Excessive physical activity and trauma affecting the replaced joint have been implicated in premature failure by loosening, fracture, and/or wear of the prosthetic implants. The patient should be cautioned to govern his/her activities accordingly as the risk of implant failure increases with weight and activity levels of the patient.

**12. MRI Safety Information**

AltıVate Shoulder Systems

	
<b>MRI Safety Information</b>	
<b>A person with DJO's AltıVate System may be safely scanned at 1.5T or 3.0T under the following conditions. Failure to follow these conditions may result in injury.</b>	
<b>Parameter</b>	<b>Condition</b>
Device Name	DJO's AltıVate System
Static Magnetic Field Strength (B0)	1.5T and 3T
MR Scanner Type	Cylindrical
B0 Field Orientation	Horizontal
Maximum Spatial Field Gradient	29.2 T/m (2,920 G/cm)
RF Excitation	Circularly Polarized (CP)

RF Transmit Coil Type	Integrated Whole Body Transmit Coil
Operating Mode	Normal Operating Mode
RF Conditions	Maximum Whole-body SAR: 2 W/kg
Scan Duration	Scan duration as detailed in the MRI Restricted Zone Summary with a cooling period of 10 minutes between scans.
Scan Regions	Landmark restrictions as detailed in the MRI Restricted Zone Summary.
Image Artifact	The presence of DJO's Altivate System may produce an image artifact of 7.0 cm. Some manipulation of scan parameters may be needed to compensate for the artifact.

### 13. Adverse Effects

- Accelerated wear of the polyethylene articulating surfaces have been reported following total shoulder replacement. Such wear may be initiated by particles of cement, metal, or other debris which can cause abrasion of the articulating surfaces.
- Accelerated wear shortens the useful life of the prosthesis and leads to early revision surgery to replace the worn prosthetic components.
- Metallosis and osteolysis may be implicated from generation of wear debris associated with the use of orthopedic implants.
- Peripheral neuropathies have been reported following total joint surgery. Subclinical nerve damage occurs more frequently, possibly the result of surgical trauma.
- Metal sensitivity reactions in patients following joint replacement have been rarely reported. Implantation of foreign material in tissues can result in histological reactions involving macrophages and fibroblasts. The clinical significance of this effect is uncertain, as similar changes may occur as a precursor to, or during the healing process. In some cases, wear debris can initiate the process of histiocytic granuloma formation and consequent osteolysis and loosening of the implant.
- Dislocation and subluxation of implant components can result from improper positioning of the components. Muscle and fibrous tissue laxity can also contribute to these conditions.
- Implants can loosen, migrate, or disassociate due to trauma or loss of fixation.
- Infection can lead to failure of the joint replacement.
- While rare, fatigue fracture of the implant can occur as a result of strenuous activity, improper alignment, or duration of service.
- Fracture of the humerus can occur while press-fitting (seating) the humeral stem into the prepared humeral canal.
- Allergic reactions.
- Soft tissue impingement.
- Scapular notching.

Intraoperative and early postoperative complications can include:

- humeral perforation, or fracture;
- humeral fracture can occur while seating the device;
- damage to blood vessels;
- temporary or permanent nerve damage resulting in pain or numbness of the affected limb;
- undesirable shortening or lengthening of the limb;
- improper implant fit; joint too loose or too tight;
- traumatic arthrosis of the shoulder from intraoperative positioning of the extremity;
- cardiovascular disorders including venous thrombosis, pulmonary embolism, or myocardial infarction;
- hematoma;
- delayed wound healing; and
- infection.

Late postoperative complications can include:

- avulsion as a result of excess muscular weakening;
- non-union due to inadequate reattachment and/or early weight bearing;
- aggravated problems of other joints of the affected limb or muscle deficiencies;
- humeral fracture by trauma or excessive loading, particularly in the presence of poor bone stock;
- fatigue fracture of the acromion or the scapular spine;
- periarticular calcification or ossification, with or without impediment to joint mobility; and
- inadequate or limited range of motion due to improper selection or positioning of components, by impingement, and calcification;
- Lack of load transfer.

Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the relevant Competent Authority as defined in EU 2017/745.

### 14. Lifetime of Device

While the expected lifetime of a reverse total shoulder replacement may be subjected to a multitude of factors, such as patient characteristics, surgeon experience, and surgical technique; evaluating bench data, published literature, and post market surveillance data of the devices under review can give an indication of the expected lifetime of the device. Based on the in-vitro longevity testing, published literature, and registry data, it can be assumed that the estimated lifetime at 10 years for the RSP® Monoblock and Altivate Reverse® is 90%. As product performance data continues to be collected through post-market surveillance activities and incorporated into the clinical evaluation report, the lifetime estimate will be re-evaluated and adjusted if required.

### 15. Sterilization

Unless opened or damaged, DJO Surgical® implants are supplied sterile in multiple pouches or barrier blister trays. Upon receipt, check all packaging for punctures or other damage. If during inspection, packaging is found opened or damaged, contact manufacturer or manufacturer's representative for instructions.

Sterilization of implants, other than Humeral Socket Inserts manufactured from Highly Cross-Linked Polyethylene with Vitamin E (e+™), is by gamma radiation at the minimum dose of 25kGy to achieve a Sterility Assurance Level (SAL) of 10<sup>-6</sup>.

Sterilization of Humeral Socket Inserts manufactured from Highly Cross-Linked Polyethylene with Vitamin E (e+™), is performed by hydrogen peroxide gas plasma to achieve a Sterility Assurance Level (SAL) of 10<sup>-6</sup>.

Implants are single-use devices.

Trials and other instruments are used to determine sizing before the sterile package needs to be opened. Should the integrity of the original sterile package be lost by being opened, punctured, or torn before implantation in the surgical field, contact manufacturer or manufacturer's representative for instructions. These products are single-use devices and CANNOT be resterilized by a healthcare facility. Contact manufacturer or manufacturer's representative for instructions.

Do not resterilize an implant or component that has been opened outside of the surgical field or in contact with or contaminated by blood or other substances. Do not try to clean an implant since standard procedures cannot be relied upon to remove contamination from the implant or component and storage of the opened implant or component should be avoided.

Instruments are provided nonsterile and should be stored in their original packaging until cleaned and sterilized according to the recommended guidelines found in the DJO Surgical® Instrumentation Instructions for Use.

**WARNING: DO NOT resterilize any sterile product distributed by DJO Surgical® (Encore Medical, L.P.) if sterile packaging is opened or damaged upon receipt. Return the implant with respective packaging to DJO Surgical® for inspection and disposition.**

**WARNING: Protect all porous coated, polished surfaces. Standard cleaning procedures cannot be relied upon to remove contamination from porous coating.**

**WARNING: DO NOT re sterilize UHMWPE (ultra-high molecular weight polyethylene), Highly Cross-Linked Polyethylene Vitamin E (e+™), or HA (hydroxyapatite) coated implants.**

DJO Surgical® has validated sterilization cycle data on file.

**NOTE: DJO Surgical® does not recommend Flash or Chemical Sterilization.**

For further information regarding the use of DJO Surgical® Reverse® Shoulder Systems, contact your DJO Surgical® representative or distributor.

DJO Surgical® Reverse® Shoulder Systems are manufactured by ENCORE MEDICAL, L.P.

9800 Metric Blvd., Austin, TX 78758 USA (Made in the USA)

**Warning:** Protect all porous coated and polished surfaces. Standard cleaning procedures cannot be relied upon to remove contamination from porous coating.

**16. \*Trademarks, Patents, and other Intellectual Property**

DJO Surgical® is a registered trademark of Encore Medical L.P. in Australia, Canada, China, European Union, United Kingdom, India, Japan, South Korea, and the United States.

DJO® is a registered trademark of DJO, LLC in Chile, China, and the United States.

AltiVate® is a registered trademark of Encore Medical L.P. in Australia, China, European Union, United Kingdom, Japan, and the United States.

e+™ and P2™ are trademarks of Encore Medical L.P. in the United States.

Reverse® is a registered trademark of Encore Medical L.P. in Germany and the United States.

RSP® is a registered trademark of Encore Medical L.P. in the United States.

Turon® is a registered trademark of Encore Medical L.P. in the United States.

DJO® reserves all rights in its intellectual property. For more information, please refer to our website.

Patented: USPNs 6,679,916; 6,790,234

© 2021 Encore Medical, L.P. and DJO, LLC All Rights Reserved.

**17. Summary of Safety and Clinical Performance (SSCP)**

The SSCP for healthcare professionals and lay persons is available in the European database on medical devices (EUDAMED; <https://ec.europa.eu/tools/eudamed>).

**18. Implant Material Composition**

Component	Composition	
<b>RSP® Monoblock Humeral Stem (Primary and Revision)</b>	Nitrogen, max: 0.05% Carbon, max: 0.08% Hydrogen, max: 0.015% Iron, max: 0.30%	Oxygen, max: 0.20% Aluminum: 5.5-6.75% Vanadium: 3.5-4.5% Yttrium, max: 0.005% Titanium: Balance
	Coating: Nitrogen, max: 0.05% Carbon, max: 0.08% Hydrogen, max: 0.015%	Iron, max: 0.50% Oxygen, max: 0.40% Titanium: Balance
<b>AltiVate Reverse® Standard Shell Humeral Stem (Primary and Revision)</b>	Nitrogen, max: 0.05% Carbon, max: 0.08% Hydrogen, max: 0.015% Iron, max: 0.30%	Oxygen, max: 0.20% Aluminum: 5.5-6.75% Vanadium: 3.5-4.5% Yttrium, max: 0.005% Titanium: Balance
	Porous Coating: Nitrogen, max: 0.05% Carbon, max: 0.08% Hydrogen, max: 0.015%	Iron, max: 0.50% Oxygen, max: 0.40% Titanium: Balance
<b>AltiVate Reverse® Small Shell Humeral Stem (Primary and Revision)</b>	Nitrogen, max: 0.05% Carbon, max: 0.08% Hydrogen, max: 0.015% Iron, max: 0.30%	Oxygen, max: 0.20% Aluminum: 5.5-6.75% Vanadium: 3.5-4.5% Yttrium, max: 0.005% Titanium: Balance
	Porous coating: Nitrogen, max: 0.05% Carbon, max: 0.08% Hydrogen, max: 0.015%	Iron, max: 0.50% Oxygen, max: 0.40% Titanium: Balance
<b>RSP® Humeral Socket Insert (Standard and Semi-Constrained)</b>	Ash, max: 125mg/kg Extraneous matter, No. Particles, max: 3 Titanium, max: 40mg/kg	Aluminum, max: 20mg/kg Calcium, max: 5mg/kg Chlorine, max: 30mg/kg
	Ultra-high molecular weight polyethylene  Ash, max: 125mg/kg Extraneous matter, No. Particles, max: 3 Titanium, max: 40mg/kg	Aluminum, max: 20mg/kg Calcium, max: 5mg/kg Chlorine, max: 30mg/kg Alpha-tocopherol: addition of 0.1% by weight

Component	Composition	
<b>AltiVate Reverse® Small Socket Insert (Neutral and Semi-Constrained)</b>	Ultra-high molecular weight polyethylene Ash, max: 125mg/kg Extraneous matter, No. Particles, max: 3 Titanium, max: 40mg/kg Aluminum, max: 20mg/kg Calcium, max: 5mg/kg Chlorine, max: 30mg/kg	
<b>AltiVate Reverse® Small Socket Insert e+™ (Neutral and Semi-Constrained)</b>	Ultra-high molecular weight polyethylene Ash, max: 125mg/kg Extraneous matter, No. Particles, max: 3 Titanium, max: 40mg/kg Aluminum, max: 20mg/kg Calcium, max: 5mg/kg Chlorine, max: 30mg/kg Alpha-tocopherol: addition of 0.1% by weight	
<b>RSP® Monoblock Hemi-Adapter with Retaining Screw</b>	Adapter: Carbon: 0.00-0.35% Aluminum: 0-1.0% Lanthanum: 0-0.20% Chromium: 26.0-30.0% Molybdenum: 5.0-7.0% Nickel, max: 1.0% Iron, max: 0.75% Silicon, max: 1.0% Manganese, max: 1.0% Nitrogen, max: 1.0% Cobalt: 56.70-69.00% Retaining Screw: Nitrogen, max: 0.05% Carbon, max: 0.08% Hydrogen, max: 0.015% Iron, max: 0.30% Oxygen, max: 0.20% Aluminum: 5.5-6.75% Vanadium: 3.5-4.5% Yttrium, max: 0.005% Titanium: Balance	
<b>RSP® Monoblock Hemi-Adapter Retaining Screw</b>	Nitrogen, max: 0.05% Carbon, max: 0.08% Hydrogen, max: 0.015% Iron, max: 0.30% Oxygen, max: 0.20% Aluminum: 5.5-6.75% Vanadium: 3.5-4.5% Yttrium, max: 0.005% Titanium: Balance	
<b>RSP® Monoblock Spacer with Retaining Screw</b>	Nitrogen, max: 0.05% Carbon, max: 0.08% Hydrogen, max: 0.015% Iron, max: 0.30% Oxygen, max: 0.20% Aluminum: 5.5-6.75% Vanadium: 3.5-4.5% Yttrium, max: 0.005% Titanium: Balance	
<b>RSP® Monoblock Spacer Retaining Screw</b>	Nitrogen, max: 0.05% Carbon, max: 0.08% Hydrogen, max: 0.015% Iron, max: 0.30% Oxygen, max: 0.20% Aluminum: 5.5-6.75% Vanadium: 3.5-4.5% Yttrium, max: 0.005% Titanium: Balance	
<b>AltiVate Reverse® Small Spacer with Retaining Screw</b>	Nitrogen, max: 0.05% Carbon, max: 0.08% Hydrogen, max: 0.015% Iron, max: 0.30% Oxygen, max: 0.20% Aluminum: 5.5-6.75% Vanadium: 3.5-4.5% Yttrium, max: 0.005% Titanium: Balance	
<b>AltiVate Reverse® Small Hemi-Adapter with Retaining Screw</b>	Adapter: Carbon: 0.00-0.35% Aluminum: 0-1.0% Lanthanum: 0-0.20% Chromium: 26.0-30.0% Molybdenum: 5.0-7.0% Nickel, max: 1.0% Iron, max: 0.75% Silicon, max: 1.0% Manganese, max: 1.0% Nitrogen, max: 1.0% Cobalt: 56.70-69.00% Screw: Nitrogen, max: 0.05% Carbon, max: 0.08% Hydrogen, max: 0.015% Oxygen, max: 0.20% Aluminum: 5.5-6.75% Vanadium: 3.5-4.5%	

Component	Composition	
	Iron, max: 0.30%	Yttrium, max: 0.005% Titanium: Balance
AltiVate® Anatomic to Reverse® Conversion Module	Nitrogen, max: 0.05% Carbon, max: 0.08% Hydrogen, max: 0.015% Iron, max: 0.30%	Oxygen, max: 0.20% Aluminum: 5.5-6.75% Vanadium: 3.5-4.5% Yttrium, max: 0.005% Titanium: Balance
Turon® Shoulder Humeral Head (Neutral and Offset)	Carbon: 0.00-0.35% Aluminum: 0-1.0% Lanthanum: 0-0.20% Chromium: 26.0-30.0% Molybdenum: 5.0-7.0%	Nickel, max: 1.0% Iron, max; 0.75% Silicon, max: 1.0% Manganese, max: 1.0% Nitrogen, max: 1.0% Cobalt: 56.70-69.00%
AltiVate® Anatomic Humeral Head (Neutral and Offset)	Carbon: 0.00-0.35% Aluminum: 0-1.0% Lanthanum: 0-0.20% Chromium: 26.0-30.0% Molybdenum: 5.0-7.0%	Nickel, max: 1.0% Iron, max; 0.75% Silicon, max: 1.0% Manganese, max: 1.0% Nitrogen, max: 1.0% Cobalt: 56.70-69.00%
RSP® Glenoid Head with Retaining Screw	Head: Carbon: 0.00-0.35% Aluminum: 0-1.0% Lanthanum: 0-0.20% Chromium: 26.0-30.0% Molybdenum: 5.0-7.0%  Screw: Nitrogen, max: 0.05% Carbon, max: 0.08% Hydrogen, max: 0.015% Iron, max: 0.30%	Nickel, max: 1.0% Iron, max; 0.75% Silicon, max: 1.0% Manganese, max: 1.0% Nitrogen, max: 1.0% Cobalt: 56.70-69.00%  Oxygen, max: 0.20% Aluminum: 5.5-6.75% Vanadium: 3.5-4.5% Yttrium, max: 0.005% Titanium: Balance
RSP® Glenoid Head Retaining Screw	Screw: Nitrogen, max: 0.05% Carbon, max: 0.08% Hydrogen, max: 0.015% Iron, max: 0.30%	Oxygen, max: 0.20% Aluminum: 5.5-6.75% Vanadium: 3.5-4.5% Yttrium, max: 0.005% Titanium: Balance
RSP® Glenoid Baseplate with HA Coating	Nitrogen, max: 0.05% Carbon, max: 0.08% Hydrogen, max: 0.015% Iron, max: 0.30%  Cp Ti Coating: Nitrogen, max: 0.05% Carbon, max: 0.08% Hydrogen, max: 0.015%  HA Coating: Arsenic, max: 3ppm Cadmium, max: 5ppm Mercury, max: 5ppm Lead, max: 30ppm	Oxygen, max: 0.20% Aluminum: 5.5-6.75% Vanadium: 3.5-4.5% Yttrium, max: 0.005% Titanium: Balance  Iron, max: 0.50% Oxygen, max: 0.40% Titanium: Balance
RSP® Glenoid Baseplate Screws (Locking and Non-Locking)	Nitrogen, max: 0.05%	Oxygen, max: 0.20%

Component	Composition	
	Carbon, max: 0.08%	Aluminum: 5.5-6.75%
	Hydrogen, max: 0.015%	Vanadium: 3.5-4.5%
	Iron, max: 0.30%	Yttrium, max: 0.005%
		Titanium: Balance


**19. Implant Card Instructions** (Applicable where implant card is available)

Implant cards are supplied as one foldable implant card and one sticker for each implant component. Each foldable implant card contains one fillable side and five blank sides. Reference the image below for each fillable field of the fillable side of the implant card. Text written in fillable fields of the implant card must be legible and at least 2mm high.

- A. Patient Name or Patient ID (to be filled by healthcare institution/provider)
- B. Date of Implantation (to be filled by healthcare institution/provider)
- C. Name and Address of Healthcare Institution (to be filled by healthcare institution/provider)


Place an implant card component sticker for each implant component used (except for screws) in Sections D of the foldable implant card. Each foldable implant card contains five Sections D.

Foldable Implant Card:

**Implant Card Component Sticker:**

P/N 509-02-432 Implant Component Sticker Shown. A sticker will be provided for each Implant.

If more than five components are used, another implant card must be completed for the remaining components. The fillable fields A-C must be completed for each foldable implant card used.

**Icon Key:**

Instructions: The table below is intended as a definition library of all symbols potentially applicable to DJO® Devices; not all symbols in this key are applicable to the devices covered by these Instructions for Use. Please refer to the product package labeling for specific symbols that are applicable to the devices covered by these Instructions for Use.

 ISO 15223-1 5.4.2	<b>Single use – do not reuse</b>
 ISO 15223-1 5.1.4	<b>Expiration Date</b>

 ISO 15223-1 5.3.4	<p style="text-align: center;"><i>Keep Dry</i></p>
 ISO 15223-1 5.1.5	<p style="text-align: center;"><i>Lot number/Batch Code</i></p>
 ISO 15223-1 5.2.1	<p style="text-align: center;"><i>Sterile</i></p>
 ISO 15223-1 5.2.4	<p style="text-align: center;"><i>Sterility symbol: R: Sterile Using Irradiation</i></p>
 	<p style="text-align: center;"><i>Sterile symbol: H<sub>2</sub>O<sub>2</sub>: Sterilized Using Hydrogen Peroxide Gas Plasma</i></p>
	<p style="text-align: center;"><i>Double Sterile Barrier</i></p>
 ISO 15223-1 5.2.7	<p style="text-align: center;"><i>Non-sterile</i></p>
	<p style="text-align: center;"><i>See "Instructions for Use"</i></p>
 ISO 15223-1 5.1.1	<p style="text-align: center;"><i>Manufacturer</i></p>
	<p style="text-align: center;"><i>Quantity of items in package</i></p>

 ISO 15223-1 5.1.2	Authorized Representative in European Community
 ISO 15223-1 5.1.6	Catalog Number
 ISO 15223-1 5.2.6	Do not resterilize
 ISO 15223-1 5.2.8	Do not use if package is damaged
 ASTM F2503:2013	MR Safe
 ASTM F2503:2013	MR Conditional
 ASTM F2503:2013	MRI Unsafe
Rx 21 CFR 801.109	Federal Law (USA) restricts this device to sale by or on the order of a physician.
	Importer
	Medical Device

	<p><i>Contains substance (cobalt) that can be carcinogenic, mutagenic or reprotoxic</i></p>
	<p><i>Country Code of Manufacturer – US</i></p>
	<p><i>Patient Name or Patient ID</i></p>
	<p><i>Date of Implantation</i></p>
	<p><i>Name and Address of the Implanting Healthcare Institution/Provider</i></p>
	<p><i>Information Website for Patients</i></p>

Bone Cement Usage – The following legends are displayed on the product labeling to indicate bone cement usage:

Usage	Legend
Implants intended to be used with bone cement	CEMENTED
Implants intended to be used without bone cement	CEMENTLESS
Implants intended to be used optionally	NO LEGEND

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