DJO Surgical Hip Systems
 DJO Surgical Hüftsysteme
 Systèmes de hanche DJO Surgical
 Sistemas de cadera DJO Surgical
 Sistemi per anca DJO Surgical
 Συστήματα ισχίου DJO Surgical
 DJO Surgical Kalça Sistemleri

Encore® Medical, L.P.
 9800 Metric Blvd.
 Austin, TX 78758-5445 USA

EC REP MDSS GmbH
 Schiffgraben 41
 Hanover, Germany

Distributed in Turkey by:
EndoProtez Ortopedi Ltd.Şti.
 Halk Sokak ¼ Sıhhiye Ankara Türkiye
 Telefon + 90 (312) 435 74 46

0400-0104 REV. U 2017-07

EN

1. Product Handling

Implants are provided sterile and should always be stored unopened in their respective protective containers. Prior to use, inspect package for damage that may compromise sterility. If packaging has been opened or damaged upon receipt, contact the manufacturer's representative. Also inspect the labeling to verify that the expiration date has not passed. If the product is expired, contact Customer Service and do not use the implant. When unpacking the implant, verify the labeling for correct Reference No. and size. When removing the implant from its packaging, the relevant aseptic instructions must be observed. Protect prosthesis from contact with objects that may damage the surface finish. Inspect each implant prior to use for visual damage. This implant is part of a system and should be used only in combination with other original DJO Surgical product belonging to the same hip system, unless otherwise specified.

2. Product Description and Implant Materials

Hip Stem	Fixation Method	Material	Applicable ASTM Standard	Applicable ISO Standard
Linear®	Cementless	Ti-6Al-4V alloy CP Ti Porous Coating	ASTM F136 / ASTM F1472 ASTM F67	ISO 5832-3 ISO 5832-2
Foundation® Press-Fit	Cementless	Ti-6Al-4V alloy CP Ti Porous Coating	ASTM F136 / ASTM F1472 ASTM F67	ISO 5832-3 ISO 5832-2
Foundation® Modular Hip	Cementless	Ti-6Al-4V alloy CP Ti Porous Coating Hydroxyapatite Coating	ASTM F136 / ASTM F1472 ASTM F67 ASTM F1185	ISO 5832-3 ISO 5832-2 ISO 13779-1
Revelation®	Cementless	Ti6Al4V alloy CP Ti Porous Coating	ASTM F136 / ASTM F1472 ASTM F67	ISO 5832-3 ISO 5832-2
Revelation® V2	Cementless	Ti6Al4V alloy CP Ti Porous Coating	ASTM F136 / ASTM F1472 ASTM F67	ISO 5832-3 ISO 5832-2
Revelation® Short Stem (Standard & V2)	Cementless	Ti6Al4V alloy CP Ti Porous Coating	ASTM F136 / ASTM F1472 ASTM F67	ISO 5832-3 ISO 5832-2
Vitality®	Cemented	CoCr	ASTM F799	ISO 5832-4
Foundation® Cemented	Cemented	CoCr	ASTM F799	ISO 5832-4
Foundation® Fracture	Cementless	Ti6Al4V alloy	ASTM F136 / ASTM F1472	ISO 5832-3
Stamina®	Cementless	Ti6Al4V alloy	ASTM F136	ISO 5832-3
CLP™ I CLP™ Offset	Cementless	Ti6Al7Nb Niobium Alloy	ASTM F1295	ISO 5832-11
R120™ / R120PC™	Cemented or Cementless	Cast CoCr	ASTM F75	ISO 5832-4
ALFA™ II	Cemented or Cementless	Cast CoCr	ASTM F75	ISO 5832-4
CLP-R™	Cemented or Cementless	Ti6Al7Nb Niobium Alloy	ASTM F1295	ISO 5832-11
Expr® Precision System: Revision Hip	Cementless	Ti6Al4V alloy	ASTM F1472	ISO 5832-3

DJO Surgical hip stems can be used with any DJO Surgical femoral heads for total joint replacement. DJO Surgical Hip Systems are for total hip replacement except for Bipolar and Unipolar which are for hemi arthroplasty applications. The CLP Offset, CLP-R, Vitality and Stamina hip systems are for either total or hemi applications.

Component	Fixation Method	Material	Applicable ASTM Standard	Applicable ISO Standard
All Poly Acetabular Cup	Cemented	Medical grade UltraHigh Molecular Weight Polyethylene Ti6Al4V alloy PMMA	ASTM F648 ASTM F67	ISO 5834-1 ISO 5832-2
Femoral Heads	Cementless Cementless ¹ Cementless ²	CoCr Ceramic Al ₂ O ₃ ¹ BioloX® <i>delta</i> Ceramic ²	ASTM F799	ISO 5832-4 ISO 6474 ISO 6474
BioloX® <i>delta</i> Option Titanium Sleeve	Cementless ²	Ti6Al4V alloy	ASTM F136 / ASTM F1472	ISO 5832-3
Modular Femoral Neck	Cemented	CoCrMo	ASTM F1537	ISO 5832-4
Acetabular Shells	Cementless	Ti6Al4V alloy CP Ti porous coating	ASTM F136 / ASTM F1472 ASTM F67	ISO 5832-3 ISO 5832-2
Polyethylene Acetabular Liners	Cementless	Medical grade UltraHigh Molecular Weight Polyethylene	ASTM F648	ISO 5834-1 & 2
Constrained Acetabular Liners	Cemented or Cementless	Medical grade UltraHigh Molecular Weight Polyethylene Ti6Al4V alloy	ASTM F648 ASTM F136 / ASTM F1472	ISO 5834-1 & 2 ISO 5832-3
X-alt™ Highly Cross-Linked Polyethylene Acetabular Liner	Cementless	Medical grade UltraHigh Molecular Weight Polyethylene (Highly Cross-Linked)	ASTM F648 ASTM F2565	ISO 5834-1 & 2
Metal/Metal Acetabular Liners	Cementless	Medical grade UltraHigh Molecular Weight Polyethylene CoCr	ASTM F648 ASTM F799	ISO 5834-1 & 2 ISO 5832-4
Bipolar/Unipolar/Sleeves	Cemented or Cementless	CoCr Ti6Al4V alloy Medical grade UltraHigh Molecular Weight Polyethylene	ASTM F799 / F-1537-07 ASTM F136 / ASTM F1472 ASTM F648	ISO 5832-4 ISO 5832-3 ISO 5834-1 & 2
Distal Centralizers	N/A	PMMA		
Screws	N/A	Ti6Al4V alloy	ASTM F136 / ASTM F1472	ISO 5832-3

¹For Ceramic-on-Ceramic applications please refer to the Keramos® Total Hip System IFU. NOTE: Keramos® is an Encore Medical Registered Trademark for CeramTec's BioloX forte™ Alumina Ceramic femoral heads.

²BioloX® *delta* and BioloX® *delta* Option Ceramic Femoral Heads are not approved for use with the Ceramic-on-Ceramic applications.

3. Indications

Joint replacement is indicated for patients suffering from disability due to:

- noninflammatory degenerative joint disease including osteoarthritis and avascular necrosis of the natural femoral head;
- rheumatoid arthritis;
- correction of functional deformity;
- femoral fracture.

This device may also be indicated in the salvage of previously failed surgical attempts.

The constrained acetabular component is indicated for primary or revision patients at high risk of hip dislocation due to a history of prior dislocation, bone loss, soft tissue laxity, neuromuscular disease, or intra-operative instability and for who all other options to constrained acetabular components have been considered.

Expri™ Revision Hip System is indicated for patients whose bone stock is of poor quality or inadequate for other reconstruction techniques as indicated by deficiencies of the femoral head, neck or portions of the proximal femur. It is intended for cementless revision hip arthroplasty on both uncemented and cemented femoral implants.

4. Intended Use

DJO Surgical hip devices are intended for treatment of patients who are candidates for total hip arthroplasty per the indications for use. While hip replacements are not intended to withstand activity levels and loads of normal healthy bone, they are a means of restoring mobility and reducing pain for many patients.

^{*}For Ceramic-on-Ceramic applications, please refer to the Keramos® Total Hip System IFU.

5. Contraindications*

Joint replacement is contraindicated where there is:

- infection or sepsis;
- insufficient bone quality which may affect the stability of the implant;
- muscular, neurological or vascular deficiencies, which compromise the affected extremity;
- skeletally immature patients and cases where there is a loss of abductor musculature, poor bone stock, poor skin coverage around hip joint which would make the procedure unjustifiable;
- osteomyelitis;
- rapid joint destruction or bone absorption apparent on roentgenogram;
- pathological conditions of the acetabulum, which would prevent achieving proper range of motion, appropriate head stability, and/or a well-seated and supported smooth articulation of the head within the acetabulum;
- alcoholism or other addictions;
- materials sensitivity;
- loss of ligamentous structures;
- high levels of physical activity (e.g. competitive sports, heavy physical labor);
- pregnancy (contraindicated for Metal on Metal applications only).

^{*}For Ceramic-on-Ceramic applications, please refer to the Keramos® Total Hip System IFU.

6. Precautions and Warnings*

An implant should never be reused. Although the implant may appear undamaged, previous stresses could create imperfections that may lead to mechanical failure. It is advised to utilize new prostheses of current design.

Familiarity with, and attention to the surgical technique recommended for this device is imperative for best results. The correct selection as well as the correct seating/placement of the prosthetic implant is extremely important. Malposition may predispose the device to

excess wear and early failure. Use of the largest stem possible is recommended. Only DJO Surgical Hip System implants, instruments, and trial prostheses should be used.

Care must be taken to protect mating surfaces (i.e. tapers) and polished bearing surfaces from nicks and scratches which could become the focal point for failure. Contouring or bending of the implant may reduce its service life and may cause immediate or eventual failure under load. An implant must not be tampered with, as tampering will adversely affect the performance of the implant. DJO Surgical (Hip, Knee, and Shoulder) systems have not been evaluated for safety and compatibility in the Magnetic resonance environment. The (Hip, Knee, shoulder) stems have not been tested for heating or migration in the Magnetic resonance environment.

1. Do not implant HA (Hydroxyapatite) coated implants with bone cement.
2. To determine the use of the hip stem with the correct femoral head (CoCr or Ceramic): please refer to Section 2 "Product Description and Implant Materials".
3. Ceramic femoral heads are only available in sizes 28 and 32 and are to be used with hip stems having a Morse Type Taper. See Section 2 "Product Description and Implant Materials" for a listing of hip stems that are cleared for use with the ceramic femoral head.
4. Ceramic femoral heads are only indicated for use with stems during total hip replacement.

Precautions and Warnings Specific to the Constrained Acetabular Liner

Precautions:

1. In order to minimize the risks of dislocation and loosening of the shell-acetabular bone or shell-bone cement interface that may occur when using a metallic shell intended for biological fixation or cemented use only, surgeons should consider providing immediate resistance to tensile forces between the metallic shell and the acetabular bone or bone cement interface through the use of orthopedic bone fixation devices such as bone screws, spikes, screw threads, fins, or other bone fixation devices.
2. To correctly position the metallic locking ring, surgeons should consult the manufacturer's instructions for appropriate device assembly.
3. Physicians should consider component malposition, component placement, and the effect on range of motion when using modular heads (with sleeves or skirts) and extended liners.
4. Regarding component malposition above, recommendation is to caution physician regarding the malpositioned acetabular components cup and the potential for impingement, premature dislocation and revision.

Warnings

1. Closed reduction of this device is not possible. Patients should be made aware that treatment of device dislocation will require additional surgery.
2. There may be a failure of the retaining ring.
3. A retaining ring that is placed incorrectly may have a reduced life.
4. Retaining ring failure, which may be due to impingement, fatigue, and/or wear, increases the probability of dislocation.
5. Failure or migration of the retaining ring may require additional surgery.

Precautions Specific to the Revelation Short Stem

1. Fracture in smaller sized stems is most likely to occur in patients who are young, physically active and/or obese.

*For Ceramic-on-Ceramic applications, please refer to the Keramos® Total Hip System IFU.

DJO Surgical (Hip, Knee, and Shoulder) systems have not been evaluated for safety and compatibility in the Magnetic resonance environment. The (Hip, Knee, Shoulder) systems have not been tested for heating or migration in the Magnetic resonance environment.

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

7. Preoperative Planning and Postoperative Care

Preoperative planning provides essential information regarding the appropriate prosthesis and likely combinations of components. Use instrument trial components for fit verification (where applicable) and extra implant components for backup. X-ray templates for all sizes of the DJO Surgical Hip system are available upon request.

Accepted surgical practices should be followed for postoperative care. The patient should be made aware of the limitation of total joint reconstruction. Excessive physical activity and trauma affecting the replaced joint have been implicated in premature failure by loosening, fracture, and/or wear of the prosthetic implants. The patient should be cautioned to govern his/her activities accordingly as the risk of implant failure increases with weight and activity levels of the patient.

Patient Counseling Information Specific to the Constrained Acetabular Liner

1. The prosthesis will not restore function to the level expected with a normal healthy joint, and the patient should be instructed as to the limitations of the device.
2. The range of motion achievable with a constrained liner is less than the range of motion of a normal joint, and less than with a semi-constrained prosthesis.
3. The patient should be told that, although the constrained hip liner provides resistance to dislocation, it can dislocate if subjected to excessive loading.
4. Once dislocated, additional surgery will be required to reduce the joint.
5. Patients should be instructed that significant reduction in the range of motion is inherent to the design characteristic of a constrained acetabular liner, and activities that may force the joint to exceed those range of motion limits should be avoided.

8. Adverse Effects*

1. Accelerated wear of the polyethylene articulating surfaces have been reported following total hip replacement. Such wear may be initiated by particles of cement, metal, or other debris which can cause abrasion of the articulating surfaces. Accelerated wear shortens the useful life of the prosthesis, and leads to early revision surgery to replace the worn prosthetic components.
2. Metallosis and osteolysis may be implicated from wear debris associated with the use of orthopedic implants.
3. Peripheral neuropathies have been reported following total joint surgery. Subclinical nerve damage occurs more frequently, possibly the result of surgical trauma.
4. Metal sensitivity reactions in patients following joint replacement have been rarely reported. Implantation of foreign material in tissues can result in histological reactions involving macrophages and fibroblasts. The clinical significance of this effect is uncertain, as similar changes may occur as a precursor to, or during the healing process. In some cases, wear debris can initiate the process of histiocytic granuloma formation and consequent osteolysis and loosening of the implant.
5. Dislocation and subluxation of implant components can result from improper positioning of the components. Muscle and fibrous tissue laxity can also contribute to these conditions.
6. Ring fracture could lead to increased risk of dislocation.
7. Implants can loosen or migrate due to trauma or loss of fixation.
8. Infection can lead to failure of the joint replacement.
9. While rare, fatigue fracture of the implant can occur as a result of strenuous activity, improper alignment, or duration of service.
10. Fracture of the femur can occur while press-fitting (seating) the femoral stem into the prepared femoral canal.
11. Allergic reactions.

Intraoperative and early postoperative complications can include:

1. acetabular perforation, or fracture;
2. femoral fracture can occur while seating the device;
3. damage to blood vessels;
4. temporary or permanent nerve damage resulting in pain or numbness of the affected limb;
5. undesirable shortening or lengthening of the limb;
6. traumatic arthrosis of the hip from intraoperative positioning of the extremity;
7. cardiovascular disorders including venous thrombosis, pulmonary embolism, or myocardial infarction;
8. hematoma;
9. delayed wound healing; and,
10. infection.

Late postoperative complications can include:

1. avulsion as a result of excess muscular weakening;
2. non-union due to inadequate reattachment and/or early weight bearing;
3. aggravated problems of other joints of the affected limb or muscle deficiencies;
4. femoral fracture by trauma or excessive loading, particularly in the presence of poor bone stock;
5. periarticular calcification or ossification, with or without impediment to joint mobility;
6. inadequate range of motion due to improper selection or positioning of components, by impingement, and calcification.

*For Ceramic-on-Ceramic applications, please refer to the Keramos® Total Hip System IFU.

9. Sterilization

Unless opened or damaged, DJO Surgical implants are supplied sterile in multiple pouches or barrier blister trays. Upon receipt, check all packaging for punctures or other damage. If during inspection, packaging is found opened or damaged, contact manufacturer or manufacturer's representative for instructions.

Sterilization of implants other than the Highly Cross-Linked Polyethylene Acetabular Liner is performed by gamma radiation at the minimum dose of 25 kGy to achieve a Sterility Assurance Level (SAL) of 10⁻⁶. Implants are single-use devices. Trials and other instruments are used to determine sizing before the sterile package needs to be opened. Should the integrity of the original sterile package be lost by being opened, punctured, or torn before implantation in the surgical field, the metallic devices may be resterilized prior to use following the guidelines listed below.

Sterilization of the Highly Cross-Linked Polyethylene Acetabular Liner is performed by hydrogen peroxide gas plasma to achieve a Sterility Assurance Level (SAL) of 10⁻⁶. These liners are single-use devices and CANNOT be resterilized by a healthcare facility. Liner trials and other instruments are used to determine sizing before the sterile package needs to be opened.

Do not resterilize an implant or component that has been opened outside of the surgical field or in contact with or contaminated by blood or other substances. Do not try to clean an implant since standard procedures cannot be relied upon to remove contamination from the implant or component and storage of the implant or component should be avoided.

Instruments are provided nonsterile and should be stored in their original packaging until cleaned and sterilized according to the recommended guidelines found in the DJO Surgical Instrumentation Instructions for Use.

User Resterilization Guidelines for Implants:

Porous coated or non-coated metallic implants or components, which have not been previously used or implanted, or contaminated with blood or other bodily substances, can be resterilized using the following steam sterilization cycle. All devices should be disassembled prior to resterilization and care should be taken to protect implant or component from mechanical damage.

Steam Resterilization for metallic implants and components ONLY!	<p>Remove from supplied packaging and wrap in protective sterilization wrap according to AAMI / AORN guidelines or place into appropriate case configuration.</p> <p>Pre-Vacuum Autoclave (HI-VAC): 270-272° F (132-134° C), 6-minute exposure time, with 4 pulses and a 5-minute dry time shall achieve a SAL of 10⁻⁶.</p> <p>Gravity Displacement Autoclave: 270-272° F (132-134° C), 30-minute exposure time, with a 5-minute dry time shall achieve a SAL of 10⁻⁶.</p>
---	--

WARNING: DO NOT resterilize any hip prosthesis distributed by DJO Surgical (Encore Medical, L.P.) if sterile packaging is opened or damaged. Return the implant with respective packaging to DJO Surgical for inspection and disposition.

WARNING: Protect all porous coated and polished surfaces. Standard cleaning procedures cannot be relied upon to remove contamination from the implant or component.

WARNING: DO NOT resterilize UHMWPE (ultra-high molecular weight polyethylene) implants, PMMA (polymethylmethacrylate) spacers, HA (Hydroxyapatite) coated implants, and ceramic implants.

WARNING: Highly Cross-Linked Polyethylene CANNOT be resterilized by a healthcare facility.

DJO Surgical has validated the above steam sterilization cycles and has data on file. Other sterilization cycles may also be suitable; however individuals or hospitals not using the recommended method are advised to validate any alternate method using appropriate laboratory techniques. Proper validation of the autoclave is essential to ensure proper sterilization temperatures and cycle times.

NOTE: DJO Surgical does not recommend Flash or Chemical Sterilization.

For further information regarding the use of the DJO Surgical Hip Systems, contact your DJO Surgical representative or distributor.

DJO Surgical Hip Systems are manufactured by ENCORE MEDICAL, L.P. ■
9800 Metric Blvd., Austin, TX 78758 USA (Made in the USA)

10. Trademarks and patents

FOUNDATION®, LINEAR®, KEYSTONE®, VITALITY®, Exprt® are registered trademarks of Encore Medical, L.P., Austin, TX 78758 USA or its affiliates.

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