
DJO Surgical Knee Systems  
 DJO Surgical Knie-Systeme  
 Systèmes de genou DJO Surgical  
 Sistemas de rodilla DJO Surgical  
 Sistemi protesici per ginocchio DJO Surgical  
 Συστήματα γονάτου DJO Surgical  
 DJO Surgical Diz Sistemleri

Encore Medical, L.P.  
 9800 Metric Blvd.  
 Austin, TX 78758

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EN

**1. Product Handling**

Implants are provided sterile and should always be stored unopened in their respective protective containers. Prior to use, inspect package for damage that may compromise sterility. If packaging has been opened or damaged upon receipt, contact the manufacturer's representative. Also inspect the labeling to verify that the expiration date has not passed. If the product is expired, contact Customer Service and do not use the implant. When unpacking the implant, verify the labeling for correct Reference No. and size. When removing the implant from its packaging, the relevant aseptic instructions must be observed. Protect prosthesis from contact with objects that may damage the surface finish. Inspect each implant prior to use for visual damage. This implant is part of a system and should be used only in combination with other original DJO Surgical product belonging to the same knee system, unless otherwise specified.

**2. Product Description and Implant Materials**

Component	Fixation Method	Material	Applicable ASTM Standard	Applicable ISO Standard
*FOUNDATION® Femur -Primary, Posterior Stabilized, Revision	Cemented	CoCrMo	ASTM F75	ISO 5832-4
*FOUNDATION Femur -Primary, Posterior Stabilized	Cemented	CoCrMo TiNbn	ASTM F75	ISO 5832-4
*FOUNDATION® Tibia -Stemmed, Modular Revision, Resurfacing	Cemented	CoCrMo	ASTM F75	ISO 5832-4
		Ti6Al4V alloy	ASTM F136 / ASTM F1472	ISO 5832-3
*FOUNDATION® Porous Coated Tibia	Cemented	Ti6Al4V alloy	ASTM F136 / ASTM F1472	ISO 5832-3
		CP Ti Porous Coating	ASTM F67	ISO 5832-2
*FOUNDATION® All Poly Tibia	Cemented	Medical grade UltraHigh Molecular Weight Polyethylene	ASTM F648	ISO 5834-1 & 2
*FOUNDATION® Tibial Inserts -Primary, Ultra-congruent, Posterior Stabilized, Lateral Pivot, Constrained	Cemented	Medical grade UltraHigh Molecular Weight Polyethylene	ASTM F648	ISO 5834-1 & 2
*FOUNDATION® Stem Extensions	Cemented	Ti6Al4V alloy	ASTM F136 / ASTM F1472	ISO 5832-3
FK/3DKNEE™ Downsize Tibia	Cemented	Ti6Al4V alloy	ASTM F136 / ASTM F1472	ISO 5832-3
3DKNEE™ Femur	Cemented	CoCrMo TiNbn	ASTM F75	ISO 5832-4
3DKNEE™ Porous Coated Femur	Cemented	CoCrMo	ASTM F75	ISO 5832-4
		CP Ti Porous Coating	ASTM F67	ISO 5832-2
3DKNEE™ Tibia	Cemented	CoCrMo	ASTM F75	ISO 5832-4
		Ti6Al4V alloy	ASTM F136 / ASTM F1472	ISO 5832-3
3DKNEE™ Porous Coated Tibia	Cementless	Ti6Al4V alloy	ASTM F136 / ASTM F1472	ISO 5832-3
		CP Ti Porous Coating	ASTM F67	ISO 5832-2
3DKNEE™ Standard Insert	Cemented	Medical grade UltraHigh Molecular Weight Polyethylene	ASTM F648	ISO 5834-1 & 2
3DKNEE™ Insert - e+ EPIK® Uni-System Inlay - e+	Cemented	Highly Cross-Linked Vitamin E UHMWPE	ASTM F2695 ASTM F2565	

Component	Fixation Method	Material	Applicable ASTM Standard	Applicable ISO Standard
EPIK® Uni-System Onlay - e+				
EMPOWR 3DKNEE™ Femur	Cemented	CoCrMo	ASTM F75	ISO 5832-4
DJO EMPOWR KNEE™ Tibia	Cemented	CoCrMo	ASTM F75	ISO 5832-4
EMPOWR 3DKNEE™ Tibial Insert - e+	Cemented	Highly Cross-Linked Vitamin E UHMWPE	ASTM F2695 ASTM F2565	
EMPOWR PS KNEE™ Femur	Cemented	CoCrMo	ASTM F75	ISO 5832-4
EMPOWR PS KNEE™ Tibial Insert - e+	Cemented	Highly Cross-Linked Vitamin E UHMWPE	ASTM F2695 ASTM F2565	
EMPOWR CR KNEE™ Tibial Insert - e+	Cemented	Highly Cross-Linked Vitamin E UHMWPE	ASTM F2695 ASTM F2565	
Exprt™ Revision Knee Femur	Cemented	CoCrMo	ASTM F75	ISO 5832-4
Exprt™ Revision Knee Tibia	Cemented	CoCrMo	ASTM F75	ISO 5832-4
Exprt™ Revision Knee Insert Standard	Cemented	Medical grade UltraHigh Molecular Weight Polyethylene	ASTM F648	ISO 5834-1 & 2
Exprt™ Revision Knee Insert e+	Cemented	Highly Cross-Linked Vitamin E UHMWPE	ASTM F2695 ASTM F2565	
Exprt™ Revision Knee Stem Extension	Cemented	CoCrMo	ASTM 1537	
Exprt™ Revision Tibial Augment	Cemented	CoCrMo	ASTM F75	ISO 5832-4
Exprt™ Revision Femoral Augment	Cemented	Ti6Al4V alloy	ASTM F136 / ASTM F1472	ISO 5832-3
EPIK® Uni Femur	Cemented	CoCrMo	ASTM F75	ISO 5832-4
EPIK® Uni Tibia	Cemented	CoCrMo	ASTM F75	ISO 5832-4
EPIK® Uni Tibial Insert	Cemented	Medical grade UltraHigh Molecular Weight Polyethylene	ASTM F648	ISO 5834-1 & 2
Patella – All Poly Domed & Sombrero	Cemented	Medical grade UltraHigh Molecular Weight Polyethylene	ASTM F648	ISO 5834-1 & 2
Patella – All Poly e+ Domed & Sombrero	Cemented	Highly Cross-Linked Vitamin E UHMWPE	ASTM F2695 ASTM F2565	
Patella – Metal Backed	Cemented	Medical grade UltraHigh Molecular Weight Polyethylene	ASTM F648	ISO 5834-1 & 2
		Ti6Al4V alloy	ASTM F136 / ASTM F1472	ISO 5832-3
		CP Ti Porous Coating	ASTM F67	ISO 5832-2
Modular Spacer Blocks	Cemented	Ti6Al4V alloy	ASTM F136 / ASTM F1472	ISO 5832-3
		CoCrMo	ASTM F75	ISO 5832-4
Attachment Screws	Cementless	Ti6Al4V alloy	ASTM F136 / ASTM F1472	ISO 5832-3
Baseplate Plugs	Cementless	Medical grade UltraHigh Molecular Weight Polyethylene	ASTM F648	ISO 5834-1 & 2

\*Not Available in the EU

Depending on the DJO Surgical Knee System, femoral prostheses may be available in left and right configurations. Femoral augmentation blocks are available in corresponding distal and posterior sizes and femoral extension stems and tibial extension stems are available in various lengths and diameters. Depending on the DJO Surgical Knee System, the stemmed baseplate is available in left and right configurations. Tibial augmentation blocks are available in corresponding medial and lateral sizes of tibial baseplates.

**Note:** Size interchangeability between the Posterior-Stabilized or Revision femur and Posterior-Stabilized or Constrained Tibial insert is limited to one size larger and one size smaller.

**Note:** Size interchangeability between the FOUNDATION® Primary femur and Lateral pivot tibial inserts is limited to one size larger and one size smaller.

**Note:** Size interchangeability between the FOUNDATION® femoral components and Ultra-congruent tibial inserts is limited to one size larger and one size smaller.

**Note:** 3DKNEE™ femurs are to be used with 3DKNEE™ tibial inserts and FK/3DKNEE™, FOUNDATION® or 3DKNEE™ tibial baseplates and patellae. If using FK/3DKNEE™ or FOUNDATION® tibial baseplates with 3DKNEE™ femurs, use baseplates that are one size smaller.

**Note:** Size 3, 5, and 7 3DKNEE™ femurs and inserts are to be used with 3DKNEE™ Tibial Components that are one size larger or FOUNDATION® or FK/3DKNEE™ downsize components that are one size smaller.

**Note:** The 3DKNEE™ femur and the 3DKNEE™ insert size must match. The insert will mate with the same size 3DKNEE™ baseplate or one size smaller FOUNDATION® or FK/3DKNEE™ baseplate.

**Note:** The 3DKNEE™ tibial insert - e+ is compatible with the same sized 3DKNEE™ femur.

**Note:** Size interchangeability between the EMPOWR 3DKNEE™ femur and EMPOWR 3DKNEE™ tibial inserts is compatible with the same size component, or one size smaller femur on a larger insert.

**Note:** Size interchangeability between the EMPOWR PS KNEE™ femur and EMPOWR PS KNEE™ tibial inserts is compatible with the same size component, one size smaller femur on a larger insert, or one size larger femur on a smaller insert.

**Note:** Size interchangeability between the EMPOWR 3D KNEE™ femur and EMPOWR CR KNEE™ tibial inserts is compatible with the same size component, up to two sizes smaller femur on a larger insert, or up to two sizes larger femur on a smaller insert.

### 3. Indications

Joint replacement is indicated for patients suffering from disability due to:

- degenerative, post-traumatic or rheumatoid arthritis;
- avascular necrosis of the femoral condyle;
- post-traumatic loss of joint configuration, particularly when there is patellofemoral erosion, dysfunction or prior patellectomy;
- moderate valgus, varus or flexion deformities;
- treatment of fractures that are unmanageable using other techniques.

This device may also be indicated in the salvage of previously failed surgical attempts. All devices are intended for cemented applications except for the 3DKNEE™ Porous Coated Femur which is intended for cementless applications.

While knee replacements are not intended to withstand activity levels and loads of normal healthy bone, they are a means of restoring mobility and reducing pain for many patients.

#### Indications for EPIK® Uni Knee:

Intended for partial replacement of the articulating surfaces of the knee when only one side of the joint is affected due to compartmental primary degenerative or post-traumatic degenerative disease, previous tibial condyle or plateau fractures, deformity or revision of previous arthroplasty. The device is a single use implant intended for implantation with bone cement.

### 4. Intended Use

DJO Surgical knee devices are intended for treatment of patients who are candidates for knee arthroplasty per the Indication for use. While total knee replacements are not intended to withstand activity levels and loads of normal healthy bone, they are a means of restoring mobility and reducing pain for many patients.

### 5. Contraindications

Joint replacement is contraindicated where there is:

- infection (or a history of infection), acute or chronic, local or systemic;
- insufficient bone quality which may affect the stability of the implant;
- muscular, neurological or vascular deficiencies, which compromise the affected extremity;
- obesity;
- alcoholism or other addictions;
- materials sensitivity;
- loss of ligamentous structures;
- high levels of physical activity (e.g. competitive sports, heavy physical labor).
- The EMPOWR 3DKNEE™ and EMPOWR CR KNEE™ are also contraindicated for patients without sufficient soft tissue integrity to provide adequate stability.

### 6. Precautions and Warnings

An implant should never be reused. Although the implant may appear undamaged, previous stresses could create imperfections that may lead to mechanical failure. It is advised to utilize new prostheses of current design.

Familiarity with, and attention to the surgical technique recommended for this device is imperative for best results. The correct selection as well as the correct seating/placement of the prosthetic implant is extremely important. Only DJO Surgical Knee System instruments and trial prostheses should be used.

Care must be taken to protect mating surfaces (i.e. tapers) and polished bearing surfaces from nicks and scratches which could become the focal point for failure. Contouring or bending of the implant may reduce its service life and may cause immediate or eventual failure under load. An implant must not be tampered with, as tampering will adversely affect the performance of the implant.

DJO Surgical (Hip, Knee, and Shoulder) systems have not been evaluated for safety and compatibility in the Magnetic resonance environment. The (Hip, Knee, Shoulder) systems have not been tested for heating or migration in the Magnetic resonance environment.

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

### 7. Preoperative Planning and Postoperative Care

Preoperative planning provides essential information regarding the appropriate prosthesis and likely combinations of components. Use instrument trial components for fit verification (where applicable) and extra implant components for backup. X-ray templates for all sizes of the DJO Surgical Knee systems are available upon request.

Accepted surgical practices should be followed for postoperative care. The patient should be made aware of the limitation of total joint reconstruction. Excessive physical activity and trauma affecting the replaced joint have been implicated in premature failure by loosening, fracture, and/or wear of the prosthetic implants. The patient should be cautioned to govern his activities accordingly as the risk of implant failure increases with weight and activity levels of the patient.

### 8. Adverse Effects

Some of the adverse effects that could occur related to total knee arthroplasty are:

- fracture of the tibia or femur;
- transient peroneal palsy secondary to surgical manipulation;
- patellar subluxation or dislocation;
- patellofemoral impingement;
- instability, changes in position, or loosening of components;
- ligamentous laxity;
- dissociation of components;
- infection;
- poor range of motion;
- shortening of limbs;
- metal sensitivity reactions.

### 9. Sterilization

Unless opened or damaged, DJO Surgical implants are supplied sterile in multiple pouches or barrier blister trays. Upon receipt, check all packaging for punctures or other damage. If packaging is opened or damaged, contact the manufacturer or manufacturer's representative for instructions.

Sterilization of implants other than the 3DKNEE™ Tibial Inserts - e+, Exprt™ Revision Knee e+ Tibial Inserts, EPIK® Onlay - e+, EPIK® Inlay - e+, and e+ Patellas are performed by gamma radiation at the minimum dose of 25 kGy to achieve a Sterility Assurance Level (SAL) of 10<sup>-6</sup>. Implants are single-use devices. Trials and other instruments are used to determine sizing before the sterile package needs to be opened. Should the original sterile package be inadvertently opened or compromised before implantation, the device cannot be implanted. Contact manufacturer or manufacturer's representative for instructions. Do not resterilize an implant or component that has been in contact with or contaminated by blood or other substances. Do not try to clean an implant since standard procedures cannot be relied upon to remove contamination from porous coating and storage of any opened implant or component should be avoided.

Sterilization of the 3DKNEE™ Tibial Inserts - e+, Exprt™ Revision Knee e+ Tibial Inserts, EPIK® Onlay - e+, EPIK® Inlay - e+, and e+ Patellas are performed by hydrogen peroxide gas plasma to achieve a Sterility Assurance Level (SAL) of 10<sup>-6</sup>. These liners are single-use devices and CANNOT be resterilized by a healthcare facility. Trials and other instruments are used to determine sizing before the sterile package needs to be opened. Should the original sterile package be inadvertently opened or compromised before implantation, the device cannot be implanted. Contact manufacturer or manufacturer's representative for instructions. Do not resterilize an implant or component that has been in contact with or contaminated by blood or other substances.

Instruments are provided nonsterile and should be stored in their original packaging until cleaned and sterilized according to the recommended guidelines found in the DJO Surgical Instrumentation Instructions for Use.

**WARNING: DO NOT resterilize any knee prosthesis distributed by DJO Surgical (Encore Medical, L.P.) if sterile packaging is opened or damaged upon receipt. Return the implant with respective packaging to DJO Surgical for inspection and disposition.**

**WARNING: Protect all porous coated and polished surfaces. Standard cleaning procedures cannot be relied upon to remove contamination from porous coating.**

**WARNING: DO NOT resterilize UHMWPE (ultra-high molecular weight polyethylene).**

DJO Surgical has validated sterilization cycle data on file.

**NOTE:** DJO Surgical does not recommend Flash or Chemical Sterilization.

For further information regarding the use of the DJO Surgical Knee Systems, contact your DJO Surgical representative or distributor.

DJO Surgical Knee Systems are manufactured by ENCORE MEDICAL, L.P.  
9800 Metric Blvd., Austin, TX 78758 USA (Made in the USA)

**10. Trademarks and patents**

FOUNDATION®, 3DKNEE™, EPIK®, is a trademark of DJO Surgical.  
U.S. patents: 5,413,604

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