


# AirHeel<sup>™</sup> w/Stabilizers

# PATIENT APPLICATION





# BEFORE USING THIS DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

#### **INTENDED USER PROFILE**

THE USER SHOULD BE ABLE TO:

- READ AND UNDERSTAND THE DIRECTIONS, WARNINGS AND CAUTIONS.

### **INTENDED USE/INDICATIONS**

Intended for plantar fasciitis and Achilles tendonitis. Note: For optimal comfort and support, wear an absorbent sock and a laced athletic shoe when wearing the brace.

### **CONTRAINDICATIONS**

None

# **APPLICATION INSTRUCTIONS**

1) Insert foot into AirHeel. Position heel in opening in back of wrap. Position Achilles aircell over tendon. The height of the Achilles aircell may be repositioned by detaching and firmly reattaching in place.

2) Place foot flat on floor and secure ankle strap. Make sure that foot aircell is located under arch of foot. To adjust AirHeel fit, tighten/loosen ankle strap.

#### WARNINGS AND PRECAUTIONS

If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, consult your medical professional immediately.

## **CLEANING INSTRUCTIONS**

Hand wash in warm water using mild soap, rinse thoroughly. AIR DRY. Note: If not rinsed thoroughly, residual soap may cause irritation and deteriorate material.

## INTENDED FOR USE ON A SINGLE PATIENT

#### NOT MADE WITH NATURAL RUBBER LATEX

#### WARRANTY:

DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

**NOTICE:** WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.



DJO, LLC T 800.336.6569 D 760.727.1280 F 800.936.6569 2900 Lake Vista Drive I Lewisville, TX 75067 I U.S.A. DJOglobal.com

COPYRIGHT © DJO, LLC 2020

MKTIS00-9866 Rev A

Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.