

REACTION™ Knee Brace

PATIENT APPLICATION

BEFORE USING THIS DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USE/INDICATIONS:

The REACTION™ knee brace is designed to help reduce pain by shock absorption and by returning kinetic energy to the knee. Indications for use are mild to moderate patellofemoral pain or maltracking and mild to moderate osteoarthritis of the knee.

CONTRAINDICATIONS:

- None

WARNINGS AND PRECAUTIONS:

If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, consult your medical professional immediately.

CLEANING INSTRUCTIONS:

Both the plastic parts and pads can be hand wash with mild soap in lukewarm water. Rinse both thoroughly. AIR DRY separately.

FOR SINGLE PATIENT USE

NOT MADE WITH NATURAL RUBBER LATEX.

WARRANTY:

DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

NOTICE: NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE OF THE ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THIS DEVICE IS ONLY ONE ELEMENT IN THE OVERALL TREATMENT PROGRAM ADMINISTERED BY A MEDICAL PROFESSIONAL. THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.



Application Information



1. Pull on sleeve until comfortable and patella is centered. (Wearing sleeve is optional.)
- 2a. Lay web over top of knee joint, centering web by aligning opening with patella (figure 2b).
3. Thread tabs on one side of back pad through inside openings on web and affix hook portions to loop face of back pad. (Note: Spacer material side of back pad should be placed toward leg and loop side away from leg. Outside of web has the logo, inside of web is dimpled.)
4. Thread loose tabs of back pad through available openings on web and then affix hook tabs to loop surface on back pad.
- 5a. Once all tabs are affixed to assure proper tightness, lift and pull both of lower tabs until snug.
- 5b. Repeat for top pair of tabs -brace should be snug and comfortable, do not over tighten.



DJO, LLC

T 800.336.6569 D 760.727.1280 F 800.936.6569
2900 Lake Vista Drive | Lewisville, TX 75067 | U.S.A.
DJOglobal.com

COPYRIGHT © DJO, LLC 2020
MKTIS00-9850 REV A

Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.