


1 Initial Preparation



- a. Pull out of the bag
- b. Open straps



- c. Roll out the brace

2 Position Brace



- a. Position under the leg.
- b. Shorter upright is applied to the thigh

3 Align first upright



- a. Centered midline to the leg and midline to the patella

4 Foam liners



- a. Begin with upper calf liner
- b. Wrap additional liners snugly around the leg

5 Position additional upright



- a. Position and press into place
- b. Be sure hinges are at equal height

6 Strapping



- a. Pull all slack out of straps
- b. Loop through the D-rings and fasten the 2 straps closest to the knee
 - i. The upper calf being first
- c. Secure additional straps

Hinge

Adjustments to be made per physician protocol



- Adjust flexion and extension by pressing in the corresponding buttons and positioning to the desired Range of Motion setting
- Quick Lock feature allows for quick and easy locked settings of 0°, 15° or 30°

Uprights



- Telescope to desired length by depressing the slider button until the desired length is reached
- Ensure that both sides are equal in length

Straps



- "Y" hook and loop strap ends allow for easy length adjustments
- Cut to desired length

Tips

- Perform preliminary fitting prior to surgery whenever possible.
- Fit to patient with the brace in a locked out position.



DJO, LLC | A DJO Global Company

T 800.336.6569 **D** 760.727.1280 **F** 800.936.6569

1430 Decision Street | Vista, CA 92081-8553 | U.S.A.

DJOglobal.com