Restoring
Life's
Vitality!

VITALSTIM™

User Manual

TABLE OF CONTENTS

VitalStim™ Therapy Unit

FOREWORD	1
ELECTROTHERAPY	2-4
Precautionary Instructions	2
Caution	2-3
Warning	2-3
Danger	3
Indications, Contraindications & Adverse Effects	4
NOMENCLATURE	5
SPECIFICATIONS	6
SETUP	7-8
Battery Installation	7
Lead Wire Installation	7
Unit Preparation	8
Amplitude Lock Feature	8
ACCESSORIES	9
Standard Accessories	9
MAINTENANCE	10
Cleaning the Unit	10
Long Term Storage	10
Factory Service	10
WARRANTY	11
NOTES	12-13

©2005 Encore Medical Corporation or its affiliates, Austin, Texas, USA. Any use of editorial, pictorial or layout composition of this publication without express written consent from the Chattanooga Group of Encore Medical, L.P. is strictly prohibited. This publication was written, illustrated and prepared for print by the Chattanooga Group of Encore Medical, L.P.

U.S. Patents: 5,725,564, 5,891,185, 5,987,359, 6,104,958, and 6,198,970.

FOREWORD

VitalStim™ Therapy Unit

Thank you for purchasing the Chattanooga Group, VitalStim Therapy Unit.

This manual contains general safety, operating, maintenance and care instructions for the owners and operators of the VitalStim Therapy Unit.

At the time of publication the information contained herein was current and up to date. However, due to continual technological improvements and increased clinical knowledge in the field of VitalStim Therapy, as well as Chattanooga Group's policy of continual improvement, Chattanooga Group reserves the right to make periodic changes and improvements to our equipment and documentation without any obligation on the part of Chattanooga Group.

Read, understand and follow the information contained in this manual.

Stay current with the latest clinical developments in the field of VitalStim Therapy and observe all applicable precautionary measures for treatment.

Keep informed on appropriate indications and contraindications for the use of VitalStim Therapy.

Federal Law restricts this device to sale on or by the order of a physician or licensed practitioner.

This equipment is to be used only by a licensed practitioner certified in the use of the VitalStim Therapy System.

Precautionary Instructions

The precautionary instructions found in this section and throughout this manual are indicated by specific symbols. Understand these symbols and their definitions before operating this equipment. The definition of these symbols are as follows:

 =CAUTION- Text with a “CAUTION” indicator will explain possible Safety infractions that could have the potential to cause minor to moderate injury or damage to equipment.

 =WARNING- Text with a “WARNING” indicator will explain possible Safety infractions that will potentially cause serious injury and equipment damage.

 =DANGER- Text with a “DANGER” indicator will explain possible Safety infractions that are imminently hazardous situations that would result in death or serious injury.

NOTE: Throughout this manual “NOTE” may be found. These notes are helpful information to aid in the particular area or function being described.

CAUTION

- Read, understand and practice the precautionary and operating instructions found in this manual. Know the limitations and hazards associated with using any electrical device. Observe the precautionary and operational decals placed on the unit.
- DO NOT operate when connected to any other unit.
- Use only fingers to operate button controls on the control panel(s). Use of sharp objects such as pencils or pens will result in damage to the unit.
- Turn unit OFF before positioning electrodes on a patient, removing electrodes from a patient or replacing batteries.
- Place the patient in a comfortable position during the VitalStim therapy session.
- Proper storage and transport temperatures for the VitalStim Therapy Unit are 40°F - 158°F (-40°C - 70°C). Relative Humidity is 10% - 100%.
- Safe use of electrotherapy during pregnancy has not been established.
- Caution should be used for patients with suspected or diagnosed heart problems.
- Caution should be used for patients with suspected or diagnosed epilepsy.
- Caution should be used in the presence of a tendency to hemorrhage following acute trauma or fracture, following recent surgical procedures when muscle contraction may disrupt the healing process, over a menstruating or pregnant uterus and over areas of skin which lack normal sensation.
- Some patients may experience skin irritation or hypersensitivity due to electrical stimulation or electrical conductive medium. The irritation can usually be reduced using an alternate conductive medium or electrode placement.

CAUTION

- Electrode placement and stimulation settings should be based on the guidance of the prescribing practitioner or other licensed health professional, certified in VitalStim Therapy.
- Powered muscle stimulators should be kept out of the reach of children.
- Powered muscle stimulators should be used only with the leads and electrodes recommended for use by the manufacturer.
- Portable powered muscle stimulators should not be used while driving, operating machinery, or during any activity in which involuntary muscle contractions may put the user at undue risk of injury.
- Always check the stimulation controls before treating a patient.
- Always turn the unit off before removing electrodes from the skin.
- Always turn the unit off before replacing the battery.
- For long term storage of the VitalStim Unit, the batteries should be removed.
- Stimulation should not be applied over, or in proximity to, cancerous lesions.

WARNING

- Federal law restricts this device to sale by, or on the order of, a physician or licensed practitioner. This device should be used only by a physician or licensed practitioner certified in VitalStim Therapy.
- Extra care should be taken when this device is used with children.
- Care must be taken when operating this equipment around other equipment. Potential electromagnetic or other interference could occur to this or to the other equipment. Try to minimize this interference by not using other equipment in conjunction with it.

WARNING

- Before administering any treatment to a patient, the operator must be certified in VitalStim Therapy.
- The long term effects of chronic electrical stimulation are unknown.
- Stimulation should not be applied transthoracically because the introduction of electrical current into the heart may cause cardiac arrhythmias.
- Stimulation should not be applied transcerebrally.
- Stimulation should not be applied over swollen, infected, or inflamed areas or skin eruptions, e.g., phlebitis, thrombophlebitis, varicose veins, etc.
- Dispose of all used accessories and spent batteries according to local, state and federal codes and regulations.

DANGER

- Stimulation should not be applied over the carotid sinus nerve particularly in patients with a known sensitivity to the carotid sinus reflex as death or serious injury is likely to occur.
- Use only electrodes and accessories designed specifically for use with the VitalStim Therapy Unit. Use of other accessories and/or techniques not approved under the VitalStim certification training may result in death, injury, or adverse effects to patient or undesirable and ineffective results.

Indications, Contraindications & Adverse Effects

Indications

The VitalStim Therapy Unit is intended for muscle re-education and for the creation of functional muscle use patterns in the throat muscles necessary for pharyngeal contraction. The device is indicated for treating from any etiology other than mechanical causes requiring surgery.

The VitalStim Therapy Unit is a prescription device intended for use by or on the order of a physician or other licensed professional, certified in the use of the VitalStim Therapy System.

Contraindications

This device should be used with caution on patients with cardiac demand pacemakers.

Use of this device is contraindicated in patients who are severely demented and exhibit non-stop verbalization. Constant verbalization could result in aspiration during trials of oral intake.

Use of this device is contraindicated in patients with significant reflux due to use of a feeding tube. Such patients are prone to repeated cases of aspiration pneumonia, and the device has not been studied in this population.


Use of this device is contraindicated in patients with dysphagia due to drug toxicity. Patients suffering from drug toxicity could aspirate during trials of oral intake.

Adverse Effects

Skin irritation and burns beneath the electrodes have been reported with the use of powered muscle stimulators. Follow the VitalStim Electrode Instruction Card, contained within the electrode packet, for recommended use.

NOMENCLATURE

VitalStim™ Therapy Unit

1. Channels 1 and 2

Lead wires connect in the respective ports.

2. Belt Clip

Attach to patient during therapy session. The Belt Clip also swings away from unit for battery installation.

3. Operator LCD Panel (See Items 10 and 11)

4. Battery Compartment

Battery compartment is on back of unit and houses two "AA" 1.5 volt batteries (included).

5. Channel 2 Intensity Up and Down Controls

Used to increase or decrease Channel 2 intensity.

6. Power On/Off Button

Used to turn unit On or Off.

7. Channel 1 Intensity Up and Down Controls.

Used to increase or decrease Channel 1 Intensity.

8. Amplitude Lock

Locks therapy amplitude after 30 seconds. Prevents accidental increase of prescribed amplitude during therapy session.

9. Timer Icon and Timer

Therapy session timer. Fixed at 60 minutes.

10. Amplitude Channel 2 Indicator

Shows Amplitude of Channel 2 during therapy session. Max. setting is 25mA, adjustable in 0.5mA increments.

11. Amplitude Channel 1 Indicator


Shows Amplitude of Channel 1 during therapy session. Max. setting is 25mA, adjustable in 0.5mA increments.

SPECIFICATIONS

- Output Configuration** Dual Channel, electrically isolated with miniature safety connectors.
- Output Waveforms** AC Mode: Rectangular symmetrical biphasic with zero net DC.
- Patient Voltage** 100 Volts Maximum, No Load
- Intensity Control** Dual intensity potentiometers: 0-25 mA peak current output, adjustable current. Constant current regulation from 0 to 4000 ohms load minimum.
- Pulse Rate** Fixed, 80 Hz
- Pulse Duration** Fixed, 300 μ Sec
- Pulse Charge** Maximum of 8 μ Coulombs under normal operation
- Output Protection** Under 15 μ Coulombs charge per pulse with any single component failure.
- Time Base** Crystal controlled oscillator operating at 3.58 MHz
- Power Source** Two AA (1.5 volt) Alkaline Batteries. ANSI Type NEDA1604A. Nickel-Metal-Hydrate or Nickel-Cadmium rechargeable batteries.
- Electrical Type** Type BF 
- Battery Life** 24 hrs or until Low Battery indicator blinks
- Circuit Design** Microcomputer with no internal adjustments
- Calibration** No Calibration Required or Needed

Safety Recognition

- United States: Conforms to UL Std 2601-1
Canada: Certified to CAN/CSA C22.2 No. 601.1
Meets 93/42/EEC Medical Device Directive
IEC/EN 60601-1


SETUP

VitalStim™ Therapy Unit

Battery Installation


Open the Belt Clip and swing away from unit back.

With a flat blade screwdriver, remove the retaining screw from the Battery Cover.

With thumb, press down on battery cover and toward unit bottom to expose battery cavity.

Install Batteries as illustrated inside unit.

Replace cover in reverse order. Install retaining screw. Snap the belt clip into position.

Lead Wire Installation

Position lead wire connector over unit connector and push down until seated.

NOTE:

Refer to the VitalStim Electrode Instruction Card contained within the electrode packet for electrode preparation and lead wire attachment to electrodes.

Attach Lead Wire Clip to Patient. Make certain there is enough slack in the lead wire to allow head movement.

Unit Preparation

After electrodes are placed on the patient per the VitalStim Certification Training Materials, turn the unit on by depressing the On/Off button once.

Set each Channel's intensity by pressing the intensity Up or Down arrow until the prescribed intensity is visible in the LCD display.

NOTE:

Electrodes must be secured to patient for amplitude to be set above 10.0 mA.

Place the unit in the Anti-Shock pouch and close securely with the hook and loop straps.


AMPLITUDE LOCK FEATURE

The Amplitude Lock Feature will automatically turn on after the Amplitude has been set and left undisturbed for approximately 30 seconds. This prevents accidental changes in the prescribed Amplitude Setting during a therapy session.

To turn off the Amplitude Lock feature, press the down arrow for either channel, then re-adjust amplitude as desired.

STANDARD ACCESSORIES

The VitalStim System can be ordered in two different kits;

Kit Number 5900- Includes VitalStim Dual Channel Unit, Lead wires, User Manual and 1 case of VitalStim Adult Size Electrodes.

Replacement Adult Electrodes Part Number: 59000

Kit Number 5905- Includes VitalStim Dual Channel Unit, Lead wires, User Manual and 1 case of VitalStim Child Size Electrodes.

Replacement Child Electrodes Part Number: 59005

Replacement Lead Wires Part Number- 59006

1. VitalStim Electrode Kits

Each pouch of Adult Electrodes and Child Electrodes contains four electrodes, Clean-Cote™ Skin Dressing Wipe, Electrode Adhesive Tape Patches and Aloe Vera Gel. The electrodes are for single use only and should be disposed of properly after each therapy session.

2. VitalStim Unit

3. Lead Wires with Clips

Lead wires connect the patient electrodes to the unit for therapy treatments. The Clips are used to secure lead wires to patient.

Anti-Shock Pouch (Not Shown)

Place unit in the Anti-Shock pouch to protect unit from damage or accidental intensity changes.

Storage Pouch (Not Shown)

Used to store and transport the VitalStim Unit and Accessories.

Cleaning the Unit

Clean the unit with a clean, lint free cloth moistened with water and mild antibacterial soap. If a more sterile cleaning is needed, use a cloth moistened with an antimicrobial cleaner.

Do not submerge the unit in water. Should the unit accidentally become submersed, remove the batteries, wrap the unit in a towel and drain as much water from the unit as possible. Set the unit upright on the towel and allow to air dry for at least 12 hrs in an area with low humidity.

Do not attempt any use of the unit until it is completely dry of all internal moisture. Any attempted use of a unit that may be wet inside could result in damage to the internal components of the unit.

Long Term Storage

If unit is to be stored or not used for 30 days or longer, remove the batteries.

Factory Service

When the VitalStim Therapy unit requires factory service, contact the selling dealer or Chattanooga Group Service Department.

All units returned to the factory for service must include the following:

WARRANTY REPAIR/OUT OF WARRANTY REPAIR

1. Written statement containing the following information:
 - RA Number- Obtain from Factory
 - Unit Model Number
 - Unit Serial Number
 - Contact person with Phone and Fax Numbers
 - Billing Address (for Out of Warranty Repair)
 - Shipping Address (Where to Ship Unit after Repair)
 - Detailed Description of Problem or Symptoms
2. Copy of original invoice issued at purchase of the unit.
3. Ship unit to Factory in the original container with all accessories and information as required in item 1 above to:

Chattanooga Group
4717 Adams Road
Hixson, TN 37343

Phone: USA: 1-866-864-0598

Canada: 1-800-361-3661

Outside USA: +1-423-870-7200

FAX: 1-423-870-7407

Web Address: www.chattgroup.com

Service to these units should be performed only by Service Technicians Certified by Chattanooga Group.

WARRANTY

VitalStim™ Therapy Unit

Chattanooga Group ("Company") warrants that the VitalStim Therapy units ("Product") are free of defects in material and workmanship. This warranty shall remain in effect for one year (12 months) from the date of original consumer purchase. If this Product fails to function during the one year warranty period due to a defect in material or workmanship, Company or the selling dealer will repair or replace this Product without charge within a period of thirty (30) days from the date on which the Product is returned to the Company or the dealer.

All repairs to the Product must be performed by a service center authorized by the Company. Any modifications or repairs performed by unauthorized centers or groups will void this warranty.

The warranty period for accessories is 90 days. Accessories consist of Lead Wires and Anti-Shock Pouch.

To participate in warranty coverage, this Product's warranty registration card (included with Product) must be filled out and returned to the Company by the original owner within ten (10) business days of purchase.

This Warranty Does Not Cover:

Replacement parts or labor furnished by anyone other than the Company, the selling dealer or a certified Company service technician.

Defects or damage caused by labor furnished by someone other than Company, the selling dealer or a certified Company service technician.

Any malfunction or failure in the Product caused by product misuse, including, but not limited to, the failure to provide reasonable and required maintenance or any use that is inconsistent with the Product User's Manual.

COMPANY SHALL NOT BE LIABLE IN ANY EVENT FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES.

Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

To Obtain Service From Company or the selling dealer under this warranty:

1. A written claim must be made within the warranty period to the Company or the selling dealer. Written claims made to the Company should be sent to:

4717 Adams Road
P.O. Box 489
Hixson, TN 37343 USA
Telephone: 1-423-870-2281 or 1-800-592-7329
Outside USA: +1-423-870-7200
Facsimile: 1-423-875-5497

and

2. The Product must be returned to the Company or the selling dealer by the owner.

This warranty gives you specific legal rights and you may also have other rights which vary from location to location.

The Company does not authorize any person or representative to create for it any other obligation or liability in connection with the sale of the Product. Any representation or agreement not contained in the warranty shall be void and of no effect.

THE FOREGOING WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OR MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

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Hixson, TN 37343 U.S.A.
1-800-592-7329 U.S.A.
+ 1-423-870-2046 OUTSIDE U.S.A. FAX