

EF1

CIRCULAR FRAME

A versatile fixation device with multiple sizing options to accommodate various indications and patient anatomies.

EXTERNAL FIXATION

VERSATILE OPTIONS

The EF1 Circular Frame is an external fixation device designed for fractures and fusions in the foot & ankle. Ring sizing is available in 120mm, 140mm, 160mm, 180mm, 200mm, 220mm, and 240mm to accommodate different patient sizes.



FEATURES & BENEFITS



EASE OF USE

ANGULAR CORRECTION CLAMPS

Angular correction clamps allow the surgeon to angulate, rotate, compress, and distract in all directions. Other products require four sizes of struts to accomplish the same thing that just one of the EF1 Circular Frame System's angular correction clamps can, potentially saving surgeons time in the OR.

The angular correction clamp technology was formulated from the EF1 Mini Rail System, which works together with the Circular Frame and allows both systems to work interchangeably.

STRONG GRIP

ROUGH COATING

EF1 Circular Frame rings have a rough, sandpaper-like coating that allows wires, nuts, bolts, etc. to grip the ring.

enovis.

T 800.495.2919 F 877.778.3864

Trilliant Surgical, LLC
727 North Shepherd Drive | Suite 100 | Houston, TX 77007 | U.S.A.
enovis.com/foot-and-ankle

Copyright © 2023 Enovis Foot and Ankle
MK-10400 Rev00

Individual results may vary. Neither Trilliant Surgical, LLC nor any of the Enovis companies dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.