



Force Line

New Premium Elastic Knitted Supports

- Smart Design
- Optimal Compression & Support
- Anatomical Fit
- Enhanced Comfort





Force Line

Common Features & Benefits

- 3D contoured knit for effective compression
- Anatomically contoured design with targeted compression & support
- Wrinkle-free, non-slip fit for comfortable wear
- Breathable material

German Engineering & Knitting Technology



MANUFORCE™

Elastic knitted wrist support provides compression & stabilization of wrist joint

EPIFORCE™

Elastic knit elbow sleeve with epicondyle silicone inserts offer compressive support

GENUFORCE™

Knitted knee sleeve provides compressive support & comfort

ACHILLOFORCE AIR™

Elastic knitted ankle support with pneumatic aircells providing intermittent compression on the achilles tendon


MALLEOFORCE™

Knitted ankle sleeve with anatomically aligned silicone inserts for compressive support around the ankle joint

MalleoForce™

INDICATIONS

- Post-operative or post-traumatic use
- Chronic soft tissue inflammation
- Chronic instability
- Sprains, moderate ligament injuries
- Recurrent swelling, load-related pain in the foot and ankle



Smart design

Correct ankle flexion angle is integrated in the knit design for an ideal fit

Protection & comfort

Anatomically shaped silicone inserts protect the malleoli & provide a massaging effect



Focused flexibility

The softness of the yarn and the knitting pattern do not hinder ankle flexion

Anatomical shape

The heel profile provides support and correct fit

ManuForce™

INDICATIONS

- Chronic, post-traumatic or postoperative inflammation of wrist joint and/or soft tissues
- Wrist sprains, Tendinitis (tenosynovitis)
- Carpal tunnel syndrome



Comfortable silicone insert

Helps reduce pressure on bony prominence and nerves



Malleable palmar aluminum stay for stabilization

Can be removed following treatment of acute phase



Removable strap

To adjust compression level

GenuForce™

INDICATIONS

- Post-operative or post-traumatic use
- Chronic soft tissue inflammation
- Recurrent swelling or load-related pain
- Chronic instability, mild ligament injuries, meniscus injuries



Brace stability

Silicon lined stays ensure non-slip grip & targeted compression zones for support where you need it most

Popliteal comfort

The softness of the yarn and the knitted pattern of the popliteal area help ensure comfort and wrinkle-free movement




Patellar stabilization

Flexible silicone buttress maintains the patella in position and provides a massaging effect

AchilloForce Air™

INDICATIONS

- Chronic or post-traumatic inflammation of the achilles tendon area
- Achilles tendonitis
- Plantar fasciitis



Interconnected pneumatic aircells for intermittent compression

The aircells provide intermittent compression and a massaging effect on the achilles tendon area while walking

Removable elastic strap

To adjust compression level



Optimal effect

Best worn in a sports shoe

EpiForce™

INDICATIONS

- Chronic, post-traumatic or postoperative inflammation of the elbow joint
- Tennis elbow (lateral epicondylitis)
- Golfer's elbow (medial epicondylitis)
- Bursitis



Removable strap

To adjust compression level


Comfort zone in elbow crease

For comfort of wear


Silicone inserts

Help relieve pressure on the elbow muscles


DONJOY® GENUFORCE™			
PART NUMBER	SIZE	THIGH CIRCUMFERENCE—6" ABOVE PATELLA (A)	CALF CIRCUMFERENCE—2" UNDER PATELLA (B)
		A (INCHES)	B (INCHES)
11-0027-1	1 = XS	14½ - 15¼"	9¾ - 11"
11-0027-2	2 = S	15¾ - 17"	11 - 11¾"
11-0027-3	3 = M	17 - 18"	11¾ - 13½"
11-0027-4	4 = L	18 - 19¼"	13½ - 14½"
11-0027-5	5 = XL	19¼ - 20½"	14½ - 15¼"
11-0027-6	6 = XXL	20½ - 21¾"	15¾ - 17"
11-0027-7	7 = XXXL	21¾ - 22¾"	17 - 18"




DONJOY® MALLEOFORCE™		
PART NUMBER	SIZE	ANKLE CIRCUMFERENCE
11-0028-1	1 = XS	6¾ - 7½"
11-0028-2	2 = S	7½ - 8¼"
11-0028-3	3 = M	8¼ - 9"
11-0028-4	4 = L	9 - 9¾"
11-0028-5	5 = XL	9¾ - 10½"
11-0028-6	6 = XXL	10½ - 11½"




DONJOY® ACHILLOFORCE AIR™		
PART NUMBER	SIZE	ANKLE CIRCUMFERENCE
11-0251-1	1 = XS	6.75 - 7½"
11-0251-2	2 = S	7.5 - 8¼"
11-0251-3	3 = M	8.25 - 9"
11-0251-4	4 = L	9 - 9¾"
11-0251-5	5 = XL	9.75 - 10½"
11-0251-6	6 = XXL	10.5 - 11½"



DONJOY® EPIFORCE™		
PART NUMBER	SIZE	FOREARM CIRCUMFERENCE
11-0253-1	1 = XS	7½ - 8¾"
11-0253-2	2 = S	8¾ - 9¾"
11-0253-3	3 = M	9¾ - 11"
11-0253-4	4 = L	11 - 12¼"
11-0253-5	5 = XL	12¼ - 13½"
11-0253-6	6 = XXL	13½ - 14½"



DONJOY® MANUFORCE™			
PART NUMBER RIGHT	PART NUMBER LEFT	SIZE	WRIST CIRCUMFERENCE
11-0254-1	11-0255-1	1 = XS	5 - 5¾"
11-0254-2	11-0255-2	2 = S	5¾ - 6¼"
11-0254-3	11-0255-3	3 = M	6¼ - 7"
11-0254-4	11-0255-4	4 = L	7 - 7½"
11-0254-5	11-0255-5	5 = XL	7½ - 8"
11-0254-6	11-0255-6	6 = XXL	8 - 8¾"



Copyright © 2016 by DJO, LLC
All rights reserved. No part of this publication may be reproduced, distributed, or transmitted in any form or by any means.



MOTION IS MEDICINE™

Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this brochure do not constitute medical, legal, or any other type of professional advice. Information related to various health, medical, and fitness conditions and their treatment is not meant to be a substitute for the advice provided by a physician or other medical professional. You should not use the information contained herein for diagnosing a health or fitness problem or disease. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you. Individual results may vary and results are not warranted or guaranteed.

DJO, LLC | A DJO Global Company
T 800.336.6569 **D** 760.727.1280 **F** 800.936.6569
 1430 Decision Street | Vista, CA 92081-8553 | U.S.A.
DJOglobal.com

DJO Global proudly sponsors:

