

CRYOKNEE™

The CryoKnee™ brace combines comfortable softgoods with four-way stretch fabric for compression, support and ability to easily integrate cold/heat therapy. The rigid and malleable Range of Motion (ROM) hinges provide protected range of motion, medial/lateral stability and a customizable fit. Included therapy pad pouch attaches to the brace for optional cold or heat therapy (cold or hot pad sold separately).

INDICATIONS:

- · Post-operative support for arthroscopic knee procedures
- · Post-operative support for total knee replacement
- · Mild to moderate ACL, PCL, MCL and/or LCL injuries and surgeries
- · Mild to moderate knee sprains and strains
- · Instability or laxity of the knee



COLD/HEAT THERAPY

Integrated cold or heat therapy pad pouch attaches easily to the brace and accommodates most cold/ hot packs (packs sold separately). Breathable material reduces condensation and is a built-in barrier to help prevent skin irritation.



RIGID CUSTOMIZABLE ROM HINGES

Removable, lightweight rigid ROM hinges are customizable to contour to the anatomy. Adjustable flexion and extension settings to limit knee motion. Extension can be limited from 0°-90° and flexion can be limited from 0°-150°. The brace can also be locked out to immobilize the knee.



4-WAY STRETCH PATELLA FABRIC

Flexible four-way stretch fabric accommodates all patella sizes and provides light compression on patella to help reduce window edema.



TRIMMABLE

Customizable softgood can be trimmed and custom-fit to the patient's leg. Trim calf strap, thigh strap and/or open up popliteal area.

Universal sizing fits both right and left knee across a wide range of patients.



enovis...

T 800.336.6569 F 800.457.4221

DJO, LLC 5919 Sea Otter Place, Suite 200 | I Carlsbad, CA 92010 | I U.S.A. enovis.com/donjoy

Copyright © 2023 DJO, LLC MKT00-12748 Rev A

WARNING!

THIS DEVICE CAN BE COLD ENOUGH TO CAUSE SERIOUS INJURY.

Please read and understand all warnings and instructions for use before using the device.

DO NOT use this device without a prescription and direction from a physician.

Individual results may vary. Neither DJO, LLC nor any of the Enovis companies dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.