

ICEMAN® FAMILY FEATURES:

DONJOY® WRAP-ON PADS

DonJoy® offers a variety of different non-sterile wrap-on pad designs to aid in delivering cold therapy to various parts of the body. Wrap-On Pads are designed with a hook engageable material and come with a wrap to secure the cold pad around the affected area. Wrap-On Pads are designed to be used in conjunction with a moisture barrier (DonJoy Sterile Dressing).

- Universal
- McGuire Knee
- Hip
- Shoulder
- Forefoot
- Ankle

**Wrap-On Pads are not available in sterile configuration.*

DONJOY COLD PADS

DonJoy offers a variety of different sterile and non-sterile cold pads to aid in delivering cold therapy to various parts of the body. Cold Pads are designed to be used in conjunction with a moisture barrier (DonJoy Sterile Dressing).

- Universal
- Hand/Wrist
- Ankle
- Large Rectangle
- Small Rectangle

**Cold Pads are available in both sterile and non-sterile configurations.*

DONJOY STERILE DRESSINGS

DonJoy sterile dressings provide a sterile moisture barrier between the skin and the cold pad. Sterile dressings are available in 4 different styles to use together with the various IceMan Cold Pads and Wrap-On Pads.

- Regular Universal
- XL Universal
- Shoulder
- Rectangle

PROVEN PERFORMANCE

- Time tested with over 20-years of experience delivering industry leading cold therapy solutions
- Developed through direct customer interaction and feedback

TECHNOLOGY ADVANTAGE

- Semi-closed loop technology helps deliver consistent cold therapy throughout the cold pad, eliminating freezing inlet temperatures

PATIENT CONVENIENCE AND EASE-OF-USE

- Coolers designed for easy set-up
- Plug-n-play operation
- Designed for easy portability and storage
- A complete line of cold pads to help deliver cold therapy to where it is needed the most

DONJOY COLD WRAPS

DonJoy offers several cold wrap styles to use with the IceMan Cold Pads. These non-sterile Cold Wraps help secure and hold the cold pad in place, cover the affected area, and are designed to be used in conjunction with a barrier.

DONJOY®

ICEMAN CLASSIC 3™

ICEMAN® CLEAR 3



T 800.336.6569 D 760.727.1280 F 800.936.6569
2900 Lake Vista Drive | Dallas, TX 75067 U.S.A.
DJOGlobal.com

Copyright © 2019 by DJO, LLC
MKT00-8100 Rev B

Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this brochure do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.

WARNING: this device can be cold enough to cause serious injury. Do not use this device without a prescription from a physician.


DONJOY® ICEMAN® CLEAR³ & CLASSIC³ THE CONSISTENTLY COOL ADVANTAGE

DONJOY® ICEMAN® CLEAR³ AND DONJOY ICEMAN CLASSIC³™

The IceMan® CLEAR³ and IceMan CLASSIC³ cold therapy units help reduce pain and swelling, speeding up rehabilitation. The IceMan helps provide extended cold therapy to aid in a variety of indications and protocols as directed by a medical professional. Both utilize DonJoy's recirculation technology, which helps maintain more consistent and accurate temperatures than other cold therapy units.

CONSISTENTLY COOL DONJOY RECIRCULATION TECHNOLOGY

The DonJoy® IceMan features a recirculation system that allows water warmed after flowing through the pad to be preserved and remixed with cooler ice water at a constant flow rate, providing consistent cool water distribution throughout the cold pad.



ICEMAN with dual pump recirculation technology. The DonJoy® IceMan features a recirculation system that allows water warmed after flowing through the cold pad to be preserved and remixed with cooler ice water at a constant flow rate, providing consistent cool water distribution throughout the cold pad.

COMPETITORS Other versions restrict or pinch water flow to the cold pad, causing variable flow rates and inconsistent cold pad temperatures.

The IceMan® CLEAR³ cold therapy unit helps reduce pain and swelling, speeding up rehabilitation and recovery. The IceMan helps provide extended cold therapy to aid a variety of indications and protocols directed by a medical professional. The IceMan® CLEAR³ utilizes DonJoy's recirculation technology, which helps maintain more consistent and accurate temperatures.



DONJOY® ICEMAN® CLEAR ³	
Description	Part Number
IceMan® CLEAR ³ Cold Therapy Unit	11-0493
IceMan CLEAR ³ w/ Universal Cold Pad, S, RH	11-1629
IceMan CLEAR ³ w/ Universal Cold Pad, NS, RH	11-1630
IceMan CLEAR ³ w/ Universal Wrap-On Pad, NS, RH	11-1631

FEATURES

Proven Performance
 • Time tested with over 20-years of experience delivering industry leading cold therapy solutions
 • Developed through direct customer interaction and feedback

Technology Advantage
 • Recirculation technology helps deliver consistent cold therapy throughout the cold pad, eliminating freezing inlet temperatures

Patient Convenience And Ease-Of-Use
 • Easy set-up and operation
 • Plug-n-play operation
 • Designed for easy portability and storage
 • A complete line of cold pads to deliver cold therapy to where it is needed the most

DONJOY® ICEMAN® CLEAR ³	
Description	Part Number
IceMan CLEAR ³ w/ Universal Wrap-On Pad, NS, EH	11-1633
IceMan CLEAR ³ w/ Shoulder, S Wrap-On Pad, NS, EH	11-1635
IceMan CLEAR ³ w/ McGuire Knee Wrap-On Pad, NS, RH	11-1636
IceMan CLEAR ³ w/ Ankle Wrap-On Pad, NS, RH	11-1638

DonJoy® IceMan CLASSIC³™ provides a distinct advantage with the controlled delivery of cold therapy helping patients and medical practitioners to easily manage post-operative pain and swelling.



DONJOY® ICEMAN CLASSIC ³ ™	
Description	Part Number
IceMan CLASSIC ³ ™ Cold Therapy Unit	11-0494
IceMan CLASSIC ³ w/ Universal Cold Pad, S, RH	11-1421
IceMan CLASSIC ³ w/ Universal Cold Pad, NS, RH	11-1423


FEATURES

Proven Performance
 • Time tested with over 20-years of experience delivering industry leading cold therapy solutions
 • Developed through direct customer interaction and feedback

Technology Advantage
 • Recirculation technology helps deliver consistent cold therapy throughout the cold pad, eliminating freezing inlet temperatures

Patient Convenience And Ease-Of-Use
 • Easy set-up and operation
 • Plug-n-play operation
 • Designed for easy portability and storage
 • A complete line of cold pads to deliver cold therapy to where it is needed the most

DONJOY® ICEMAN CLASSIC ³ ™	
Description	Part Number
IceMan CLASSIC ³ ™ w/ Universal Wrap-On Pad, NS, RH	11-1422
IceMan CLASSIC ³ w/ Universal Wrap-On Pad, NS, EH	11-1424
IceMan CLASSIC ³ w/ Shoulder, Wrap-On Pad, NS, EH	11-9098
IceMan CLASSIC ³ w/ McGuire Knee Wrap-On Pad, NS, RH	11-9099



STERILE DRESSINGS	
Description (Non-Sterile)	Part Number
Universal Sterile Dressing	11-0822-9-00000
Universal XL Sterile Dressing	11-0823-9-00000
Rectangle XL Sterile Dressing	11-0820-9-00000
Shoulder XL Sterile Dressing	11-1360-9-00000

S = Sterile

COLD WRAPS	
Description (Non-Sterile)	Part Number
Shoulder S/M Cold Wraps	11-0611-9-00000
Shoulder L/XL Cold Wraps	11-0612-9-00000
Back S/M Cold Wraps	11-0613-9-00000
Back L/XL Cold Wraps	11-0614-9-00000
Knee/Thigh Cold Wraps	11-0610-9-00000

NS = Non-Sterile

ACCESSORIES	
Description (Non-Sterile)	Part Number
Bed Hanger	11-0657-9-00000
Dual Cold Pad Connector	11-0674-9-02000
Power Supply 8.5V (Domestic)	25-4882
Power Supply 12V (Domestic)	25-4041

RH = Regular Hose (6.5')

EH = Extended Hose (15')

WRAP-ON PADS	
Description (Non-Sterile)	Part Number
Universal Wrap-On Pad, NS, RH (10.85" x 11.34")	11-0679-9-00000
Universal Wrap-On Pad, NS, EH (10.85" x 11.34")	11-0680-9-00000
Universal Wrap-On XL, NS, RH (12" x 11.34")	11-0681-9-00000
Universal Wrap-On Pad XL, NS, EH (12" x 11.34")	11-0682-9-00000
McGuire Knee Wrap-On Pad, NS, RH (10.85" x 11.34")	11-1268-9-00000
Shoulder Wrap-On Pad S, NS, EH (12.75" x 19" x 12.25")	11-1397-9-00000
Shoulder Wrap-On Pad, NS, EH (14" x 23" x 12.5")	11-1334-9-00000
Ankle Wrap-On Pad, NS, RH	11-1522
Fore-Foot Wrap-On Pad, NS, RH	11-1521
Hip-Left Wrap-On Pad, NS, RH	11-9112
Hip-Right Wrap-On Pad, NS, RH	11-9111

COLD PADS	
Description (Sterile)	Part Number
Universal Cold Pad, S, RH (10.85" x 11.34")	11-0801-9-00000
Universal Cold Pad, S, EH (10.85" x 11.34")	11-0802-9-00000
Universal Cold Pad XL, S, EH (12" x 11.34")	11-0804-9-00000
Rectangle Cold Pad S, S, RH (4.81" x 10.75")	11-0817-9-00000
Rectangle Cold Pad XL, L, S, EH (9.76" x 13.75")	11-0813-9-00000
Wrist/Hand Cold Pad, S, RH (2.9" x 4.2")	11-1269-9-00000
Ankle Cold Pad, S, RH (12.26" x 8.28")	11-0815-9-00000
Universal Cold Pad, NS, RH (10.85" x 11.34")	11-0675-9-00000
Universal Cold Pad, NS, EH (10.85" x 11.34")	11-0676-9-00000
Universal Cold Pad XL, NS, EH (12" x 11.34")	11-0678-9-00000