


DYNACLIP®

BONE FIXATION THAT EMPOWERS
ACTIVE, ADAPTIVE HEALING

The DynaClip's versatile design, multiple configurations, and numerous size offerings allow for dependable durability and sustained dynamic compression across a variety of surgical applications.¹


STAPLE SYSTEM

DYNACLIP®




A x B
1 x 2 x 3

DYNACLIP FORTE™




A x B x C
1 x 2 x 3

DYNACLIP DELTA™








A x B x C x D x E
1 x 2 x 3 x 4

DYNACLIP QUATTRO™



A x B x C
1 x 2 x 3 x 4

Dimensions in millimeters (mm)



DYNACLIP DELTA™



24 x 16 x 14 x 16 x 18
10.5 x 1.8 x 5.5 x 2.5



28 x 16 x 14 x 20 x 22
10.5 x 1.9 x 5.5 x 2.5

DYNACLIP QUATTRO™



18 x 16 x 16
10.5 x 1.5 x 5.5 x 2.5




22 x 16 x 20
10.5 x 1.5 x 5.5 x 2.5



22 x 20 x 20
10.5 x 1.5 x 5.5 x 2.5

PROCEDURE GUIDE

INDICATION		BRIDGE LENGTH x LEG LENGTH
1	HALLUX IPJ FUSION	8x8, 10x10, 12x 12
2	LESSER PIPJ FUSION	8x8, 10x10
3	AKIN OSTEOTOMY	8x8, 10x10
4	MTPJ FUSION	10x10, 12x12, 14x14, 18x18 DELTA: 24x16x14x12, 28x16x14x12
5	CHEVRON OSTEOTOMY	14x14, 18x18
6	GENERAL MT OSTEOTOMY	12x12, 14x14, 18x18
7	BASE WEDGE OSTEOTOMY	10x10, 12x12, 14x14, 18x18
8	1ST TMT (LAPIDUS) FUSION	12x12, 14x14, 18x18, 20x18 FORTE: 26x20, 30x20, 26x16, 30x16
9	COTTON OSTEOTOMY	18x18, 18x20, 20x18, 20x22
10	2 ND /3 RD TMT (LISFRANC) FUSION	18x18, 18x20, 20x18, 20x22 FORTE: 26x20, 30x20, 26x16, 30x16
11	NAVICULOCUNEIFORM FUSION	14x14, 18x18, 18x20 FORTE: 26x20, 30x20
12	TALONAVICULAR FUSION	20x18, 20x22, 26x20 FORTE: 26x20, 30x20 QUATTRO: 18x16, 22x16, 22x20
13	CALCANEOCUBOID FUSION	20x18, 20x22, 26x20 FORTE: 26x20, 30x20 QUATTRO: 18x16, 22x16, 22x20
14	EVANS OSTEOTOMY	20x18, 20x22, 26x20
15	DYER OSTEOTOMY	20x18, 20x22, 26x20



enovis

T 800.495.2919 F 877778.3864

Medshape, Inc.
1575 Northside Drive NW | Suite 440 | Atlanta, GA 30318 | U.S.A.
enovis.com/foot-and-ankle

Copyright © 2023 MedShape, Inc.
MK-10171 Rev 06

1. Data on file. MedShape, Inc. 2018.

Individual results may vary. Neither Medshape, Inc. nor any of the Enovis companies dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.