


MATRIXCELLECT® 100 DBM PUTTY

MatrixCelect 100 is a human bone derived, flowable demineralized bone matrix putty with no additional carriers, indicated for homologous use for the treatment of surgically created or traumatic skeletal defects.

BIOLOGICS

- Extremely irrigation resistant properties intended to aid in securing final allograft placement
- In-vivo test results demonstrate all 5 bone forming elements present (chondrocytes, osteocytes, cartilage, bone marrow cells, and new bone)
- Every lot is verified for osteoinductivity post-sterilization
- Provided sterile via low dose, low temperature, gamma sterilization with a sterility assurance level (SAL) of 10⁻⁶
- Packaged with a shelf life of 2 years and may be stored in ambient temperatures
- Available in 1.0, 2.5, and 5.0cc syringes



| Part Number | Volume |
|-------------|---------------|
| 452010-TS | 1.0cc syringe |
| 452025-TS | 2.5cc syringe |
| 452050-TS | 5.0cc syringe |

Trilliant is a distributor for CellRight Technologies.