



MOVING REHABILITATION FORWARD™

Hot/Cold Therapy & Clinic Supplies | 2013



## INTRODUCTION



Chattanooga is the world's largest manufacturer of rehabilitation equipment for treating musculoskeletal, neurological and soft tissue disorders. For over six decades, Chattanooga has set a benchmark for leadership, reliability and excellence. Our products contribute to better treatment outcomes in hospitals, clinics and home settings worldwide. We lead the rehabilitation industry by example, and through continuous innovation, provide real world solutions for clinicians and their patients.

As the world's leading brand in rehabilitation, Chattanooga is sold in over 80 countries worldwide. Chattanooga distributors and customers are supported by a highly trained and seasoned team of personnel who oversee new product development, product training and education, order placement, product delivery, and after sales service in addition to providing dedicated sales and marketing support.

Chattanooga is a brand of DJO Global Inc. DJO Global provides solutions for musculoskeletal health, vascular health and pain management. Our products help prevent injuries or rehabilitate after surgery, injury or from degenerative disease, enabling patients to regain or maintain their natural motion.

Visit [www.DJOGlobal.com/chattanooga](http://www.DJOGlobal.com/chattanooga)



AIRCAST

chattanooga

CMFV

Compex

djosurgical

DONJOY

Dr. Comfort

Empi

PROCARE


## HOT/COLD THERAPY AND CLINIC SUPPLIES

A world wide leader in therapeutic products, Chattanooga offers a diverse selection of high quality hot/cold therapy and clinical rehabilitation supplies.

Chattanooga is the proud inventor of the original Hydrocollator HotPacs and heating unit, an all natural product that was first to market more than 60 years ago, and is still the first choice of physical therapists and other medical professionals when using heat therapy. The HotPac generates deep moist heat which delivers increased advantages over dry heat to relieve chronic pain, acute muscle strains and general musculoskeletal discomfort.

Chattanooga also sets the standard in cold therapy with the original ColPac, which come in a variety of shapes and sizes. Chilling units are also available for convenient storage and efficient cooling.

For precise measurements; hand dynamometers, goniometers, pinch gauge, and calipers are all available. For athletic training rooms and physical therapy clinics; carts, stools, lotion, and electrodes can all be ordered. Combining industry leading hot and cold therapy products with clinic supplies, Chattanooga can meet any clinics needs.





## CONTENTS

<b>Introduction</b>	<b>1</b>
<b>Hydrocollator Mobile Heating Units</b>	<b>3</b>
<b>Hydrocollator Stationary Heating Units</b>	<b>4</b>
<b>Hydrocollator HotPacs</b>	<b>5</b>
<b>Terry Covers</b>	<b>6</b>
<b>ColPac Chilling Unit</b>	<b>7</b>
<b>ColPac, boo-boo pac, Sport-Pac</b>	<b>8</b>
<b>Theratherm Digital Moist Heat</b>	<b>9</b>
<b>Fluidotherapy Dry Heat Therapy</b>	<b>10</b>
<b>Pedal Exercisers</b>	<b>11</b>
<b>Saunders Total Back</b>	<b>12</b>
<b>Stabilizer</b>	<b>12</b>
<b>Chattanooga Evaluation Devices</b>	<b>13</b>
<b>Baseline Evaluation Devices</b>	<b>15</b>
<b>Sully Shoulder Support</b>	<b>17</b>
<b>Carts &amp; Stools</b>	<b>18</b>
<b>Gels &amp; Lotions</b>	<b>19</b>
<b>Education Manuals</b>	<b>20</b>
<b>Electrodes &amp; Accessories</b>	<b>21</b>

## ORDERING INFORMATION

## Part Number Description

**M-2 Mobile Hydrocollator**

2402	M2 Mobile - Includes 12 Standard (1006) HotPacs
2402-2	M2 Mobile - Includes 3 Standard (1006), 3 Oversize (1004), and 3 Cervical (1002) HotPacs
2402-3	M2 Mobile - Includes 4 Standard (1006), 2 Oversize (1004), and 4 Cervical (1002) HotPacs
2402-4	M2 Mobile - Includes 4 Standard (1006), and 4 Oversize (1004)

**M-4 Mobile Hydrocollator**

2502	M4 Mobile - Includes 24 Standard (1006) HotPacs
2502-2	M4 Mobile - Includes 8 Oversize (1004), and 8 Cervical (1002) HotPacs
2502-3	M4 Mobile - Includes 6 Standard (1006), 6 Oversize (1004), and 6 Cervical (1002) HotPacs

**SS-2 Mobile Hydrocollator**

2302	SS2 Mobile - Includes 8 Standard (1006) HotPacs
2302-2	SS2 Mobile - Includes 4 Standard (1006), and 4 Cervical (1002), HotPacs
2302-3	SS2 Mobile - Includes 2 Standard (1006), 2 Oversize (1004), and 2 Cervical (1002) HotPacs

## TECHNICAL SPECIFICATIONS

**M-2 Mobile Hydrocollator**

Mains Power:	110-120 V, 50/60 Hz
Power Consumption:	1000 W
Weight:	48 lb (22 kg) /Shipping - 93 lb (42 kg)
Dimensions <sup>o</sup> :	27" x 16" x 33" (69 x 41 x 84 cm)
Electrical Safety Class:	Class 1, Type B
Safety Tests:	Conforms to UL 60601-1, certified to Can/CSA C222 No. 601.1

**M-4 Mobile Hydrocollator**

Mains Power:	110-120 V, 50/60 Hz
Power Consumption:	1500 W
Weight:	132 lb (60 kg) /Shipping - 160 lb (72 kg)
Dimensions <sup>o</sup> :	35" x 20" x 33" (89 x 51 x 84 cm)
Electrical Safety Class:	Class 1, Type B
Safety Tests:	Conforms to UL 60601-1, certified to Can/CSA C222 No. 601.1

**SS-2 Mobile Hydrocollator**

Mains Power:	110-120 V, 50/60 Hz
Power Consumption:	1000 W
Weight:	66 lb (30 kg) /Shipping - 80 lb (36 kg)
Dimensions <sup>o</sup> :	21" x 16" x 33" (53 x 41 x 84 cm)
Electrical Safety Class:	Class 1, Type B
Safety Tests:	Conforms to UL 60601-1, certified to Can/CSA C222 No. 601.1

<sup>o</sup> Dimensions are listed as L x W x H

## HYDROCOLLATOR® MOBILE HEATING UNITS



Our benchmark product is the standard all others can only hope to match. Durable and easy to maintain, these high-quality stainless steel units give you a constant supply of temperature-consistent HotPacs.



## FEATURES

## Hydrocollator Mobile Heating Units

- High-quality stainless steel
- Thermostatically controlled temperature
- 3" (8 cm) swivel-type rubber casters for silent, friction-free movement
- Easy maintenance. Simple to fill and drain. No plumbing required
- One-year warranty

\* Approximate with HotPacs

	M-2	M-4	SS-2	E-2	E-1/E-1 MHT
Tank Capacity	14 gal (52 L)	27 gal (102 L)	10 gal (37 L)	8 gal (37 L)	3 gal (37 L)
Temperature Range	160° - 165°F (71° - 74°C)	160° - 165°F (71° - 74°C)	160° - 165°F (71° - 74°C)	160° - 165°F (71° - 74°C)	160° - 165°F (71° - 74°C)
Thermal Cut-Out Temp.	180° - 185°F (82° - 85°C)	180° - 185°F (82° - 85°C)	180° - 185°F (82° - 85°C)	180° - 185°F (82° - 85°C)	180° - 185°F (82° - 85°C)
Temperature Accuracy	+/- 10%	+/- 10%	+/- 10%	+/- 10%	+/- 10%
Heat up Time to 160° F (to 70°C)	6 Hours	8 Hours	4 Hours	4 Hours	2 Hours
Cool Down Time from 160° F (from 70°C)	3 Hours	4 Hours	4 Hours	2 Hours	1 Hour
Fiberglass Insulation	Yes	Yes	Yes	No	No

## HYDROCOLLATOR® STATIONARY HEATING UNITS



E-1

E-2

### ORDERING INFORMATION

Part Number	Description
<b>E-1 Stationary Hydrocollator</b>	
2102	E1 Stationary - Includes 4 Standard (1006) HotPacs
2102-2	E1 Stationary - Includes 3 Standard (1006) HotPacs and 1 Cervical (1002) HotPac
2102-2	E1 Stationary - Includes 2 Standard (1006) HotPacs and 2 Cervical (1002) HotPacs
<b>E-2 Stationary Hydrocollator</b>	
2802	E1 Stationary - Includes 2 Oversize (1004), 3 Standard (1006) HotPacs and 1 Cervical (1002) HotPac
<b>Hydrocollator Accessories</b>	
4016	Wall Mounted Towel Rack
4010	Side Table Rack
4230	Extra Shelf for Side Table Rack (not shown)
22100	Mini Non-Electric Lotion Warmer
4012	TM-1 Electric Warmer
4228	Dial Thermometer (not shown)

### HYDROCOLLATOR ACCESSORIES

#### Wall Mounted Towel Rack



Six hooks extend 13" (33 cm) from wall. Measures 2" x 33" (5 cm x 84 cm).

#### Side Table Rack



Stainless steel work surface and two wire shelves add efficiency to units. Use with M-4, M-2 and SS-2.

#### Mini Non-Electric Lotion Warmer



Hangs on the front of the Hydrocollator E-1 or E-2 heating unit and uses the unit's warmth to keep lotion or gel at a comfortable temperature. Bottles are not included.

#### TM-1 Electric Lotion Warmer



Easily attaches to any treatment table or wall to keep massage lotion or ultrasound gel at a comfortable temperature and within reach. Bottles are not included.

### TECHNICAL SPECIFICATIONS

#### E-1 Stationary Hydrocollator

Mains Power:	110-120 V, 50/60 Hz
Power Consumption:	1000 W
Weight:	15 lb (7 kg) / Shipping - 21 lb (10 kg)
Dimensions*:	13" x 8" x 16" (33 x 20 x 41 cm)
Electrical Safety Class:	Class 1, Type B
Safety Tests:	Conforms to UL 60601-1, certified to Can/CSA C222 No. 601.1

#### E-2 Stationary Hydrocollator

Mains Power:	110-120 V, 50/60 Hz
Power Consumption:	1000 W
Weight:	29 lb (13 kg) / Shipping - 46 lb (21 kg)
Dimensions*:	15" x 13" x 20" (38 x 33 x 51 cm)
Electrical Safety Class:	Class 1, Type B
Safety Tests:	Conforms to UL 60601-1, certified to Can/CSA C222 No. 601.1

\* Dimensions are listed as L x W x H

## ORDERING INFORMATION

## Part Number Description

## Hydrocollator HotPac - Moist Heat Hot Pack

A 1008	10" x 18" (25 cm x 46 cm)
B 1010	10" x 24" (25 cm x 61 cm)
C 1004	Oversize 15" x 24" (38 cm x 61 cm)
D 1002	Neck Contour 24" long (61 cm)
E 1014	Half Size 5" x 12" (13 cm x 30 cm)
F 1006	Standard 10" x 12" (25 cm x 30 cm)
G 1012	Knee-Shoulder 10" x 20" (25 cm x 50 cm)

## HYDROCOLLATOR® HOTPACS™



The first choice for physical therapists, physicians, athletic trainers, and chiropractors for over 60 years.



## FEATURES

## The Original Hydrocollator® Moist Heat HotPac

Chattanooga is the inventor of the original Hydrocollator HotPac and proud to still be producing this all natural product today. The benefits of moist heat therapy over dry heat are well established for pain relief due to arthritis, myofibrosis, back pain, muscle strains and spasms.

## Hydrocollator HotPacs

- Seven HotPac sizes designed to fit every body contour
- Provides up to 30 minutes of deep, soothing, moist heat
- Brown tabs on most HotPacs for easy rotation in busy clinics
- One HotPac provides hundreds of cost effective treatments
- One-year warranty



## HYDROCOLLATOR® TERRY COVERS



### FEATURES

#### All-Terry Covers

- For high-volume institutional use or therapy at home
- Hook and loop closures to prevent unwrapping or slipping
- Washable and reusable
- One-year warranty

#### Foam-Filled Terry Covers

- For home or light clinic use
- Polyurethane foam layer between two layers of terry
- Available with or without pockets
- Washable and reusable
- One-year warranty

### ORDERING INFORMATION

Part Number	Description
<b>All-Terry Covers</b>	
1118	Standard 19" x 27" (48 cm x 69 cm)
1126	Oversize 24" x 36" (61 cm x 91 cm)
1104	Neck Contour 25" x 18" (64 cm x 46 cm)
<b>Foam-Filled Terry Covers</b>	
1102	Standard 19" x 27" (48 cm x 69 cm)
1124	Oversize 24" x 36" (61 cm x 91 cm)
1120	Neck Contour 25" x 16" (64 cm x 41 cm)
<b>Foam-Filled Pocket Terry Covers</b>	
1108	Standard 17" x 27" (43 cm x 69 cm)
1110	Oversize 24" x 36" (61 cm x 91 cm)
<b>Hydrocollator HotPac Set in Retail Packaging</b>	
1064	Standard Set in Retail Package includes one Standard HotPac and one Foam Filled Terry Cover
1062	Neck Set in Retail Package includes one Neck HotPac and one Foam Filled Terry Cover



**ORDERING INFORMATION****Part Number Description****ColPac Chilling Unit Model C-5**

3102	C-5 ColPac Chilling Unit Includes 6 Standard (1500) and 6 Half Size (1506) ColPacs
------	--

**Standard Accessories**

4010	Side Table Rack <i>See page 4 for Description and image</i>
4230	Extra Shelf for Side Table Rack ( <i>not shown</i> )

**ColPac Freezer**

90910	ColPac Freezer 21 1/4" x 26" x 33 1/2" (53 cm x 66 cm x 86 cm) Includes 12 Standard ColPacs (1500) Shipping weight 150 lbs (68 kg)
-------	---

**TECHNICAL SPECIFICATIONS****ColPac Chilling Unit Model C-5**

Mains Power:	115 VAC ~ 60 Hz
Current:	4A
Weight:	165 lb (74 kg) Shipping Weight - 190 lb (86 kg)
Dimensions*:	27" x 16" x 33" (69 x 41 x 84 cm)
Electrical Safety Class:	Type B
Safety Tests:	UL 544 (U.S. Only)

\* Dimensions are listed as L x W x H

**COLPAC® CHILLING UNIT**

Large capacity and fast chilling time means your patients will never have to wait for a ColPac. The ColPac Chilling Units help you keep pace throughout a busy day. No plumbing required, simply plug into a standard electrical outlet.

**FEATURES****ColPac Chilling Unit**

- Closed-cell foam insulation for efficiency
- Stainless steel, welded construction for long life
- Heavy-duty compressor for more efficient cooling
- Specially designed interior chilling coils cool packs faster
- 3-inch (8 cm) swivel-type rubber casters for silent, friction-free movement of mobile units
- Drain valve for easy cleaning and defrosting
- One-year warranty

**ColPac Freezer**

- Five cubic feet of spacious capacity
- Four durable stainless steel shelves
- Comes complete with trays and cup holders
- Shipped with 12 standard size ColPacs
- UL listed and cUL listed
- One-year warranty





## COLPAC®, BOO-BOO PAC™, FLEXIPAC®, SPORT-PAC™



Seven different sizes of blue vinyl and four sizes of black polyurethane ensure effective cold therapy for any body area or shape.



### FEATURES

#### The Original ColPac – There's a ColPac® for Every Application!

- Designed for heavy clinical use.
- Effective pain relief
- Filled with a non-toxic silica gel that remains pliable throughout treatment
- Delivers up to 30 minutes of soothing relief from acute pain, swelling and fever
- Latex free
- One-year warranty

#### boo-boo pac

- Fun, bear-shaped design especially made for kids
- Ideal for pediatric applications
- Filled with non-toxic silica gel
- Delivers up to 30 minutes of soothing relief from acute pain, swelling and fever
- Shipped in a display box containing 12 packs

#### FlexiPac

- Freezable to provide penetrating cold therapy
- Microwavable for soothing heat
- Flexible gel conforms to any body contour at any temperature

#### Sport-Pac

- Soccer ball design is ideal for youth sport applications
- Easily transferred from the freezer to the team cooler for sideline use
- Delivers up to 30 minutes of effective cold therapy for acute pain due to strains, sprains and bruises
- Shipped in a display box containing 10 packs

### ORDERING INFORMATION

Part Number	Description
<b>Blue Vinyl ColPac</b>	
A 1512	Override 11" x 21" (28 cm x 53 cm)
B 1500	Standard Size 11" x 14" (28 cm x 36 cm)
C 1506	Half Size 7½" x 11" (19 cm x 28 cm)
D 1508	Neck Contour 23" long (58 cm)
E 1502	3" x 11" (8 cm x 28 cm)
F 1504	Quarter Size 5½" x 7½" (14 cm x 19 cm)
G 1510	Eye (Specially shaped to cover closed eyelids)
<b>Black Polyurethane ColPac</b>	
H 1554	Neck Contour 21" long (53 cm)
I 1556	Override 12½" x 18½" (32 cm x 47 cm)
J 1552	Standard Size 10" x 13½" (25 cm x 34 cm)
K 1562	Half Size 6½" x 11" (17 cm x 28 cm)
<b>L boo-boo pac</b>	
Shipped in a display box containing 12 packs. NOTE: Must be purchased in quantities of 12 each	
1534	Royal Blue boo-boo pac (12 per case)
<b>M FlexiPac Hot &amp; Cold Compresses</b>	
4026	5" x 6" (13 cm x 15 cm) 48/case
4020	5" x 10" (13 cm x 25 cm) 24/case
4029	8" x 14" (20 cm x 35 cm) 12/case
4032	5" x 10" (13 cm x 25 cm) 25/case Chiropractically Yours (not shown)
<b>N Sport-Pac</b>	
Shipped in a display box containing 10 packs. NOTE: Must be purchased in quantities of 10 each	
1517	Sport-Pac (10 per case)

## ORDERING INFORMATION

## Part Number Description

## Theratherm Digital Moist Heat Pack

1030	Small 7" x 15" (18 cm x 38 cm)
1031	Medium 14" x 14" (36 cm x 36 cm)
1032	Standard 14" x 27" (36 cm x 68 cm)
1033	Shoulder/Neck 23" x 20" (58 cm x 50 cm)

## Replacement Flannel Covers

13330	Small 7" x 15" (18 cm x 38 cm)
13331	Medium 14" x 14" (36 cm x 36 cm)
13332	Standard 14" x 27" (36 cm x 68 cm)
13343	Shoulder/Neck 23" x 20" (58 cm x 50 cm)

## THERATHERM® DIGITAL MOIST HEAT PACK



Theratherm Digital Moist Heat Pack provides intense, moist heat therapy for temporary relief from pain associated with arthritis, back pain caused by muscle spasm and inflammation caused by strain and tension. Its flannel cover draws moisture from the air and as the pack's temperature rises, the moisture is released onto the skin, speeding relief to your patients. Its unique, easy-to-use digital control is convenient and safe.



## FEATURES

## Theratherm Digital Moist Heat Pack

- Convenient preset time and temperature
- Digital hand control for operator to program and monitor temperature and/or treatment time
- Lockout mode protects against inconsistent setting changes throughout treatment sessions
- Moist heat therapy temperatures range from 88° to 166° F (32° to 74° C)
- Treatment sessions up to 60 minutes
- Specially designed flannel cover for efficient moisture retention
- UL and cUL classified
- One-year warranty



## FLUIDOTHERAPY® – DRY HEAT THERAPY



Dry Heat Therapy relieves localized pain, increases localized blood circulation, increases range of motion and relieves minor pain and stiffness associated with nonrheumatoid arthritis.



### FEATURES

#### Fluidotherapy – the only name in dry heat therapy

##### Key effects of Fluidotherapy:

- **Dry Heat** increases circulation and assists in pain relief.
- **Massaging Action** provides warmth, increased local circulation and pain diminishing effects to help increase joint range of motion.
- **Limb Buoyancy** allows patients to freely perform a variety of resistive exercises with increased mobility and decreased pain. The fluidized Cellex particles act like a low viscosity fluid, allowing limbs to be suspended much like in a liquid state.

Fluidotherapy FLU110D



Single extremity unit for treating the hand, wrist, elbow, foot, and ankle. Includes thirty pounds (13.6 kg) of Cellex Dry Heat Media.

Fluidotherapy FLU115D



Double extremity unit for treating up to two hands or two feet independently or simultaneously. Includes 40 lb (18.1 kg) of Cellex Dry Heat Media.

### ORDERING INFORMATION

Part Number	Description
<b>Single Extremity Fluidotherapy</b>	
FLU110D	Fluido 110 V
<b>Double Extremity Fluidotherapy</b>	
FLU115D	Fluido 110 V
<b>Standard Accessories</b>	
MED0001	Cellex Dry Heat Media 10 lb (4.5 kg) container

### TECHNICAL SPECIFICATIONS

#### Single and Double Extremity Fluidotherapy

Mains Power:	120V, 50/60 Hz, 11A
Variable Adjustments:	Time, Temperature and Air Speed
Modes of Operation:	Continuous and Pulsed
Pulse Mode:	OFF to 6 seconds; ON/OFF
Treatment Time:	1-99 minutes and Continuous
Operating Temperature:	Dependent on air speed selected. 88° F to 125° F (31° C to 52° C)

Optimal Performance: 115° to 118° F at 50% air speed

Air Speed: 0% - 100% (5% increments)

Safety Class: Class 1, Type B

Safety Tests: UL 60601-1

#### Single Extremity Fluidotherapy

Weight: 70 lbs (31.7 kg) without  
Cellex Dry Heat Media

Shipping Weight: 100 lb (45.4 kg)

Dimensions\*: 34" x 11.5" x 33"  
(86.4 cm x 29.2 cm x 83.8 cm)

#### Double Extremity Fluidotherapy

Weight: 60 lbs (27.2 kg) without  
Cellex Dry Heat Media

Shipping Weight: 100 lb (45.4 kg)

Dimensions\*: 34" x 18.5" x 33"  
(86.4 cm x 47 cm x 83.8 cm)

\* Dimensions are listed as L x W x H

NOTE: Dedicated circuit is recommended.



#### Cellex® Dry Heat Media

Cellex Dry Heat Media is designed for use in all Chattanooga Fluidotherapy units. Made of natural cellulose, Cellex is environmentally safe and easily disposable. Cellex Dry Heat Media is packaged in a 10-lb (4.5 kg) container making it easy to store and handle.

NOTE: Problems will occur in a Fluidotherapy system if a non-processed media is used.

**ORDERING INFORMATION****Part Number Description****Chattanooga Standard Exerciser™**

18000 Retail Package

18020 Standard Package

**Chattanooga Deluxe Exerciser™**

18010 Chattanooga Deluxe Exerciser

**CHATTANOOGA STANDARD EXERCISER™****FEATURES****Chattanooga Standard Exerciser™**

- Perfect for upper or lower body exercises
- Extremely lightweight: 5 lbs (2.25 kg)
- Wrap-around foot straps
- Comes completely assembled
- Retail packaging available
- Shipping dimensions: 9" x 16" x 45"
- One-year warranty

**CHATTANOOGA DELUXE EXERCISER™****FEATURES****Chattanooga Deluxe Exerciser™**

- Use for upper or lower body exercises at home, the office or while traveling
- Automatic, multi-functional electronic monitor with LED display
- Quick dial with variable resistance
- Shipping dimensions: 18" x 9.5" x 12.5"
- One-year warranty



## SAUNDERS® TOTAL BACK EXERCISE DEVICE



### FEATURES

#### Saunders® Total Back Exercise Device

Promotes strength, flexibility and coordination through controlled functional exercise

- Pelvic stabilization allows isolation of lumbar extensors in any range of motion
- Adjusting the table angle allows either mild exercise in a pain free range, or intense conditioning
- Angle of recline reduces the load on the spine by 20% to 30% and allows for spinal positions in flexion, extension or side lying through a range of 110°
- "Side sit-ups" are great for the obliques, transversus abdominus, multifundus, quadratus lumborum and erector spinae muscles
- Includes a full-size wall poster feature 28 numbered exercises



### ORDERING INFORMATION

Part Number	Description
-------------	-------------

#### Saunders® Total Back Exercise Device

640243-002	Saunders Total Back Black
640243-003	Saunders Total Back Imperial Blue
640243-005	Saunders Total Back Burgundy

## STABILIZER™

### FEATURES

#### Stabilizer™

The Stabilizer Pressure Biofeedback unit gives valuable information to ensure quality and precision in exercise performance and muscle testing. It allows the clinician and patient to determine if the patient is able to selectively isolate and maintain contractions of the cervical or lumbopelvic core stabilization muscles. The measuring range is 0-200 mmHg analog pressure with an accuracy of +/- 3 mmHg pressure.

- Use with exercise to protect and stabilize joints
- Prevent and treat low back and neck pain
- Improve core stabilization exercises
- Enhance home exercise programs



### ORDERING INFORMATION

Part Number	Description
-------------	-------------

#### Stabilizer™

9296	Stabilizer Unit Comes complete with display gauge, pressure cell and instruction booklet.
------	--

## ORDERING INFORMATION

## Part Number Description

## Chattanooga Hydraulic Hand Dynamometer

43050 200 lbs (91 kg) Dial Gauge

## Hydraulic Pinch Gauge

43051 50 lbs (23 kg) Dial Gauge

## Skinfold Caliper

43063 Skinfold Caliper

## CHATTANOOGA HYDRAULIC HAND DYNAMOMETER



## FEATURES

## Hydraulic Hand Dynamometer

- Gives accurate grip strength readings without being able to “feel” the handle move
- Dynamometer accurately measures grip strength up to 200 lbs (91 kg)
- Adjustable handle with five positions to accommodate any hand size
- Indicator remains at maximum reading until reset
- Includes carrying case
- One-year warranty



## CHATTANOOGA HYDRAULIC PINCH GAUGE

## FEATURES

## Hydraulic Pinch Gauge

- Reliably and accurately register up to 50 lbs (23 kg)
- Measures all pinch tests (tip, key and palmer)
- Includes carrying case
- One-year warranty



## CHATTANOOGA SKINFOLD CALIPER

## FEATURES

## Skinfold Caliper

- Constructed of lightweight aluminum
- Large 60 mm scale is easy to read
- Measurement arms and “floating tips” assure accurate and reproducible results
- Includes carrying case, booklet and body fat percentage tables
- One-year warranty



## CHATTANOOGA 3-PIECE HAND EVALUATION SET



### FEATURES

#### 3-Piece Hand Evaluation Set

- All the instruments necessary to evaluate strength and range of motion for hand and finger
- Hydraulic Hand Dynamometer
- Hydraulic Pinch Gauge (50 lbs/23 kg)
- Stainless steel finger goniometer
- Includes carrying case
- One-year warranty



## CHATTANOOGA 7-PIECE HAND EVALUATION SET

### FEATURES

#### 7-Piece Hand Evaluation Set

- All instruments needed for a full hand evaluation
- Hydraulic Hand Dynamometer
- 2-point discriminator with 3rd point
- Finger circumference gauge and finger motion gauge
- Wartenburg pinwheel
- Hydraulic pinch gauge registering up to 50 lbs (23 kg)
- Stainless steel goniometer
- Digit-O-Meter
- Includes carrying case
- One-year warranty




## CHATTANOOGA GONIOMETERS

### FEATURES

#### Goniometers

- The plastic 360° scale ISOM (International Standards of Measurement) goniometers (A, B, C) have a head with three scales calibrated for use with the ISOM system. The clear plastic permits observation of the joint's axis of motion and range of motion
- The stainless steel finger and small joint goniometer (E) measures the range of motion of finger joints (metacarpophalangeal and interphalangeal) and other small joints including toes. The head has two opposing 180° scales, which are marked in 5° increments
- All goniometers have a linear scale in inches and centimeters
- One-year warranty



### ORDERING INFORMATION

Part Number	Description
<b>3-Piece Hand Evaluation Set</b>	
43055	200 lbs (91 kg) Dial Gauge Dynamometer
<b>7-Piece Hand Evaluation Set</b>	
43054	200 lbs (91 kg) Dial Gauge Dynamometer 30 lbs (14 kg) Pinch Gauge
<b>Goniometers</b>	
A 43058	Plastic 360° ISOM 12" (30 cm)
B 43059	Plastic 360° ISOM 8" (20 cm)
C 43060	Plastic 360° ISOM 6" (15 cm)
D 43057	Stainless Steel Finger and Small Joint 6" (15 cm)

**ORDERING INFORMATION****Part Number Description****Digital Hand Dynamometer**

43100 300 lbs (136 kg) Digital Gauge

**Digital Pinch Gauge**

43101 100 lbs (45 kg) Digital Gauge

**Inclinometer**43062 Digital Inclinometer  
Includes carrying case.

43061 Bubble Inclinometer

**Goniometers**

A 43119 Stainless Steel 180° Digit 4" (10 cm)

B 43118 Plastic 180° Digit 4" (10 cm)

**BASELINE® DIGITAL HAND DYNAMOMETER****FEATURES****Digital Hand Dynamometer**

- Patient will not feel the handle move
- Dynamometer accurately measures grip strength up to 300 lbs (136 kg)
- Handle adjusts to accommodate any hand size
- Indicator remains at maximum reading until reset
- Includes carrying case
- One-year warranty

**BASELINE® DIGITAL PINCH GAUGE****FEATURES****Digital Pinch Gauge**


- Reliably and accurately register up to 100 lbs (45 kg)
- Measures all pinch tests (tip, key and palmer)
- Includes carrying case
- One-year warranty

**BASELINE® INCLINOMETER****FEATURES****Inclinometer**

- Allows total range of motion to be accurately measured
- Digital model has convenient LED display and carrying case
- Our most popular model, the Bubble Inclinometer is simple and easy to use
- One-year warranty

**BASELINE® GONIOMETERS****FEATURES****Goniometers**

- The plastic and stainless steel 180° scale digit goniometers (F and G) require the use of only one hand to measure 110° flexion and 40° hyperextension in 5° increments
- All goniometers have a linear scale in inches and centimeters
- One-year warranty





## 3-PIECE BASELINE® DIGITAL HAND EVALUATION SET



### FEATURES

#### 3-Piece Digital Hand Evaluation Set

- All the instruments necessary to evaluate strength and range of motion for hand and finger
- Digital Hand Dynamometer,
- Digital Pinch Gauge (50 lbs/23 kg)
- Stainless steel finger goniometer
- Includes carrying case
- One-year warranty



## 7-PIECE BASELINE® DIGITAL HAND EVALUATION SET

### FEATURES

#### 7-Piece Digital Hand Evaluation Set

- All instruments needed for a full hand evaluation
- Digital Hand Dynamometer,
- 2-point discriminator with 3rd point
- Finger circumference gauge and finger motion gauge
- Wartenburg pinwheel
- Mechanical pinch gauge registering up to 60 lbs (27 kg)
- Stainless steel goniometer
- Digit-O-Meter
- One-year warranty



## BASELINE® ORTHOPEDIC AND SPORTS MEDICINE MECHANICAL PINCH GAUGE

### FEATURES

#### Orthopedic and Sports Medicine Mechanical Pinch Gauge

- Measures tip, key and palmer pinch strength of healthy hands
- Results are consistent with the published Mathiowetz studies
- Indicator remains at the subject's maximum reading until reset
- Strength readings are in pounds and kilograms up to 60 lbs (27 kg)
- Lightweight anodized aluminum with wrist strap
- Includes carrying case
- One-year warranty



### ORDERING INFORMATION

Part Number	Description
<b>3-Piece Baseline Digital Hand Evaluation Set</b>	
43102	300 lbs (136 kg) Digital Gauge Dynamometer
<b>7-Piece Baseline Digital Hand Evaluation Set</b>	
43103	300 lbs (136 kg) Digital Gauge Dynamometer 60 lbs (27 kg) Pinch Gauge
<b>Baseline Orthopedic and Sports Medicine Mechanical Pinch Gauge</b>	
43052	30 lbs (14 kg) Dial Gauge (Blue)
43053	60 lbs (27 kg) Dial Gauge (Red)

## ORDERING INFORMATION

Part Number	Size	Chest	Bicep
<b>Sully® Shoulder Support</b>			
660087-001	X-Small	34"-40"	9"-10½"
660087-002	Small	34"-40"	11"-13"
660087-003	Medium	38"-42"	13"-15"
660087-004	Large	42"-46"	14"-17"
660087-005	X-Large	45"-51"	15"-17"

## SULLY® SHOULDER SUPPORT




## FEATURES

**Sully® Shoulder Support**


On the outside, elastic straps attach with hook to the neoprene at any point, in any direction, and with any amount of force you choose. You can functionally stabilize, assist or restrict movement according to the specific needs of each athlete.

On the inside, the perforated, breathable, rubber-like neoprene grips the skin so the support moves with the athlete. It's as if the stabilization straps are attached directly to the body.


The Sully comes complete with a highly adaptable, bifurcated neoprene strap, and two additional 4" neoprene straps. Worn by professional and amateur athletes in all sports! Measure biceps in a non-flexed position.

**Multiple Applications**


Posterior Instabilities




Muscle Strains (pectoral etc.)



Rotator Cuff Deceleration



First Strap




Second Strap

Anterior Instabilities




First Strap




Second Strap

Inferior Instabilities



First Strap



Second Strap

Multi-Directional Instabilities

## STAINLESS STEEL CARTS



### FEATURES

#### Stainless Steel Carts

Versatile, functional and mobile, these sturdy carts go wherever you need them.



### ORDERING INFORMATION

Part Number	Description
-------------	-------------

#### Stainless Steel Carts

<b>A 4002</b>	Model MB Specially designed to hold the SS, E-2 and E-1 heating units 20" H x 19" W x 17" D (51 cm x 48 cm x 43 cm)
<b>B 4004</b>	Model MB-T Accommodates stationary heating, ultrasound and electrotherapy units. 30" H x 19" W x 17" D (76 cm x 48 cm x 43 cm)
<b>C 4018</b>	Model MB-TD Accommodates stationary heating and electrotherapy units. Drawer provides convenient storage for clinic supplies. 30" H x 19" W x 17" D (76 cm x 48 cm x 43 cm)

#### Pneumatic Stools with back

77064	Imperial Blue
77065	Burgundy
77066	Dove Gray
77067	Black

#### Pneumatic Stools without back

77060	Imperial Blue
77061	Burgundy
77062	Dove Gray
77063	Black

## PNEUMATIC THERAPY STOOLS

### FEATURES

#### Pneumatic Therapy Stools

A high-quality, low-priced therapy stool with round, 3" (8 cm) thick foam upholstered swivel seat. Available with or without back in the four colors listed below. A convenient gas spring makes raising and lowering the stool simple 18" - 22" (46 cm - 56 cm). Maximum weight 350 lbs (158 kg).

#### Available Colors

	Black
	Burgundy
	Dove
	Imperial Blue

*This color chart is produced within the limitations of lithographic color printing and is intended as a general guide only.*



## ORDERING INFORMATION

## Part Number Description

**A Ultra-Myossage® Lotion**

4262 1 Gallon (3.8 liter) Bottle

4260 8.5 oz. (250 ml) Bottle

**B Myossage® Lotion**

4210 1 Gallon (3.8 liter) Bottle

4208 8.5 oz. (250 ml) Bottle (12/case)

**C Utility Straps**

631374-001 Utility Strap 48"

631325-001 Utility Strap 72"

630483-001 Utility Strap 96"

**D Nylatex Assortment (1 pack of 9 rolls)**

1200 Nylatex Assortment

**D Nylatex Individual (3 rolls per package)**

1202 \*2.5" x 18" (6 cm x 46 cm)

1264 2.5" x 24" (6 cm x 61 cm)

1204 \*2.5" x 36" (6 cm x 91 cm)

1206 \*2.5" x 48" (6 cm x 122 cm)

1208 \*4" x 18" (10 cm x 46 cm)

1210 \*4" x 36" (10 cm x 91 cm)

1212 \*4" x 48" (10 cm x 122 cm)

1213 4" x 60" (10 cm x 152 cm)

1214 \*6" x 18" (15 cm x 46 cm)

1216 \*6" x 36" (15 cm x 91 cm)

1218 \*6" x 48" (15 cm x 122 cm)

1219 6" x 60" (15 cm x 152 cm)

\* Included in Nylatex Assortment Pack (1200)

## GELS AND LOTIONS



## FEATURES

**A Ultra-Myossage® Lotion**

Formulated with aloe vera, this top quality lotion is ideal for massage or as an ultrasound coupling agent. No menthol, hexachlorophene or polysorbates. Hypoallergenic, water-soluble and greaseless.

**B Myossage® Lotion**

Soothing, effective alcohol-free, non-greasy and non-staining lotion.

**C Utility Straps**

A versatile black webbed strap, available in three lengths. Use as a transfer belt, an extension for the EZ-ON traction harness belts, or as shown to assist manual cervical traction.

**D Nylatex Wraps**

- Ideal for holding cold packs in place
- Several sizes available
- Reusable

## EDUCATION MANUALS



Manuals to help educate about back pain and daily activities for preventative measures.



### FEATURES

#### Managing Back Pain

By H. Duane Saunders, MS PT & Michael S. Melnik, MS, OTR

A photo-filled guide to daily activities for back pain sufferers. Helps your patients in acute pain to safely participate in essential activities such as dressing, child care, kitchen activities, and more. 31 pages.

#### Your Healthy Back

By Richard A Pollock, CSP, H. Duane Saunders, MS PT and Michael S. Melnik, MS OTR

The complete guide for back injury prevention, body mechanics, work station design and claims management for supervisors. 24 pages.

#### For Your Back

By H. Duane Saunders, MS PT

The most comprehensive of our manuals. It thoroughly covers back facts, simple anatomy, and causes and treatments for the most common types of back problems. 48 pages.

#### For Your Neck

By H. Duane Saunders, MS PT

Discusses common causes of neck injuries and explains preventive strategies. 48 pages.

#### Save Your Back

By H. Duane Saunders, MS PT and Michael S Melnik, MS OTR

The photographic guide to patient handling techniques. Essential for health care provider training. 48 pages.

### ORDERING INFORMATION

Part Number	Description
-------------	-------------

#### Managing Back Pain

630309-001	Managing Back Pain
------------	--------------------

#### Your Healthy Back

621060-001	Your Healthy Back
------------	-------------------

#### For Your Back

630308-001	For Your Back
------------	---------------

#### For Your Neck

630310-001	For Your Neck
------------	---------------

#### Save Your Back

621059-001	Save Your Back
------------	----------------

## ORDERING INFORMATION

## Part Number Description

**DURA-STICK® Premium Self-Adhesive Electrodes**

42170	Dura-Stick Premium 1.25" Round
42171	Dura-Stick Premium 2" Round
42172	Dura-Stick Premium 1.5" x 2.5" Oval
42173	Dura-Stick Premium 3" x 5" Oval (20/case) (10 packs of 2)
42174	Dura-Stick Premium 2" Square
42175	Dura-Stick Premium 2" x 3.5" Rectangle
42176	Dura-Stick Premium 2" Square <b>Blue Gel</b>
42177	Dura-Stick Premium 1.5" x 3.5" Rectangle <b>Blue Gel</b>

**DURA-STICK Self-Adhesive Electrodes**

42160	Dura-Stick 2" x 4" Oval
42178	Dura-Stick 2" Square
42179	Dura-Stick 2" x 3.5" Rectangle

**DURA-STICK Disposable Electrodes**

42180	Dura-Stick 2.25" x 2.5" Disposable 40/Package
-------	--

NOTE: ALL QUANTITIES ARE 10 PACKS OF 4 EACH,  
40 PER CASE, UNLESS OTHERWISE STATED

**DURA-STICK® PREMIUM SELF-ADHESIVE ELECTRODES**

## FEATURES

**DURA-STICK® Premium Self-Adhesive Electrodes**

- 10 packs of 4 each, 40/case unless otherwise noted
- Cloth Backing
- Stainless steel mesh
- Choice wire and molded connectors
- Provides optimal current dispersion with low impedance

**DURA-STICK SELF-ADHESIVE ELECTRODES**

## FEATURES

**DURA-STICK Self-Adhesive Electrodes**

- 10 packs of 4 each, 40/case
- Carbon silver film layer incorporated into a low profile design
- Cloth backing for increased "stretchability" and memory
- Innovative polymer for improved adhesion and durability

**DURA-STICK DISPOSABLE ELECTRODES**

## FEATURES

**DURA-STICK Disposable Electrodes**

- 40/Package
- Multi-day" electrode for extended use
- May be used for 2-5 days
- Solid gel center surrounded by adhesive
- Beige mesh backing, looks like a bandage
- Dispose of when finished



## DURA-STICK® PLUS SELF-ADHESIVE ELECTRODES



### FEATURES

#### Dura-Stick® Plus Self-Adhesive Electrodes

- 10 packs of 4 each, 40/case
- Our most economical electrode.
- Permagem™ Hydrogel for superior conformability to body contours
- Blue foam backing
- Carbon conductive film incorporated into a flexible spun lace backing
- Flexible and sturdy design for durability and ease of use



## DURA-STICK sEMG ELECTRODE

### FEATURES

#### DURA-STICK sEMG Electrode

- Round white
- Foam Backing
- 100 pk/individual



## CONDUCTIVE RUBBER ELECTRODES

### FEATURES

#### Conductive Rubber Electrodes

- Conductive rubber electrodes
- 3" (8 cm) round
- Available in black or red



### ORDERING INFORMATION

Part Number	Description
<b>Dura-Stick® Plus Self-Adhesive Electrodes</b>	
42045	Dura-Stick Plus 2.75" x 5" Rectangle
42159	Dura-Stick Plus 1.5" x 2.5" Oval
42181	Dura-Stick Plus 2" x 3.5" Rectangle
42182	Dura-Stick Plus 2" Round
42183	Dura-Stick Plus 2" Square
42160	Dura-Stick Plus 2" x 4" Oval
<b>DURA-STICK sEMG Electrode</b>	
42185	Dura-Stick sEMG Round White (100 pk/individual)
<b>Conductive Rubber Electrodes</b>	
72852	Black Rubber Carbon 3" (8 cm) Round (2X)
72853	Red Rubber Carbon 3" (8 cm) Round (2X)

NOTE: ALL QUANTITIES ARE 10 PACKS OF 4 EACH, 40 PER CASE, UNLESS OTHERWISE STATED



**DJO, LLC | A DJO Global Company**  
**T** 800.592.7329   **D** 760.727.1280   **F** 800.936.6569  
1430 Decision Street | Vista, CA 92081-8553 | U.S.A.  
[DJOglobal.com/chattanooga](http://DJOglobal.com/chattanooga)

*Together in Motion™*