

MOVING REHABILITATION FORWARD**

Physical Agent Modalities | 2013





INTRODUCTION



Chattanooga is the world's largest manufacturer of rehabilitation equipment for treating musculoskeletal, neurological and soft tissue disorders. For over six decades, Chattanooga has set a benchmark for leadership, reliability and excellence. Our products contribute to better treatment outcomes in hospitals, clinics and home settings worldwide. We lead the rehabilitation industry by example, and through continuous innovation, provide real world solutions for clinicians and their patients.

As the world's leading brand in rehabilitation, Chattanooga is sold in over 80 countries worldwide. Chattanooga distributors and customers are supported by a highly trained and seasoned team of personnel who oversee new product development, product training and education, order placement, product delivery, and after sales service in addition to providing dedicated sales and marketing support.

Chattanooga is a brand of DJO Global Inc. DJO Global provides solutions for musculoskeletal health, vascular health and pain management. Our products help prevent injuries or rehabilitate after surgery, injury or from degenerative disease, enabling patients to regain or maintain their natural motion.

Visit www.DJOglobal.com/chattanooga



















PHYSICAL AGENT MODALITIES

For over 25 years, Chattanooga has been the world leader in electrotherapy devices and has gradually expanded its product offering to include other essential physical agent modalities, such as ultrasound, shortwave and laser therapy.

Chattanooga developed these treatment options under the well established Intelect product range, appreciated by therapists for its ease of use and functional design.

Chattanooga is committed to continually working with customers to develop smart solutions that empower therapists to provide the best treatment to their patients.



chattanooga

Salling Salling
Total ferential of onnuru Meroparrank India ferential of
The state of the s
ET finding For Company Control of
00

CONTENTS Introduction Intelect Legend XT 3 Intelect TranSport® Ultrasound 5 Intelect TranSport® Electrotherapy 6 Intelect TranSport® Combo Intelect® Legend 8 Intelect® Portable Electrotherapy 9 **Gels and Lotions** 10 **Electrodes & Accessories** 11

Part Number

Description

ORDERING INFORMATION

Intelect® Legend XT System			
2763	2-Channel Electrotherapy		

2786 4-Channel Electrotherapy
2760 2-Channel Electrotherapy

2788 4-Channel Electrotherapy

Intelect Legend XT System with Cart

2793 2-Channel Electrotherapy
2797 4-Channel Electrotherapy

2791 2-Channel Combination

2795 4-Channel Combination

Standard Accessories

27383	5 cm ² Sound Head Applicator
27312	Electrotherapy Leadwire Channel 1
27313	Electrotherapy Leadwire Channel 2
27314	Electrotherapy Leadwire Channel 3
27215	Electrotherapy Leadwire Channel A

2/313	Electronierapy Leadwire Charmer 4	
Optional Accessories		
2781	Channel 3 and 4 Electrotherapy Module	
2780ASY	Therapy Cart (assembled)	
2767	Battery Module (Self-contained Nickel Metal Hydride battery automatically and continually charges when unit is plugged into mains power supply.)	
27469	Patient Interrupt Switch Channels 1 and 2	
27470	Patient Interrupt Switch Channels 3 and 4	
72852	Black Rubber Carbon 3" (8 cm) Electrodes (2x)	
72853	Red Rubber Carbon 3" (8 cm) Electrodes (2x)	
10648	Nylatex® Wrap, 2.5" x 24" (2x)	
27465	Patient Data Cards (25/pack)	
27508	Operator Remote Channels 1 and 2	
27079	Operator Remote Channels 3 and 4	
27381	1 cm² Sound Head Applicator	
27382	2 cm² Sound Head Applicator	
27384	10 cm ² Sound Head Applicator	
27133	Intelect Legend XT Carry Bag, Large	

INTELECT® LEGEND XT



The Intelect® Legend XT is one of the most advanced therapy systems available. Beneath its award-winning, ergonomic design is a sophisticated, versatile, user-friendly system. The Intelect Legend XT lets you offer a comprehensive range of therapies in a single, state-of-the-art unit.



FEATURES

Intelect® Legend XT

- Patented award-winning, ergonomic modular design
- Easy-to-add modularity allows the addition of two additional channels of electrotherapy or a battery module
- High contrast 5" (13 cm) FSTN LCD monochromatic user interface
- Two independent electrotherapy channels, expandable to four independent channels
- Multiple waveform electrotherapy-six clinical waveforms: Interferential, Premodulated, High Volt, Russian, Symmetrical Biphasic and Microcurrent
- Documentation of treatment data with Patient Data Cards
- User defined protocols for your specific needs

- Dual frequency ultrasound at 1 or 3.3 MHz
- Pulsed or Continuous Duty Cycles (10%, 20%, 50% & 100%)
- Four ergonomically designed soundhead applicator sizes: 1 cm², 2 cm², 5 cm² & 10 cm²
- Watertight soundheads for use in underwater therapy
- Therapy System Cart maximizes storage of clinical supplies and organizes cables
- Optional Operator Remotes allow treatment stop, start and pause, as well as change in amplitude at any time during a treatment session. Each remote is color-coded for designation between channels
- Two-year warranty for all clinical electrotherapy



Intelect Legend XT Carry Bag

Ruggedly constructed to easily transport the Legend XT and all accessories required for therapy on-the-go.



Modularity

Flexible modular design easily accommodates either two additional channels or a battery module.

chattanooga





User Defined Protocols

Customize waveforms to your specific needs and organize alphabetically by patient name for easy recall.



Patient Documentation

Patient Data Cards allow clinicians to capture and record accurate patient information using programs such as pain map profile, modality records and electrode placement diagrams.



Multiple Waveforms

Six essential clinical waveforms to manage a broad spectrum of clinical conditions: Interferential, Premodulated, High Volt, Russian, Microcurrent and Symmetrical Biphasic.

TECHNICAL SPECIFICATIONS

Intelect Legend XT	
Mains Power:	100-240 VAC~50/60 Hz, 1.0 A
Power Consumption:	1000 W
Weight:	7 lb (3.2 kg) /Shipping - 16 lb (7.3 kg)
Dimensions*:	12.75" x 11.4" x 8.75" (32.4 x 28.9 x 22.2 cm)
Electrical Safety Class:	Class 1, Electrotherapy Type BF, Ultrasound Type B
Safety Tests:	UL/IEC/EN 60601-1, IEC/ EN 60601-1-2, IEC 60601-2-5, IEC 60601-2-10
* Dimensions are listed	as L x W x H

WAVEFORMS



INTERFERENTIAL

Continuous

• Manual Vector Scan



PREMODULATED

Continuous

• Cycle Timed



HIGH VOLT

• Twin peak



RUSSIAN

• Continuous

• Cycle Timed

MM MM

MICROCURRENT

Polarity Based



SYMMETRICAL BIPHASIC

• Continuous

Cycle Timed

Patented Technology

The Intelect Legend XT Therapy System is protected by one or more of the following patents:

D503,807 D511,008 D497,462 D513,615 D506,260 D504,593 D505,726 D508,998

Additional Patents Pending

ORDERING INFORMATION

OKDEKING	INFORMATION
Part Number	Description
Intelect TranS _I	oort® Ultrasound
2782	Intelect Transport Ultrasound
Standard Acces	ssories
27383	5 cm² Sound Head Applicator
4248	ConductorTM Ultrasound Gel 8.5 oz (250 ml) Bottle (1x)
Optional Acces	sories
27478	Intelect Transport Nickel Metal Hydride Battery Pack
27381	1 cm² Sound Head Applicator
27382	2 cm² Sound Head Applicator
27384	10 cm² Sound Head Applicator
27467	Intelect Transport Carry Bag
2780ASY	Therapy Cart (assembled)
2884	Cart Adapter (Gray) Allows the Intelect Transport unit to attach securely to the Therapy Cart.

Note: The Cart Adapter, part number 2884, is necessary when ordering a cart for a Transport Unit.

TECHNICAL SPECIFICATIONS

Intelect Legend XT	
Mains Power:	100-240 VAC~50/60 Hz, 1.0 A
Weight:	5 lb (2.3 kg) /Shipping - 9 lb (4 kg)
Dimensions*:	12.9" x 11.3" x 6.4" (32.8 x 28.8 x 16.3 cm)
Electrical Safety Class:	Class 1, Type B
Safety Tests:	UL/IEC/EN 60601-1, IEC/EN 60601-1-2, IEC 60601-2-5

^{*} Dimensions are listed as L x W x H

INTELECT TRANSPORT® ULTRASOUND





FEATURES

Intelect TranSport® Ultrasound

- Fully functional 1 and 3.3 MHz frequencies
- Ten user-defined memory positions for User Protocols
- Pulsed and continuous therapy operation (10%, 20%, 50% and 100%)
- Ergonomically designed Sound Head Applicators in variety of sizes (1 cm², 2 cm², 5 cm² and 10 cm²)
- Sound heads are watertight for use in underwater therapy
- All Sound Head Applicators are interchangeable due to Electronic Signature™ residing on all transducers
- Head-warming feature on all Sound Head Applicators
- Unit uniquely designed for table mount, therapy cart or mobile use
- 110V or battery powered option
- Optional Therapy Cart integrates with the unit to maximize storage of clinical supplies
- Two-year warranty



Cart Adapter



Optional Battery Pack



INTELECT TRANSPORT® ELECTROTHERAPY





FEATURES

Intelect TranSport® Electrotherapy

- Two channels of electrical stimulation output
- Logical control system and a large easy-to-read graphical LCD display
- Four standard waveforms on board: Interferential, Premodulated, High Volt, and Russian
- Ten user-defined memory positions for User Protocols
- Independent intensity and parameter controls for each channel
- Unit uniquely designed for tabletop, wall mount, therapy cart or mobile use
- 110V or battery-powered option
- Optional Therapy Cart integrates with the unit to maximize storage of clinical supplies
- Two-year warranty





The custom designed carrying bag (optional) is lightweight. Protect and transport the unit and accessories for therapy on the go.

ORDERING INFORMATION

Intelect TranSport® Electrotherapy

2783 Intelect Transport®
2-Channel Electrotherapy

Standard Accessories

Part Number

27312 Electrotherapy Leadwire Channel 1
27313 Electrotherapy Leadwire Channel 2

Description

Optional Accessories

27478	Intelect Transport Nickel Metal Hydride Battery Pack
72852	Black Rubber Carbon 3" (8 cm) Electrodes (2x)
72853	Red Rubber Carbon 3" (8 cm) Electrodes (2x)
10648	Nylatex® Wrap 2.5" x 24" (2x)
27467	Intelect Transport Carry Bag
2780ASY	Therapy Cart (assembled)
2884	Cart Adapter (Gray) Allows the Intelect Transport unit to attach securely to the Therapy Cart.

Note: The Cart Adapter, part number 2884, is necessary when ordering a cart for a Transport Unit.

WAVEFORMS



INTERFERENTIAL

- Continuous
- Manual Vector Scan



PREMODULATED

- Continuous
- Cycle Timed



HIGH VOLT

• Twin peak



RUSSIAN

ContinuousCycle Timed

TECHNICAL SPECIFICATIONS

Intelect Legend XT

Mains Power: 100-240 VAC-50/60 Hz, 1.0 A

Weight: 5 lb (2.3 kg) /Shipping - 9 lb (4 kg)

Dimensions*: 12.9" x 11.3" x 6.4" (32.8 x 28.8 x 16.3 cm)

Electrical Safety Class: Class 1, Type BF

Safety Tests: UL/IEC/EN 60601-1, IEC/EN 60601-1-2, IEC 60601-2-10

* Dimensions are listed as L x W x H

ORDERING INFORMATION

Part Number	Description	
Intelect TranSport® Combo		
2738	Intelect Transport Combo	
Standard Acce	ssories	
27383	5 cm ² Sound Head Applicator	
4248	ConductorTM Ultrasound Gel 8.5 oz (250 ml) Bottle (1x)	
27312	Electrotherapy Leadwire Channel 1	
27313	Electrotherapy Leadwire Channel 2	
Optional Acces	ssories	
27478	Intelect Transport Nickel Metal Hydride Battery Pack	

Optional Accessories	
27478	Intelect Transport Nickel Metal Hydride Battery Pack
72852	Black Rubber Carbon 3" (8 cm) Electrodes (2x)
72853	Red Rubber Carbon 3" (8 cm) Electrodes (2x)
10648	Nylatex® Wrap 2.5" x 24" (2x)
27381	1 cm² Sound Head Applicator
27382	2 cm ² Sound Head Applicator
27384	10 cm² Sound Head Applicator
27467	Intelect Transport Carry Bag
2780ASY	Therapy Cart (assembled)
2884	Cart Adapter (Gray) Allows the Intelect Transport unit to

Note: The Cart Adapter, part number 2884, is necessary when ordering a cart for a Transport Unit.

WAVEFORMS



INTERFERENTIAL

attach securely to the Therapy Cart.

- Continuous
- Manual Vector Scan



PREMODULATED

- Continuous
- Cycle Timed



• Twin peak



RUSSIAN

- Continuous
- Cycle Timed

TECHNICAL SPECIFICATIONS

Intelect Legend XT	
Mains Power:	100-240 VAC~50/60 Hz
Weight:	5 lb (2.5 kg) /Shipping - 16 lb (7.3 kg)
Dimensions*:	12.9" x 11.3" x 10.1" (32.8 x 28.8 x 25.7 cm)
Electrical Safety Class:	Class 1, Electrotherapy Type BF, Ultrasound Type B
Safety Tests:	UL/IEC/EN 60601-1, IEC/EN 60601-1-2, IEC 60601-2-5, IEC 60601-2-10
* Dimensions are listed	as L x W x H

INTELECT TRANSPORT® COMBO





FEATURES

Intelect TranSport® Combo

- Fully functional 1 and 3.3 MHz ultrasound frequencies
- Two channels of electrical stimulation output
- Ten user-defined memory positions for User Protocols
- Four standard waveforms on board: Interferential, Premodulated, High Volt, and Russian
- Pulsed and continuous therapy operation (10%, 20%, 50% and 100%)
- Independent intensity and parameter controls for each channel

 Ergonomically designed Sound Head Applicators in a variety of sizes (1 cm², 2 cm², 5 cm² and 10 cm²)

• All Sound Head Applicators are interchangeable with Electronic Signature™ residing on all transducers

- Watertight sound heads for use in underwater therapy
- Head-warming feature on all Sound Head Applicators
- Unit uniquely designed for tabletop, wall mount, therapy cart or mobile use
- 110V or battery powered option
- Optional Therapy Cart integrates with the unit to maximize storage of clinical supplies
- Two-year warranty



INTELECT® LEGEND





FEATURES

Intelect® Legend Combination, Ultrasound and Stimulator Units

Intelect® is a name the clinical rehabilitation community relies on. The Intelect Legend Stimulator, Combination and Ultrasound units are well known for their exceptional clinical performance, durability, and ease of use. Units carry a two-year warranty.

Ultrasound Applicators

Waterproof and precalibrated applicators operate with 1 or 3.3 MHz frequencies. Available in 2 cm², 5 cm² and 10 cm².





Stim Lead Wires

Designed with a rip stop feature and colorcoding on the connectors to easily distinguish between channels.



EMG Lead Wires

Designed with a rip stop feature and colorcoding on the connectors to easily distinguish between channels.

ORDERING INFORMATION

A Intelect L	egend Combination	
INTCB	Two Channel Combo	

Four Channel Combo

B Intelect Legend Ultrasound

7550

Part Number Description

INT001	Dual Frequency System

C Intelect Legend Stimulator

INTO02	Two Channel Stimulator	
7560	Four Channel Stimulator	

Ultrasound Applicators

78048	2 cm ² Applicator	
78047	5 cm ² Applicator	
78046	10 cm ² Applicator	

Stim Lead Wires

27312	Channel 1	
27313	Channel 2	
27314	Channel 3	
27315	Channel 4	

EMG Lead Wires

27321	Channel 1
27322	Channel 2

WAVEFORMS



INTERFERENTIAL

- Continuous
- Manual Vector Scan



PREMODULATED

- Continuous
- Cycle Timed



HIGH VOLT

• Twin peak



RUSSIAN

ContinuousCycle Timed

77595

ORDERING INFORMATION Part Number Description Intelect TENS **A** 77710 Intelect TENS Basic **B** 77600 Intelect TENS Standard **C** 77712 Intelect TENS Digital **Intelect NMES D** 77715 Intelect NMES Standard E 77717 Intelect NMES Digital Intelect IFC **F** 77719 Intelect IFC Portable Batteries (not shown) Rechargeable NiCd 9-volt Battery **Battery Recharger** Lead Wires (not shown) Lead Wire, TENS 360° Pivot, 77593 Unshielded Connector (Pair) Lead Wire, TENS 360° Pivot, 77594 Shielded Connector (Pair)

Lead Wire, TENS Non-Pivot,

Shielded Connector (Pair)

IENS Economy, T	ENS Standard & TENS Digital
Channel:	Dual, isolated
Waveform:	Asymmetric biphasic square pulse
Pulse Amplitude:	0-80 mA each channel adjustable (500 ohm load)
Pulse Width:	60-250µ sec. adjustable (Economy & Standard)
Pulse Width:	50-300µ sec. adjustable (Digital)
Pulse Frequency:	50-300µ sec. adjustable (Digital)
NMES Portable - S	Standard and Digital
Channel:	Dual, isolated
Waveform:	Asymmetric biphasic square pulse
Pulse Amplitude:	0-80 mA each channel adjustable (500 ohm load)
Pulse Width:	250μ sec. (Standard)
Pulse Width:	50-300µ sec. adjustable (Digital)
Pulse Frequency:	5, 30, 100 Hz (Standard)
Pulse Frequency:	2-120 Hz, adjustable (Digital)
IFC	
Channel:	Single for 2 pole or Dual for 4 pole
Waveform:	Symmetric biphasic square
Pulse Amplitude:	0~32 mA, adjustable (at 500 ohm load)
Pulse Width:	125μ sec./phase
Pulse Frequency:	4000+/-10% Hz fixed

INTELECT® PORTABLE ELECTROTHERAPY



Long known for innovative clinical electrotherapy, Chattanooga offers you a full line of portable electrotherapy products with Intelect TENS units, Intelect NMES (Neuromuscular Electrical Stimulation) and Intelect IFC (Interferential Current) units. Treat multiple patient needs from pain management to the prevention or retardation of disuse muscle atrophy. All devices come with lead wires, 9-volt battery, DURA-STICK® Self-Adhesive Electrodes and a convenient carrying case.













FEATURES

A-C Intelect TENS

- Dual Independent Channels
- Adjustable Timer on Standard and Digital units
- Adjustable Frequency and Pulse Duration parameters
- Basic and Standard units feature three modes: Burst, Normal and Modulation
- Digital unit features five modes: Burst, Normal, Modulated Rate and Width, Strength Duration and Bi-Pulse
- Digital unit has easy-to-read LCD display and convenient digital controls

D-E Intelect NMES

- Dual Independent Channels
- Adjustable Pulse Amplitude, Pulse Frequency and Contraction/Relaxation Time parameters
- Digital unit features
- Adjustable Pulse Width and Timer parameters
- Three stim modes: Synchronous, Constant and Alternate
- Easy-to-read LCD display and convenient digital controls

F Intelect IFC

- Dual Channel, Manual Selection for two or four Electrodes
- Adjustable Pulse Amplitude, Interference Frequency and Difference Frequency parameters
- Four stimulation modes (frequency shifts): Continuous, 1/1, 8/8 and 10/10 Ramped
- Includes 120V AC line cord adapter
- Sliding cover for controls and battery

GELS AND LOTIONS













FEATURES

A Intelect® Ultrasound Gel

Non-greasy, non-staining and hypoallergenic gel enhances conductivity.

B Ultrasound Transmission Gel

Water soluable, non-greasy ultrasound and coupling gel.

C Ultra-Myossage® Lotion

Formulated with aloe vera, this top quality lotion is ideal for massage or as an ultrasound coupling agent. No menthol, hexachlorophene or polysorbates. Hypoallergenic, water-soluble and greaseless.

D Myossage® Lotion

Soothing, effective alcohol-free, non-greasy and non-staining lotion.

E Nylatex Wraps

- Ideal for holding cold packs in place
- Several sizes available
- Reusable

ORDERING INFORMATION

Part Number	Description
A Intelect® Ult	trasound Gel
4266	1.3 Gallon (5 Liter) Container Includes one 8.5 oz. container for refill
B Conductor™	Transmission Gel
4238	1.3 Gallon (5 Liter) Container Includes one 8.5 oz. container for refill
4248	8.5 oz. (250 ml) Bottle (24/case)
C Ultra-Myoss	age® Lotion
4262	1 Gallon (3.8 liter) Bottle
4260	8.5 oz. (250 ml) Bottle
D Myossage® L	Lotion
4210	1 Gallon (3.8 liter) Bottle
4208	8.5 oz. (250 ml) Bottle (12/case)
E Nylatex Asso	ortment (1 pack of 9 rolls)
1200	Nylatex Assortment
E Nylatex Indi	vidual (3 rolls per package)
1202	*2.5" x 18" (6 cm x 46 cm)
1264	2.5" x 24" (6 cm x 61 cm)
1204	*2.5" x 36" (6 cm x 91 cm)
1206	*2.5" x 48" (6 cm x 122 cm)
1208	*4" x 18" (10 cm x 46 cm)
1210	*4" x 36" (10 cm x 91 cm)
1212	*4" x 48" (10 cm x 122 cm)
1213	4" x 60" (10 cm x 152 cm)
1214	*6" x 18" (15 cm x 46 cm)
1216	*6" x 36" (15 cm x 91 cm)
1218	*v6" x 48" (15 cm x 122 cm)

^{*} Included in Nylatex Assortment Pack (1200)

6" x 60" (15 cm x 152 cm)

1219

ORDERING INFORMATION

Part Number	Description	
DURA-STICK® Premium Self-Adhesive Electrodes		
42170	Dura-Stick Premium 1.25" Round	
42171	Dura-Stick Premium 2" Round	
42172	Dura-Stick Premium 1.5" x 2.5" Oval	
42173	Dura-Stick Premium 3" x 5" Oval (20/case) (10 packs of 2)	
42174	Dura-Stick Premium 2" Square	
42175	Dura-Stick Premium 2" x 3.5" Rectangle	
42176	Dura-Stick Premium 2" Square Blue Gel	
42177	Dura-Stick Premium 1.5" x 3.5" Rectangle Blue Gel	

DURA-STICK Self-Adhesive Electrodes

42160	Dura-Stick 2" x 4" Oval
42178	Dura-Stick 2" Square
42179	Dura-Stick 2" x 3.5" Rectangle

DURA-STICK Disposable Electrodes

42180	Dura-Stick 2.25" x 2.5" Disposable
42100	40/Package

NOTE: ALL QUANTITIES ARE 10 PACKS OF 4 EACH, 40 PER CASE, UNLESS OTHERWISE STATED

DURA-STICK® PREMIUM SELF-ADHESIVE ELECTRODES



FEATURES

DURA-STICK® Premium Self-Adhesive Electrodes

- 10 packs of 4 each, 40/case unless otherwise noted
- Cloth Backing
- Stainless steel mesh
- Choice wire and molded connectors
- Provides optimal current dispersion with low impedance





DURA-STICK SELF-ADHESIVE ELECTRODES

FEATURES

DURA-STICK Self-Adhesive Electrodes

- 10 packs of 4 each, 40/case
- Carbon silver film layer incorporated into a low profile design
- Cloth backing for increased "stretchability" and memory
- Innovative polymer for improved adhesion and durability





DURA-STICK DISPOSABLE ELECTRODES

FEATURES

DURA-STICK Disposable Electrodes

- 40/Package
- Multi-day" electrode for extended use
- May be used for 2-5 days
- Solid gel center surrounded by adhesive
- Beige mesh backing, looks like a bandage
- Dispose of when finished





DURA-STICK® PLUS SELF-ADHESIVE ELECTRODES



FEATURES

Dura-Stick® Plus Self-Adhesive Electrodes

- 10 packs of 4 each, 40/case
- Our most economical electrode.
- Permagel™ Hydrogel for superior conformability to body contours
- Blue foam backing
- Carbon conductive film incorporated into a flexible spun lace backing
- Flexible and sturdy design for durability and ease of use





DURA-STICK SEMG ELECTRODE

FEATURES

DURA-STICK sEMG Electrode

- Round white
- Foam Backing
- 100 pk/individual



CONDUCTIVE RUBBER ELECTRODES

FEATURES

Conductive Rubber Electrodes

- Conductive rubber electrodes
- 3" (8 cm) round
- Available in black or red



ORDERING INFORMATION

Part Number	Description
Dura-Stick® Plus Self-Adhesive Electrodes	
42045	Dura-Stick Plus 2.75" x 5" Rectangle
42159	Dura-Stick Plus 1.5" x 2.5" Oval
42181	Dura-Stick Plus 2" x 3.5" Rectangle
42182	Dura-Stick Plus 2" Round
42183	Dura-Stick Plus 2" Square
42160	Dura-Stick Plus 2" x 4" Oval

DURA-STICK SEMG Electrode

42185	Dura-Stick sEMG Round White
	(100 pk/individual)

Conductive Rubber Electrodes

72852	Black Rubber Carbon 3" (8 cm) Round (2X)
72853	Red Rubber Carbon 3" (8 cm) Round (2X)

NOTE: ALL QUANTITIES ARE 10 PACKS OF 4 EACH, 40 PER CASE, UNLESS OTHERWISE STATED





DJO, LLC | A DJO Global Company
T 800.592.7329 | D 760.727.1280 | F 800.936.6569 1430 Decision Street | Vista, CA 92081-8553 | U.S.A. DJOglobal.com/chattanooga

Together in Motion.

©2013 DJO, LLC 00-2456 Rev. A