

The Kelvi® Pro System is the safest product on the market as it is the only device safety rated for 24 hours of continuous use.







For every step of your healing journey there's Kelvi



Back Pro



Ankle Pro




Knee Pro



Shoulder Pro



For more information please visit us online



enovis...

Copyright © 2023 DJO, LLC | MKT00-13088 Rev A

Kelvi® and KelviTec® are registered trademarks of Hypothermia Devices, Inc. dba Kelvi

Individual results may vary. Neither DJO, LLC nor any of the Enovis companies dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.

WARNING: This device can be cold or hot enough to cause serious injury. Serious adverse reactions & safety hazards may occur when using this device. Read and understand all warnings and Instructions for Use, including contraindications, before using this device. Additional warnings appear in the Instructions for Use. Patients can only use this device with a prescription from a healthcare professional. The prescription must state how long and how often the device should be used and the length of breaks between uses.



SAFE. PRECISE. SMART. Cooling, Heating & Contrast Therapy



COLD, HEAT, AND CONTRAST TREATMENT IN ONE CONVENIENT, DIGITAL DEVICE.



KelviTec® technology measures the temperature of your skin, not the temperature of the device using digital thermometers in the wrap.

Continuous feedback from the wrap to the console so you get a therapeutic treatment throughout, even when you sleep.

No ice, water, or solution to be adjusted making it easy for you or your caregiver to use.



TREATMENT OPTIONS

Easy application in clinic, on the road, or at home.

CRYOTHERAPY

39°-59°F

THERMOTHERAPY

95°-108°F

CONTRAST THERAPY

Rapid comfortable transitions between cooling and heating

39°-108°F

COLD THERAPY

- · Decreases pain
- · May increase patient satisfaction during recovery
- · Clinically proven to decrease opioid consumption post-operatively

HEAT THERAPY

- · Increase blood flow and range of motion
- · Decrease pain

CONTRAST THERAPY

- Rapid change of temperature between cold and heat therapy provides a unique treatment option not available in other at-home devices.
- Contrast therapy has shown clinical improvements in pain and motion.

| Device | Cold Therapy | Heat Therapy | Rapid Contrast Therapy | Pressure or Compression | Technology | Software based with digital interface |
|-----------------------|--------------|--------------|---------------------------|----------------------------|---|---|
| Kelvi® | 39F-59F | 95F-108F | 39F-59F & 95F-108F | Up to 75mmHG | Wrap measures temperature of the skin | Yes; Pre-set option on the console with app for customizable treatments |
| Game Ready Pro 2.1 | 35F-50F | No | No | Up to 75mmHG | Ice and water | No |
| NICE | 42F-58F | No | No | Up to 50mmHG | Glycol and distilled water solution | No |
| Vascutherm 4 | 43F-50F | 105F | No | Up to 50mmHG | Alcohol and distilled water solution | No |