PROCARE



XCEL**T**RAX[™]
Walking Braces



An affordable and comprehensive family of Walking Braces, the XCELTRAX provides ideal support and comfort to help promote healing for patients with mild to severe foot and ankle injuries.



Superior Fit

The lightweight design and unique pre-shaped ergonomic frame helps modify the level of pressure on the injury at both the malleoli and lower leg.



Ideal Comfort and Sizing

The anatomical design of the struts accommodate a wide variety of leg sizes.

Integrated pump design to make the pneumatic liner simple to inflate and adjust.



Available in

Standard & Pneumatic

Tall and Ankle, XSmall to XLarge



Natural and Stable Gait

The XcelTrax walking brace edgeless sole design adjusts to each patient's unique gait pattern. This unique sole design, wider foot bed and low profile rocker sole help promote a more stable and natural gait.





PROCARE®

	SHOES	SHOE SIZE		
P/N	MENS	WOMENS	SIZES	REPLACEMENT LINER
79-95512	2-4	3.5-5.5	XSmall	79-95532
79-95513	4.5-7	6-8	Small	79-95533
79-95515	7.5-10.5	8.5-11.5	Medium	79-95535
79-95517	10.5-12.5	11.5-13.5	Large	79-95537
79-95518	12.5+	13.5+	XLarge	79-95538
RAX Air Ankle				
79-95522	2-4	3.5-5.5	XSmall	79-95542
79-95523	4.5-7	6-8	Small	79-95543
79-95525	7.5-10.5	8.5-11.5	Medium	79-95545
79-95527	10.5-12.5	11.5-13.5	Large	79-95547
79-95528	12.5+	13.5+	XLarge	79-95548
RAX Standard Tall				
79-95492	2-4	3.5-5.5	XSmall	79-95552
79-95493	4.5-7	6-8	Small	79-95553
79-95495	7.5-10.5	8.5-11.5	Medium	79-95555
79-95497	10.5-12.5	11.5-13.5	Large	79-95557
79-95498	12.5+	13.5+	XLarge	79-95558
RAX Standard Ankle				
79-95502	2-4	3.5-5.5	XSmall	79-95562
79-95503	4.5-7	6-8	Small	79-95563
79-95505	7.5-10.5	8.5-11.5	Medium	79-95565
79-95507	10.5-12.5	11.5-13.5	Large	79-95567
79-95508	12.5+	13.5+	XLarge	79-95568

Accessories		
P/N	PART DESCRIPTION	
79-95093	Weather Cover, Small	
79-95095	Weather Cover, Medium	
79-95097	Weather Cover, Large	
79-95210	Walking Brace Liner Extension Pad	



DJO, LLC | A DJO Global Company

T 800.793.6065 **D** 760.727.1280 **F** 800.580.9625 1430 Decision Street I Vista, CA 92081-8553 I U.S.A. procaresoftgoods.com

Together in Motion...

DJO Global proudly sponsors:









