



## DONJOY BACK BRACE II TLSO

The DonJoy Back Brace II TLSO accommodates the appropriate motion restriction throughout the healing process. The lightweight design includes a patented Mechanical Advantage Pulley System, a quick-release shoulder strap buckle, and an easy to use thoracic adjustment. All these features allow for the patient to comfortably control and relieve acute pain often associated with anomalies of the spine such as osteoporosis, compression fractures, spinal stenosis, strain, and excessive kyphosis.

Padded shoulder straps  
enhance patient comfort.

Midline cutout  
avoids contact with  
surgical incision.

Rigid overlapping  
side panels for  
lateral stability.

Rigid anterior plate  
helps to resist flexion.






Patented Mechanical Advantage Pulley System easily adjusts trunk compression.




Quick release shoulder buckles simplify donning and doffing.



The sternal bar easily adjusts without tools for patient height.

Part #	Description	Size	Waist Circumference	Back Panel Height
11-1585-2	DONJOY BACK BRACE II - TLSO	S	30"–35"	17"
11-1585-3	DONJOY BACK BRACE II - TLSO	M	35"–40"	17"
11-1585-4	DONJOY BACK BRACE II - TLSO	L	40"–45"	19"
11-1585-5	DONJOY BACK BRACE II - TLSO	XL	45"–50"	19"
11-1585-6	DONJOY BACK BRACE II - TLSO	XXL	50"–55"	19"
25-1585-2	STRAP KIT	S		
25-1585-3	STRAP KIT	M,L,XL,XXL		



DONJOY BACK BRACE PRODUCT LINE	SI-S5	LI-L5	LI-SI	LI-SI	T9-SI	T9-SI	T2/T3-SI
	SI BELT	LSO LITE	LSO LOW PROFILE	LSO REGULAR	LSO PREMIUM PLUS	LSO LSO W/CHAIRBACK	TLSO
INDICATIONS FOR USE							
Sacroiliac Pain	✓						
Low Back Pain	✓	✓	✓	✓	✓	✓	✓
Low Back Sprains/Strains	✓	✓	✓	✓	✓	✓	✓
Lumbar Disc Displacement		✓	✓	✓	✓	✓	✓
Osteoporosis			✓	✓	✓	✓	✓
Disc Herniation & Degeneration			✓	✓	✓	✓	✓
Post-Operative Laminectomy			✓	✓	✓	✓	✓
Post-Operative Disectomy			✓	✓	✓	✓	✓
IDET Procedure			✓	✓	✓	✓	✓
Spondylolisthesis			✓	✓	✓	✓	✓
Post-Operative Fusion					✓	✓	✓
Spinal Stenosis					✓	✓	✓
Compression Fractures							✓
Kyphosis							✓
Scoliosis							✓



**DJO, LLC | A DJO Global Company**  
**T** 800.336.6569 **D** 760.727.1280 **F** 800.936.6569  
 1430 Decision Street | Vista, CA 92081-8553 | U.S.A.  
[DJOglobal.com](http://DJOglobal.com)



DJO Global proudly sponsors:

