



# REGENERATION BONE GROWTH STIMULATION

## DEVICE MEASUREMENT GUIDE



FRACTURE SITE		DUAL COIL PART# 01-203-0001 L 9" x W 6 1/2" x H 7" 1.85 LBS.	SIZE 1 PART# 01-211-0001 L 4" x W 4 1/4" x H 1 1/4" 0.9 LBS.	SIZE 2 PART# 01-211-0002 L 9" x W 8 1/2" x H 4 1/2" 1.9 LBS.	SIZE 3 PART# 01-211-0003 L 10 1/2" x W 11 1/2" x H 5 3/4" 2.45 LBS.	SIZE 4 PART# 01-211-0004 L 13" x W 14" x H 9" 4.95 LBS.	SPINE PART# 01-207-0007 L 10 1/2" x W 9" x H 1" 1.85 LBS.
Clavicle	Proximal		•				
	Mid		•				
	Distal		•				
Humerus	Proximal			•			
	Mid			•			
	Distal	•		•			
Radius/Ulna	Proximal	•		•			
	Mid	•		•			
	Distal	•					
Hand/Wrist	Phalanges	•	•				
	Metacarpals	•	•				
Femur	Proximal			•	•	•	
	Mid			•	•	•	
	Distal			•	•		
Tibia/Fibula	Proximal			•	•		
	Mid	•		•			
	Distal	•		•			
Foot/Ankle	Phalanges	•	•				
	Metatarsals	•	•				
	Ankle	•		•			
Scaphoid	•	•					
Lumbar/Spine						•	

Individual results may vary. This therapy is not for everyone. Please consult your physician. A prescription is required. For more information, please call DJO® at 800-336-6569. Prior to use, refer to the Instructions for Use supplied with these devices for indications, contraindications, side effects, suggested procedure, warnings and precautions.