



DONJOY®



THE CLASSIC ADVANTAGE The IceMan CLASSIC³™ cold therapy unit helps reduce pain and swelling, speeding up rehabilitation. The IceMan® helps provide extended cold therapy to aid in a variety of indications and protocols as directed by a medical professional. The IceMan CLASSIC³ utilizes DonJoy's patented recirculation technology.



ICEMAN with dual pump recirculation technology. The DonJoy® IceMan® features a patented recirculation system that allows water warmed after flowing through the cold pad to be preserved and remixed with cooler ice water at a constant flow rate, providing consistent cool water distribution throughout the cold pad.



COMPETITORS Other versions restrict or pinch water flow to the cold pad, causing variable flow rates and inconsistent cold pad temperatures.



DonJoy® IceMan CLASSIC3™ provides a distinct advantage with the controlled delivery of cold therapy helping patients and medical practitioners to easily manage post-operative pain and swelling.

FEATURES

Proven Performance

- Time tested with over 20-years of experience delivering industry leading cold therapy solutions
- Developed through direct customer interaction and feedback

Technology Advantage

- Patented recirculation technology helps deliver consistent cold therapy throughout the cold pad, eliminating freezing inlet temperatures

Patient Convenience and Ease-Of-Use

- Cooler designed for easy set-up
- Plug-n-play operation
- Designed for easy portability and storage
- A complete line of cold pads to help deliver cold therapy to where it is needed the most

DONJOY® ICEMAN CLASSIC3™

Description	Part Number
IceMan CLASSIC3™ Cold Therapy Unit	11-0494
IceMan CLASSIC3™ w/ Universal Cold Pad, S, RH	11-1421
IceMan CLASSIC3™ w/ Universal Cold Pad, NS, RH	11-1422
IceMan CLASSIC3™ w/ Universal Wrap-On Pad, NS, RH	11-1423
IceMan CLASSIC3™ w/ Universal Wrap-On Pad, NS, EH	11-1424
IceMan CLASSIC3™ w/ Shoulder, S Wrap-On Pad, NS, EH	11-9098
IceMan CLASSIC3™ w/ McGuire Knee Wrap-On Pad, NS, RH	11-9099

WRAP-ON PADS

Description (Non-Sterile)	Part Number
Universal Wrap-On Pad, NS, RH (10.85" x 11.34")	11-0679-9-00000
Universal Wrap-On Pad, NS, EH (10.85" x 11.34")	11-0680-9-00000
Universal Wrap-On XL, NS, RH (12"x11.34")	11-0681-9-00000
Universal Wrap-On Pad XL, NS, EH (12"x11.34")	11-0682-9-00000
McGuire Knee Wrap-On Pad, NS, RH (10.85"x11.34")	11-1268-9-00000
Shoulder Wrap-On Pad S, NS, EH (12.75"x19"x12.25")	11-1397-9-00000
Shoulder Wrap-On Pad, NS, EH (14"x23"x12.5")	11-1334-9-00000
Ankle Wrap-On Pad, NS, RH	11-1522
Fore-Foot Wrap-On Pad, NS, RH	11-1521
Hip-Left Wrap-On Pad, NS, RH	11-9112
Hip-Right Wrap-On Pad, NS, RH	11-9111

COLD PADS

Description (Sterile)	Part Number
Universal Cold Pad, S, RH (10.85"x11.34")	11-0801-9-00000
Universal Cold Pad, S, EH (10.85"x11.34")	11-0802-9-00000
Universal Cold Pad XL, S, EH (12"x11.34")	11-0804-9-00000
Rectangle Cold Pad S, S, RH (4.81"x10.75")	11-0817-9-00000
Rectangle Cold Pad XL, L, S, EH (9.76"x13.75")	11-0813-9-00000
Wrist/Hand Cold Pad, S, RH (2.9"x4.2")	11-1269-9-00000
Ankle Cold Pad, S, RH (12.26"x8.28")	11-0815-9-00000

Description (Sterile)	Part Number
Universal Cold Pad, NS, RH (10.85"x11.34")	11-0675-9-00000
Universal Cold Pad, NS, EH (10.85"x11.34")	11-0676-9-00000
Universal Cold Pad XL, NS, EH (12"x11.34")	11-0678-9-00000

STERILE DRESSINGS

Description (Non-Sterile)	Part Number
Universal Sterile Dressing	11-0822-9-00000
Universal XL Sterile Dressing	11-0823-9-00000
Rectangle XL Sterile Dressing	11-0820-9-00000
Shoulder XL Sterile Dressing	11-1360-9-00000

S = Sterile
NS = Non-Sterile
RH = Regular Hose (6.5")
EH = Extended Hose (15")

WARNING: THIS DEVICE CAN BE COLD ENOUGH TO CAUSE SERIOUS INJURY. DO NOT USE THIS DEVICE WITHOUT A PRESCRIPTION FROM A PHYSICIAN.



T 800.336.6569 D 760.727.1280 F 800.936.6569
 2900 Lake Vista Drive | Dallas, TX 75067 | U.S.A.
 DJOglobal.com/donjoy

Copyright © 2019 by DJO, LLC
 MKT00-8244 Rev A

Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this sheet do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.