

DONJOY®



Revolutionary Cold Therapy Compression Sleeve

The Freeze Sleeve™ MD provides 360° of cold therapy. This 100% coverage helps reduce inflammation to aching muscles and joints, pain, and post operative procedures.

SIMPLE

The easy slip on process eliminates the need for messy ice bags, cumbersome wraps and multi-step products on the market today.

RELIEF & RECOVERY

Helps provide relief and recovery for aching muscles and joints.

COVERAGE

360 degrees of even coverage, and also allows full range of motion to the user.

TREATMENT

The Freeze Sleeve™ MD helps provide the recommended 20 minutes of cold therapy treatment.

COMFORT

The compression design allows the Freeze Sleeve™ MD to be used comfortably on multiple areas of the arms and legs.

MOBILITY

Allows the user to apply daily cold therapy treatment on the go whenever and wherever you are!

CLEAN

Specially formulated with antimicrobial gel that eliminates odors and is nontoxic.



ELBOW



WRIST



ANKLE



HAMSTRING



CALF



KNEE

HOW TO MEASURE FOR FIT

The Freeze Sleeve™ MD should fit snug around the treated area. Its “firm fit” is not intended for support, but to hold the sleeve in place aiding in mobile recovery. If patients are between sizes, it is recommended they order the larger size to help ensure an adequate fit.



WRIST

Measure circumference of wrist 2 inches above where the hand meets the wrist.



SHIN & CALF

Measure circumference of area around the largest portion of the calf.



ELBOW

With arm bent at 90 degrees, measure from outer point and around inner crease of elbow.



KNEE

While standing, measure circumference of knee 3 inches above where the knee bends.



ANKLE

Measure circumference of ankle from the heel to 2 inches above the bridge of the foot.



HAMSTRING & QUAD

While standing, measure circumference of area around the largest portion of the quad & hamstring.

ORDERING INFORMATION

PART NUMBER	DESCRIPTION	SIZE
11-3340-2	8-10 in (20-26 cm)	SMALL
11-3340-3	10-13 in (26-33 cm)	MEDIUM
11-3340-4	13-16 in (33-42 cm)	LARGE
11-3340-5	16-20 in (42-52 cm)	X-LARGE
11-3340-6	20-25 in (52-63 cm)	XX-LARGE
11-3340-7	25-30 in (63-76 cm)	XXX-LARGE

Not made with natural rubber latex.



T 800.336.6569 D 760.727.1280 F 800.936.6569
2900 Lake Vista Drive | Lewisville, TX 75067 | U.S.A.
DJOGlobal.com/donjoy

Copyright © 2019 by DJO, LLC
MKT00-9263 Rev A

 and Freeze Sleeve™ are trademarks of Freeze Sleeve, LLC

Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this sheet do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.