


EVOLVE34™

LAPIDUS CORRECTION SYSTEM

Your Lapidus, your way. A comprehensive assistant allowing multiple fixation options to support your preferred approach.

BUNION SYSTEMS



The **Evolve34™ Lapidus Correction System** is designed to provide powerful correction, predictable results, and procedural flexibility. The design allows for 2-in-1 correction of both the intermetatarsal angle and the frontal plane while also giving surgeons the ability to joint prep and fixate according to their preferred method.

FEATURES & BENEFITS



FLEXIBILITY

Allows for any joint prep desired and fixation choices based on surgeon preference, empowering the surgeon to perform the procedure the way *they* want.

ANATOMICAL CORRECTION REFERENCE

The jig is anchored on the 2nd metatarsal, providing the surgeon with an anatomical reference for correction during the procedure. The optional cut guide is designed to minimize shortening of the 1st metatarsal and allows for cuts perpendicular to the 2nd metatarsal.

COMPREHENSIVE ASSISTANT

Evolve34™ is an all-in-one, hands-free tool that helps the surgeon with correction, distraction, and compression. No need for multiple pieces.

EFFICIENT INSTRUMENTATION

The Evolve34 Lapidus Correction System comes in a small sterilization tray with the reusable jig and other instruments, along with a procedure pack for disposables.



FLEXIBLE FIXATION

Your Lapidus, your way. Evolve34 is compatible with multiple fixation options.

DYNACLIP®
BONE
FIXATION
SYSTEM



ARSENAL
FOOT
PLATING
SYSTEM™



TIGER
CANNULATED
SCREWS



enovis.

T 800.495.2919 F 877.778.3864

Medshape, Inc.
1575 Northside Drive NW | Suite 440 | Atlanta, GA 30318 | U.S.A.
enovis.com/foot-and-ankle

Copyright © 2023 MedShape, Inc.
MK-10363 Rev 00
CRO-2023-0123

EVOLVE34 SINGLE-USE INSTRUMENTS

PART #	DESCRIPTION
3400-04-000	LAPIDUS SMOOTH FIXATION AND SAWBLADE KIT

Individual results may vary. Neither Medshape, Inc. nor any of the Enovis companies dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.