

# Exos Form<sup>™</sup> II Spine Bracing System

Reshaping Form, Fit and Function










### V-STX™

A patented, thermoformable material used with the posterior, anterior and lordotic inserts, allows for true customization and ideal patient fit.



# **ComfortCORE™ Foam**

Conforms to the unique contours of a patient's body, adjusts instantly to individual movements and delivers improved patient comfort.



### Diamond<sup>2</sup> Grid™

Uniquely designed construct, delivers patients' a more durable, longer-wearing brace.



## Boa® Technology

Independent superior and inferior compression delivered through an innovative closure system.

## Adjustable Belt Wings

Modifiable belt accommodates varying patient body structures for optimal fit.

# Modular Design

Step-up / Step-down design allows for single brace use through rehabilitation and recovery.

## Semi-Universal Sizing

Removable circumference tape allows for accurate patient sizing.



# Exos Form<sup>™</sup> II Spine Bracing System

Exos Form II is designed with patented thermoformable VST-X™ material for true customization, ComfortCORE™ Foam that contours to the patient's body, Diamond² Grid™ construct for a more durable and longerwearing brace and Boa® Technology's advanced closure system. The end result is a spine brace that delivers superior decompression and support, a quick and easily customizable fit, and provides patient comfort and satisfaction by reshaping form, fit and function.



# Exos FORM™ II 627

The Exos FORM™ II 627 (with anterior panel) provides relief from a wide range of indications from mild to acute lower back pain, as well as post-operative support. The Exos FORM™ II 627 provides superior sagittal control and support from L-1 to L-5.

PART NO.	DESCRIPTION
305627-40	EXOS FORM™ II 627, S/M
305627-60	EXOS FORM™ II 627, L/XL

PDAC Assigned Code: Lo627 (CST) / Lo642 (OTS)



# Exos FORM™II 637

The Exos FORM™ II 637 provides a higher degree of spinal support and relief from a wide range of indicationsfrom acute lower back pain to post-operative support. The Exos FORM™ II 637 provides superior sagittal and coronal control and support from Tg-S1.

PART NO.	DESCRIPTION
305637-40	EXOS FORM™ II 637, S/M
305637-60	EXOS FORM™ II 637, L/XL

PDAC Assigned Code: Lo637 (CST) / Lo650 (OTS)



## Exos FORM™II 626

The Exos FORM™ II 626 provides relief from mild lower back pain. The Exos FORM™ II 626 provides superior sagittal control and support from L-1 to L-5.

PART NO.	DESCRIPTION
305626-40	EXOS FORM™ II 626, S/M
305626-60	EXOS FORM™ II 626, L/XL

PDAC Assigned Code: Lo626 (CST) / Lo641 (OTS)



# Exos FORM™ II 631

The Exos FORM™ II 631 provides a higher degree of spinal support and relief from a wide range of indications from acute lower back pain to post-operative support. The Exos FORM™ II 631 provides superior sagittal control and support from T9-S1.

PART NO.	DESCRIPTION
305631-40	EXOS FORM™ II 631, S/M
305631-60	EXOS FORM™ II 631, L/XL

PDAC Assigned Code: Lo631 (CST) / Lo648 (OTS)



# Exos FORM™ II 621

The Exos FORM™ II 621 provides superior stabilization of the sacroiliac (SI) joint through direct compression. Silicone grip strips prevent the brace from migrating.


PART NO.	DESCRIPTION
305621-40	EXOS FORM™ II 621, S/M
305621-60	EXOS FORM™ II 621, L/XL

PDAC Assigned Code: Lo621

#### **COMMON INDICATIONS**

	Sacrum	L1 - L5	L1 - L5	T9 - S1	T9 - S1
INDICATIONS	Exos FORM™ II 621	Exos FORM™ II 626	Exos FORM™ II 627	Exos FORM™ II 631	Exos FORM™ II 637
Sacroiliac Pain	•				
Acute and Chronic Lower Back Pain		•	•	•	•
Lower Back Sprains / Strains		•	•	•	•
Lumbar Disc Displacement				•	
Osteoporosis		•	•	•	•
Disc Herniation and Degeneration			•		•
Post-Operative Laminectomy		•	•	•	•
Post-Operative Disectomy					•
Spondylolisthesis		•	•	•	•
Post-Operative Fusion				•	
Spinal Stenosis				•	•
Stable / Non-Displaced Fractures					

#### **RANGE OF SPINAL SUPPORT**





# ORDERING INFORMATION

#### SIZING GUIDE WAIST CIRCUMFERENCE

DESCRIPTION	Exos FORM™ II 637	Exos FORM™ II 631	Exos FORM™ II 627	Exos FORM™ II 626	Exos FORM™ II 621
S/M	25" - 49" (63cm - 125cm)	less than 35" (89 cm)			
L/XL	49" - 60" (125cm - 152cm)	36" - 55" (91 cm - 140 cm)			

#### EXOS FORM™ II SPINE BRACING SYSTEM

PART NO.	DESCRIPTION	SIZE
305621-40	EXOS FORM II 621	S/M
305621-60	EXOS FORM II 621	L/XL
305626-40	EXOS FORM II 626	S/M
305626-60	EXOS FORM II 626	L/XL
305627-40	EXOS FORM II 627	S/M
305627-60	EXOS FORM II 627	L/XL
305631-40	EXOS FORM II 631	S/M
305631-60	EXOS FORM II 631	L/XL
305637-40	EXOS FORM II 637	S/M
305637-60	EXOS FORM II 637	L/XL

#### **EXOS FORM™ II SPINE BRACING ACCESSORIES**

PART NO.	DESCRIPTION	SIZE			
905-631-40	POSTERIOR T9 PANEL	S/M			
905-631-60	POSTERIOR T9 PANEL	L/XL			
905-637-40	LATERAL EXTENSIONS	S/M			
905-637-60	LATERAL EXTENSIONS	L/XL			
905-AP-40	ANTERIOR PANEL	S/M			
905-AP-60	ANTERIOR PANEL	L/XL			
905-1-WINGS-40	BELT WING SET	S/M			
905-1-WINGS-60	BELT WING SET	L/XL			
905-1-BWING-40	BRIDGE WING	S/M			
905-1-BWING-60	BRIDGE WING	L/XL			
905-1-INSERT-40	LORDOTIC INSERTS	S/M			
905-1-INSERT-60	LORDOTIC INSERTS	L/XL			
905-JITT	SI BELT PAD SET	UNIVERSAL			
905-POSTPAD-40	POSTERIOR T9 PANEL PADDING	S/M			
905-POSTPAD-60	POSTERIOR T9 PANEL PADDING	L/XL			
905-ANTPAD-40	ANTERIOR PANEL PADDING	S/M			
905-ANTPAD-60	ANTERIOR PANEL PADDING	L/XL			
905-HOTCOLD-40	HOT/COLD PACK ACCESSORY	S/M			
905-HOTCOLD-60	HOT/COLD PACK ACCESSORY	L/XL			

# Embrace the Reform



DJO, LLC | A DJO Global Company

**T** 800.336.6569 **D** 760.727.1280 **F** 800.936.6569 1430 Decision Street | Vista, CA 92081-8553 | U.S.A.



