

A uniquely designed, versatile mini external fixation device offering flexible angular ability.

#### **EXTERNAL FIXATION**

## FLEXIBLE ANGULAR ABILITY

Other mini rail frames limit surgeons to a rail or bar with two pin cars connected to it. Each car holds two pins that are parallel, completely limiting the surgeon's ability to place a pin. The EF1 Mini Rail offers flexible angular abilities since it is not restricted to parallel pin placement





# **FEATURES & BENEFITS**



The sector of th

Pulais and a state of the state

#### 360° CONE OF ANGULATION

#### PIN, COLLET, & CLAMP

The EF1 Mini Rail pin, collet, and clamp designs allow the surgeon to select the optimal pin position first, then connect to the clamp for flexible angulation ability.



### VARIABLE PIN SIZES

Color-coded pin collets enable the surgeon to easily choose a 3mm, 4mm, or 5mm pin per their preference.



# enovis.

T 800.495.2919 F 877.778.3864

Trilliant Surgical, LLC 727 North Shepherd Drive | Suite 100 | Houston, TX 77007 | U.S.A. enovis.com/foot-and-ankle

Copyright S 2023 Enovis Foot and Ankle MK-10401 Rev00

Individual results may vary. Neither Trilliant Surgical, LLC nor any of the Enovis companies dispense medical advice. The contents of this document do not constitute medical, (egal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.