



c€ XMS01038-1 Ascaplasty. Adjustable positioning based on the point of entry. LOCKING DIAL Locks the targeting arc in position.

ACCURATE.

The "Sure-Shot" targeting arc enhances the surgeon's ability, like a marksman, to position Peca Implants. The design ensures a safe distance between where the proximal Peca implant exits the lateral cortex from the osteotomy.






EXACT-T® RECESSAllows Exact driver positioning and provides optimal torque.

CONTROLLED.

The freedom paddle offers full tri-planar control to the surgeon. With the device positioned medially, the surgeon has the freedom to manipulate the transverse, frontal, and sagittal planes, in addition to the DMAA, which would otherwise be limited by a transfix pin into the head.



CONTROLLED TRI-PLANAR CORRECTION





LEARN MORE

enovis

T 800.495.2919 F 877.778.3864

Trilliant Surgical, LLC 727 North Shepherd Drive | Suite 100 | Houston, TX 7707 | U.S.A. enovis.com/foot-and-ankle

Copyright © 2023 Enovis Foot & Ankle MK-10483 Rev 00

Individual results may vary. Neither Trilliant Surgical, LLC nor any of the Enovis companies dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.