

# STERISHIELD® II DUAL LAYER AMNION

SteriShield II is a dual layer amnion patch intended to create an environment conducive for regenerating healthy tissue and improving clinical outcomes.

## BIOLOGICS

- Published benefits of amniotic membranes may include reduced inflammation, reduced fibrosis/scarring, decreased pain, and anti-adhesion properties\*
- A dual layer amnion preparation allows for reversible application at the surgical site and eliminates site specific placement requirements
- Provided sterile via low dose, low temperature, gamma sterilization with a sterility assurance level (SAL) of 10<sup>-6</sup>
- Packaged with a shelf life of 3 years and may be stored in ambient temperatures
- Available as a patch in a variety of dimensions



Part Number	Size Dimension (L x W)
830-TS	1 x 1cm patch
831-TS	2 x 2cm patch
832-TS	2 x 3cm patch

Part Number	Size Dimension (L x W)
833-TS	4 x 4cm patch
834-TS	4 x 6cm patch
835-TS	4 x 8cm patch

\*Please refer to [djoglobal.com](http://djoglobal.com) for references (not specific to SteriShield II). All biologic products are minimally manipulated and intended for homologous use only. Trilliant is a distributor for BoneBank® Allografts.