

### UDI Cross Reference Tool for Trilliant Surgical Products

This catalog contains applicable GTIN information for Trilliant Surgical products in order to assist in deriving required UDI information, along with other pertinent product information as specified on respective labeling (i.e., materials, symbols, etc.).

All product shipped from Trilliant Surgical in Houston, TX contains all the necessary UDI information on the respective product's package labeling; however, for non-sterile implants and instruments that are cleaned and sterilized prior to use, they may be separated from their original label and packaging in order to undergo cleaning and sterilization wherein the original labeling may be discarded.

When product becomes separated from their original labeling, UDI information can be derived by combining available product identifying information, Product Identifiers (PI), on devices with a cross referencing tool containing associated Device Identifiers (DI).

*Note: Trilliant Surgical utilizes Lot Numbers as Product Identifiers for non-sterile product.*

Within the GTIN catalog, for each product, there is an applicable GTIN, which is inclusive of the respective product's Device Identifier (DI), the portion of the UDI that identifies the specific version or model a device and the device labeler.

For example, the DI for part 101-00-001, (DI=00812926020075), can identify the part as a 3S Hemi Implant, Small, Non-Sterile and Trilliant Surgical as the device labeler for the part on global scale.

Together, the UDI information can be generated;  $DI + PI = UDI$

The format to generate a UDI (for Trilliant Surgical product) is as follows:

**(01)<Device Identifier>(10)<Product Identifier>**

For example, if the part is 101-00-001 (DI=00812926020075) and the associated Lot is TSL001234, the UDI is:

**(01)00812926020075(10)TSL001234**

#### **Important to Note:**

- Copies of the most current reference catalog are available on the website, <https://www.djoglobal.com/udi-gtin-information>
- Obtain the most recent version of this document before use to ensure information is current and accurate for each represented product
- This catalog may include active and inactive product; all products may not currently be offered and/or available at the time of use

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
101-00-001	3S Hemi Implant, Small, Non-Sterile	00812926020075	Yes	Cobalt Chromium (CoCrMo)
101-00-002	3S Hemi Implant, Medium, Non-Sterile	00812926020082	Yes	Cobalt Chromium (CoCrMo)
101-00-003	3S Hemi Implant, Large, Non-Sterile	00812926020099	Yes	Cobalt Chromium (CoCrMo)
101-00-004	3S Hemi Implant, XLarge, Non-Sterile	00812926020105	Yes	Cobalt Chromium (CoCrMo)
102-10-007	7mm Twist Subtalar Implant	00812926020235	Yes	Titanium Alloy (Ti6Al4V)
102-10-008	8mm Twist Subtalar Implant	00812926020242	Yes	Titanium Alloy (Ti6Al4V)
102-10-009	9mm Twist Subtalar Implant	00812926020259	Yes	Titanium Alloy (Ti6Al4V)
102-10-010	10mm Twist Subtalar Implant	00812926020266	Yes	Titanium Alloy (Ti6Al4V)
102-10-011	11mm Twist Subtalar Implant	00812926020273	Yes	Titanium Alloy (Ti6Al4V)
102-10-012	12mm Twist Subtalar Implant	00812926020280	Yes	Titanium Alloy (Ti6Al4V)
102-20-007	7mm Disco Subtalar Implant	00812926020297	Yes	Titanium Alloy (Ti6Al4V)
102-20-008	8mm Disco Subtalar Implant	00812926020303	Yes	Titanium Alloy (Ti6Al4V)
102-20-009	9mm Disco Subtalar Implant	00812926020310	Yes	Titanium Alloy (Ti6Al4V)
102-20-010	10mm Disco Subtalar Implant	00812926020327	Yes	Titanium Alloy (Ti6Al4V)
102-20-011	11mm Disco Subtalar Implant	00812926020334	Yes	Titanium Alloy (Ti6Al4V)
102-20-012	12mm Disco Subtalar Implant	00812926020341	Yes	Titanium Alloy (Ti6Al4V)
111-00-001	3S Hemi Non-Sterile Sterilization Tray	00812926020112	No	Polypropylene, Polyphenylsulfone
111-00-002	3S Hemi Handle w/ Impactor Cap	00812926020129	No	Silicone, Stainless Steel, Aluminum Alloy
111-00-003	3S Hemi Sizer S/M, Non-Sterile	00812926020136	No	Titanium Alloy (Ti6Al4V)
111-00-004	3S Hemi Sizer L/XL, Non-Sterile	00812926020143	No	Titanium Alloy (Ti6Al4V)
111-00-005	3S Hemi, Non-Sterile, Awl	00812926020150	No	Stainless Steel
111-00-006	3S Hemi S/M Broach, Non-Sterile	00812926020167	No	Stainless Steel
111-00-007	3S Hemi L/XL Broach, Non-Sterile	00812926020174	No	Stainless Steel
111-00-008	3S Hemi Trial, Small, Non-Sterile	00812926020181	No	Titanium Alloy (Ti6Al4V)
111-00-009	3S Hemi Trial, Medium, Non-Sterile	00812926020198	No	Titanium Alloy (Ti6Al4V)
111-00-010	3S Hemi Trial, Large, Non-Sterile	00812926020204	No	Titanium Alloy (Ti6Al4V)
111-00-011	3S Hemi Trial, XLarge, Non-Sterile	00812926020211	No	Titanium Alloy (Ti6Al4V)
111-00-012	3S Hemi Impactor Bit, Non-Sterile	00812926020228	No	Stainless Steel
111-00-013	3S Hemi Impactor Head, Non-Sterile	00812926028033	No	Polyphenylsulfone
111-99-002	3S Hemi, Non-Sterile, Tray Lid	00812926027852	No	Aluminum
112-00-001	Subtalar Implant System Sterilization Tray Base Assembly	00812926020358	No	Polypropylene, Aluminum
112-00-002	Probe	00812926020365	No	Stainless Steel
112-00-003	Insertion Tool	00812926020372	No	Stainless Steel
112-00-004	Implant Pin	00812926020389	Yes	Stainless Steel
112-00-005	Removal Tool	00812926020396	No	Stainless Steel
112-00-006	Subtalar Tray Shipping Insert	00812926020402	No	Polypropylene
112-10-007	7mm Twist Subtalar Trial	00812926020419	No	Titanium Alloy (Ti6Al4V)
112-10-008	8mm Twist Subtalar Trial	00812926020426	No	Titanium Alloy (Ti6Al4V)
112-10-009	9mm Twist Subtalar Trial	00812926020433	No	Titanium Alloy (Ti6Al4V)
112-10-010	10mm Twist Subtalar Trial	00812926020440	No	Titanium Alloy (Ti6Al4V)
112-10-011	11mm Twist Subtalar Trial	00812926020457	No	Titanium Alloy (Ti6Al4V)
112-10-012	12mm Twist Subtalar Trial	00812926020464	No	Titanium Alloy (Ti6Al4V)
112-20-007	7mm Disco Subtalar Trial	00812926020471	No	Titanium Alloy (Ti6Al4V)
112-20-008	8mm Disco Subtalar Trial	00812926020488	No	Titanium Alloy (Ti6Al4V)
112-20-009	9mm Disco Subtalar Trial	00812926020495	No	Titanium Alloy (Ti6Al4V)

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
112-20-010	10mm Disco Subtalar Trial	00812926020501	No	Titanium Alloy (Ti6Al4V)
112-20-011	11mm Disco Subtalar Trial	00812926020518	No	Titanium Alloy (Ti6Al4V)
112-20-012	12mm Disco Subtalar Trial	00812926020525	No	Titanium Alloy (Ti6Al4V)
112-99-011	Subtalar Implant System Sterilization Tray Lid	00812926027869	No	Aluminum
200-20-008	2.0mm x 8mm Tiger Cannulated Screw	00812926020532	Yes	Titanium Alloy (Ti6Al4V)
200-20-010	2.0mm x 10mm Tiger Cannulated Screw	00812926020549	Yes	Titanium Alloy (Ti6Al4V)
200-20-012	2.0mm x 12mm Tiger Cannulated Screw	00812926020556	Yes	Titanium Alloy (Ti6Al4V)
200-20-014	2.0mm x 14mm Tiger Cannulated Screw	00812926020563	Yes	Titanium Alloy (Ti6Al4V)
200-20-016	2.0mm x 16mm Tiger Cannulated Screw	00812926020570	Yes	Titanium Alloy (Ti6Al4V)
200-20-018	2.0mm x 18mm Tiger Cannulated Screw	00812926020587	Yes	Titanium Alloy (Ti6Al4V)
200-20-020	2.0mm x 20mm Tiger Cannulated Screw	00812926020594	Yes	Titanium Alloy (Ti6Al4V)
200-20-022	2.0mm x 22mm Tiger Cannulated Screw	00812926020600	Yes	Titanium Alloy (Ti6Al4V)
200-20-024	2.0mm x 24mm Tiger Cannulated Screw	00812926020617	Yes	Titanium Alloy (Ti6Al4V)
200-20-026	2.0mm x 26mm Tiger Cannulated Screw	00812926020624	Yes	Titanium Alloy (Ti6Al4V)
200-20-028	2.0mm x 28mm Tiger Cannulated Screw	00812926020631	Yes	Titanium Alloy (Ti6Al4V)
200-20-030	2.0mm x 30mm Tiger Cannulated Screw	00812926020648	Yes	Titanium Alloy (Ti6Al4V)
200-20-032	2.0mm x 32mm Tiger Cannulated Screw	00812926020655	Yes	Titanium Alloy (Ti6Al4V)
200-20-034	2.0mm x 34mm Tiger Cannulated Screw	00812926020662	Yes	Titanium Alloy (Ti6Al4V)
200-20-036	2.0mm x 36mm Tiger Cannulated Screw	00812926020679	Yes	Titanium Alloy (Ti6Al4V)
200-20-038	2.0mm x 38mm Tiger Cannulated Screw	00812926020686	Yes	Titanium Alloy (Ti6Al4V)
200-20-040	2.0mm x 40mm Tiger Cannulated Screw	00812926020693	Yes	Titanium Alloy (Ti6Al4V)
200-20-042	2.0mm x 42mm Tiger Cannulated Screw	00812926020709	Yes	Titanium Alloy (Ti6Al4V)
200-24-008	2.4mm x 8mm Tiger Cannulated Screw	00812926020716	Yes	Titanium Alloy (Ti6Al4V)
200-24-010	2.4mm x 10mm Tiger Cannulated Screw	00812926020723	Yes	Titanium Alloy (Ti6Al4V)
200-24-012	2.4mm x 12mm Tiger Cannulated Screw	00812926020730	Yes	Titanium Alloy (Ti6Al4V)
200-24-014	2.4mm x 14mm Tiger Cannulated Screw	00812926020747	Yes	Titanium Alloy (Ti6Al4V)
200-24-016	2.4mm x 16mm Tiger Cannulated Screw	00812926020754	Yes	Titanium Alloy (Ti6Al4V)
200-24-018	2.4mm x 18mm Tiger Cannulated Screw	00812926020761	Yes	Titanium Alloy (Ti6Al4V)
200-24-020	2.4mm x 20mm Tiger Cannulated Screw	00812926020778	Yes	Titanium Alloy (Ti6Al4V)
200-24-022	2.4mm x 22mm Tiger Cannulated Screw	00812926020785	Yes	Titanium Alloy (Ti6Al4V)
200-24-024	2.4mm x 24mm Tiger Cannulated Screw	00812926020792	Yes	Titanium Alloy (Ti6Al4V)
200-24-026	2.4mm x 26mm Tiger Cannulated Screw	00812926020808	Yes	Titanium Alloy (Ti6Al4V)
200-24-028	2.4mm x 28mm Tiger Cannulated Screw	00812926020815	Yes	Titanium Alloy (Ti6Al4V)
200-24-030	2.4mm x 30mm Tiger Cannulated Screw	00812926020822	Yes	Titanium Alloy (Ti6Al4V)
200-24-032	2.4mm x 32mm Tiger Cannulated Screw	00812926020839	Yes	Titanium Alloy (Ti6Al4V)
200-24-034	2.4mm x 34mm Tiger Cannulated Screw	00812926020846	Yes	Titanium Alloy (Ti6Al4V)
200-24-036	2.4mm x 36mm Tiger Cannulated Screw	00812926020853	Yes	Titanium Alloy (Ti6Al4V)
200-24-038	2.4mm x 38mm Tiger Cannulated Screw	00812926020860	Yes	Titanium Alloy (Ti6Al4V)
200-24-040	2.4mm x 40mm Tiger Cannulated Screw	00812926020877	Yes	Titanium Alloy (Ti6Al4V)
200-24-042	2.4mm x 42mm Tiger Cannulated Screw	00812926020884	Yes	Titanium Alloy (Ti6Al4V)
200-24-044	2.4mm x 44mm Tiger Cannulated Screw	00812926020891	Yes	Titanium Alloy (Ti6Al4V)
200-24-046	2.4mm x 46mm Tiger Cannulated Screw	00812926020907	Yes	Titanium Alloy (Ti6Al4V)
200-24-048	2.4mm x 48mm Tiger Cannulated Screw	00812926020914	Yes	Titanium Alloy (Ti6Al4V)
200-30-010	3.0mm x 10mm Tiger Cannulated Screw	00812926020921	Yes	Titanium Alloy (Ti6Al4V)
200-30-012	3.0mm x 12mm Tiger Cannulated Screw	00812926020938	Yes	Titanium Alloy (Ti6Al4V)
200-30-014	3.0mm x 14mm Tiger Cannulated Screw	00812926020945	Yes	Titanium Alloy (Ti6Al4V)

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
200-30-016	3.0mm x 16mm Tiger Cannulated Screw	00812926020952	Yes	Titanium Alloy (Ti6Al4V)
200-30-018	3.0mm x 18mm Tiger Cannulated Screw	00812926020969	Yes	Titanium Alloy (Ti6Al4V)
200-30-020	3.0mm x 20mm Tiger Cannulated Screw	00812926020976	Yes	Titanium Alloy (Ti6Al4V)
200-30-022	3.0mm x 22mm Tiger Cannulated Screw	00812926020983	Yes	Titanium Alloy (Ti6Al4V)
200-30-024	3.0mm x 24mm Tiger Cannulated Screw	00812926020990	Yes	Titanium Alloy (Ti6Al4V)
200-30-026	3.0mm x 26mm Tiger Cannulated Screw	00812926021003	Yes	Titanium Alloy (Ti6Al4V)
200-30-028	3.0mm x 28mm Tiger Cannulated Screw	00812926021010	Yes	Titanium Alloy (Ti6Al4V)
200-30-030	3.0mm x 30mm Tiger Cannulated Screw	00812926021027	Yes	Titanium Alloy (Ti6Al4V)
200-30-032	3.0mm x 32mm Tiger Cannulated Screw	00812926021034	Yes	Titanium Alloy (Ti6Al4V)
200-30-034	3.0mm x 34mm Tiger Cannulated Screw	00812926021041	Yes	Titanium Alloy (Ti6Al4V)
200-30-036	3.0mm x 36mm Tiger Cannulated Screw	00812926021058	Yes	Titanium Alloy (Ti6Al4V)
200-30-038	3.0mm x 38mm Tiger Cannulated Screw	00812926021065	Yes	Titanium Alloy (Ti6Al4V)
200-30-040	3.0mm x 40mm Tiger Cannulated Screw	00812926021072	Yes	Titanium Alloy (Ti6Al4V)
200-40-010	4.0mm x 10mm Tiger Cannulated Screw	00812926021089	Yes	Titanium Alloy (Ti6Al4V)
200-40-012	4.0mm x 12mm Tiger Cannulated Screw	00812926021096	Yes	Titanium Alloy (Ti6Al4V)
200-40-014	4.0mm x 14mm Tiger Cannulated Screw	00812926021102	Yes	Titanium Alloy (Ti6Al4V)
200-40-016	4.0mm x 16mm Tiger Cannulated Screw	00812926021119	Yes	Titanium Alloy (Ti6Al4V)
200-40-018	4.0mm x 18mm Tiger Cannulated Screw	00812926021126	Yes	Titanium Alloy (Ti6Al4V)
200-40-020	4.0mm x 20mm Tiger Cannulated Screw	00812926021133	Yes	Titanium Alloy (Ti6Al4V)
200-40-022	4.0mm x 22mm Tiger Cannulated Screw	00812926021140	Yes	Titanium Alloy (Ti6Al4V)
200-40-024	4.0mm x 24mm Tiger Cannulated Screw	00812926021157	Yes	Titanium Alloy (Ti6Al4V)
200-40-026	4.0mm x 26mm Tiger Cannulated Screw	00812926021164	Yes	Titanium Alloy (Ti6Al4V)
200-40-028	4.0mm x 28mm Tiger Cannulated Screw	00812926021171	Yes	Titanium Alloy (Ti6Al4V)
200-40-030	4.0mm x 30mm Tiger Cannulated Screw	00812926021188	Yes	Titanium Alloy (Ti6Al4V)
200-40-032	4.0mm x 32mm Tiger Cannulated Screw	00812926021195	Yes	Titanium Alloy (Ti6Al4V)
200-40-034	4.0mm x 34mm Tiger Cannulated Screw	00812926021201	Yes	Titanium Alloy (Ti6Al4V)
200-40-036	4.0mm x 36mm Tiger Cannulated Screw	00812926021218	Yes	Titanium Alloy (Ti6Al4V)
200-40-038	4.0mm x 38mm Tiger Cannulated Screw	00812926021225	Yes	Titanium Alloy (Ti6Al4V)
200-40-040	4.0mm x 40mm Tiger Cannulated Screw	00812926021232	Yes	Titanium Alloy (Ti6Al4V)
200-40-042	4.0mm x 42mm Tiger Cannulated Screw	00812926021249	Yes	Titanium Alloy (Ti6Al4V)
200-40-044	4.0mm x 44mm Tiger Cannulated Screw	00812926021256	Yes	Titanium Alloy (Ti6Al4V)
200-40-046	4.0mm x 46mm Tiger Cannulated Screw	00812926021263	Yes	Titanium Alloy (Ti6Al4V)
200-40-048	4.0mm x 48mm Tiger Cannulated Screw	00812926021270	Yes	Titanium Alloy (Ti6Al4V)
200-40-050	4.0mm x 50mm Tiger Cannulated Screw	00812926021287	Yes	Titanium Alloy (Ti6Al4V)
200-40-052	4.0mm x 52mm Tiger Cannulated Screw	00812926021294	Yes	Titanium Alloy (Ti6Al4V)
200-40-054	4.0mm x 54mm Tiger Cannulated Screw	00812926021300	Yes	Titanium Alloy (Ti6Al4V)
200-40-056	4.0mm x 56mm Tiger Cannulated Screw	00812926021317	Yes	Titanium Alloy (Ti6Al4V)
200-45-001	4.5mm Smooth Washer	00812926029580	Yes	Titanium Alloy (Ti6Al4V)
200-45-002	4.5mm Spiked Washer	00812926029573	Yes	Titanium Alloy (Ti6Al4V)
200-45-003	Ridged Washer	00810028392120	Yes	Titanium Alloy (Ti6Al4V)
200-45-020	4.5mm x 20mm STTCS	00812926028811	Yes	Titanium Alloy (Ti6Al4V)
200-45-022	4.5mm x 22mm STTCS	00812926028828	Yes	Titanium Alloy (Ti6Al4V)
200-45-024	4.5mm x 24mm STTCS	00812926028835	Yes	Titanium Alloy (Ti6Al4V)
200-45-026	4.5mm x 26mm STTCS	00812926028842	Yes	Titanium Alloy (Ti6Al4V)
200-45-028	4.5mm x 28mm STTCS	00812926028859	Yes	Titanium Alloy (Ti6Al4V)
200-45-030	4.5mm x 30mm STTCS	00812926028866	Yes	Titanium Alloy (Ti6Al4V)

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
200-45-032	4.5mm x 32mm STTCS	00812926028873	Yes	Titanium Alloy (Ti6Al4V)
200-45-034	4.5mm x 34mm STTCS	00812926028880	Yes	Titanium Alloy (Ti6Al4V)
200-45-036	4.5mm x 36mm STTCS	00812926028897	Yes	Titanium Alloy (Ti6Al4V)
200-45-038	4.5mm x 38mm STTCS	00812926028903	Yes	Titanium Alloy (Ti6Al4V)
200-45-040	4.5mm x 40mm STTCS	00812926028910	Yes	Titanium Alloy (Ti6Al4V)
200-45-045	4.5mm x 45mm STTCS	00812926028927	Yes	Titanium Alloy (Ti6Al4V)
200-45-050	4.5mm x 50mm STTCS	00812926028934	Yes	Titanium Alloy (Ti6Al4V)
200-45-055	4.5mm x 55mm STTCS	00812926028941	Yes	Titanium Alloy (Ti6Al4V)
200-45-060	4.5mm x 60mm STTCS	00812926028958	Yes	Titanium Alloy (Ti6Al4V)
200-45-065	4.5mm x 65mm STTCS	00812926028965	Yes	Titanium Alloy (Ti6Al4V)
200-45-070	4.5mm x 70mm STTCS	00812926028972	Yes	Titanium Alloy (Ti6Al4V)
200-45-140	4.5mm x 40mm LTTCS	00812926028989	Yes	Titanium Alloy (Ti6Al4V)
200-45-145	4.5mm x 45mm LTTCS	00812926028996	Yes	Titanium Alloy (Ti6Al4V)
200-45-150	4.5mm x 50mm LTTCS	00812926029009	Yes	Titanium Alloy (Ti6Al4V)
200-45-155	4.5mm x 55mm LTTCS	00812926029016	Yes	Titanium Alloy (Ti6Al4V)
200-45-160	4.5mm x 60mm LTTCS	00812926029023	Yes	Titanium Alloy (Ti6Al4V)
200-45-165	4.5mm x 65mm LTTCS	00812926029030	Yes	Titanium Alloy (Ti6Al4V)
200-45-170	4.5mm x 70mm LTTCS	00812926029047	Yes	Titanium Alloy (Ti6Al4V)
201-40-020	4.0mm x 20mm Cannulated Screw	00810028392755	Yes	Titanium Alloy (Ti6Al4V)
201-40-022	4.0mm x 22mm Cannulated Screw	00810028392762	Yes	Titanium Alloy (Ti6Al4V)
201-40-024	4.0mm x 24mm Cannulated Screw	00810028392779	Yes	Titanium Alloy (Ti6Al4V)
201-40-026	4.0mm x 26mm Cannulated Screw	00810028392786	Yes	Titanium Alloy (Ti6Al4V)
201-40-028	4.0mm x 28mm Cannulated Screw	00810028392793	Yes	Titanium Alloy (Ti6Al4V)
201-40-030	4.0mm x 30mm Cannulated Screw	00810028392809	Yes	Titanium Alloy (Ti6Al4V)
201-40-032	4.0mm x 32mm Cannulated Screw	00810028392816	Yes	Titanium Alloy (Ti6Al4V)
201-40-034	4.0mm x 34mm Cannulated Screw	00810028392823	Yes	Titanium Alloy (Ti6Al4V)
201-40-036	4.0mm x 36mm Cannulated Screw	00810028392830	Yes	Titanium Alloy (Ti6Al4V)
201-40-038	4.0mm x 38mm Cannulated Screw	00810028392847	Yes	Titanium Alloy (Ti6Al4V)
201-40-040	4.0mm x 40mm Cannulated Screw	00810028392854	Yes	Titanium Alloy (Ti6Al4V)
201-40-042	4.0mm x 42mm Cannulated Screw	00810028392861	Yes	Titanium Alloy (Ti6Al4V)
201-40-044	4.0mm x 44mm Cannulated Screw	00810028392878	Yes	Titanium Alloy (Ti6Al4V)
201-40-046	4.0mm x 46mm Cannulated Screw	00810028392885	Yes	Titanium Alloy (Ti6Al4V)
201-40-048	4.0mm x 48mm Cannulated Screw	00810028392892	Yes	Titanium Alloy (Ti6Al4V)
201-40-050	4.0mm x 50mm Cannulated Screw	00810028392908	Yes	Titanium Alloy (Ti6Al4V)
201-40-055	4.0mm x 55mm Cannulated Screw	00810028392915	Yes	Titanium Alloy (Ti6Al4V)
201-40-060	4.0mm x 60mm Cannulated Screw	00810028392922	Yes	Titanium Alloy (Ti6Al4V)
202-20-010	2.0 x 10mm Tiger Cann. Headless Screw	00812926021492	Yes	Titanium Alloy (Ti6Al4V)
202-20-012	2.0 x 12mm Tiger Cann. Headless Screw	00812926021508	Yes	Titanium Alloy (Ti6Al4V)
202-20-014	2.0 x 14mm Tiger Cann. Headless Screw	00812926021515	Yes	Titanium Alloy (Ti6Al4V)
202-20-016	2.0 x 16mm Tiger Cann. Headless Screw	00812926021522	Yes	Titanium Alloy (Ti6Al4V)
202-20-018	2.0 x 18mm Tiger Cann. Headless Screw	00812926021539	Yes	Titanium Alloy (Ti6Al4V)
202-20-020	2.0 x 20mm Tiger Cann. Headless Screw	00812926021546	Yes	Titanium Alloy (Ti6Al4V)
202-20-022	2.0 x 22mm Tiger Cann. Headless Screw	00812926021553	Yes	Titanium Alloy (Ti6Al4V)
202-20-024	2.0 x 24mm Tiger Cann Headless Screw	00812926021560	Yes	Titanium Alloy (Ti6Al4V)
202-24-010	2.4 x 10mm Tiger Cann. Headless Screw	00812926021577	Yes	Titanium Alloy (Ti6Al4V)
202-24-012	2.4 x 12mm Tiger Cann. Headless Screw	00812926021584	Yes	Titanium Alloy (Ti6Al4V)

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
202-24-014	2.4 x 14mm Tiger Cann. Headless Screw	00812926021591	Yes	Titanium Alloy (Ti6Al4V)
202-24-016	2.4 x 16mm Tiger Cann. Headless Screw	00812926021607	Yes	Titanium Alloy (Ti6Al4V)
202-24-018	2.4 x 18mm Tiger Cann. Headless Screw	00812926021614	Yes	Titanium Alloy (Ti6Al4V)
202-24-020	2.4 x 20mm Tiger Cann. Headless Screw	00812926021621	Yes	Titanium Alloy (Ti6Al4V)
202-24-022	2.4 x 22mm Tiger Cann. Headless Screw	00812926021638	Yes	Titanium Alloy (Ti6Al4V)
202-24-024	2.4 x 24mm Tiger Cann. Headless Screw	00812926021645	Yes	Titanium Alloy (Ti6Al4V)
202-24-026	2.4 x 26mm Tiger Cann. Headless Screw	00812926021652	Yes	Titanium Alloy (Ti6Al4V)
202-24-028	2.4 x 28mm Tiger Cann. Headless Screw	00812926021669	Yes	Titanium Alloy (Ti6Al4V)
202-24-030	2.4 x 30mm Tiger Cann. Headless Screw	00812926021676	Yes	Titanium Alloy (Ti6Al4V)
202-30-010	3.0 x 10mm Tiger Cann. Headless Screw	00812926021683	Yes	Titanium Alloy (Ti6Al4V)
202-30-012	3.0 x 12mm Tiger Cann. Headless Screw	00812926021690	Yes	Titanium Alloy (Ti6Al4V)
202-30-014	3.0 x 14mm Tiger Cann. Headless Screw	00812926021706	Yes	Titanium Alloy (Ti6Al4V)
202-30-016	3.0 x 16mm Tiger Cann. Headless Screw	00812926021713	Yes	Titanium Alloy (Ti6Al4V)
202-30-018	3.0 x 18mm Tiger Cann. Headless Screw	00812926021720	Yes	Titanium Alloy (Ti6Al4V)
202-30-020	3.0 x 20mm Tiger Cann. Headless Screw	00812926021737	Yes	Titanium Alloy (Ti6Al4V)
202-30-022	3.0 x 22mm Tiger Cann. Headless Screw	00812926021744	Yes	Titanium Alloy (Ti6Al4V)
202-30-024	3.0 x 24mm Tiger Cann. Headless Screw	00812926021751	Yes	Titanium Alloy (Ti6Al4V)
202-30-026	3.0 x 26mm Tiger Cann. Headless Screw	00812926021768	Yes	Titanium Alloy (Ti6Al4V)
202-30-028	3.0 x 28mm Tiger Cann. Headless Screw	00812926021775	Yes	Titanium Alloy (Ti6Al4V)
202-30-030	3.0 x 30mm Tiger Cann. Headless Screw	00812926021782	Yes	Titanium Alloy (Ti6Al4V)
202-30-032	3.0 x 32mm Tiger Cann. Headless Screw	00812926021799	Yes	Titanium Alloy (Ti6Al4V)
202-30-034	3.0 x 34mm Tiger Cann. Headless Screw	00812926021805	Yes	Titanium Alloy (Ti6Al4V)
202-45-020	4.5mm x 20mm Headless, STTHCS	00812926029221	Yes	Titanium Alloy (Ti6Al4V)
202-45-022	4.5mm x 22mm Headless, STTHCS	00812926029238	Yes	Titanium Alloy (Ti6Al4V)
202-45-024	4.5mm x 24mm Headless, STTHCS	00812926029245	Yes	Titanium Alloy (Ti6Al4V)
202-45-026	4.5mm x 26mm Headless, STTHCS	00812926029252	Yes	Titanium Alloy (Ti6Al4V)
202-45-028	4.5mm x 28mm Headless, STTHCS	00812926029269	Yes	Titanium Alloy (Ti6Al4V)
202-45-030	4.5mm x 30mm Headless, STTHCS	00812926029276	Yes	Titanium Alloy (Ti6Al4V)
202-45-032	4.5mm x 32mm Headless, STTHCS	00812926029283	Yes	Titanium Alloy (Ti6Al4V)
202-45-034	4.5mm x 34mm Headless, STTHCS	00812926029290	Yes	Titanium Alloy (Ti6Al4V)
202-45-036	4.5mm x 36mm Headless, STTHCS	00812926029306	Yes	Titanium Alloy (Ti6Al4V)
202-45-038	4.5mm x 38mm Headless, STTHCS	00812926029313	Yes	Titanium Alloy (Ti6Al4V)
202-45-040	4.5mm x 40mm Headless, STTHCS	00812926029320	Yes	Titanium Alloy (Ti6Al4V)
202-45-045	4.5mm x 45mm Headless, STTHCS	00812926029337	Yes	Titanium Alloy (Ti6Al4V)
202-45-050	4.5mm x 50mm Headless, STTHCS	00812926029344	Yes	Titanium Alloy (Ti6Al4V)
202-45-055	4.5mm x 55mm Headless, STTHCS	00812926029351	Yes	Titanium Alloy (Ti6Al4V)
202-45-060	4.5mm x 60mm Headless, STTHCS	00812926029368	Yes	Titanium Alloy (Ti6Al4V)
202-45-065	4.5mm x 65mm Headless, STTHCS	00812926029375	Yes	Titanium Alloy (Ti6Al4V)
202-45-070	4.5mm x 70mm Headless, STTHCS	00812926029382	Yes	Titanium Alloy (Ti6Al4V)
204-24-006	2.6 x 2.4 x 6mm Hammertoe Implant	00812926028699	Yes	Titanium Alloy (Ti6Al4V)
204-24-008	2.6 x 2.4 x 8mm Hammertoe Implant	00812926028705	Yes	Titanium Alloy (Ti6Al4V)
204-24-010	2.6 x 2.4 x 10mm Hammertoe Implant	00812926028712	Yes	Titanium Alloy (Ti6Al4V)
204-24-012	2.6 x 2.4 x 12mm Hammertoe Implant	00812926028729	Yes	Titanium Alloy (Ti6Al4V)
204-24-014	2.6 x 2.4 x 14mm Hammertoe Implant	00812926028736	Yes	Titanium Alloy (Ti6Al4V)
204-24-016	2.6 x 2.4 x 16mm Hammertoe Implant	00812926028743	Yes	Titanium Alloy (Ti6Al4V)
204-24-018	2.6 x 2.4 x 18mm Hammertoe Implant	00812926028750	Yes	Titanium Alloy (Ti6Al4V)

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
204-30-006	3.4 x 3.0 x 6mm Hammertoe Implant	00812926021935	Yes	Titanium Alloy (Ti6Al4V)
204-30-008	3.4 x 3.0 x 8mm Hammertoe Implant	00812926021942	Yes	Titanium Alloy (Ti6Al4V)
204-30-010	3.4 x 3.0 x 10mm Hammertoe Implant	00812926021959	Yes	Titanium Alloy (Ti6Al4V)
204-30-012	3.4 x 3.0 x 12mm Hammertoe Implant	00812926021966	Yes	Titanium Alloy (Ti6Al4V)
204-30-014	3.4 x 3.0 x 14mm Hammertoe Implant	00812926021973	Yes	Titanium Alloy (Ti6Al4V)
204-30-016	3.4 x 3.0 x 16mm Hammertoe Implant	00812926021980	Yes	Titanium Alloy (Ti6Al4V)
204-30-018	3.4 x 3.0 x 18mm Hammertoe Implant	00812926021997	Yes	Titanium Alloy (Ti6Al4V)
204-30-020	3.4 x 3.0 x 20mm Hammertoe Implant	00812926022000	Yes	Titanium Alloy (Ti6Al4V)
204-30-022	3.4 x 3.0 x 22mm Hammertoe Implant	00812926022017	Yes	Titanium Alloy (Ti6Al4V)
204-30-024	3.4 x 3.0 x 24mm Hammertoe Implant	00812926022024	Yes	Titanium Alloy (Ti6Al4V)
204-30-026	3.4 x 3.0 x 26mm Hammertoe Implant	00812926022031	Yes	Titanium Alloy (Ti6Al4V)
204-30-028	3.4 x 3.0 x 28mm Hammertoe Implant	00812926022048	Yes	Titanium Alloy (Ti6Al4V)
204-30-030	3.4 x 3.0 x 30mm Hammertoe Implant	00812926022055	Yes	Titanium Alloy (Ti6Al4V)
204-30-106	4.5 x 3.0 x 6mm Hammertoe Implant	00812926022062	Yes	Titanium Alloy (Ti6Al4V)
204-30-108	4.5 x 3.0 x 8mm Hammertoe Implant	00812926022079	Yes	Titanium Alloy (Ti6Al4V)
204-30-110	4.5 x 3.0 x 10mm Hammertoe Implant	00812926022086	Yes	Titanium Alloy (Ti6Al4V)
204-30-112	4.5 x 3.0 x 12mm Hammertoe Implant	00812926022093	Yes	Titanium Alloy (Ti6Al4V)
204-30-114	4.5 x 3.0 x 14mm Hammertoe Implant	00812926022109	Yes	Titanium Alloy (Ti6Al4V)
204-30-116	4.5 x 3.0 x 16mm Hammertoe Implant	00812926022116	Yes	Titanium Alloy (Ti6Al4V)
204-30-118	4.5 x 3.0 x 18mm Hammertoe Implant	00812926022123	Yes	Titanium Alloy (Ti6Al4V)
204-30-120	4.5 x 3.0 x 20mm Hammertoe Implant	00812926022130	Yes	Titanium Alloy (Ti6Al4V)
204-30-122	4.5 x 3.0 x 22mm Hammertoe Implant	00812926022147	Yes	Titanium Alloy (Ti6Al4V)
204-30-124	4.5 x 3.0 x 24mm Hammertoe Implant	00812926022154	Yes	Titanium Alloy (Ti6Al4V)
204-30-126	4.5 x 3.0 x 26mm Hammertoe Implant	00812926022161	Yes	Titanium Alloy (Ti6Al4V)
204-30-128	4.5 x 3.0 x 28mm Hammertoe Implant	00812926022178	Yes	Titanium Alloy (Ti6Al4V)
204-30-130	4.5 x 3.0 x 30mm Hammertoe Implant	00812926022185	Yes	Titanium Alloy (Ti6Al4V)
204-40-206	5.5mm x 4.0mm x 6mm Hammertoe Implant	00812926022192	Yes	Titanium Alloy (Ti6Al4V)
204-40-208	5.5mm x 4.0mm x 8mm Hammertoe Implant	00812926022208	Yes	Titanium Alloy (Ti6Al4V)
204-40-210	5.5mm x 4.0 x 10mm Hammertoe Implant	00812926022215	Yes	Titanium Alloy (Ti6Al4V)
205-55-001	5.5mm Smooth Washer	00812926025070	Yes	Titanium Alloy (Ti6Al4V)
205-55-030	5.5mm x 30mm STTCS	00812926024073	Yes	Titanium Alloy (Ti6Al4V)
205-55-035	5.5mm x 35mm STTCS	00812926024080	Yes	Titanium Alloy (Ti6Al4V)
205-55-040	5.5mm x 40mm STTCS	00812926024097	Yes	Titanium Alloy (Ti6Al4V)
205-55-045	5.5mm x 45mm STTCS	00812926024103	Yes	Titanium Alloy (Ti6Al4V)
205-55-050	5.5mm x 50mm STTCS	00812926024110	Yes	Titanium Alloy (Ti6Al4V)
205-55-055	5.5mm x 55mm STTCS	00812926024127	Yes	Titanium Alloy (Ti6Al4V)
205-55-060	5.5mm x 60mm STTCS	00812926024134	Yes	Titanium Alloy (Ti6Al4V)
205-55-065	5.5mm x 65mm STTCS	00812926024141	Yes	Titanium Alloy (Ti6Al4V)
205-55-070	5.5mm x 70mm STTCS	00812926024158	Yes	Titanium Alloy (Ti6Al4V)
205-55-075	5.5mm x 75mm STTCS	00812926024165	Yes	Titanium Alloy (Ti6Al4V)
205-55-080	5.5mm x 80mm STTCS	00812926024172	Yes	Titanium Alloy (Ti6Al4V)
205-70-001	7.0mm Smooth Washer	00812926025087	Yes	Titanium Alloy (Ti6Al4V)
205-70-050	7.0mm x 50mm STTCS	00812926024189	Yes	Titanium Alloy (Ti6Al4V)
205-70-055	7.0mm x 55mm STTCS	00812926024196	Yes	Titanium Alloy (Ti6Al4V)
205-70-060	7.0mm x 60mm STTCS	00812926024202	Yes	Titanium Alloy (Ti6Al4V)
205-70-065	7.0mm x 65mm STTCS	00812926024219	Yes	Titanium Alloy (Ti6Al4V)

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
205-70-070	7.0mm x 70mm STTCS	00812926024226	Yes	Titanium Alloy (Ti6Al4V)
205-70-075	7.0mm x 75mm STTCS	00812926024233	Yes	Titanium Alloy (Ti6Al4V)
205-70-080	7.0mm x 80mm STTCS	00812926024240	Yes	Titanium Alloy (Ti6Al4V)
205-70-085	7.0mm x 85mm STTCS	00812926024257	Yes	Titanium Alloy (Ti6Al4V)
205-70-090	7.0mm x 90mm STTCS	00812926024264	Yes	Titanium Alloy (Ti6Al4V)
205-70-095	7.0mm x 95mm STTCS	00812926024271	Yes	Titanium Alloy (Ti6Al4V)
205-70-100	7.0mm x 100mm STTCS	00812926024288	Yes	Titanium Alloy (Ti6Al4V)
205-70-105	7.0mm x 105mm STTCS	00812926024295	Yes	Titanium Alloy (Ti6Al4V)
205-70-110	7.0mm x 110mm STTCS	00812926024301	Yes	Titanium Alloy (Ti6Al4V)
205-70-115	7.0mm x 115mm STTCS	00812926024318	Yes	Titanium Alloy (Ti6Al4V)
205-70-120	7.0mm x 120mm STTCS	00812926024325	Yes	Titanium Alloy (Ti6Al4V)
206-55-050	5.5mm x 50mm LTTCS	00812926024332	Yes	Titanium Alloy (Ti6Al4V)
206-55-055	5.5mm x 55mm LTTCS	00812926024349	Yes	Titanium Alloy (Ti6Al4V)
206-55-060	5.5mm x 60mm LTTCS	00812926024356	Yes	Titanium Alloy (Ti6Al4V)
206-55-065	5.5mm x 65mm LTTCS	00812926024363	Yes	Titanium Alloy (Ti6Al4V)
206-55-070	5.5mm x 70mm LTTCS	00812926024370	Yes	Titanium Alloy (Ti6Al4V)
206-55-075	5.5mm x 75mm LTTCS	00812926024387	Yes	Titanium Alloy (Ti6Al4V)
206-55-080	5.5mm x 80mm LTTCS	00812926024394	Yes	Titanium Alloy (Ti6Al4V)
206-55-085	5.5mm x 85mm LTTCS	00812926024400	Yes	Titanium Alloy (Ti6Al4V)
206-55-090	5.5mm x 90mm LTTCS	00812926024417	Yes	Titanium Alloy (Ti6Al4V)
206-70-050	7.0mm x 50mm LTTCS	00812926024424	Yes	Titanium Alloy (Ti6Al4V)
206-70-055	7.0mm x 55mm LTTCS	00812926024431	Yes	Titanium Alloy (Ti6Al4V)
206-70-060	7.0mm x 60mm LTTCS	00812926024448	Yes	Titanium Alloy (Ti6Al4V)
206-70-065	7.0mm x 65mm LTTCS	00812926024455	Yes	Titanium Alloy (Ti6Al4V)
206-70-070	7.0mm x 70mm LTTCS	00812926024462	Yes	Titanium Alloy (Ti6Al4V)
206-70-075	7.0mm x 75mm LTTCS	00812926024479	Yes	Titanium Alloy (Ti6Al4V)
206-70-080	7.0mm x 80mm LTTCS	00812926024486	Yes	Titanium Alloy (Ti6Al4V)
206-70-085	7.0mm x 85mm LTTCS	00812926024493	Yes	Titanium Alloy (Ti6Al4V)
206-70-090	7.0mm x 90mm LTTCS	00812926024509	Yes	Titanium Alloy (Ti6Al4V)
206-70-095	7.0mm x 95mm LTTCS	00812926024516	Yes	Titanium Alloy (Ti6Al4V)
206-70-100	7.0mm x 100mm LTTCS	00812926024523	Yes	Titanium Alloy (Ti6Al4V)
206-70-105	7.0mm x 105mm LTTCS	00812926024530	Yes	Titanium Alloy (Ti6Al4V)
206-70-110	7.0mm x 110mm LTTCS	00812926024547	Yes	Titanium Alloy (Ti6Al4V)
206-70-115	7.0mm x 115mm LTTCS	00812926024554	Yes	Titanium Alloy (Ti6Al4V)
206-70-120	7.0mm x 120mm LTTCS	00812926024561	Yes	Titanium Alloy (Ti6Al4V)
207-55-030	5.5 x 30mm STHCS	00812926024578	Yes	Titanium Alloy (Ti6Al4V)
207-55-035	5.5 x 35mm STHCS	00812926024585	Yes	Titanium Alloy (Ti6Al4V)
207-55-040	5.5 x 40mm STHCS	00812926024592	Yes	Titanium Alloy (Ti6Al4V)
207-55-045	5.5 x 45mm STHCS	00812926024608	Yes	Titanium Alloy (Ti6Al4V)
207-55-050	5.5 x 50mm STHCS	00812926024615	Yes	Titanium Alloy (Ti6Al4V)
207-55-055	5.5 x 55mm STHCS	00812926024622	Yes	Titanium Alloy (Ti6Al4V)
207-55-060	5.5 x 60mm STHCS	00812926024639	Yes	Titanium Alloy (Ti6Al4V)
207-55-065	5.5 x 65mm STHCS	00812926024646	Yes	Titanium Alloy (Ti6Al4V)
207-55-070	5.5 x 70mm STHCS	00812926024653	Yes	Titanium Alloy (Ti6Al4V)
207-55-075	5.5 x 75mm STHCS	00812926024660	Yes	Titanium Alloy (Ti6Al4V)
207-55-080	5.5 x 80mm STHCS	00812926024677	Yes	Titanium Alloy (Ti6Al4V)



## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
207-70-045	7.0 x 45mm STTHCS	00812926027494	Yes	Titanium Alloy (Ti6Al4V)
207-70-050	7.0 x 50mm STTHCS	00812926024684	Yes	Titanium Alloy (Ti6Al4V)
207-70-055	7.0 x 55mm STTHCS	00812926024691	Yes	Titanium Alloy (Ti6Al4V)
207-70-060	7.0 x 60mm STTHCS	00812926024707	Yes	Titanium Alloy (Ti6Al4V)
207-70-065	7.0 x 65mm STTHCS	00812926024714	Yes	Titanium Alloy (Ti6Al4V)
207-70-070	7.0 x 70mm STTHCS	00812926024721	Yes	Titanium Alloy (Ti6Al4V)
207-70-075	7.0 x 75mm STTHCS	00812926024738	Yes	Titanium Alloy (Ti6Al4V)
207-70-080	7.0 x 80mm STTHCS	00812926024745	Yes	Titanium Alloy (Ti6Al4V)
207-70-085	7.0 x 85mm STTHCS	00812926024752	Yes	Titanium Alloy (Ti6Al4V)
207-70-090	7.0 x 90mm STTHCS	00812926024769	Yes	Titanium Alloy (Ti6Al4V)
207-70-095	7.0 x 95mm STTHCS	00812926024776	Yes	Titanium Alloy (Ti6Al4V)
207-70-100	7.0 x 100mm STTHCS	00812926024783	Yes	Titanium Alloy (Ti6Al4V)
207-70-105	7.0 x 105mm STTHCS	00812926024790	Yes	Titanium Alloy (Ti6Al4V)
207-70-110	7.0 x 110mm STTHCS	00812926024806	Yes	Titanium Alloy (Ti6Al4V)
207-70-115	7.0 x 115mm STTHCS	00812926024813	Yes	Titanium Alloy (Ti6Al4V)
207-70-120	7.0 x 120mm STTHCS	00812926024820	Yes	Titanium Alloy (Ti6Al4V)
208-55-050	5.5 x 50mm LTHCS	00812926024837	Yes	Titanium Alloy (Ti6Al4V)
208-55-055	5.5 x 55mm LTHCS	00812926024844	Yes	Titanium Alloy (Ti6Al4V)
208-55-060	5.5 x 60mm LTHCS	00812926024851	Yes	Titanium Alloy (Ti6Al4V)
208-55-065	5.5 x 65mm LTHCS	00812926024868	Yes	Titanium Alloy (Ti6Al4V)
208-55-070	5.5 x 70mm LTHCS	00812926024875	Yes	Titanium Alloy (Ti6Al4V)
208-55-075	5.5 x 75mm LTHCS	00812926024882	Yes	Titanium Alloy (Ti6Al4V)
208-55-080	5.5 x 80mm LTHCS	00812926024899	Yes	Titanium Alloy (Ti6Al4V)
208-55-085	5.5 x 85mm LTHCS	00812926024905	Yes	Titanium Alloy (Ti6Al4V)
208-55-090	5.5 x 90mm LTHCS	00812926024912	Yes	Titanium Alloy (Ti6Al4V)
208-70-050	7.0 x 50mm LTHCS	00812926024929	Yes	Titanium Alloy (Ti6Al4V)
208-70-055	7.0 x 55mm LTHCS	00812926024936	Yes	Titanium Alloy (Ti6Al4V)
208-70-060	7.0 x 60mm LTHCS	00812926024943	Yes	Titanium Alloy (Ti6Al4V)
208-70-065	7.0 x 65mm LTHCS	00812926024950	Yes	Titanium Alloy (Ti6Al4V)
208-70-070	7.0 x 70mm LTHCS	00812926024967	Yes	Titanium Alloy (Ti6Al4V)
208-70-075	7.0 x 75mm LTHCS	00812926024974	Yes	Titanium Alloy (Ti6Al4V)
208-70-080	7.0 x 80mm LTHCS	00812926024981	Yes	Titanium Alloy (Ti6Al4V)
208-70-085	7.0 x 85mm LTHCS	00812926024998	Yes	Titanium Alloy (Ti6Al4V)
208-70-090	7.0 x 90mm LTHCS	00812926025001	Yes	Titanium Alloy (Ti6Al4V)
208-70-095	7.0 x 95mm LTHCS	00812926025018	Yes	Titanium Alloy (Ti6Al4V)
208-70-100	7.0 x 100mm LTHCS	00812926025025	Yes	Titanium Alloy (Ti6Al4V)
208-70-105	7.0 x 105mm LTHCS	00812926025032	Yes	Titanium Alloy (Ti6Al4V)
208-70-110	7.0 x 110mm LTHCS	00812926025049	Yes	Titanium Alloy (Ti6Al4V)
208-70-115	7.0 x 115mm LTHCS	00812926025056	Yes	Titanium Alloy (Ti6Al4V)
208-70-120	7.0 x 120mm LTHCS	00812926025063	Yes	Titanium Alloy (Ti6Al4V)
210-00-001	Tiger Screw Sterilization Tray Base Assembly	00812926021324	No	Polypropylene, Polyphenylsulfone
210-00-003	Cannulated Driver Handle	00812926021331	No	Silicone, Stainless Steel, Aluminum Alloy
210-00-004	Ratcheting Cannulated Driver Handle	00812926021348	No	Silicone, Stainless Steel, Aluminum Alloy
210-00-005	Cannulated Depth Gauge	00812926021355	No	Aluminum

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
210-00-006	Screw Remover	00812926021362	No	Stainless Steel
210-00-007	Osteotomy Guide	00812926027012	No	Stainless Steel
210-00-008	4.5mm Screw Remover	00812926029504	No	Stainless Steel
210-00-101	Combined Tiger / Tiger Headless Screw Tray	00812926027029	No	Polypropylene, Aluminum
210-00-102	Compact Tiger Cannulated Screw Tray Base Assembly	00812926027036	No	Polypropylene, Polyphenylsulfone
210-00-103	4.5mm Tiger Cannulated Screw System Tray	00812926029535	No	Aluminum, Stainless Steel, Silicone
210-24-001	2.0/2.4mm Cannulated Drill Bit	00812926021379	No	Stainless Steel
210-24-002	2.0/2.4mm Cannulated Countersink	00812926021386	No	Stainless Steel
210-24-003	2.0/2.4mm Cannulated Screw Driver Bit	00812926021393	No	Stainless Steel
210-24-004	0.035" x 5" K-Wire Standard	00812926021409	Yes	Stainless Steel
210-24-005	0.035" x 5" K-Wire Partially Threaded	00812926021416	Yes	Stainless Steel
210-24-006	0.035" x 5" K-Wire Double Ended	00812926021423	Yes	Stainless Steel
210-24-007	2.0/2.4mm Tiger Solid Core Slip Fit Driver	00812926029399	No	Stainless Steel
210-24-008	2.0/2.4mm Tiger Cannulated Slip Fit Driver	00812926029405	No	Stainless Steel
210-24-014	.035" Marked K-wire	00812926027821	Yes	Stainless Steel
210-30-002	Combo Depth Gauge Countersink	00812926027814	No	Stainless Steel
210-40-001	3.0/4.0mm Cannulated Drill Bit	00812926021430	No	Stainless Steel
210-40-002	3.0/4.0mm Cann. Screw Countersink	00812926021447	No	Stainless Steel
210-40-003	3.0/4.0mm Cannulated Screw Driver Bit	00812926021454	No	Stainless Steel
210-40-004	0.045" x 6" K-Wire Standard	00812926021461	Yes	Stainless Steel
210-40-005	0.045" x 6" K-Wire Partially Threaded	00812926021478	Yes	Stainless Steel
210-40-006	0.045" x 6" K-Wire Double Ended	00812926021485	Yes	Stainless Steel
210-40-007	3.0/4.0mm Tiger Solid Core Slip Fit Driver	00812926029412	No	Stainless Steel
210-40-008	3.0/4.0mm Tiger Cannulated Slip Fit Driver	00812926029429	No	Stainless Steel
210-40-014	0.045" x 6" Marked K-Wire Standard	00812926027838	Yes	Stainless Steel
210-45-001	4.5mm Cannulated Drill Bit	00812926029436	No	Stainless Steel
210-45-002	4.5mm Cannulated Countersink	00812926029450	No	Stainless Steel
210-45-003	4.5mm Cannulated Driver	00812926029467	No	Stainless Steel
210-45-005	4.5mm Cannulated Depth Gauge	00812926029498	No	Aluminum
210-45-007	4.5mm Solid-Core Screw Remover	00812926029511	No	Stainless Steel
210-54-006	0.054" x 6" K-Wire Double Ended	00812926022314	Yes	Stainless Steel
210-60-004	0.062" x 6" K-Wire Standard	00812926023946	Yes	Stainless Steel
210-60-005	0.062" x 6" K-Wire Partially Threaded	00812926023953	Yes	Stainless Steel
210-60-006	0.062" x 6" K-Wire Double Ended	00812926023960	Yes	Stainless Steel
210-99-011	Tiger Cannulated Screw Tray Lid	00812926027876	No	Polyphenylsulfone
210-99-101	Combined Tiger / Tiger Headless Screw Tray Lid	00812926027883	No	Aluminum
210-99-103	0.062" K-Wire Tube	00810028391369	No	Aluminum
210-99-112	Tiger Screw Compact Lid	00812926027890	No	Polyphenylsulfone
211-40-001	4.0mm Cannulated Drill Bit	00810028391703	Yes	Stainless Steel
211-40-002	4.0mm Cannulated Countersink	00810028391710	No	Stainless Steel
211-40-004	Cannulated Depth Gauge	00810028391758	No	Aluminum
212-00-001	Tiger Headless Screw Sterilization Tray Base Assembly	00812926021812	No	Polypropylene, Aluminum
212-00-005	Cannulated Depth Gauge: Headless	00812926021829	No	Aluminum
212-20-002	2.0mm Cann. Prox. Drill	00812926021836	No	Stainless Steel
212-20-003	2.0mm Cann. Screw Driver Bit: Headless	00812926021843	No	Stainless Steel

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
212-20-007	2.0mm Removal Bit: Headless	00812926021850	No	Stainless Steel
212-24-002	2.4mm Cann. Prox. Drill	00812926021867	No	Stainless Steel
212-24-003	2.4mm Cann. Screw Driver Bit: Headless	00812926021874	No	Stainless Steel
212-24-007	2.4mm Removal Bit: Headless	00812926021881	No	Stainless Steel
212-30-002	3.0mm Cann. Prox. Drill	00812926021898	No	Stainless Steel
212-40-003	3.0/4.0mm Cann. Screw Driver Bit: Headless	00812926021911	No	Stainless Steel
212-40-007	3.0mm Removal Bit: Headless	00812926021928	No	Stainless Steel
212-45-002	4.5mm Cannulated Proximal Drill	00812926029474	No	Stainless Steel
212-99-011	Tiger Headless Cannulated Screw Tray Lid	00812926027906	No	Polyphenylsulfone
214-00-001	Hammer Toe Implant Tray Base Assembly	00812926022222	No	Polypropylene, Polyphenylsulfone, Aluminum
214-00-006	3.4mm Spade Pilot Drill	00812926022239	No	Stainless Steel
214-00-007	4.5mm Spade Pilot Drill	00812926022246	No	Stainless Steel
214-00-008	5.5mm Spade Pilot Drill	00812926022253	No	Stainless Steel
214-00-009	3.4mm Spade Driver	00812926022260	No	Stainless Steel
214-00-010	5.5mm Spade Driver	00812926022277	No	Stainless Steel
214-00-011	4.5mm Spade Driver	00812926022284	No	Stainless Steel
214-00-012	Cannulated Depth Gauge	00812926022291	No	Aluminum
214-00-013	3.0mm Thread Pilot Drill	00812926022307	No	Stainless Steel
214-00-014	Two-Step Hammer Toe Impactor	00812926027043	No	Stainless Steel
214-00-015	2.6mm Spade Pilot Drill	00812926028767	No	Stainless Steel
214-00-016	2.6mm Spade Driver	00812926028774	No	Stainless Steel
214-00-017	2.4mm Thread Pilot Drill	00812926028781	No	Stainless Steel
214-99-004	Hammer Toe Implant Tray Lid	00812926027913	No	Polyphenylsulfone
214-99-008	2.6mm x 2.4mm Two Step Pocket Insert	00812926029931	No	Radel
215-00-001	Tiger Large Cannulated Screw Tray Assembly	00812926025094	No	Aluminum, Stainless Steel, Polypropylene, Silicone, Polyphenylsulfone
215-00-003	Large Cannulated Fixed Driver Handle	00812926025117	No	Silicone, Stainless Steel, Aluminum Alloy
215-00-004	Large Cannulated Ratcheting Driver Handle	00812926025124	No	Silicone, Stainless Steel, Aluminum Alloy
215-00-005	Large Cannulated Depth Gauge	00812926025131	No	Aluminum
215-00-006	Large Cannulated Screw Sleeve	00812926025148	No	Stainless Steel
215-00-010	Large Cannulated Parallel Wire Guide	00812926025216	No	Stainless Steel
215-00-011	Large Cannulated Screw Extractor	00812926025223	No	Stainless Steel
215-00-012	1/4 SQ / Jacobs Adapter	00812926027050	No	Stainless Steel
215-55-001	5.5mm Large Cannulated Pilot Drill	00812926025230	No	Stainless Steel
215-55-002	5.5mm Large Cann. Countersink	00812926025247	No	Stainless Steel
215-55-003	5.5mm Large Cannulated Driver	00812926025254	No	Stainless Steel
215-55-004	0.094" x 12" K-Wire Standard	00812926027067	Yes	Stainless Steel
215-55-005	0.094" x 12" K-Wire Partially Threaded	00812926025315	Yes	Stainless Steel
215-55-007	5.5mm Screw Drill Sleeve	00812926025155	No	Stainless Steel
215-55-008	5.5mm Guide Pin Sleeve	00812926025162	No	Stainless Steel
215-55-009	5.5mm Screw Trocar	00812926025179	No	Stainless Steel
215-70-001	7.0mm Large Cannulated Pilot Drill	00812926025278	No	Stainless Steel
215-70-002	7.0mm Large Cann. Countersink	00812926025285	No	Stainless Steel
215-70-003	7.0mm Large Cannulated Driver	00812926025292	No	Stainless Steel

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
215-70-004	0.110" x 12" K-Wire Standard	00812926027074	Yes	Stainless Steel
215-70-005	0.110" x 12" K-Wire Partially Threaded	00812926025322	Yes	Stainless Steel
215-70-007	7.0mm Screw Drill Sleeve	00812926025186	No	Stainless Steel
215-70-008	7.0mm Guide Pin Sleeve	00812926025193	No	Stainless Steel
215-70-009	7.0mm Screw Trocar	00812926025209	No	Stainless Steel
217-55-002	5.5mm Cannulated Proximal Drill	00812926025261	No	Stainless Steel
217-55-003	5.5mm Headless Cannulated Driver	00812926027470	No	Stainless Steel
217-70-002	7.0mm Cannulated Proximal Drill	00812926025308	No	Stainless Steel
217-70-003	7.0mm Headless Cannulated Driver	00812926027487	No	Stainless Steel
220-08-001	8mm Concave Reamer	00812926022321	Yes	Stainless Steel
220-08-002	8mm Convex Reamer	00812926022338	Yes	Stainless Steel
220-08-003	8mm Sterile HTR Kit	00812926022345	Yes	Stainless Steel
220-08-004	8mm HTR Implant System, Non-Sterile	00812926022352	Yes	Stainless Steel
220-10-001	10mm CONCAVE REAMER	00812926022376	Yes	Stainless Steel
220-10-002	10mm CONVEX REAMER	00812926022383	Yes	Stainless Steel
220-10-003	10mm Sterile HTR Kit	00812926022390	Yes	Stainless Steel
220-10-004	HTR Implant, Non-Sterile	00812926022406	Yes	Stainless Steel
220-12-001	12mm Concave Reamer	00812926027098	Yes	Stainless Steel
220-12-002	12mm Convex Reamer	00812926027104	Yes	Stainless Steel
220-12-004	12mm HTR Implant System, Non-Sterile	00812926027128	Yes	Stainless Steel
300-10-001	Arsenal Ankle 8 Hole Anterior Tibia Plate	00810028393226	Yes	Titanium Alloy (Ti6Al4V)
300-10-002	Arsenal Ankle 10 Hole Anterior Tibia Plate	00810028393233	Yes	Titanium Alloy (Ti6Al4V)
300-11-001	Arsenal Ankle 7 Hole Anterolateral Tibia Plate, L	00810028393240	Yes	Titanium Alloy (Ti6Al4V)
300-11-002	Arsenal Ankle 7 Hole Anterolateral Tibia Plate, R	00810028393257	Yes	Titanium Alloy (Ti6Al4V)
300-11-003	Arsenal Ankle 9 Hole Anterolateral Tibia Plate, L	00810028393264	Yes	Titanium Alloy (Ti6Al4V)
300-11-004	Arsenal Ankle 9 Hole Anterolateral Tibia Plate, R	00810028393271	Yes	Titanium Alloy (Ti6Al4V)
300-11-005	Arsenal Ankle 11 Hole Anterolateral Tibia Plate, L	00810028393288	Yes	Titanium Alloy (Ti6Al4V)
300-11-006	Arsenal Ankle 11 Hole Anterolateral Tibia Plate, R	00810028393295	Yes	Titanium Alloy (Ti6Al4V)
300-11-007	Arsenal Ankle 13 Hole Anterolateral Tibia Plate, L	00810028393301	Yes	Titanium Alloy (Ti6Al4V)
300-11-008	Arsenal Ankle 13 Hole Anterolateral Tibia Plate, R	00810028393318	Yes	Titanium Alloy (Ti6Al4V)
300-11-009	Arsenal Ankle 16 Hole Anterolateral Tibia Plate, L	00810028393325	Yes	Titanium Alloy (Ti6Al4V)
300-11-010	Arsenal Ankle 16 Hole Anterolateral Tibia Plate, R	00810028393332	Yes	Titanium Alloy (Ti6Al4V)
300-12-001	Arsenal Ankle 4 Hole Posterior Malleolar Plate	00810028393202	Yes	Titanium Alloy (Ti6Al4V)
300-12-002	Arsenal Ankle 6 Hole Posterior Malleolar Plate	00810028393219	Yes	Titanium Alloy (Ti6Al4V)
300-50-001	2 Hole VLC Gridlock Plate	00812926022420	Yes	Titanium
300-50-002	3 Hole VLC Gridlock Plate	00812926022437	Yes	Titanium
300-50-003	4 Hole VLC Gridlock Plate	00812926022444	Yes	Titanium
300-50-004	Custom 4 Hole VLC Gridlock Plate	00812926022451	Yes	Titanium
300-50-005	5 Hole VLC Gridlock Plate	00812926022468	Yes	Titanium
300-50-006	6 Hole VLC Gridlock Plate	00812926022475	Yes	Titanium
300-50-007	7 Hole VLC Gridlock Plate	00812926025339	Yes	Titanium
300-50-008	Custom 5 Hole VLC Gridlock Plate	00812926025346	Yes	Titanium
300-51-001	Small MPX Gridlock Plate	00812926022482	Yes	Titanium
300-51-002	Medium MPX Gridlock Plate	00812926022499	Yes	Titanium
300-51-003	Large MPX Gridlock Plate	00812926022505	Yes	Titanium
300-52-001	Stepped MPJ VLC Gridlock Plate, Right	00812926022512	Yes	Titanium

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
300-52-002	Stepped MPJ VLC Gridlock Plate, Left	00812926022529	Yes	Titanium
300-52-003	Medium MPJ VLC Gridlock Plate	00812926022536	Yes	Titanium
300-52-004	Petite MPJ VLC Gridlock Plate, Right	00812926025353	Yes	Titanium
300-52-005	Petite MPJ VLC Gridlock Plate, Left	00812926025360	Yes	Titanium
300-53-001	Lapidus VLC Gridlock Plate, Right	00812926022543	Yes	Titanium
300-53-002	Lapidus VLC Gridlock Plate, Left	00812926022550	Yes	Titanium
300-53-003	Lapidus VLC Gridlock Plate	00812926022567	Yes	Titanium
300-54-002	2mm Wedge Gridlock Plate	00812926022574	Yes	Titanium
300-54-004	4mm Wedge Gridlock Plate	00812926022581	Yes	Titanium
300-54-006	6mm Wedge Gridlock Plate	00812926022598	Yes	Titanium
300-54-102	2mm Pre-Bent Wedge VL Gridlock Plate	00812926022604	Yes	Titanium
300-54-104	4mm Pre-Bent Wedge VL Gridlock Plate	00812926022611	Yes	Titanium
300-54-106	6mm Pre-Bent Wedge VL Gridlock Plate	00812926022628	Yes	Titanium
300-54-208	8mm Evans Plate	00812926027197	Yes	Titanium
300-54-210	10mm Evans Plate	00812926027203	Yes	Titanium
300-54-212	12mm Evans Plate	00812926027210	Yes	Titanium
300-55-003	4 Hole VL Gridlock Plate	00812926027227	Yes	Titanium
300-56-002	Small VLC "L" Plate, Right	00812926025377	Yes	Titanium
300-56-003	Small VLC "L" Plate, Left	00812926025384	Yes	Titanium
300-56-004	Large VLC "L" Plate, Right	00812926025391	Yes	Titanium
300-56-005	Large VLC "L" Plate, Left	00812926025407	Yes	Titanium
300-60-001	7 Hole VL Gridlock Fibula Plate	00812926025438	Yes	Titanium
300-60-002	9 Hole VL Gridlock Fibula Plate	00812926025445	Yes	Titanium
300-60-003	11 Hole VL Gridlock Fibula Plate	00812926025452	Yes	Titanium
300-61-001	7 Hole Contoured VL Gridlock Fibula Plate	00812926025469	Yes	Titanium
300-61-002	9 Hole Contoured VL Gridlock Fibula Plate	00812926025476	Yes	Titanium
300-61-003	11 Hole Contoured VL Gridlock Fibula Plate	00812926025483	Yes	Titanium
300-61-004	7 Hole VL Gridlock Distal Fibula Plate	00812926027562	Yes	Titanium
300-61-005	8 Hole Contoured VL Gridlock Fibula Plate	00812926027579	Yes	Titanium
300-61-006	10 Hole Contoured VL Gridlock Fibula Plate	00812926027586	Yes	Titanium
300-61-007	12 Hole Contoured Fibula Plate	00812926027593	Yes	Titanium
300-62-001	7 Hole Contour Tibia Plate	00812926025490	Yes	Titanium
300-62-002	9 Hole Contour Tibia Plate	00812926025506	Yes	Titanium
300-63-001	6 Hole 1/3 Tubular Plate	00812926027609	Yes	Titanium
300-63-002	9 Hole 1/3 Tubular Plate	00812926027616	Yes	Titanium
300-63-003	12 Hole 1/3 Tubular Plate	00812926027623	Yes	Titanium
300-64-001	6 Hole Malleolar Gridlock Hook Plate	00812926027647	Yes	Titanium
300-64-002	Buttress Plate	00812926027630	Yes	Titanium
300-65-001	8 Hole Gridlock Helical Plate, Left	00812926028484	Yes	Titanium
300-65-002	10 Hole Gridlock Helical Plate, Left	00812926028491	Yes	Titanium
300-65-003	8 Hole Gridlock Helical Plate, Right	00812926028507	Yes	Titanium
300-65-004	10 Hole Gridlock Helical Plate, Right	00812926028514	Yes	Titanium
300-70-001	Minimally Invasive Bunion Plate, Right	00812926027708	Yes	Titanium
300-70-002	Minimally Invasive Bunion Plate, Left	00812926027715	Yes	Titanium
300-80-001	Arsenal Custom 4 Hole VLC Plate	00812926029948	Yes	Titanium
300-80-002	Arsenal 2 Hole VLC Plate	00812926029955	Yes	Titanium

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
300-80-003	Arsenal 3 Hole VLC Plate	00812926029962	Yes	Titanium
300-80-004	Arsenal 4 Hole VLC Plate	00812926029979	Yes	Titanium
300-80-005	Arsenal 5 Hole VLC Plate	00812926029986	Yes	Titanium
300-80-006	Arsenal 6 Hole VLC Plate	00812926029993	Yes	Titanium
300-81-001	Arsenal Mini Condylar Plate, Short	00810028390003	Yes	Titanium
300-81-002	Arsenal Mini Condylar Plate, Long	00810028390010	Yes	Titanium
300-81-003	Arsenal 3 Hole Mini VLC Plate	00810028390027	Yes	Titanium
300-81-004	Arsenal 4 Hole Mini VLC Plate	00810028390034	Yes	Titanium
300-81-005	Arsenal 5 Hole Mini VLC Plate	00810028390041	Yes	Titanium
300-81-006	Arsenal 6 Hole Mini VLC Plate	00810028390058	Yes	Titanium
300-81-007	Arsenal Mini Z Plate, Short	00810028390065	Yes	Titanium
300-81-008	Arsenal Mini Z Plate, Long	00810028390072	Yes	Titanium
300-82-001	Arsenal Petite MPJ Plate, L	00810028390089	Yes	Titanium
300-82-002	Arsenal Petite MPJ Plate, R	00810028390096	Yes	Titanium
300-82-003	Arsenal Standard MPJ Plate, L	00810028390102	Yes	Titanium
300-82-004	Arsenal Standard MPJ Plate, R	00810028390119	Yes	Titanium
300-82-005	Arsenal Revision MPJ Plate, L	00810028390126	Yes	Titanium
300-82-006	Arsenal Revision MPJ Plate, R	00810028390133	Yes	Titanium
300-82-007	Arsenal Petite MPJ Plate, L, 0°	00810028391512	Yes	Titanium
300-82-008	Arsenal Petite MPJ Plate, R, 0°	00810028391529	Yes	Titanium
300-82-009	Arsenal Petite MPJ Plate, L, 5°	00810028391536	Yes	Titanium
300-82-010	Arsenal Petite MPJ Plate, R, 5°	00810028391543	Yes	Titanium
300-83-001	Arsenal Petite Lapidus Plate, L	00810028390140	Yes	Titanium
300-83-002	Arsenal Petite Lapidus Plate, R	00810028390157	Yes	Titanium
300-83-003	Arsenal Standard Lapidus Plate, L	00810028390164	Yes	Titanium
300-83-004	Arsenal Standard Lapidus Plate, R	00810028390171	Yes	Titanium
300-83-005	Arsenal Large Lapidus Plate, L	00810028390188	Yes	Titanium
300-83-006	Arsenal Large Lapidus Plate, R	00810028390195	Yes	Titanium
300-84-001	Arsenal 5th Met Hook Plate, L	00810028390201	Yes	Titanium
300-84-002	Arsenal 5th Met Hook Plate, R	00810028390218	Yes	Titanium
300-84-003	Arsenal Short 5th Met Hook Plate, L	00810028391550	Yes	Titanium
300-84-004	Arsenal Short 5th Met Hook Plate, R	00810028391567	Yes	Titanium
300-85-001	Arsenal Lisfranc Plate	00810028390225	Yes	Titanium
300-85-002	Arsenal TN Plate	00810028390232	Yes	Titanium
300-85-003	Arsenal NC Plate, Left	00810028390249	Yes	Titanium
300-85-004	Arsenal NC Plate, Right	00810028391345	Yes	Titanium
300-85-005	Arsenal U-Plate, Small, L	00810028391574	Yes	Titanium
300-85-006	Arsenal U-Plate, Small, R	00810028391581	Yes	Titanium
300-85-007	Arsenal U-Plate, Standard, L	00810028391598	Yes	Titanium
300-85-008	Arsenal U-Plate, Standard, R	00810028391604	Yes	Titanium
300-86-001	Arsenal Evans Plate	00810028390256	Yes	Titanium
300-86-002	Arsenal Universal Calc Slide Plate	00810028390263	Yes	Titanium
300-86-003	Arsenal Calc Perimeter Plate, Left	00812926029825	Yes	Titanium
300-86-004	Arsenal Calc Perimeter Plate, Right	00812926029832	Yes	Titanium
300-86-005	Arsenal Sinus Tarsi Calc Plate, L	00812926029849	Yes	Titanium
300-86-006	Arsenal Sinus Tarsi Calc Plate, R	00812926029856	Yes	Titanium

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
300-86-007	Arsenal Sinus Support Plate	00812926029863	Yes	Titanium
300-87-001	Arsenal Short Medial Column Plate, Left	00812926029870	Yes	Titanium
300-87-002	Arsenal Short Medial Column Plate, Right	00812926029887	Yes	Titanium
300-87-003	Arsenal Long Medial Column Plate, Left	00812926029894	Yes	Titanium
300-87-004	Arsenal Long Medial Column Plate, Right	00812926029900	Yes	Titanium
300-88-001	Arsenal L Plate, Left	00810028390270	Yes	Titanium
300-88-002	Arsenal L Plate, Right	00810028390287	Yes	Titanium
300-90-001	Arsenal Ankle 6 Hole 1/3 Tubular Plate	00810028392939	Yes	Titanium Alloy (Ti6Al4V)
300-90-002	Arsenal Ankle 8 Hole 1/3 Tubular Plate	00810028392946	Yes	Titanium Alloy (Ti6Al4V)
300-90-003	Arsenal Ankle 10 Hole 1/3 Tubular Plate	00810028392953	Yes	Titanium Alloy (Ti6Al4V)
300-90-004	Arsenal Ankle 12 Hole 1/3 Tubular Plate	00810028392960	Yes	Titanium Alloy (Ti6Al4V)
300-91-001	Arsenal Ankle Buttress Plate	00810028392977	Yes	Titanium Alloy (Ti6Al4V)
300-91-002	Arsenal Ankle 6 Hole Straight Plate	00810028392984	Yes	Titanium Alloy (Ti6Al4V)
300-91-003	Arsenal Ankle 8 Hole Straight Plate	00810028392991	Yes	Titanium Alloy (Ti6Al4V)
300-91-004	Arsenal Ankle 10 Hole Straight Plate	00810028393004	Yes	Titanium Alloy (Ti6Al4V)
300-91-005	Arsenal Ankle 12 Hole Straight Plate	00810028393011	Yes	Titanium Alloy (Ti6Al4V)
300-92-001	Arsenal Ankle 8 Hole Lateral Fibula Plate, Symmetric	00810028393028	Yes	Titanium Alloy (Ti6Al4V)
300-92-002	Arsenal Ankle 10 Hole Lateral Fibula Plate, Symmetric	00810028393035	Yes	Titanium Alloy (Ti6Al4V)
300-92-003	Arsenal 12 Hole Lateral Fibula Plate, Symmetric	00810028393042	Yes	Titanium Alloy (Ti6Al4V)
300-93-001	Arsenal Ankle 6 Hole Lateral Fibula Plate, L	00810028393059	Yes	Titanium Alloy (Ti6Al4V)
300-93-002	Arsenal Ankle 6 Hole Lateral Fibula Plate, R	00810028393066	Yes	Titanium Alloy (Ti6Al4V)
300-93-003	Arsenal Ankle 8 Hole Lateral Fibula Plate, L	00810028393073	Yes	Titanium Alloy (Ti6Al4V)
300-93-004	Arsenal Ankle 8 Hole Lateral Fibula Plate, R	00810028393080	Yes	Titanium Alloy (Ti6Al4V)
300-93-005	Arsenal Ankle 10 Hole Lateral Fibula Plate, L	00810028393097	Yes	Titanium Alloy (Ti6Al4V)
300-93-006	Arsenal Ankle 10 Hole Lateral Fibula Plate, R	00810028393103	Yes	Titanium Alloy (Ti6Al4V)
300-94-001	Arsenal Ankle 6 Hole Posterior Fibula Plate, L	00810028393110	Yes	Titanium Alloy (Ti6Al4V)
300-94-002	Arsenal Ankle 6 Hole Posterior Fibula Plate, R	00810028393127	Yes	Titanium Alloy (Ti6Al4V)
300-94-003	Arsenal Ankle 9 Hole Posterior Fibula Plate, L	00810028393134	Yes	Titanium Alloy (Ti6Al4V)
300-94-004	Arsenal Ankle 9 Hole Posterior Fibula Plate, R	00810028393141	Yes	Titanium Alloy (Ti6Al4V)
300-95-001	Arsenal Ankle 7 Hole Medial Malleolar Plate	00810028393158	Yes	Titanium Alloy (Ti6Al4V)
300-95-002	Arsenal Ankle 9 Hole Medial Malleolar Plate	00810028393165	Yes	Titanium Alloy (Ti6Al4V)
300-96-001	Arsenal Ankle Tension Band Plate	00810028393172	Yes	Titanium Alloy (Ti6Al4V)
300-97-001	Arsenal Ankle 4 Hole Hook Plate	00810028393189	Yes	Titanium Alloy (Ti6Al4V)
300-97-002	Arsenal Ankle 6 Hole Hook Plate	00810028393196	Yes	Titanium Alloy (Ti6Al4V)
301-24-010	2.4mm x 10mm GL Screw	00812926022635	Yes	Titanium Alloy (Ti6Al4V)
301-24-012	2.4mm x 12mm GL Screw	00812926022642	Yes	Titanium Alloy (Ti6Al4V)
301-24-014	2.4mm x 14mm GL Screw	00812926022659	Yes	Titanium Alloy (Ti6Al4V)
301-24-016	2.4mm x 16mm GL Screw	00812926022666	Yes	Titanium Alloy (Ti6Al4V)
301-24-018	2.4mm x 18mm GL Screw	00812926022673	Yes	Titanium Alloy (Ti6Al4V)
301-24-020	2.4mm x 20mm GL Screw	00812926022680	Yes	Titanium Alloy (Ti6Al4V)
301-24-022	2.4mm x 22mm GL Screw	00812926022697	Yes	Titanium Alloy (Ti6Al4V)
301-24-024	2.4mm x 24mm GL Screw	00812926022703	Yes	Titanium Alloy (Ti6Al4V)
301-24-026	2.4mm x 26mm GL Screw	00812926022710	Yes	Titanium Alloy (Ti6Al4V)
301-24-028	2.4mm x 28mm GL Screw	00812926022727	Yes	Titanium Alloy (Ti6Al4V)
301-24-030	2.4mm x 30mm GL Screw	00812926022734	Yes	Titanium Alloy (Ti6Al4V)
301-24-032	2.4mm x 32mm GL Screw	00812926022741	Yes	Titanium Alloy (Ti6Al4V)

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
301-24-034	2.4mm x 34mm GL Screw	00812926022758	Yes	Titanium Alloy (Ti6Al4V)
301-24-036	2.4mm x 36mm GL Screw	00812926022765	Yes	Titanium Alloy (Ti6Al4V)
301-24-038	2.4mm x 38mm GL Screw	00812926022772	Yes	Titanium Alloy (Ti6Al4V)
301-24-040	2.4mm x 40mm GL Screw	00812926022789	Yes	Titanium Alloy (Ti6Al4V)
301-30-010	3.0mm x 10mm GL Screw	00812926022796	Yes	Titanium Alloy (Ti6Al4V)
301-30-012	3.0mm x 12mm GL Screw	00812926022802	Yes	Titanium Alloy (Ti6Al4V)
301-30-014	3.0mm x 14mm GL Screw	00812926022819	Yes	Titanium Alloy (Ti6Al4V)
301-30-016	3.0mm x 16mm GL Screw	00812926022826	Yes	Titanium Alloy (Ti6Al4V)
301-30-018	3.0mm x 18mm GL Screw	00812926022833	Yes	Titanium Alloy (Ti6Al4V)
301-30-020	3.0mm x 20mm GL Screw	00812926022840	Yes	Titanium Alloy (Ti6Al4V)
301-30-022	3.0mm x 22mm GL Screw	00812926022857	Yes	Titanium Alloy (Ti6Al4V)
301-30-024	3.0mm x 24mm GL Screw	00812926022864	Yes	Titanium Alloy (Ti6Al4V)
301-30-026	3.0mm x 26mm GL Screw	00812926022871	Yes	Titanium Alloy (Ti6Al4V)
301-30-028	3.0mm x 28mm GL Screw	00812926022888	Yes	Titanium Alloy (Ti6Al4V)
301-30-030	3.0mm x 30mm GL Screw	00812926022895	Yes	Titanium Alloy (Ti6Al4V)
301-30-032	3.0mm x 32mm GL Screw	00812926022901	Yes	Titanium Alloy (Ti6Al4V)
301-30-034	3.0mm x 34mm GL Screw	00812926022918	Yes	Titanium Alloy (Ti6Al4V)
301-30-036	3.0mm x 36mm GL Screw	00812926022925	Yes	Titanium Alloy (Ti6Al4V)
301-30-038	3.0mm x 38mm GL Screw	00812926022932	Yes	Titanium Alloy (Ti6Al4V)
301-30-040	3.0mm x 40mm GL Screw	00812926022949	Yes	Titanium Alloy (Ti6Al4V)
301-40-010	4.0mm x 10mm GL Screw	00812926022956	Yes	Titanium Alloy (Ti6Al4V)
301-40-012	4.0mm x 12mm GL Screw	00812926022963	Yes	Titanium Alloy (Ti6Al4V)
301-40-014	4.0mm x 14mm GL Screw	00812926022970	Yes	Titanium Alloy (Ti6Al4V)
301-40-016	4.0mm x 16mm GL Screw	00812926022987	Yes	Titanium Alloy (Ti6Al4V)
301-40-018	4.0mm x 18mm GL Screw	00812926022994	Yes	Titanium Alloy (Ti6Al4V)
301-40-020	4.0mm x 20mm GL Screw	00812926023007	Yes	Titanium Alloy (Ti6Al4V)
301-40-022	4.0mm x 22mm GL Screw	00812926023014	Yes	Titanium Alloy (Ti6Al4V)
301-40-024	4.0mm x 24mm GL Screw	00812926023021	Yes	Titanium Alloy (Ti6Al4V)
301-40-026	4.0mm x 26mm GL Screw	00812926023038	Yes	Titanium Alloy (Ti6Al4V)
301-40-028	4.0mm x 28mm GL Screw	00812926023045	Yes	Titanium Alloy (Ti6Al4V)
301-40-030	4.0mm x 30mm GL Screw	00812926023052	Yes	Titanium Alloy (Ti6Al4V)
301-40-032	4.0mm x 32mm GL Screw	00812926023069	Yes	Titanium Alloy (Ti6Al4V)
301-40-034	4.0mm x 34mm GL Screw	00812926023076	Yes	Titanium Alloy (Ti6Al4V)
301-40-036	4.0mm x 36mm GL Screw	00812926023083	Yes	Titanium Alloy (Ti6Al4V)
301-40-038	4.0mm x 38mm GL Screw	00812926023090	Yes	Titanium Alloy (Ti6Al4V)
301-40-040	4.0mm x 40mm GL Screw	00812926023106	Yes	Titanium Alloy (Ti6Al4V)
302-24-010	2.4mm x 10mm GL Locking Screw	00812926023113	Yes	Titanium Alloy (Ti6Al4V)
302-24-012	2.4mm x 12mm GL Locking Screw	00812926023120	Yes	Titanium Alloy (Ti6Al4V)
302-24-014	2.4mm x 14mm GL Locking Screw	00812926023137	Yes	Titanium Alloy (Ti6Al4V)
302-24-016	2.4mm x 16mm GL Locking Screw	00812926023144	Yes	Titanium Alloy (Ti6Al4V)
302-24-018	2.4mm x 18mm GL Locking Screw	00812926023151	Yes	Titanium Alloy (Ti6Al4V)
302-24-020	2.4mm x 20mm GL Locking Screw	00812926023168	Yes	Titanium Alloy (Ti6Al4V)
302-24-022	2.4mm x 22mm GL Locking Screw	00812926023175	Yes	Titanium Alloy (Ti6Al4V)
302-24-024	2.4mm x 24mm GL Locking Screw	00812926023182	Yes	Titanium Alloy (Ti6Al4V)
302-24-026	2.4mm x 26mm GL Locking Screw	00812926023199	Yes	Titanium Alloy (Ti6Al4V)
302-24-028	2.4mm x 28mm GL Locking Screw	00812926023205	Yes	Titanium Alloy (Ti6Al4V)



## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
302-24-030	2.4mm x 30mm GL Locking Screw	00812926023212	Yes	Titanium Alloy (Ti6Al4V)
302-24-032	2.4mm x 32mm GL Locking Screw	00812926023229	Yes	Titanium Alloy (Ti6Al4V)
302-24-034	2.4mm x 34mm GL Locking Screw	00812926023236	Yes	Titanium Alloy (Ti6Al4V)
302-24-036	2.4mm x 36mm GL Locking Screw	00812926023243	Yes	Titanium Alloy (Ti6Al4V)
302-24-038	2.4mm x 38mm GL Locking Screw	00812926023250	Yes	Titanium Alloy (Ti6Al4V)
302-24-040	2.4mm x 40mm GL Locking Screw	00812926023267	Yes	Titanium Alloy (Ti6Al4V)
302-30-010	3.0mm x 10mm GL Locking Screw	00812926023274	Yes	Titanium Alloy (Ti6Al4V)
302-30-012	3.0mm x 12mm GL Locking Screw	00812926023281	Yes	Titanium Alloy (Ti6Al4V)
302-30-014	3.0mm x 14mm GL Locking Screw	00812926023298	Yes	Titanium Alloy (Ti6Al4V)
302-30-016	3.0mm x 16mm GL Locking Screw	00812926023304	Yes	Titanium Alloy (Ti6Al4V)
302-30-018	3.0mm x 18mm GL Locking Screw	00812926023311	Yes	Titanium Alloy (Ti6Al4V)
302-30-020	3.0mm x 20mm GL Locking Screw	00812926023328	Yes	Titanium Alloy (Ti6Al4V)
302-30-022	3.0mm x 22mm GL Locking Screw	00812926023335	Yes	Titanium Alloy (Ti6Al4V)
302-30-024	3.0mm x 24mm GL Locking Screw	00812926023342	Yes	Titanium Alloy (Ti6Al4V)
302-30-026	3.0mm x 26mm GL Locking Screw	00812926023359	Yes	Titanium Alloy (Ti6Al4V)
302-30-028	3.0mm x 28mm GL Locking Screw	00812926023366	Yes	Titanium Alloy (Ti6Al4V)
302-30-030	3.0mm x 30mm GL Locking Screw	00812926023373	Yes	Titanium Alloy (Ti6Al4V)
302-30-032	3.0mm x 32mm GL Locking Screw	00812926023380	Yes	Titanium Alloy (Ti6Al4V)
302-30-034	3.0mm x 34mm GL Locking Screw	00812926023397	Yes	Titanium Alloy (Ti6Al4V)
302-30-036	3.0mm x 36mm GL Locking Screw	00812926023403	Yes	Titanium Alloy (Ti6Al4V)
302-30-038	3.0mm x 38mm GL Locking Screw	00812926023410	Yes	Titanium Alloy (Ti6Al4V)
302-30-040	3.0mm x 40mm GL Locking Screw	00812926023427	Yes	Titanium Alloy (Ti6Al4V)
302-40-010	4.0mm x 10mm GL Locking Screw	00812926023434	Yes	Titanium Alloy (Ti6Al4V)
302-40-012	4.0mm x 12mm GL Locking Screw	00812926023441	Yes	Titanium Alloy (Ti6Al4V)
302-40-014	4.0mm x 14mm GL Locking Screw	00812926023458	Yes	Titanium Alloy (Ti6Al4V)
302-40-016	4.0mm x 16mm GL Locking Screw	00812926023465	Yes	Titanium Alloy (Ti6Al4V)
302-40-018	4.0mm x 18mm GL Locking Screw	00812926023472	Yes	Titanium Alloy (Ti6Al4V)
302-40-020	4.0mm x 20mm GL Locking Screw	00812926023489	Yes	Titanium Alloy (Ti6Al4V)
302-40-022	4.0mm x 22mm GL Locking Screw	00812926023496	Yes	Titanium Alloy (Ti6Al4V)
302-40-024	4.0mm x 24mm GL Locking Screw	00812926023502	Yes	Titanium Alloy (Ti6Al4V)
302-40-026	4.0mm x 26mm GL Locking Screw	00812926023519	Yes	Titanium Alloy (Ti6Al4V)
302-40-028	4.0mm x 28mm GL Locking Screw	00812926023526	Yes	Titanium Alloy (Ti6Al4V)
302-40-030	4.0mm x 30mm GL Locking Screw	00812926023533	Yes	Titanium Alloy (Ti6Al4V)
302-40-032	4.0mm x 32mm GL Locking Screw	00812926023540	Yes	Titanium Alloy (Ti6Al4V)
302-40-034	4.0mm x 34mm GL Locking Screw	00812926023557	Yes	Titanium Alloy (Ti6Al4V)
302-40-036	4.0mm x 36mm GL Locking Screw	00812926023564	Yes	Titanium Alloy (Ti6Al4V)
302-40-038	4.0mm x 38mm GL Locking Screw	00812926023571	Yes	Titanium Alloy (Ti6Al4V)
302-40-040	4.0mm x 40mm GL Locking Screw	00812926023588	Yes	Titanium Alloy (Ti6Al4V)
303-50-001	Smooth GL Screw Washer	00812926023595	Yes	Titanium Alloy (Ti6Al4V)
303-50-002	Spiked GL Screw Washer	00812926023601	Yes	Titanium Alloy (Ti6Al4V)
303-50-003	Large Smooth GL Screw Washer	00812926023618	Yes	Titanium Alloy (Ti6Al4V)
303-50-004	12mm GL Spiked Washer	00812926023625	Yes	Titanium Alloy (Ti6Al4V)
303-60-001	Smooth Gridlock Ankle Screw Washer	00812926026633	Yes	Titanium Alloy (Ti6Al4V)
303-60-002	Spiked Gridlock Ankle Screw Washer	00812926026640	Yes	Titanium Alloy (Ti6Al4V)
303-90-001	Arsenal 2.7/3.5mm/4.0mm Screw Washer, Smooth	00810028392106	Yes	Titanium Alloy (Ti6Al4V)
303-90-002	Arsenal 2.7/3.5mm/4.0mm Screw Washer, Ridged	00810028392113	Yes	Titanium Alloy (Ti6Al4V)

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
304-27-010	2.7mm x 10mm Gridlock Ankle Screw	00812926025513	Yes	Titanium Alloy (Ti6Al4V)
304-27-012	2.7mm x 12mm Gridlock Ankle Screw	00812926025520	Yes	Titanium Alloy (Ti6Al4V)
304-27-014	2.7mm x 14mm Gridlock Ankle Screw	00812926025537	Yes	Titanium Alloy (Ti6Al4V)
304-27-016	2.7mm x 16mm Gridlock Ankle Screw	00812926025544	Yes	Titanium Alloy (Ti6Al4V)
304-27-018	2.7mm x 18mm Gridlock Ankle Screw	00812926025551	Yes	Titanium Alloy (Ti6Al4V)
304-27-020	2.7mm x 20mm Gridlock Ankle Screw	00812926025568	Yes	Titanium Alloy (Ti6Al4V)
304-27-022	2.7mm x 22mm Gridlock Ankle Screw	00812926025575	Yes	Titanium Alloy (Ti6Al4V)
304-27-024	2.7mm x 24mm Gridlock Ankle Screw	00812926025582	Yes	Titanium Alloy (Ti6Al4V)
304-35-010	3.5mm x 10mm Gridlock Ankle Screw	00812926025599	Yes	Titanium Alloy (Ti6Al4V)
304-35-012	3.5mm x 12mm Gridlock Ankle Screw	00812926025605	Yes	Titanium Alloy (Ti6Al4V)
304-35-014	3.5mm x 14mm Gridlock Ankle Screw	00812926025612	Yes	Titanium Alloy (Ti6Al4V)
304-35-016	3.5mm x 16mm Gridlock Ankle Screw	00812926025629	Yes	Titanium Alloy (Ti6Al4V)
304-35-018	3.5mm x 18mm Gridlock Ankle Screw	00812926025636	Yes	Titanium Alloy (Ti6Al4V)
304-35-020	3.5mm x 20mm Gridlock Ankle Screw	00812926025643	Yes	Titanium Alloy (Ti6Al4V)
304-35-022	3.5mm x 22mm Gridlock Ankle Screw	00812926025650	Yes	Titanium Alloy (Ti6Al4V)
304-35-024	3.5mm x 24mm Gridlock Ankle Screw	00812926025667	Yes	Titanium Alloy (Ti6Al4V)
304-35-026	3.5mm x 26mm Gridlock Ankle Screw	00812926025674	Yes	Titanium Alloy (Ti6Al4V)
304-35-028	3.5mm x 28mm Gridlock Ankle Screw	00812926025681	Yes	Titanium Alloy (Ti6Al4V)
304-35-030	3.5mm x 30mm Gridlock Ankle Screw	00812926025698	Yes	Titanium Alloy (Ti6Al4V)
304-35-032	3.5mm x 32mm Gridlock Ankle Screw	00812926025704	Yes	Titanium Alloy (Ti6Al4V)
304-35-034	3.5mm x 34mm Gridlock Ankle Screw	00812926025711	Yes	Titanium Alloy (Ti6Al4V)
304-35-036	3.5mm x 36mm Gridlock Ankle Screw	00812926025728	Yes	Titanium Alloy (Ti6Al4V)
304-35-038	3.5mm x 38mm Gridlock Ankle Screw	00812926025735	Yes	Titanium Alloy (Ti6Al4V)
304-35-040	3.5mm x 40mm Gridlock Ankle Screw	00812926025742	Yes	Titanium Alloy (Ti6Al4V)
304-35-042	3.5mm x 42mm Gridlock Ankle Screw	00812926025759	Yes	Titanium Alloy (Ti6Al4V)
304-35-044	3.5mm x 44mm Gridlock Ankle Screw	00812926025766	Yes	Titanium Alloy (Ti6Al4V)
304-35-046	3.5mm x 46mm Gridlock Ankle Screw	00812926025773	Yes	Titanium Alloy (Ti6Al4V)
304-35-048	3.5mm x 48mm Gridlock Ankle Screw	00812926025780	Yes	Titanium Alloy (Ti6Al4V)
304-35-050	3.5mm x 50mm Gridlock Ankle Screw	00812926025797	Yes	Titanium Alloy (Ti6Al4V)
304-40-010	4.0mm x 10mm Gridlock Ankle Screw	00812926027692	Yes	Titanium Alloy (Ti6Al4V)
304-40-012	4.0mm x 12mm Gridlock Ankle Screw	00812926025803	Yes	Titanium Alloy (Ti6Al4V)
304-40-014	4.0mm x 14mm Gridlock Ankle Screw	00812926025810	Yes	Titanium Alloy (Ti6Al4V)
304-40-016	4.0mm x 16mm Gridlock Ankle Screw	00812926025827	Yes	Titanium Alloy (Ti6Al4V)
304-40-018	4.0mm x 18mm Gridlock Ankle Screw	00812926025834	Yes	Titanium Alloy (Ti6Al4V)
304-40-020	4.0mm x 20mm Gridlock Ankle Screw	00812926025841	Yes	Titanium Alloy (Ti6Al4V)
304-40-022	4.0mm x 22mm Gridlock Ankle Screw	00812926025858	Yes	Titanium Alloy (Ti6Al4V)
304-40-024	4.0mm x 24mm Gridlock Ankle Screw	00812926025865	Yes	Titanium Alloy (Ti6Al4V)
304-40-026	4.0mm x 26mm Gridlock Ankle Screw	00812926025872	Yes	Titanium Alloy (Ti6Al4V)
304-40-028	4.0mm x 28mm Gridlock Ankle Screw	00812926025889	Yes	Titanium Alloy (Ti6Al4V)
304-40-030	4.0mm x 30mm Gridlock Ankle Screw	00812926025896	Yes	Titanium Alloy (Ti6Al4V)
304-40-032	4.0mm x 32mm Gridlock Ankle Screw	00812926025902	Yes	Titanium Alloy (Ti6Al4V)
304-40-034	4.0mm x 34mm Gridlock Ankle Screw	00812926025919	Yes	Titanium Alloy (Ti6Al4V)
304-40-036	4.0mm x 36mm Gridlock Ankle Screw	00812926025926	Yes	Titanium Alloy (Ti6Al4V)
304-40-038	4.0mm x 38mm Gridlock Ankle Screw	00812926025933	Yes	Titanium Alloy (Ti6Al4V)
304-40-040	4.0mm x 40mm Gridlock Ankle Screw	00812926025940	Yes	Titanium Alloy (Ti6Al4V)
304-40-042	4.0mm x 42mm Gridlock Ankle Screw	00812926025957	Yes	Titanium Alloy (Ti6Al4V)

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
304-40-044	4.0mm x 44mm Gridlock Ankle Screw	00812926025964	Yes	Titanium Alloy (Ti6Al4V)
304-40-046	4.0mm x 46mm Gridlock Ankle Screw	00812926025971	Yes	Titanium Alloy (Ti6Al4V)
304-40-048	4.0mm x 48mm Gridlock Ankle Screw	00812926025988	Yes	Titanium Alloy (Ti6Al4V)
304-40-050	4.0mm x 50mm Gridlock Ankle Screw	00812926025995	Yes	Titanium Alloy (Ti6Al4V)
305-27-010	2.7mm x 10mm Gridlock Ankle Locking Screw	00812926026008	Yes	Titanium Alloy (Ti6Al4V)
305-27-012	2.7mm x 12mm Gridlock Ankle Locking Screw	00812926026015	Yes	Titanium Alloy (Ti6Al4V)
305-27-014	2.7mm x 14mm Gridlock Ankle Locking Screw	00812926026022	Yes	Titanium Alloy (Ti6Al4V)
305-27-016	2.7mm x 16mm Gridlock Ankle Locking Screw	00812926026039	Yes	Titanium Alloy (Ti6Al4V)
305-27-018	2.7mm x 18mm Gridlock Ankle Locking Screw	00812926026046	Yes	Titanium Alloy (Ti6Al4V)
305-27-020	2.7mm x 20mm Gridlock Ankle Locking Screw	00812926026053	Yes	Titanium Alloy (Ti6Al4V)
305-27-022	2.7mm x 22mm Gridlock Ankle Locking Screw	00812926026060	Yes	Titanium Alloy (Ti6Al4V)
305-27-024	2.7mm x 24mm Gridlock Ankle Locking Screw	00812926026077	Yes	Titanium Alloy (Ti6Al4V)
305-35-010	3.5mm x 10mm Gridlock Ankle Locking Screw	00812926026084	Yes	Titanium Alloy (Ti6Al4V)
305-35-012	3.5mm x 12mm Gridlock Ankle Locking Screw	00812926026091	Yes	Titanium Alloy (Ti6Al4V)
305-35-014	3.5mm x 14mm Gridlock Ankle Locking Screw	00812926026107	Yes	Titanium Alloy (Ti6Al4V)
305-35-016	3.5mm x 16mm Gridlock Ankle Locking Screw	00812926026114	Yes	Titanium Alloy (Ti6Al4V)
305-35-018	3.5mm x 18mm Gridlock Ankle Locking Screw	00812926026121	Yes	Titanium Alloy (Ti6Al4V)
305-35-020	3.5mm x 20mm Gridlock Ankle Locking Screw	00812926026138	Yes	Titanium Alloy (Ti6Al4V)
305-35-022	3.5mm x 22mm Gridlock Ankle Locking Screw	00812926026145	Yes	Titanium Alloy (Ti6Al4V)
305-35-024	3.5mm x 24mm Gridlock Ankle Locking Screw	00812926026152	Yes	Titanium Alloy (Ti6Al4V)
305-35-026	3.5mm x 26mm Gridlock Ankle Locking Screw	00812926026169	Yes	Titanium Alloy (Ti6Al4V)
305-35-028	3.5mm x 28mm Gridlock Ankle Locking Screw	00812926026176	Yes	Titanium Alloy (Ti6Al4V)
305-35-030	3.5mm x 30mm Gridlock Ankle Locking Screw	00812926026183	Yes	Titanium Alloy (Ti6Al4V)
305-35-032	3.5mm x 32mm Gridlock Ankle Locking Screw	00812926026190	Yes	Titanium Alloy (Ti6Al4V)
305-35-034	3.5mm x 34mm Gridlock Ankle Locking Screw	00812926026206	Yes	Titanium Alloy (Ti6Al4V)
305-35-036	3.5mm x 36mm Gridlock Ankle Locking Screw	00812926026213	Yes	Titanium Alloy (Ti6Al4V)
305-35-038	3.5mm x 38mm Gridlock Ankle Locking Screw	00812926026220	Yes	Titanium Alloy (Ti6Al4V)
305-35-040	3.5mm x 40mm Gridlock Ankle Locking Screw	00812926026237	Yes	Titanium Alloy (Ti6Al4V)
305-35-042	3.5mm x 42mm Gridlock Ankle Locking Screw	00812926026244	Yes	Titanium Alloy (Ti6Al4V)
305-35-044	3.5mm x 44mm Gridlock Ankle Locking Screw	00812926026251	Yes	Titanium Alloy (Ti6Al4V)
305-35-046	3.5mm x 46mm Gridlock Ankle Locking Screw	00812926026268	Yes	Titanium Alloy (Ti6Al4V)
305-35-048	3.5mm x 48mm Gridlock Ankle Locking Screw	00812926026275	Yes	Titanium Alloy (Ti6Al4V)
305-35-050	3.5mm x 50mm Gridlock Ankle Locking Screw	00812926026282	Yes	Titanium Alloy (Ti6Al4V)
305-40-010	4.0mm x 10mm Gridlock Ankle Locking Screw	00812926026299	Yes	Titanium Alloy (Ti6Al4V)
305-40-012	4.0mm x 12mm Gridlock Ankle Locking Screw	00812926026305	Yes	Titanium Alloy (Ti6Al4V)
305-40-014	4.0mm x 14mm Gridlock Ankle Locking Screw	00812926026312	Yes	Titanium Alloy (Ti6Al4V)
305-40-016	4.0mm x 16mm Gridlock Ankle Locking Screw	00812926026329	Yes	Titanium Alloy (Ti6Al4V)
305-40-018	4.0mm x 18mm Gridlock Ankle Locking Screw	00812926026336	Yes	Titanium Alloy (Ti6Al4V)
305-40-020	4.0mm x 20mm Gridlock Ankle Locking Screw	00812926026343	Yes	Titanium Alloy (Ti6Al4V)
305-40-022	4.0mm x 22mm Gridlock Ankle Locking Screw	00812926026350	Yes	Titanium Alloy (Ti6Al4V)
305-40-024	4.0mm x 24mm Gridlock Ankle Locking Screw	00812926026367	Yes	Titanium Alloy (Ti6Al4V)
305-40-026	4.0mm x 26mm Gridlock Ankle Locking Screw	00812926026374	Yes	Titanium Alloy (Ti6Al4V)
305-40-028	4.0mm x 28mm Gridlock Ankle Locking Screw	00812926026381	Yes	Titanium Alloy (Ti6Al4V)
305-40-030	4.0mm x 30mm Gridlock Ankle Locking Screw	00812926026398	Yes	Titanium Alloy (Ti6Al4V)
305-40-032	4.0mm x 32mm Gridlock Ankle Locking Screw	00812926026404	Yes	Titanium Alloy (Ti6Al4V)
305-40-034	4.0mm x 34mm Gridlock Ankle Locking Screw	00812926026411	Yes	Titanium Alloy (Ti6Al4V)

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
305-40-036	4.0mm x 36mm Gridlock Ankle Locking Screw	00812926026428	Yes	Titanium Alloy (Ti6Al4V)
305-40-038	4.0mm x 38mm Gridlock Ankle Locking Screw	00812926026435	Yes	Titanium Alloy (Ti6Al4V)
305-40-040	4.0mm x 40mm Gridlock Ankle Locking Screw	00812926026442	Yes	Titanium Alloy (Ti6Al4V)
305-40-042	4.0mm x 42mm Gridlock Ankle Locking Screw	00812926026459	Yes	Titanium Alloy (Ti6Al4V)
305-40-044	4.0mm x 44mm Gridlock Ankle Locking Screw	00812926026466	Yes	Titanium Alloy (Ti6Al4V)
305-40-046	4.0mm x 46mm Gridlock Ankle Locking Screw	00812926026473	Yes	Titanium Alloy (Ti6Al4V)
305-40-048	4.0mm x 48mm Gridlock Ankle Locking Screw	00812926026480	Yes	Titanium Alloy (Ti6Al4V)
305-40-050	4.0mm x 50mm Gridlock Ankle Locking Screw	00812926026497	Yes	Titanium Alloy (Ti6Al4V)
306-40-040	4.0mm x 40mm Gridlock Ankle Syndesmotic Screw	00812926026503	Yes	Titanium Alloy (Ti6Al4V)
306-40-042	4.0mm x 42mm Gridlock Ankle Syndesmotic Screw	00812926026510	Yes	Titanium Alloy (Ti6Al4V)
306-40-044	4.0mm x 44mm Gridlock Ankle Syndesmotic Screw	00812926026527	Yes	Titanium Alloy (Ti6Al4V)
306-40-046	4.0mm x 46mm Gridlock Ankle Syndesmotic Screw	00812926026534	Yes	Titanium Alloy (Ti6Al4V)
306-40-048	4.0mm x 48mm Gridlock Ankle Syndesmotic Screw	00812926026541	Yes	Titanium Alloy (Ti6Al4V)
306-40-050	4.0mm x 50mm Gridlock Ankle Syndesmotic Screw	00812926026558	Yes	Titanium Alloy (Ti6Al4V)
306-40-052	4.0mm x 52mm Gridlock Ankle Syndesmotic Screw	00812926026565	Yes	Titanium Alloy (Ti6Al4V)
306-40-054	4.0mm x 54mm Gridlock Ankle Syndesmotic Screw	00812926026572	Yes	Titanium Alloy (Ti6Al4V)
306-40-056	4.0mm x 56mm Gridlock Ankle Syndesmotic Screw	00812926026589	Yes	Titanium Alloy (Ti6Al4V)
306-40-058	4.0mm x 58mm Gridlock Ankle Syndesmotic Screw	00812926026596	Yes	Titanium Alloy (Ti6Al4V)
306-40-060	4.0mm x 60mm Gridlock Ankle Syndesmotic Screw	00812926026602	Yes	Titanium Alloy (Ti6Al4V)
306-40-062	4.0mm x 62mm Gridlock Ankle Syndesmotic Screw	00812926026619	Yes	Titanium Alloy (Ti6Al4V)
306-40-064	4.0mm x 64mm Gridlock Ankle Syndesmotic Screw	00812926026626	Yes	Titanium Alloy (Ti6Al4V)
306-41-040	4.0mm x 40mm Gridlock Ankle Syndesmotic Screw LP	00812926027234	Yes	Titanium Alloy (Ti6Al4V)
306-41-042	4.0mm x 42mm Gridlock Ankle Syndesmotic Screw LP	00812926027241	Yes	Titanium Alloy (Ti6Al4V)
306-41-044	4.0mm x 44mm Gridlock Ankle Syndesmotic Screw LP	00812926027258	Yes	Titanium Alloy (Ti6Al4V)
306-41-046	4.0mm x 46mm Gridlock Ankle Syndesmotic Screw LP	00812926027265	Yes	Titanium Alloy (Ti6Al4V)
306-41-048	4.0mm x 48mm Gridlock Ankle Syndesmotic Screw LP	00812926027272	Yes	Titanium Alloy (Ti6Al4V)
306-41-050	4.0mm x 50mm Gridlock Ankle Syndesmotic Screw LP	00812926027289	Yes	Titanium Alloy (Ti6Al4V)
306-41-052	4.0mm x 52mm Gridlock Ankle Syndesmotic Screw LP	00812926027296	Yes	Titanium Alloy (Ti6Al4V)
306-41-054	4.0mm x 54mm Gridlock Ankle Syndesmotic Screw LP	00812926027302	Yes	Titanium Alloy (Ti6Al4V)
306-41-056	4.0mm x 56mm Gridlock Ankle Syndesmotic Screw LP	00812926027319	Yes	Titanium Alloy (Ti6Al4V)
306-41-058	4.0mm x 58mm Gridlock Ankle Syndesmotic Screw LP	00812926027326	Yes	Titanium Alloy (Ti6Al4V)
306-41-060	4.0mm x 60mm Gridlock Ankle Syndesmotic Screw LP	00812926027333	Yes	Titanium Alloy (Ti6Al4V)
306-41-062	4.0mm x 62mm Gridlock Ankle Syndesmotic Screw LP	00812926027340	Yes	Titanium Alloy (Ti6Al4V)
306-41-064	4.0mm x 64mm Gridlock Ankle Syndesmotic Screw LP	00812926027357	Yes	Titanium Alloy (Ti6Al4V)
307-22-006	2.2mm x 6mm Arsenal Screw	00810028390294	Yes	Titanium Alloy (Ti6Al4V)
307-22-008	2.2mm x 8mm Arsenal Screw	00810028390300	Yes	Titanium Alloy (Ti6Al4V)
307-22-010	2.2mm x 10mm Arsenal Screw	00810028390317	Yes	Titanium Alloy (Ti6Al4V)
307-22-012	2.2mm x 12mm Arsenal Screw	00810028390324	Yes	Titanium Alloy (Ti6Al4V)
307-22-014	2.2mm x 14mm Arsenal Screw	00810028390331	Yes	Titanium Alloy (Ti6Al4V)
307-22-016	2.2mm x 16mm Arsenal Screw	00810028390348	Yes	Titanium Alloy (Ti6Al4V)
307-22-018	2.2mm x 18mm Arsenal Screw	00810028390355	Yes	Titanium Alloy (Ti6Al4V)
307-22-020	2.2mm x 20mm Arsenal Screw	00810028390362	Yes	Titanium Alloy (Ti6Al4V)
307-22-022	2.2mm x 22mm Arsenal Screw	00810028390379	Yes	Titanium Alloy (Ti6Al4V)
307-22-024	2.2mm x 24mm Arsenal Screw	00810028390386	Yes	Titanium Alloy (Ti6Al4V)
307-22-026	2.2mm x 26mm Arsenal Screw	00810028390393	Yes	Titanium Alloy (Ti6Al4V)
307-22-028	2.2mm x 28mm Arsenal Screw	00810028390409	Yes	Titanium Alloy (Ti6Al4V)

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
307-22-030	2.2mm x 30mm Arsenal Screw	00810028390416	Yes	Titanium Alloy (Ti6Al4V)
307-27-008	2.7mm x 8mm Arsenal Screw	00810028390423	Yes	Titanium Alloy (Ti6Al4V)
307-27-010	2.7mm x 10mm Arsenal Screw	00810028390430	Yes	Titanium Alloy (Ti6Al4V)
307-27-012	2.7mm x 12mm Arsenal Screw	00810028390447	Yes	Titanium Alloy (Ti6Al4V)
307-27-014	2.7mm x 14mm Arsenal Screw	00810028390454	Yes	Titanium Alloy (Ti6Al4V)
307-27-016	2.7mm x 16mm Arsenal Screw	00810028390461	Yes	Titanium Alloy (Ti6Al4V)
307-27-018	2.7mm x 18mm Arsenal Screw	00810028390478	Yes	Titanium Alloy (Ti6Al4V)
307-27-020	2.7mm x 20mm Arsenal Screw	00810028390485	Yes	Titanium Alloy (Ti6Al4V)
307-27-022	2.7mm x 22mm Arsenal Screw	00810028390492	Yes	Titanium Alloy (Ti6Al4V)
307-27-024	2.7mm x 24mm Arsenal Screw	00810028390508	Yes	Titanium Alloy (Ti6Al4V)
307-27-026	2.7mm x 26mm Arsenal Screw	00810028390515	Yes	Titanium Alloy (Ti6Al4V)
307-27-028	2.7mm x 28mm Arsenal Screw	00810028390522	Yes	Titanium Alloy (Ti6Al4V)
307-27-030	2.7mm x 30mm Arsenal Screw	00810028390539	Yes	Titanium Alloy (Ti6Al4V)
307-27-032	2.7mm x 32mm Arsenal Screw	00810028390546	Yes	Titanium Alloy (Ti6Al4V)
307-27-034	2.7mm x 34mm Arsenal Screw	00810028390553	Yes	Titanium Alloy (Ti6Al4V)
307-27-036	2.7mm x 36mm Arsenal Screw	00810028390560	Yes	Titanium Alloy (Ti6Al4V)
307-27-038	2.7mm x 38mm Arsenal Screw	00810028390577	Yes	Titanium Alloy (Ti6Al4V)
307-27-040	2.7mm x 40mm Arsenal Screw	00810028390584	Yes	Titanium Alloy (Ti6Al4V)
307-27-042	2.7mm x 42mm Arsenal Screw	00810028391925	Yes	Titanium Alloy (Ti6Al4V)
307-27-044	2.7mm x 44mm Arsenal Screw	00810028391932	Yes	Titanium Alloy (Ti6Al4V)
307-27-046	2.7mm x 46mm Arsenal Screw	00810028391949	Yes	Titanium Alloy (Ti6Al4V)
307-27-048	2.7mm x 48mm Arsenal Screw	00810028391956	Yes	Titanium Alloy (Ti6Al4V)
307-27-050	2.7mm x 50mm Arsenal Screw	00810028391963	Yes	Titanium Alloy (Ti6Al4V)
307-35-008	3.5mm x 8mm Arsenal Screw	00810028390591	Yes	Titanium Alloy (Ti6Al4V)
307-35-010	3.5mm x 10mm Arsenal Screw	00810028390607	Yes	Titanium Alloy (Ti6Al4V)
307-35-012	3.5mm x 12mm Arsenal Screw	00810028390614	Yes	Titanium Alloy (Ti6Al4V)
307-35-014	3.5mm x 14mm Arsenal Screw	00810028390621	Yes	Titanium Alloy (Ti6Al4V)
307-35-016	3.5mm x 16mm Arsenal Screw	00810028390638	Yes	Titanium Alloy (Ti6Al4V)
307-35-018	3.5mm x 18mm Arsenal Screw	00810028390645	Yes	Titanium Alloy (Ti6Al4V)
307-35-020	3.5mm x 20mm Arsenal Screw	00810028390652	Yes	Titanium Alloy (Ti6Al4V)
307-35-022	3.5mm x 22mm Arsenal Screw	00810028390669	Yes	Titanium Alloy (Ti6Al4V)
307-35-024	3.5mm x 24mm Arsenal Screw	00810028390676	Yes	Titanium Alloy (Ti6Al4V)
307-35-026	3.5mm x 26mm Arsenal Screw	00810028390683	Yes	Titanium Alloy (Ti6Al4V)
307-35-028	3.5mm x 28mm Arsenal Screw	00810028390690	Yes	Titanium Alloy (Ti6Al4V)
307-35-030	3.5mm x 30mm Arsenal Screw	00810028390706	Yes	Titanium Alloy (Ti6Al4V)
307-35-032	3.5mm x 32mm Arsenal Screw	00810028390713	Yes	Titanium Alloy (Ti6Al4V)
307-35-034	3.5mm x 34mm Arsenal Screw	00810028390720	Yes	Titanium Alloy (Ti6Al4V)
307-35-036	3.5mm x 36mm Arsenal Screw	00810028390737	Yes	Titanium Alloy (Ti6Al4V)
307-35-038	3.5mm x 38mm Arsenal Screw	00810028390744	Yes	Titanium Alloy (Ti6Al4V)
307-35-040	3.5mm x 40mm Arsenal Screw	00810028390751	Yes	Titanium Alloy (Ti6Al4V)
307-35-042	3.5mm x 42mm Arsenal Screw	00810028390768	Yes	Titanium Alloy (Ti6Al4V)
307-35-044	3.5mm x 44mm Arsenal Screw	00810028390775	Yes	Titanium Alloy (Ti6Al4V)
307-35-046	3.5mm x 46mm Arsenal Screw	00810028390782	Yes	Titanium Alloy (Ti6Al4V)
307-35-048	3.5mm x 48mm Arsenal Screw	00810028390799	Yes	Titanium Alloy (Ti6Al4V)
307-35-050	3.5mm x 50mm Arsenal Screw	00810028390805	Yes	Titanium Alloy (Ti6Al4V)
307-35-055	3.5mm x 55mm Arsenal Screw	00810028391970	Yes	Titanium Alloy (Ti6Al4V)

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
307-35-060	3.5mm x 60mm Arsenal Screw	00810028391987	Yes	Titanium Alloy (Ti6Al4V)
307-35-065	3.5mm x 65mm Arsenal Screw	00810028391994	Yes	Titanium Alloy (Ti6Al4V)
307-35-070	3.5mm x 70mm Arsenal Screw	00810028392007	Yes	Titanium Alloy (Ti6Al4V)
308-22-006	2.2mm x 6mm Arsenal Locking Screw	00810028390812	Yes	Titanium Alloy (Ti6Al4V)
307-40-010	4.0mm x 10mm Arsenal Screw	00810028392137	Yes	Titanium Alloy (Ti6Al4V)
307-40-012	4.0mm x 12mm Arsenal Screw	00810028392144	Yes	Titanium Alloy (Ti6Al4V)
307-40-014	4.0mm x 14mm Arsenal Screw	00810028392151	Yes	Titanium Alloy (Ti6Al4V)
307-40-016	4.0mm x 16mm Arsenal Screw	00810028392168	Yes	Titanium Alloy (Ti6Al4V)
307-40-018	4.0mm x 18mm Arsenal Screw	00810028392175	Yes	Titanium Alloy (Ti6Al4V)
307-40-020	4.0mm x 20mm Arsenal Screw	00810028392182	Yes	Titanium Alloy (Ti6Al4V)
307-40-022	4.0mm x 22mm Arsenal Screw	00810028392199	Yes	Titanium Alloy (Ti6Al4V)
307-40-024	4.0mm x 24mm Arsenal Screw	00810028392205	Yes	Titanium Alloy (Ti6Al4V)
307-40-026	4.0mm x 26mm Arsenal Screw	00810028392212	Yes	Titanium Alloy (Ti6Al4V)
307-40-028	4.0mm x 28mm Arsenal Screw	00810028392229	Yes	Titanium Alloy (Ti6Al4V)
307-40-030	4.0mm x 30mm Arsenal Screw	00810028392236	Yes	Titanium Alloy (Ti6Al4V)
307-40-032	4.0mm x 32mm Arsenal Screw	00810028392243	Yes	Titanium Alloy (Ti6Al4V)
307-40-034	4.0mm x 34mm Arsenal Screw	00810028392250	Yes	Titanium Alloy (Ti6Al4V)
307-40-036	4.0mm x 36mm Arsenal Screw	00810028392267	Yes	Titanium Alloy (Ti6Al4V)
307-40-038	4.0mm x 38mm Arsenal Screw	00810028392274	Yes	Titanium Alloy (Ti6Al4V)
307-40-040	4.0mm x 40mm Arsenal Screw	00810028392281	Yes	Titanium Alloy (Ti6Al4V)
307-40-042	4.0mm x 42mm Arsenal Screw	00810028392298	Yes	Titanium Alloy (Ti6Al4V)
307-40-044	4.0mm x 44mm Arsenal Screw	00810028392304	Yes	Titanium Alloy (Ti6Al4V)
307-40-046	4.0mm x 46mm Arsenal Screw	00810028392311	Yes	Titanium Alloy (Ti6Al4V)
307-40-048	4.0mm x 48mm Arsenal Screw	00810028392328	Yes	Titanium Alloy (Ti6Al4V)
307-40-050	4.0mm x 50mm Arsenal Screw	00810028392335	Yes	Titanium Alloy (Ti6Al4V)
307-40-055	4.0mm x 55mm Arsenal Screw	00810028392342	Yes	Titanium Alloy (Ti6Al4V)
307-40-060	4.0mm x 60mm Arsenal Screw	00810028392359	Yes	Titanium Alloy (Ti6Al4V)
307-40-065	4.0mm x 65mm Arsenal Screw	00810028392366	Yes	Titanium Alloy (Ti6Al4V)
307-40-070	4.0mm x 70mm Arsenal Screw	00810028392373	Yes	Titanium Alloy (Ti6Al4V)
308-22-008	2.2mm x 8mm Arsenal Locking Screw	00810028390829	Yes	Titanium Alloy (Ti6Al4V)
308-22-010	2.2mm x 10mm Arsenal Locking Screw	00810028390836	Yes	Titanium Alloy (Ti6Al4V)
308-22-012	2.2mm x 12mm Arsenal Locking Screw	00810028390843	Yes	Titanium Alloy (Ti6Al4V)
308-22-014	2.2mm x 14mm Arsenal Locking Screw	00810028390850	Yes	Titanium Alloy (Ti6Al4V)
308-22-016	2.2mm x 16mm Arsenal Locking Screw	00810028390867	Yes	Titanium Alloy (Ti6Al4V)
308-22-018	2.2mm x 18mm Arsenal Locking Screw	00810028390874	Yes	Titanium Alloy (Ti6Al4V)
308-22-020	2.2mm x 20mm Arsenal Locking Screw	00810028390881	Yes	Titanium Alloy (Ti6Al4V)
308-22-022	2.2mm x 22mm Arsenal Locking Screw	00810028390898	Yes	Titanium Alloy (Ti6Al4V)
308-22-024	2.2mm x 24mm Arsenal Locking Screw	00810028390904	Yes	Titanium Alloy (Ti6Al4V)
308-22-026	2.2mm x 26mm Arsenal Locking Screw	00810028390911	Yes	Titanium Alloy (Ti6Al4V)
308-22-028	2.2mm x 28mm Arsenal Locking Screw	00810028390928	Yes	Titanium Alloy (Ti6Al4V)
308-22-030	2.2mm x 30mm Arsenal Locking Screw	00810028390935	Yes	Titanium Alloy (Ti6Al4V)
308-27-008	2.7mm x 8mm Arsenal Locking Screw	00810028390942	Yes	Titanium Alloy (Ti6Al4V)
308-27-010	2.7mm x 10mm Arsenal Locking Screw	00810028390959	Yes	Titanium Alloy (Ti6Al4V)
308-27-012	2.7mm x 12mm Arsenal Locking Screw	00810028390966	Yes	Titanium Alloy (Ti6Al4V)
308-27-014	2.7mm x 14mm Arsenal Locking Screw	00810028390973	Yes	Titanium Alloy (Ti6Al4V)
308-27-016	2.7mm x 16mm Arsenal Locking Screw	00810028390980	Yes	Titanium Alloy (Ti6Al4V)

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
308-27-018	2.7mm x 18mm Arsenal Locking Screw	00810028390997	Yes	Titanium Alloy (Ti6Al4V)
308-27-020	2.7mm x 20mm Arsenal Locking Screw	00810028391000	Yes	Titanium Alloy (Ti6Al4V)
308-27-022	2.7mm x 22mm Arsenal Locking Screw	00810028391017	Yes	Titanium Alloy (Ti6Al4V)
308-27-024	2.7mm x 24mm Arsenal Locking Screw	00810028391024	Yes	Titanium Alloy (Ti6Al4V)
308-27-026	2.7mm x 26mm Arsenal Locking Screw	00810028391031	Yes	Titanium Alloy (Ti6Al4V)
308-27-028	2.7mm x 28mm Arsenal Locking Screw	00810028391048	Yes	Titanium Alloy (Ti6Al4V)
308-27-030	2.7mm x 30mm Arsenal Locking Screw	00810028391055	Yes	Titanium Alloy (Ti6Al4V)
308-27-032	2.7mm x 32mm Arsenal Locking Screw	00810028391062	Yes	Titanium Alloy (Ti6Al4V)
308-27-034	2.7mm x 34mm Arsenal Locking Screw	00810028391079	Yes	Titanium Alloy (Ti6Al4V)
308-27-036	2.7mm x 36mm Arsenal Locking Screw	00810028391086	Yes	Titanium Alloy (Ti6Al4V)
308-27-038	2.7mm x 38mm Arsenal Locking Screw	00810028391093	Yes	Titanium Alloy (Ti6Al4V)
308-27-040	2.7mm x 40mm Arsenal Locking Screw	00810028391109	Yes	Titanium Alloy (Ti6Al4V)
308-27-042	2.7mm x 42mm Arsenal Locking Screw	00810028392014	Yes	Titanium Alloy (Ti6Al4V)
308-27-044	2.7mm x 44mm Arsenal Locking Screw	00810028392021	Yes	Titanium Alloy (Ti6Al4V)
308-27-046	2.7mm x 46mm Arsenal Locking Screw	00810028392038	Yes	Titanium Alloy (Ti6Al4V)
308-27-048	2.7mm x 48mm Arsenal Locking Screw	00810028392045	Yes	Titanium Alloy (Ti6Al4V)
308-27-050	2.7mm x 50mm Arsenal Locking Screw	00810028392052	Yes	Titanium Alloy (Ti6Al4V)
308-35-008	3.5mm x 8mm Arsenal Locking Screw	00810028391116	Yes	Titanium Alloy (Ti6Al4V)
308-35-010	3.5mm x 10mm Arsenal Locking Screw	00810028391123	Yes	Titanium Alloy (Ti6Al4V)
308-35-012	3.5mm x 12mm Arsenal Locking Screw	00810028391130	Yes	Titanium Alloy (Ti6Al4V)
308-35-014	3.5mm x 14mm Arsenal Locking Screw	00810028391147	Yes	Titanium Alloy (Ti6Al4V)
308-35-016	3.5mm x 16mm Arsenal Locking Screw	00810028391154	Yes	Titanium Alloy (Ti6Al4V)
308-35-018	3.5mm x 18mm Arsenal Locking Screw	00810028391161	Yes	Titanium Alloy (Ti6Al4V)
308-35-020	3.5mm x 20mm Arsenal Locking Screw	00810028391178	Yes	Titanium Alloy (Ti6Al4V)
308-35-022	3.5mm x 22mm Arsenal Locking Screw	00810028391185	Yes	Titanium Alloy (Ti6Al4V)
308-35-024	3.5mm x 24mm Arsenal Locking Screw	00810028391192	Yes	Titanium Alloy (Ti6Al4V)
308-35-026	3.5mm x 26mm Arsenal Locking Screw	00810028391208	Yes	Titanium Alloy (Ti6Al4V)
308-35-028	3.5mm x 28mm Arsenal Locking Screw	00810028391215	Yes	Titanium Alloy (Ti6Al4V)
308-35-030	3.5mm x 30mm Arsenal Locking Screw	00810028391222	Yes	Titanium Alloy (Ti6Al4V)
308-35-032	3.5mm x 32mm Arsenal Locking Screw	00810028391239	Yes	Titanium Alloy (Ti6Al4V)
308-35-034	3.5mm x 34mm Arsenal Locking Screw	00810028391246	Yes	Titanium Alloy (Ti6Al4V)
308-35-036	3.5mm x 36mm Arsenal Locking Screw	00810028391253	Yes	Titanium Alloy (Ti6Al4V)
308-35-038	3.5mm x 38mm Arsenal Locking Screw	00810028391260	Yes	Titanium Alloy (Ti6Al4V)
308-35-040	3.5mm x 40mm Arsenal Locking Screw	00810028391277	Yes	Titanium Alloy (Ti6Al4V)
308-35-042	3.5mm x 42mm Arsenal Locking Screw	00810028391284	Yes	Titanium Alloy (Ti6Al4V)
308-35-044	3.5mm x 44mm Arsenal Locking Screw	00810028391291	Yes	Titanium Alloy (Ti6Al4V)
308-35-046	3.5mm x 46mm Arsenal Locking Screw	00810028391307	Yes	Titanium Alloy (Ti6Al4V)
308-35-048	3.5mm x 48mm Arsenal Locking Screw	00810028391314	Yes	Titanium Alloy (Ti6Al4V)
308-35-050	3.5mm x 50mm Arsenal Locking Screw	00810028391321	Yes	Titanium Alloy (Ti6Al4V)
308-35-055	3.5mm x 55mm Arsenal Locking Screw	00810028392069	Yes	Titanium Alloy (Ti6Al4V)
308-35-060	3.5mm x 60mm Arsenal Locking Screw	00810028392076	Yes	Titanium Alloy (Ti6Al4V)
308-35-065	3.5mm x 65mm Arsenal Locking Screw	00810028392083	Yes	Titanium Alloy (Ti6Al4V)
308-35-070	3.5mm x 70mm Arsenal Locking Screw	00810028392090	Yes	Titanium Alloy (Ti6Al4V)
308-40-010	4.0mm x 10mm Arsenal Locking Screw	00810028392380	Yes	Titanium Alloy (Ti6Al4V)
308-40-012	4.0mm x 12mm Arsenal Locking Screw	00810028392397	Yes	Titanium Alloy (Ti6Al4V)
308-40-014	4.0mm x 14mm Arsenal Locking Screw	00810028392403	Yes	Titanium Alloy (Ti6Al4V)

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
308-40-016	4.0mm x 16mm Arsenal Locking Screw	00810028392410	Yes	Titanium Alloy (Ti6Al4V)
308-40-018	4.0mm x 18mm Arsenal Locking Screw	00810028392427	Yes	Titanium Alloy (Ti6Al4V)
308-40-020	4.0mm x 20mm Arsenal Locking Screw	00810028392434	Yes	Titanium Alloy (Ti6Al4V)
308-40-022	4.0mm x 22mm Arsenal Locking Screw	00810028392441	Yes	Titanium Alloy (Ti6Al4V)
308-40-024	4.0mm x 24mm Arsenal Locking Screw	00810028392458	Yes	Titanium Alloy (Ti6Al4V)
308-40-026	4.0mm x 26mm Arsenal Locking Screw	00810028392465	Yes	Titanium Alloy (Ti6Al4V)
308-40-028	4.0mm x 28mm Arsenal Locking Screw	00810028392472	Yes	Titanium Alloy (Ti6Al4V)
308-40-030	4.0mm x 30mm Arsenal Locking Screw	00810028392489	Yes	Titanium Alloy (Ti6Al4V)
308-40-032	4.0mm x 32mm Arsenal Locking Screw	00810028392496	Yes	Titanium Alloy (Ti6Al4V)
308-40-034	4.0mm x 34mm Arsenal Locking Screw	00810028392502	Yes	Titanium Alloy (Ti6Al4V)
308-40-036	4.0mm x 36mm Arsenal Locking Screw	00810028392519	Yes	Titanium Alloy (Ti6Al4V)
308-40-038	4.0mm x 38mm Arsenal Locking Screw	00810028392526	Yes	Titanium Alloy (Ti6Al4V)
308-40-040	4.0mm x 40mm Arsenal Locking Screw	00810028392533	Yes	Titanium Alloy (Ti6Al4V)
308-40-042	4.0mm x 42mm Arsenal Locking Screw	00810028392540	Yes	Titanium Alloy (Ti6Al4V)
308-40-044	4.0mm x 44mm Arsenal Locking Screw	00810028392557	Yes	Titanium Alloy (Ti6Al4V)
308-40-046	4.0mm x 46mm Arsenal Locking Screw	00810028392564	Yes	Titanium Alloy (Ti6Al4V)
308-40-048	4.0mm x 48mm Arsenal Locking Screw	00810028392571	Yes	Titanium Alloy (Ti6Al4V)
308-40-050	4.0mm x 50mm Arsenal Locking Screw	00810028392588	Yes	Titanium Alloy (Ti6Al4V)
308-40-055	4.0mm x 55mm Arsenal Locking Screw	00810028392595	Yes	Titanium Alloy (Ti6Al4V)
308-40-060	4.0mm x 60mm Arsenal Locking Screw	00810028392601	Yes	Titanium Alloy (Ti6Al4V)
308-40-065	4.0mm x 65mm Arsenal Locking Screw	00810028392618	Yes	Titanium Alloy (Ti6Al4V)
308-40-070	4.0mm x 70mm Arsenal Locking Screw	00810028392625	Yes	Titanium Alloy (Ti6Al4V)
309-45-065	65 mm Malleolus Implant	00810028392632	Yes	Titanium Alloy (Ti6Al4V)
309-45-070	70 mm Malleolus Implant	00810028392649	Yes	Titanium Alloy (Ti6Al4V)
309-45-075	75 mm Malleolus Implant	00810028392656	Yes	Titanium Alloy (Ti6Al4V)
309-45-080	80 mm Malleolus Implant	00810028392663	Yes	Titanium Alloy (Ti6Al4V)
309-45-085	85 mm Malleolus Implant	00810028392670	Yes	Titanium Alloy (Ti6Al4V)
309-45-090	90 mm Malleolus Implant	00810028392687	Yes	Titanium Alloy (Ti6Al4V)
309-45-095	95 mm Malleolus Implant	00810028392694	Yes	Titanium Alloy (Ti6Al4V)
309-45-100	100 mm Malleolus Implant	00810028392700	Yes	Titanium Alloy (Ti6Al4V)
309-45-105	105 mm Malleolus Implant	00810028392717	Yes	Titanium Alloy (Ti6Al4V)
309-45-110	110 mm Malleolus Implant	00810028392724	Yes	Titanium Alloy (Ti6Al4V)
309-45-115	115 mm Malleolus Implant	00810028392731	Yes	Titanium Alloy (Ti6Al4V)
309-45-120	120 mm Malleolus Implant	00810028392748	Yes	Titanium Alloy (Ti6Al4V)
310-00-001	Large Bone Clamp	00812926023632	No	Stainless Steel
310-00-002	Small Bone Clamp	00812926023649	No	Stainless Steel
310-00-003	5" Joint Distraction Plier	00812926023656	No	Stainless Steel
310-00-004	4" Osteotomy Distraction Plier	00812926023663	No	Stainless Steel
310-00-006	Plate Bending Irons	00812926023670	No	Stainless Steel
310-00-008	Large Bending Pliers	00812926023687	No	Stainless Steel
310-00-009	GL Calibrated Screw Sleeve	00812926028040	No	Stainless Steel
310-00-010	Plate Screw Depth Gauge	00812926023694	No	Stainless Steel, Titanium Alloy (Ti6Al4V), Polyphenylsufone
310-00-019	GL Plate Forceps	00812926024066	No	Stainless Steel
310-00-020	6" Distraction Pliers	00812926028019	No	Stainless Steel
310-00-021	8" Distraction Pliers	00812926028026	No	Stainless Steel



## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
310-01-001	Universal IM Guide Pin	00812926023755	Yes	Stainless Steel
310-01-002	Cannulated Reamer Guide	00812926023762	No	Stainless Steel
310-01-016	16mm Concave Reamer Head	00812926023779	No	Stainless Steel
310-01-018	18mm Concave Reamer Head	00812926023793	No	Stainless Steel
310-01-020	20mm Concave Reamer Head	00812926023816	No	Stainless Steel
310-01-022	22mm Concave Reamer Head	00812926023830	No	Stainless Steel
310-01-100	Gridlock Sterilization Tray Assembly	00812926023700	No	Aluminum, Stainless Steel, Silicone
310-01-101	Gridlock Plating System Screw Module Base Assembly	00812926023717	No	Polypropylene, Aluminum
310-01-102	Gridlock VLC MPX L Plating Module Base Assembly	00812926023724	No	Polypropylene, Polyphenylsulfone
310-01-103	Gridlock Instrument Module Base Assembly	00812926023731	No	Polypropylene, Aluminum
310-01-104	Gridlock Specialty Plate Module Base Assembly	00812926023748	No	Polypropylene, Polyphenylsulfone
310-02-016	16mm Convex Reamer Head	00812926023786	No	Stainless Steel
310-02-018	18mm Convex Reamer Head	00812926023809	No	Stainless Steel
310-02-020	20mm Convex Reamer Head	00812926023823	No	Stainless Steel
310-02-022	22mm Convex Reamer Head	00812926023847	No	Stainless Steel
310-24-001	2.4mm GL Screw Drill Bit	00812926023854	No	Stainless Steel
310-24-005	2.4mm Variable/Static Angle Drill Guide	00812926023984	No	Stainless Steel
310-24-006	2.4mm Drill Sleeves	00812926024011	No	Stainless Steel
310-24-007	2.4mm Compression/Neutral Drill Guide	00812926024042	No	Stainless Steel
310-24-008	2.4mm GL Static Threaded Drill Guide	00812926025414	No	Stainless Steel
310-24-009	2.4mm Variable Angle Threaded Drill Guide	00812926028293	No	Stainless Steel
310-25-001	2.4mm GL Screw Overdrill Bit	00812926023885	No	Stainless Steel
310-30-001	3.0mm GL Screw Drill Bit	00812926023861	No	Stainless Steel
310-30-002	2.4/3.0mm GL Screw Countersink	00812926023915	No	Stainless Steel
310-30-003	2.4/3.0/4.0mm GL Screw Driver Bit	00812926023977	No	Stainless Steel
310-30-005	3.0mm Variable/Static Angle Drill Guide	00812926023991	No	Stainless Steel
310-30-006	3.0mm Drill Sleeves	00812926024028	No	Stainless Steel
310-30-007	3.0mm Compression/Neutral Drill Guide	00812926024059	No	Stainless Steel
310-30-008	3.0mm GL Static Threaded Drill Guide	00812926025421	No	Stainless Steel
310-30-009	3.0mm Variable Angle Threaded Drill Guide	00812926028309	No	Stainless Steel
310-31-001	3.0mm GL Screw Overdrill Bit	00812926023892	No	Stainless Steel
310-31-003	Gridlock Driver Bit, Short	00812926027548	No	Stainless Steel
310-40-001	4.0mm GL Pilot Drill Bit	00812926023878	No	Stainless Steel
310-40-002	2.4/3.0/4.0mm GL Screw Countersink	00812926023922	No	Stainless Steel
310-40-005	4.0mm Variable/Static Angle Drill Guide	00812926024004	No	Stainless Steel
310-40-006	4.0mm Drill Sleeves	00812926024035	No	Stainless Steel
310-40-009	4.0mm Variable Angle Threaded Drill Guide	00812926028316	No	Stainless Steel
310-41-001	4.0mm GL Overdrill Bit	00812926023908	No	Stainless Steel
310-60-005	0.062" x 4.5" Olive Wire Partially Threaded	00812926023939	Yes	Stainless Steel
310-99-033	Gridlock Plating System Screw Module Lid	00812926027920	No	Aluminum
310-99-034	Gridlock Plating System Instrument Module Lid	00812926027937	No	Aluminum
310-99-105	Gridlock VLC MPX L Plating Module Lid	00812926027944	No	Polyphenylsulfone
310-99-107	Gridlock Specialty Plate Module Lid	00812926027951	No	Polyphenylsulfone
311-00-001	Small Joint Reamer Sterilization Tray Base Assembly	00812926027463	No	Polypropylene

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
311-01-016	16mm Concave Reamer	00812926027364	No	Stainless Steel
311-01-018	18mm Concave Reamer	00812926027388	No	Stainless Steel
311-01-020	20mm Concave Reamer	00812926027401	No	Stainless Steel
311-01-022	22mm Concave Reamer	00812926027425	No	Stainless Steel
311-01-024	24mm Concave Reamer	00812926027449	No	Stainless Steel
311-02-016	16mm Convex Reamer	00812926027371	No	Stainless Steel
311-02-018	18mm Convex Reamer	00812926027395	No	Stainless Steel
311-02-020	20mm Convex Reamer	00812926027418	No	Stainless Steel
311-02-022	22mm Convex Reamer	00812926027432	No	Stainless Steel
311-02-024	24mm Convex Reamer	00812926027456	No	Stainless Steel
311-99-011	Small Joint Reaming System Lid	00812926027968	No	Polyphenylsulfone
317-00-001	Minimally Invasive Bunion Plate System Sterilization Tray	00812926027722	No	Polypropylene
317-00-002	Plate Placement Guide, Right	00812926027746	No	Stainless Steel
317-00-003	Plate Placement Guide, Left	00812926027753	No	Stainless Steel
317-00-004	MIB Threaded Drill Guide	00812926027760	No	Stainless Steel
317-00-005	MIB Interfrag Drill Guide	00812926027777	No	Stainless Steel
317-00-006	Broach	00812926027784	No	Stainless Steel
317-00-007	Tissue Retractor	00812926027791	No	Stainless Steel
317-24-001	2.4mm Calibrated Drill Bit	00812926027807	No	Stainless Steel
317-99-011	Minimally Invasive Bunion Plate System Sterilization Tray Lid	00812926027739	No	Polyphenylsulfone
320-00-001	Large Serrated Bone Clamp	00812926026695	No	Stainless Steel
320-00-003	Syndesmotic Bone Clamp	00812926026701	No	Stainless Steel
320-00-005	Fixed Driver Handle Gridlock Ankle Plating System	00812926027005	No	Silicone, Stainless Steel, Aluminum Alloy
320-00-006	Gridlock Ankle Plate Bending Irons	00812926026718	No	Stainless Steel
320-00-009	GL Calibrated Screw Sleeve, Ankle	00812926028057	No	Stainless Steel
320-00-010	Gridlock Ankle Plating Depth Gauge	00812926026725	No	Stainless Steel, Titanium Alloy (Ti6Al4V), Polyphenylsulfone
320-00-011	8mm Hohmann Retractor	00812926026732	No	Stainless Steel
320-00-012	15mm Hohmann Retractor	00812926026749	No	Stainless Steel
320-00-013	Bone Pick	00812926026756	No	Stainless Steel
320-00-014	Periosteal Elevator	00812926026763	No	Aluminum, Stainless Steel
320-00-016	GL Ankle Screw Remover	00812926026770	No	Stainless Steel
320-00-017	Hook Drill Guide	00812926027654	No	Stainless Steel
320-00-019	Guide Handle	00812926028798	No	Stainless Steel
320-00-020	Hook Tamp	00812926028804	No	Stainless Steel
320-01-100	Gridlock Ankle Plate System Tray Assembly	00812926026657	No	Aluminum, Stainless Steel, Silicone
320-01-101	Ankle Screw Module Total Base Assembly	00812926026664	No	Polypropylene
320-01-102	Gridlock Ankle Plate Module Base Assembly	00812926026671	No	Polyphenylsulfone
320-01-103	Ankle Instrument Module Base Assembly	00812926026688	No	Polypropylene
320-01-104	Ankle Plate Module, Ext.	00812926028286	No	Polyphenylsulfone
320-01-999	Ankle Plate Module Insert	00812926027975	No	Polypropylene
320-27-001	2.7mm Gridlock Ankle Screw Drill Bit	00812926026794	No	Stainless Steel
320-27-005	2.7mm Ankle Variable/Static Angle Drill Guide	00812926026909	No	Stainless Steel
320-27-006	2.7mm Ankle Drill Sleeves	00812926026947	No	Stainless Steel
320-27-008	2.7mm Ankle Threaded Static Drill Guide	00812926026978	No	Stainless Steel

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
320-27-009	2.7mm Variable Angle Threaded Drill Guide	00812926028323	No	Stainless Steel
320-27-011	2.7mm Gridlock Ankle Screw Overdrill Bit	00812926026824	No	Stainless Steel
320-35-001	3.5mm Gridlock Ankle Screw Drill Bit	00812926026800	No	Stainless Steel
320-35-002	Gridlock Ankle Screw Countersink	00812926026862	No	Stainless Steel
320-35-003	Gridlock Ankle Screw Driver Bit	00812926026893	No	Stainless Steel
320-35-005	3.5mm Ankle Variable/Static Angle Drill Guide	00812926026916	No	Stainless Steel
320-35-006	3.5mm Ankle Drill Sleeves	00812926026954	No	Stainless Steel
320-35-008	3.5mm Ankle Threaded Static Drill Guide	00812926026985	No	Stainless Steel
320-35-009	3.5mm Variable Angle Threaded Drill Guide	00812926028330	No	Stainless Steel
320-35-011	3.5mm Gridlock Ankle Screw Overdrill Bit	00812926026831	No	Stainless Steel
320-36-003	Gridlock Ankle Driver Bit, Short	00812926027555	No	Stainless Steel
320-40-001	4.0mm Gridlock Ankle Screw Drill Bit	00812926026817	No	Stainless Steel
320-40-002	Gridlock Ankle Washer Countersink	00812926026879	No	Stainless Steel
320-40-005	4.0mm Ankle Variable/Static Angle Drill Guide	00812926026923	No	Stainless Steel
320-40-006	4.0mm Ankle Drill Sleeves	00812926026961	No	Stainless Steel
320-40-008	4.0mm Ankle Threaded Static Drill Guide	00812926026992	No	Stainless Steel
320-40-009	4.0mm Variable Angle Threaded Drill Guide	00812926028347	No	Stainless Steel
320-40-011	4.0mm Gridlock Ankle Screw Overdrill Bit	00812926026848	No	Stainless Steel
320-40-012	Syndesmotic Screw Drill Bit	00812926026855	No	Stainless Steel
320-40-015	Syndesmotic Screw Drill Guide	00812926026930	No	Stainless Steel
320-40-016	Syndesmotic Screw Removal Bit	00812926026787	No	Stainless Steel
320-60-005	0.062" x 4.65" GL Olive Wire Partially Threaded	00812926026886	Yes	Stainless Steel
320-99-111	Gridlock Ankle Screw Module Lid	00812926027982	No	Aluminum
320-99-112	Gridlock Ankle Plate Module Lid	00812926027999	No	Polyphenylsulfone
320-99-113	Gridlock Ankle Instrument Module Lid	00812926028002	No	Aluminum
330-00-001	Calibrated Spreader	00812926029528	No	Stainless Steel
330-00-002	Small Bone Distractor	00812926029542	No	Stainless Steel
330-00-003	Small Bone Compressor	00812926029559	No	Stainless Steel
330-00-004	McGlamry Elevator	00812926029566	No	Stainless Steel
330-00-005	Arsenal Variable Angle Drill Guide	00812926029740	No	Stainless Steel
330-00-006	Arsenal Compression Drill Guide	00812926029757	No	Stainless Steel
330-00-007	Arsenal Static Angle Drill Guide	00812926029764	No	Stainless Steel
330-00-008	2.7mm/3.5mm Arsenal Threaded Static Drill Guide	00812926029795	No	Stainless Steel
330-00-009	Arsenal Screw Remover	00810028391338	No	Stainless Steel
330-00-010	Large Bending Pliers with Springs	00810028391352	No	Stainless Steel
330-00-020	Arsenal Drill Tower Outer Sleeve	00810028391765	No	Stainless Steel
330-00-021	Arsenal Drill Tower 2.7/4.0mm Insert	00810028391772	No	Stainless Steel
330-00-022	Arsenal Drill Tower 3.5mm Insert	00810028391789	No	Stainless Steel
330-01-100	Arsenal Plating System Tray	00812926029801	No	Aluminum, Stainless Steel, Silicone, Polyphenylsulfone, Polypropylene
330-01-200	Arsenal Modular Tray Base	00812926029924	No	Aluminum, Stainless Steel, Silicone
330-01-201	Arsenal Modular Tray Lid	00810028391420	No	Aluminum, Stainless Steel
330-01-202	Arsenal Modular Tray Sublevel	00810028391437	No	Aluminum, Stainless Steel
330-01-203	Arsenal 2.7mm/3.5mm Screw Module	00810028391444	No	Aluminum, Stainless Steel, Silicone, Polyphenylsulfone (Radel)

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
330-01-204	Arsenal Modular Caddy Base	00810028391451	No	Polypropylene
330-01-205	Arsenal Modular Caddy Lid	00810028391468	No	Polyphenylsulfone (Radel)
330-01-206	Arsenal MPJ Plating Caddy	00810028391475	No	Polypropylene
330-01-207	Arsenal Reamer Caddy	00810028391482	No	Polypropylene
330-01-208	Arsenal Lapidus Plating Caddy	00810028391499	No	Polypropylene
330-01-209	Arsenal Universal Plating Caddy	00810028391505	No	Polypropylene
330-01-210	K-Wire Tube	00812926029917	No	Stainless Steel, Polyacetal
330-07-003	2.2mm Arsenal Screw Driver Bit	00812926029658	No	Stainless Steel
330-10-003	2.7mm/3.5mm/4.0mm Arsenal Screw Driver Bit	00812926029665	No	Stainless Steel
330-22-001	2.2mm Arsenal Screw Drill Bit	00812926029672	No	Stainless Steel
330-22-005	2.2mm Arsenal Screw Depth Gauge	00812926029597	No	Stainless Steel
330-22-006	2.2mm Arsenal Drill Sleeves	00812926029733	No	Stainless Steel
330-22-011	2.2mm Arsenal Screw Overdrill Bit	00812926029702	No	Stainless Steel
330-24-004	0.094" x 6" K-Wire Standard	00812926029610	Yes	Stainless Steel
330-27-001	2.7mm Arsenal Screw Drill Bit	00812926029689	No	Stainless Steel
330-27-002	2.7mm/3.5mm Arsenal Screw Countersink	00812926029641	No	Stainless Steel
330-27-005	2.7mm/3.5mm Arsenal Screw Depth Gauge	00812926029603	No	Stainless Steel
330-27-006	2.7mm Arsenal Drill Sleeves	00812926029771	No	Stainless Steel
330-27-011	2.7mm Arsenal Screw Overdrill Bit	00812926029719	No	Stainless Steel
330-35-001	3.5mm Arsenal Screw Drill Bit	00812926029696	No	Stainless Steel
330-35-006	3.5mm Arsenal Drill Sleeves	00812926029788	No	Stainless Steel
330-35-011	3.5mm Arsenal Screw Overdrill Bit	00812926029726	No	Stainless Steel
330-60-005	0.062" x 3" Olive Wire Partially Threaded	00812926029627	Yes	Stainless Steel
330-86-001	Calc Slide Plate Placement Guide	00812926029634	No	Stainless Steel
340-00-001	Arsenal Ankle Variable/Static Angle Drill Guide	00810028391796	No	Stainless Steel
340-00-002	Arsenal Ankle Hook Plate Guide	00810028391833	No	Stainless Steel
340-00-003	Arsenal Ankle Hook Plate Tamp	00810028391840	No	Stainless Steel
340-00-004	Parallel Distractor	00810028391857	No	Stainless Steel
340-00-005	2.7mm/3.5mm/4.0 Arsenal Screw Depth Gauge	00810028391734	No	Stainless Steel
340-00-006	Bending Irons	00810028391864	No	Stainless Steel
340-00-007	Pointed Bone Clamp	00810028391871	No	Stainless Steel
340-00-008	Syndesmotic Clamp	00810028391888	No	Stainless Steel
340-00-009	Arsenal Ankle Screw Remover	00810028391741	No	Stainless Steel
340-00-010	Periosteal Elevator	00810028391918	No	Stainless Steel
340-00-011	15mm Hohmann Retractor (6" length)	00810028391901	No	Stainless Steel
340-00-019	Arsenal Ankle Forceps	00810028391895	No	Stainless Steel
340-01-100	Arsenal Ankle Tray Base	00810028393349	No	Aluminum, Stainless Steel, Silicone
340-01-101	Arsenal Ankle Tray Lid	00810028393356	No	Aluminum, Stainless Steel
340-01-200	Arsenal Ankle Upper Sublevel	00810028393363	No	Aluminum, Stainless Steel, Silicone
340-01-201	Arsenal Ankle Lower Sublevel	00810028393370	No	Aluminum, Stainless Steel, Silicone, Polyphenylsulfone (Radel)
340-01-202	Arsenal Ankle Screw Caddy	00810028393387	No	Aluminum, Stainless Steel, Silicone, Polyphenylsulfone (Radel)
340-01-203	Arsenal Ankle Straight Plate Caddy	00810028393394	No	Aluminum, Stainless Steel, Silicone

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
340-01-204	Arsenal Ankle Lateral Fibula Plate Caddy	00810028393400	No	Aluminum, Stainless Steel, Silicone
340-01-205	Arsenal Ankle Malleolus and Posterior Plate Caddy	00810028393417	No	Aluminum, Stainless Steel, Silicone
340-01-206	Arsenal Ankle Pilon Plate Caddy	00810028393424	No	Aluminum, Stainless Steel, Silicone
340-15-003	Cannulated Driver Bit	00810028391666	No	Stainless Steel
340-15-007	Solid-Core Removal Tool	00810028391697	No	Stainless Steel
340-27-001	2.7mm Arsenal Ankle Screw Drill Bit	00810028391611	Yes	Stainless Steel
340-27-006	2.7mm Arsenal Ankle Drill Sleeves	00810028391802	No	Stainless Steel
340-27-011	2.7mm Arsenal Ankle Screw Overdrill Bit	00810028391642	Yes	Stainless Steel
340-35-001	3.5mm Arsenal Ankle Screw Drill Bit	00810028391628	Yes	Stainless Steel
340-35-006	3.5mm Arsenal Ankle Drill Sleeves	00810028391819	No	Stainless Steel
340-35-011	3.5mm Arsenal Ankle Screw Overdrill Bit	00810028391659	Yes	Stainless Steel
340-40-001	4.0mm Arsenal Ankle Screw Drill Bit	00810028391635	Yes	Stainless Steel
340-40-006	4.0mm Arsenal Ankle Drill Sleeves	00810028391826	No	Stainless Steel
340-45-001	Malleolus Implant Drill Bit	00810028391673	Yes	Stainless Steel
340-45-002	Malleolus Implant Countersink	00810028391680	No	Stainless Steel
340-60-005	.062" x 3" Olive Wire, Arsenal Ankle	00810028391727	Yes	Stainless Steel
500-08-001	8mm Sterile Staple Kit Assembly	00812926027500	Yes	Nitinol, Stainless Steel
500-08-101	8mm x 8mm x 8mm Staple	00812926028217	Yes	Nitinol
500-10-001	10mm Sterile Staple Kit Assembly	00812926027517	Yes	Nitinol, Stainless Steel
500-10-101	10mm x 10mm x 10mm Staple	00812926028224	Yes	Nitinol
500-12-101	12mm x 12mm x 12mm Staple	00812926028606	Yes	Nitinol
500-15-101	15mm x 15mm x 15mm Staple	00812926028231	Yes	Nitinol
500-15-102	15mm x 12mm x 12mm Staple	00812926028248	Yes	Nitinol
500-18-101	18mm x 18mm x 18mm Staple	00812926028255	Yes	Nitinol
500-18-102	18mm x 15mm x 15mm Staple	00812926028262	Yes	Nitinol
500-20-101	20mm x 20mm x 20mm Staple	00812926028279	Yes	Nitinol
500-20-201	20mm x 20mm x 20mm High Strength Staple	00812926028613	Yes	Nitinol
500-25-101	25mm x 20mm x 20mm Staple	00812926028620	Yes	Nitinol
510-00-001	Sniper Staple System Sterilization Tray Base Assembly	00812926028194	No	Polypropylene
510-00-003	8mm, 10mm Staple Spreader	00812926028149	No	Stainless Steel
510-00-004	15mm, 18mm, 20mm Staple Spreader	00812926028156	No	Stainless Steel
510-00-005	Staple Tamp	00812926028187	No	Stainless Steel
510-00-006	8-12mm Staple Spreader	00812926028637	No	Stainless Steel
510-00-007	15-25mm Staple Spreader	00812926028644	No	Stainless Steel
510-00-010	Lever Assist	00812926027845	No	Stainless Steel
510-08-101	8mm Drill Guide Assembly	00812926028064	No	Stainless Steel
510-10-002	8mm/10mm Staple Pilot Drill	00812926028118	No	Stainless Steel
510-10-101	10mm Drill Guide Assembly	00812926028071	No	Stainless Steel
510-12-001	12mm Drill Guide Assembly	00812926028668	No	Stainless Steel
510-15-001	15mm Drill Guide Assembly	00812926028088	No	Stainless Steel
510-16-001	1.6mm Post / Pilot Drill	00812926028521	No	Stainless Steel
510-16-002	1.6mm Pilot Drill	00812926028538	No	Stainless Steel
510-16-005	1.6mm Locating Pin	00812926028545	No	Stainless Steel
510-18-001	18mm Drill Guide Assembly	00812926028095	No	Stainless Steel

## Trilliant Surgical GTIN Reference Catalog

Product / Catalog Number	Product / Catalog Description	GTIN (HRI)	Single Use Only?	Material(s)
510-20-001	20mm Drill Guide Assembly	00812926028101	No	Stainless Steel
510-20-002	2.0mm Calibrated Pilot Drill	00812926028125	No	Stainless Steel
510-20-005	2.0mm Locating Pin	00812926028163	No	Stainless Steel
510-20-101	20mm, HS, Drill Guide Assembly	00812926028675	No	Stainless Steel
510-25-001	25mm Drill Guide Assembly	00812926028682	No	Stainless Steel
510-25-002	2.5mm Calibrated Pilot Drill	00812926028132	No	Stainless Steel
510-25-005	2.5mm Locating Pin	00812926028170	No	Stainless Steel
510-30-002	3.0mm Calibrated Pilot Drill	00812926028651	No	Stainless Steel
510-30-005	3.0mm Locating Pin	00812926028590	No	Stainless Steel
510-99-101	Sniper Staple System Sterilization Tray Lid	00812926028200	No	Polyphenylsulfone