



When portability  
is the priority



# Patients are leaving the hospital sooner

and DVT risk can last for 30 plus days post discharge. DJO Global introduces VenaPro, a truly portable device to help prevent post-operative DVT for moderate to high risk patients.


VenaPro is designed for optimal patient comfort and compliance.

- *No tubes or cords*
- *One button operation*
- *Weighs less than one lb*
- *Designed based on clinically proven technology*
- *Graduated and asymmetric compression*
- *Air releases to cool patients' legs*
- *USB Downloadable compliance data*
- *Re-chargeable battery (4-5 hour life)*
- *Applies 50mmhg once per minute*
- *Patients don't need to return the device!*



## Indications

- Aid in the prevention of DVT
- Can also be used as an aid in the prophylaxis for DVT by persons expecting to be stationary for long periods of time
- Enhance blood circulation
- Diminish post-operative pain and swelling
- Reduce wound healing time
- Aid in the treatment and healing of: stasis dermatitis, venous stasis ulcers, arterial and diabetic leg ulcers, chronic venous insufficiency and reduction of edema in the lower limbs





## How can I get VenaPro?

**Purchase VenaPro:** Contact your local sales representative or call DJO Directly at 800-336-6569

**DJO Insurance Billing:** DJO Global's hospital to home program (H2HDVT) is a simple program where we handle all the billing, shipping and patient education for a completely hassle free way to protect your patients at home. To set up H2HDVT in your clinic contact your local sales representative or call 800-328-2536 ext 7308

Part #	Description
30M	VenaPro System
3111-024	Compliance Data Communicator Cable
3111-020	VenaPro Sleeve Extender
3111-057	VenaPro Charger
3111-058	VenaPro Clinical Charger (3 prong)



**Distributed by DJO Global**

1430 Decision Street | Vista, CA 92081-8553 | U.S.A.

[DJOGlobal.com](http://DJOGlobal.com)