










| Prescription ONLY | Available Devices  | Primera TENS & NMES 77615  | Intelect TENS 77712 77710 77600  | Intelect NMES 77715 77717  | Intelect Digital IFC 77719  |
|--|---|--|--|--|--|
| Discontinued Devices  | Features | Intensity: 0 - 80 mA Frequency: 0 - 150 Hz Pulse Width: 0 - 200 μ s Programs: 8 TENS and 6 NMES | Intensity: 0 - 80 mA Frequency: 50 - 300 Hz Pulse Width: 60 - 250 μ s Programs: Adjustable | Intensity: 0 - 80 mA Frequency: 2 - 120 Hz Pulse Width: 50 - 300 μ s Programs: Adjustable | Intensity: 0 - 32 mA Frequency: 4000-4125 Hz (Carrier) Pulse Width: 125 μ s Programs: 4 Modes, Adjustable |
| Select (TENS) 199584  | Intensity: 0-60 mA Frequency: 0-100 Hz Pulse Width: 0-400 μ s Programs: 5 pre-set, 4 customizable | ✓ | ✓ | | |
| Active (TENS) 199611  | Intensity: 0-60 mA Frequency: 2-125 Hz Pulse Width: 48-400 μ s Program: 1 Modulated | ✓ | ✓ | | |
| Phoenix (TENS & NMES) 199696  | Intensity: 0-100 mA Frequency: 4-222 Hz Pulse Width: 0-300 μ s Programs: 4 preset - 2 NMES, 1 TENS and 1 Edema | ✓ | ✓ | ✓ | |
| Continuum (TENS & NMES) 199610  | Intensity: 0-100 mA Frequency: 2-150 Hz Pulse Width: 0-400 μ s Programs: 13 pre-programmed, 3 custom | ✓ | ✓ | ✓ | |
| Elite (NMES) 199729  | Intensity: 0-120 mA Frequency: 1-120 Hz Pulse Width: 0-380 μ s Programs: 3 programs (endurance, strength, massage) | ✓ | | ✓ | |
| Axon (IFC) 1-99720  | Intensity: 0-33 mA Frequency: 2-160 Hz (Beat), 4000-4160 Hz (Carrier) Programs: 6 IFC programs | | | | ✓ |
| IF3 Wave (IFC & NMES) 7110S  | Intensity: 0-50 mA Frequency: 1-150 Hz (Beat), 4000-4150 Hz (Carrier) Pulse Width: 0-125 μ s Programs: 3 modalities - Interferential Current, NMES & Pulsed DC | | | ✓ | ✓ |

Prescription Device Cross Reference Chart

| <p>Prescription ONLY</p> | <p>Available Devices</p>  | <p>Primera TENS & NMES 77615</p>  | <p>Intelect TENS 77712 77710 77600</p>  | <p>Intelect NMES 77715 77717</p>  | <p>Intelect Digital IFC 77719</p>  |
|--|--|--|---|--|--|
| <p>Discontinued Devices</p>  | <p>Features</p> | <p>Intensity: 0 - 80 mA Frequency: 0 - 150 Hz Pulse Width: 0 - 200 μs Programs: 8 TENS and 6 NMES</p> | <p>Intensity: 0 - 80 mA Frequency: 50 - 300 Hz Pulse Width: 60 - 250 μs Programs: Adjustable</p> | <p>Intensity: 0 - 80 mA Frequency: 2 - 120 Hz Pulse Width: 50 - 300 μs Programs: Adjustable</p> | <p>Intensity: 0 - 32 mA Frequency: 4000-4125 Hz (Carrier) Pulse Width: 125 μs Programs: 4 Modes, Adjustable</p> |
| <p>Sensus (TENS) 700740</p>  | <p>Intensity: 100 mA Frequency: 60-100 Hz Pulse Width: 230-430 μs Programs: 1 program</p> |  |  | | |